

DICOM Conformance Statement, X500

Revision Data

| Revision | Description | Submit Date |
|----------|---|-------------|
| 00 | DICOM Conformance Statement for 1.0 Release | 12/05/06 |

| | | |
|------------|---|-----------|
| 1.0 | Purpose | 5 |
| 2.0 | Scope | 5 |
| 3.0 | Definitions | 6 |
| 4.0 | Implementation Model | 8 |
| 4.0.1 | Application Data Flow Diagram | 9 |
| 4.0.2 | Verification | 10 |
| 4.0.3 | DICOM Store | 10 |
| 4.0.4 | DICOM Print | 11 |
| 4.0.5 | Patient Registration using Worklist | 12 |
| 4.0.6 | Modality Performed Procedure Step | 15 |
| 4.0.7 | Removable Media Storage | 15 |
| 4.0.8 | Storage Commitment | 15 |
| 4.1 | AE Functional Definition | 16 |
| 4.1.1 | Verification Real-World Activities | 16 |
| 4.1.2 | Store Real-World Activities | 16 |
| 4.1.3 | Storage Commitment - Push Model Real-World Activities | 16 |
| 4.1.4 | Print Real-World Activities | 17 |
| 4.2 | Modality Worklist Real-World Activities | 17 |
| 4.3 | Modality Performed Procedure Step Real-World Activities | 17 |
| 4.4 | Removable Media Storage Real-World Activities | 17 |
| 4.5 | Sequencing of Real-World Activities | 18 |
| 5.0 | AE Specifications | 19 |
| 5.1 | Acuson X500 AE Specification | 19 |
| 5.1.1 | Association Establishment Policies | 19 |
| 5.1.2 | Association Initiation by Real-World Activities | 20 |
| 6.0 | Removable Media Interchange Specifications | 59 |
| 6.1 | Supported Application Profiles | 59 |
| 6.2 | Supported SOP Classes | 59 |

| | | |
|-------------|--|------------|
| 6.2.1 | Supported SOP Classes and Transfer Syntaxes | 59 |
| 6.3 | Information Object Definition and DICOMDIR Keys | 60 |
| 6.3.1 | DICOM File Meta Information | 60 |
| 6.3.2 | Basic Directory Information Object Definitions - File-set Identification Module 61 | |
| 6.3.3 | Basic Directory Information Object Definitions - Directory Identification Module 61 | |
| 6.3.4 | Physical Storage Media and Media Formats | 62 |
| 7.0 | Communication Profiles | 63 |
| 7.1 | TCP/IP Stack Supported | 63 |
| 7.1.1 | Physical Media Supported | 63 |
| 7.1.2 | Chapter Extensions/Specializations/Privatizations | 63 |
| 8.0 | Configuration | 64 |
| 8.1 | General System Configuration | 64 |
| 8.1.1 | Hospital Name | 64 |
| 8.2 | DICOM Network Configuration | 64 |
| 8.2.1 | Local | 64 |
| 8.2.2 | Remote | 65 |
| 8.3 | External Equipment Configuration | 69 |
| 8.4 | Support of Extended Character Sets | 69 |
| 9.0 | Security | 70 |
| 9.1 | Security Profiles | 70 |
| 9.2 | Association Level Security | 70 |
| 9.3 | Application Level Security | 70 |
| 9.4 | Virus Protection | 70 |
| 10.0 | Appendix A: Mapping for Obstetric DICOM SR | 71 |
| 11.0 | Appendix B: Mapping for Cardiac DICOM SR | 108 |

1.0 Purpose

This document describes the conformance to the ACR-NEMA DICOM 3.0 Standard by the Acuson X500 ultrasound system software version 1.0 from Siemens Medical Solutions USA, Inc. Ultrasound Division. It shall establish the conformance specifications for this system only, and does not apply to other products offered by Siemens Medical Solutions USA, Inc., or its affiliates.

The Acuson X500 system is a device that generates ultrasound images that can be sent using DICOM standard protocols and definitions to other DICOM compliant devices that support SOP classes as defined in Table 2 in this document.

2.0 Scope

The DICOM standard provides a well-defined set of structures and protocols that allow inter-operability to a wide variety of medical imaging devices.

When configured with the DICOM option, the Acuson X500 system provides support for essential services related to ultrasound scanning and connectivity to DICOM compliant devices. Acuson X500 system products will not support all features supported by the DICOM standard. This document clearly states the DICOM services and data classes that are supported by the applications included with the Acuson X500. The intent of this document is to allow users and other vendors who also conform to the DICOM standard to exchange information within the specific context of those elements of the DICOM standard that Acuson X500 system supports.

This document is written with respect to the adopted portions of the DICOM standard, Revision 3. The following sections of this document follow the outline specified in the DICOM Standard NEMA publication PS3.2.¹

¹ Second part of the DICOM standard: NEMA Standards Publication PS 3.2-2003, Digital Imaging and Communications in Medicine (DICOM), Part 2: Conformance

3.0 Definitions

The following table provides a list of terms, their acronyms (if applicable), and their descriptions.

Table 1 Terms, Acronyms, and Descriptions.

| Term | Acronym | Description |
|--|-----------|--|
| American College of Radiology - National Electrical Manufacturer's Association | ACR-NEMA | The American College of Radiology (ACR) and the National Electrical Manufacturers Association (NEMA) formed a joint committee to develop a standard for Digital Imaging and Communications in Medicine (DICOM). |
| Application Entity | AE | An application that supports DICOM communication with other DICOM applications. |
| DICOM Conformance Statement | DCS | A formal statement associated with a specific implementation of the DICOM Standard. It specifies the Service Classes, Information Objects, Communications Protocols and Media Storage Application Profiles supported by the implementation. |
| DICOM Message Service Element | DIMSE | Defines an Application Service Element (both the service and protocol) used by peer DICOM Application Entities for the purpose of exchanging medical images and related information. |
| Digital Imaging and Communications in Medicine, Version 3.0 | DICOM 3.0 | A well-defined set of structures and protocols that allow inter-operability to a wide variety of medical imaging devices. |
| Ethernet | - | Network methodology devised in 1976 by Digital Equipment Corporation, Intel and Xerox which is the most common in practice today. Ethernet is the IEEE standard 802.3 |
| Information Object Definition | IOD | A data abstraction of a class of similar Real-World Objects which defines the nature and attributes relevant to the class of Real-World objects represented. |
| Integrating the Healthcare Enterprise | IHE | An initiative sponsored by the Radiological Society of North America (RSNA) to document and demonstrate standards-based methods of sharing information in support of optimal patient care. For additional information see www.rsna.org/ihe . |
| Picture Archiving and Communications Systems | PACS | A DICOM server that accepts medical images from another DICOM system and stores the images for later retrieval. |
| Protocol Data Unit | PDU | The PDUs are message formats exchanged between peer entities within a layer. A PDU shall consist of protocol control information and user data. |

Table 1 Terms, Acronyms, and Descriptions. (Continued)

| Term | Acronym | Description |
|---|----------------|--|
| Real-World Activity | RWA | That which exists in the real world which pertains to specific area of information processing within the area of interest of the DICOM Standard. Such a Real-World Activity may be represented by one or more computer information metaphors called SOP Classes. |
| Request | RQ | A request from one DICOM AE for service from another DICOM AE |
| Response | RSP | A response from one DICOM AE to the request for service from another DICOM AE |
| Service Class Provider | SCP | The role played by a DICOM Application Entity (DIMSE-Service-User) which performs operations and invokes notifications on a specific Association. |
| Service Class User | SCU | The role played by a DICOM Application Entity (DIMSE-Service-User) which invokes notifications and performs operations on a specific Association. |
| Service-Object Pairs | SOP | The union of a specific set of DIMSE Services and one related Information Object Definition which completely defines a precise context for communication. |
| Structured Report | SR | Also called Procedure Report. A DICOM object which contains measurement, calculations, diagnoses, image references and other information concerning a patient exam. |
| Unique identifier | UID | A series of digits and periods (.) used to uniquely identify an object such as an Ultrasound image in DICOM. |
| VA Hospital Information System Technology Architecture DICOM Conformance Requirements | VISTA | DICOM requirements document of the US Department of Veteran's Affairs (VA) Hospital Information System Technology Architecture. For additional information see www.va.gov/imaging . |

4.0 Implementation Model

Acuson X500 system users can store images and other data directly on the Acuson X500 system hard disk. Images and structured reports can be exported to a DICOM archive server or workstation on a network. In the following sections, Acuson X500 system Real World Activities are indicated by “Real World Activity” name while “X500 AE” indicates the invoked Application Entity. Similarly, the activities associated with service providers are indicated as “Real World Service Activity.”

4.0.1 Application Data Flow Diagram

Figure 1 illustrates the Acuson X500 system's Application Entity (AE), which is shown in the box. Relationships between user invoked activities (in the circles at the left of the AE) and the associated real-world activities provided by DICOM service providers (in the circles on the right side of the diagram) are shown.




Figure 1 Implementation Model

4.0.2 Verification

Verification is a part of the DICOM configuration located on the 'DICOM' page of the System Presets. Verification can be used to send a DICOM Verification request to a remote Application Entity (AE) and listen for a response.

When used as a diagnostic tool, Verification returns the following messages to the user:

- If the verification succeeds: "DICOM - Successfully contacted system"
- If the verification fails: "DICOM - Unable to communicate with system"




Figure 2 Verification Model.

4.0.3 DICOM Store

When requested the Acuson X500 sends images and/or structured reports to the preconfigured DICOM Storage server.

DICOM Store can be seen as two sub-operations:

- queueing images and/or structured reports for transfer
- transferring images and/or structured reports to the storage server.

Queueing images and structured reports for transfer:

Acuson X500 can be configured to automatically queue up images and structured reports for transfer as they are being created. "AutoStore to DICOM" option in DICOM presets has to be set for this.

Alternatively, user can select exams or individual images and manually queue them up from Review mode. When an exam is selected for DICOM store all images and structured reports (generally zero or one) will be queued. Structured reports can't be selected individually for store, the entire exam must be stored.

Transfer of images and structured reports to the storage server:

Further, once images and/or structured reports are queued they may be immediately transferred or delayed till the end of study using the transfer storage configuration.

Acuson X500 supports two storage configurations: “Store At End of Exam” and “Store During Exam”.

If the storage configuration is set to “Store At End of Exam” transfer attempts begin when the user selects “Close Study” or “New Patient”.

If the storage configuration is set to “Store During Exam”, transfer attempts to destination devices begin immediately after they are queued.

For both “Store At End of Exam” and “Store During Exam” settings, image and/or structured report transfer will be delayed if the Acuson X500 is busy performing another DICOM Store operation.




Figure 3 Store Model.

4.0.4 DICOM Print

Acuson X500 system is capable of grayscale (B/W) and color printing.

When requested, single frame images will be printed to a pre-configured DICOM network printer.

DICOM Print can be seen as two sub-operations:

- paging images for transfer
- transferring pages to printer

Paging images for transfer:

Acuson X500 can be configured to automatically queue up images to be printed on B/W Printer and/or Color printer as they are being created.

Alternatively, user can select exams or individual images and manually queue them up from Review mode for print.

Every image queued up is added into a page in the respective printer layout (DICOM B/W Printer Layout or DICOM Color Printer Layout).

Transfer of pages to the Printer:

Further, pages may be immediately transferred to the printer or delayed till the end of study based on the transfer configuration.

Acuson X500 supports two configurations: “Print At End of Exam” and “Print When Page Full”.

If the configuration is set to “Print At End of Exam,” transfer attempts of all pages to the destination DICOM printer begin as a batch when the user ends the exam.

If the configuration is set to “Print When Page Full”, transfer attempt of a page to the destination DICOM printer begins as soon as it becomes full.

For both “Print At End of Exam” and “Print when page full” settings, page transfer will be delayed if the Acuson X500 is busy performing another DICOM Print operation.




Figure 4 Print Model.

4.0.5 Patient Registration using Worklist

Patient registration can be automated by using the 'Worklist' Real World Activity. Pressing the 'New Patient' key on the keyboard initiates the patient data registration process and closes the previous active study. Pressing the 'Worklist' button on the patient data display screen invokes the Worklist query screen. The Worklist query screen can also be initiated from the Study screen.

Pressing the 'Search' button will attempt to find all matching patient data using the information entered on the Worklist Query screen. Patient name fields that are partially filled or empty will be treated as though an implicit wildcard was appended at the end of each field. Patient ID, Requested Procedure ID and Accession number will be exact match only. If no matches are found, a message will be presented to the operator indicating so. If more than one matching patient is found, a pick list of patient procedures will be presented to the user to select from. Each of the fields will be sortable in ascending and descending order.

The pick list of patient procedures will be limited to a number of preset entries. If more than this number of matching records are found in the query, the search will terminate and a the user will be notified. The search list criteria will contain:

- Patient name
- Patient ID
- Accession number
- Exam start date/time range
- Requested Procedure ID
- US/All modalities
- Scheduled station AE title




Figure 5 Modality Worklist Model

Once a Worklist query is initiated, a “Retrieving worklist, please wait ...” dialog will be presented to the user. The user will only have one option, “Cancel,” which will abort the query operation.

The following data fields in Modality Worklist Screen are initially populated from the New Patient Screen, if filled in, and can be used for query:

| Attribute Name | Tag |
|---------------------|-------------|
| Patient's Full Name | (0010,0010) |
| Patient ID | (0010,0020) |
| Accession Number | (0008,0050) |

The following data fields will be populated on the worklist screen for each return:

| Attribute Name | Tag |
|---------------------------------------|--------------------------|
| Patient's Full Name | (0010,0010) |
| Patient ID | (0010,0020) |
| Accession Number | (0008,0050) |
| Exam Start Date/Time | (0040,0002), (0040,0003) |
| Scheduled Procedure Step Sequence | (0040,0100)* |
| >Scheduled Procedure Step Description | (0040,0007) |
| >Scheduled Protocol Code Sequence | (0040,0008) |
| >>Code Value | (0008,0100) |
| Requested Procedure Description | (0032,1060) |
| Exam Type | (0008,1030)** |

*<code1>, ..., <codeN>: <sched1>, ..., <schedn>

where:

code<i> = Sequence item code value(0008,0100)

for a given sequence or value multiplicity

sched<i> = Scheduled procedure step(0040,0007)

for a given sequence or value multiplicity

**if a value exists for (0008,1030). Otherwise, Exam

Type is set to value of Scheduled procedure step

(0040,0007). If (0040,0007) is also empty, Exam

Type is set to Requested procedure Description

(0032,1060) if it exists.

The user will have the option to select a patient procedure step, or cancel the operation. Selection of a procedure step from the list will cause demographic information for the patient to be loaded in to the patient data fields.

The following data fields will be populated on the patient data screen:

| Attribute Name | Tag |
|----------------------------------|--------------------------|
| Patient Name (first,middle,last) | (0010,0010) |
| Patient ID | (0010,0020) |
| Accession number | (0008,0050) |
| Exam start date/time | (0040,0002), (0040,0003) |
| DOB | (0010,0030) |
| Sex | (0010,0040) |
| Weight | (0010,1030) |
| Height | (0010,1020) |

| Attribute Name | Tag |
|----------------|-------------|
| Physician | (0008,0090) |
| Indication | (0080,1080) |
| LMP | (0010,21D0) |

4.0.6 Modality Performed Procedure Step

The Acuson X500 System supports reporting of Modality Performed Procedure Step (MPPS) orders when the patient registration process utilizes the 'Worklist' Real World Activity. Procedure steps are presented to the operator after successful query of a server that supports the MPPS option. A detail window allows the operator access to individual scheduled procedure steps. Pressing the 'Procedures' push button on the Review Screen actualizes the detail window when multiple procedure steps are listed for the patient.




Figure 6 MPPS Model

4.0.7 Removable Media Storage

The Acuson X500 can perform DICOM operations to its standard on-board 120mm CD/DVD disk drive.

The Acuson X500 performs the File Set Creator and File Set Reader Roles for CD and DVD disks. The File Set Reader functionality does not support import of DICOM Structured Reports or measurements of imported images. Both limitations are overcome when DICOM and TIFF/AVI format is exported to CD/DVD. A DICOM conforming CD/DVD media is created when the user saves studies in DICOM format to the CD/DVD. A DICOM 3.0 conforming DICOMDIR file is created together with the directory structures, image files and structured reports (if any exist).

4.0.8 Storage Commitment

The user can exercise the Storage Commitment option by configuring and selecting a Storage Commitment server from the DICOM Presets menu. The Acuson X500 system requests commitment of images and structured reports (if any exist) and upon successful acknowledgment from the Storage server marks the study on the system hard drive as 'Archived'.




Figure 7 Storage Commitment Model

4.1 AE Functional Definition

4.1.1 Verification Real-World Activities

The Acuson X500 application entity performs Verification Service Class as an SCU and SCP allowing the operator to verify the ability of an application on a remote device to receive DICOM messages and allowing the operator of a remote DICOM device to verify the Acuson X500 system's ability to receive DICOM messages. (C-ECHO DIMSE)

4.1.2 Store Real-World Activities

The Acuson X500 Application Entity (AE) performs all of the functions to transmit ultrasound images, structured reports and associated data to network servers or workstations. The Acuson X500 AE supports the Ultrasound Image, Ultrasound Multi-Frame Image, Ultrasound Image (Retired), Ultrasound Multi-Frame (Retired) and Secondary Capture storage SOP classes SOP classes as an SCU.

The Acuson X500 AE also supports Structured Reports, for Obstetric and Cardiac studies, using the Comprehensive SR SOP Class as an SCU.

The Acuson X500 AE initiates an association for C-STORE Requests to store providers when the user invokes "DICOM Store". The association may be used to store multiple images, clips and structured reports and is closed when no images, clips or structured reports are available to be stored to the remote device for five seconds.

4.1.3 Storage Commitment - Push Model Real-World Activities

The Acuson X500 AE supports Storage Commitment Push Model SOP class to inform servers when all the store operations for a study have been completed. The Storage Commitment SCU uses the N-ACTION primitive to request safekeeping of a set of SOP Instances. The Storage Commitment SCU also processes the N-EVENT-REPORT primitives that are received from the SCP indicating 'successful' or 'non-successful' commitment status. The N-EVENT-REPORT information is used to mark a study as being successfully archived to a DICOM SCP.

The successful commit status and archival indication on the X500 does not ensure permanent archival of the images and Structured Reports. The operations performed by the SCP are dependent on its capabilities and configuration.

4.1.4 Print Real-World Activities

The Acuson X500 AE provides all aspects of the Print Management SCU. The Acuson X500 AE initiates an association to the printer when the user invokes "DICOM Print". The association may be used to print multiple pages and is closed when no pages are available to be printed to the remote device for five seconds.

4.2 Modality Worklist Real-World Activities

The Acuson X500 AE supports the DICOM Basic Worklist Management Service as an SCU. The AE initiates an association to the active Worklist server when a Worklist query is selected (via the "Worklist" button). The association is closed upon the completion of each query. A preset maximum number of matching results is accepted, at which point, the Acuson X500 AE issues a C-CANCEL-RQ request.

4.3 Modality Performed Procedure Step Real-World Activities

The Acuson X500 AE supports Modality Performed Procedure Step (MPPS) in the role of SCU. The Acuson X500 is capable of displaying scheduled procedure steps via the User Interface (UI) for Modality Performed Procedure Step. The operator can select a single PPS. The operator can notify the MPPS server that a MPPS is 'In Progress', 'Discontinued' or 'Completed'.

4.4 Removable Media Storage Real-World Activities

The Acuson X500 AE provides a standard implementation of DICOM Store to CD or DVD. The Acuson X500 AE selects one or more studies and exports the same to CD or DVD. Acuson X500 AE creates a DICOM File Format Image File for every image, clip and structured report in each of the selected studies.

A DICOMDIR file is created along with the files.

Measurements are not supported on imported images unless TIFF/AVI format is exported.

The DICOM SR cannot be imported from media unless the TIFF/AVI format is exported along with the DICOM SR.

4.5 Sequencing of Real-World Activities

Print, Store, Echo, Worklist, Storage Commit and MPPS commands can be transmitted simultaneously within the limits described below.

Storage Commit

The Storage Commitment (if enabled) command is sent in the following situations:

- a. On series close, when all images have previously stored successfully.
- b. The series is closed before all images are stored successfully, all previous stores have succeeded and the last image stores successfully.
- c. The series is closed before all images are stored successfully, at least one store has succeeded, at least one store has failed and the last store with non-zero retry count fails or succeeds.
- d. A series has been partially committed as in c. Later, due to "Retry Job" button press on the Store Status UI screen the store jobs are retried. Another Storage Commit is sent when at least one store has succeeded and the last store with non-zero retry count fails or succeeds.

MPPS

The MPPS (if enabled) command is sent in the following situations:

- a. N-CREATE command is sent whenever a new procedure step is selected. The state of the MPPS command is set to "In-Progress".
- b. N-SET command is sent when the Procedure Step is closed by the user pressing either the Completed or Discontinued button on the Close Procedure dialog. The state of the MPPS command is set, according to the state (Completed or Discontinued) set by the user.

5.0 AE Specifications

The following specifications apply to the Acuson X500 AE as depicted in Figure 1.

5.1 Acuson X500 AE Specification

The Acuson X500 AE provides conformance to the following DICOM Service SOP Classes as an SCU.

Table 2 Supported SOP Classes.

| Service SOP Class Name | SOP Class UID |
|---|-------------------------------|
| Verification | 1.2.840.10008.1.1 |
| Ultrasound Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.6 |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 |
| Ultrasound Multi-Frame Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.3 |
| Ultrasound Multi-Frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 |
| Storage Commitment - Push Model | 1.2.840.10008.1.20.1 |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 |
| Basic Color Print Management Meta SOP Class | 1.2.840.10008.5.1.1.18 |
| Basic Grayscale Image Box SOP Class | 1.2.840.10008.5.1.1.4 |
| Basic Color Image Box SOP Class | 1.2.840.10008.5.1.1.4.1 |
| Printer SOP Class | 1.2.840.10008.5.1.1.16 |
| Modality Worklist Information Model C- FIND | 1.2.840.10008.5.1.4.31 |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 |
| Comprehensive SR | 1.2.840.10008.5.1.4.1.1.88.33 |

5.1.1 Association Establishment Policies

5.1.1.1 General

The Acuson X500 system utilizes TCP/IP. The Maximum Length PDU negotiation is included in all association establishment requests. The maximum length PDU offered for an association initiated by Acuson X500 is:

- Maximum PDU Offered: 28672

5.1.1.2 Association Establishment Order

Acuson X500 initiates each C-Store Request one at a time, one for each transfer request being processed.

Image format on Acuson X500 can be set to one of “Automatic”, “Old Ultrasound” or “Secondary Capture”.

For the “Automatic” setting, Acuson X500 proposes Ultrasound Multi-Frame Image, Ultrasound Image, Ultrasound Multi-Frame Image (Retired), Ultrasound Image (Retired), Secondary Capture Image and Comprehensive SR sequentially.

For the “Old Ultrasound” setting, Acuson X500 proposes Ultrasound Multi-Frame Image (Retired), Ultrasound Image (Retired), Secondary Capture and Comprehensive SR Image to be negotiated sequentially.

For the “Secondary Capture” setting, Acuson X500 proposes Secondary Capture Image and Comprehensive SR to be negotiated sequentially.

5.1.1.3 Asynchronous Nature

All associations use the default synchronous mode of operation. Asynchronous Operations Window negotiations are not supported on the Acuson X500 system.

5.1.1.4 Implementation Identifying Information

- Implementation Class UID: “1.3.12.2.1107.5.5.5” (See below).
- Implementation Version Name:
“MergeCOM3_351MergeCOM3_351”

Siemens has provided registration for all Siemens Medical Solutions Groups. This unique Class UID is defined as:

“1.3.12.2.1107.5.5.product”

Where the interpretation is:

1. = International Standards Organization (ISO)

3. = International branch of ISO

12.2.1107.5. = Assigned to Siemens-UB MED

5. = Ultrasound Modality (SMS-UG)

Product = 5 - DICOM implementation for SONOLINE G20, G40, G50, G60 S, X300 and Acuson CV70, X500

5.1.2 Association Initiation by Real-World Activities

5.1.2.1 Real World Activity – Verification

The Acuson X500 is capable of supporting Verification service class as SCU or SCP. Verification can be initiated as a singular event from the Systems Presets menu to any configured SCP that supports Verification.

Proposed Presentation Contexts – Verification

The Acuson X500 will propose Presentation contexts as shown in table 3.

Table 3 Verification Presentation Context.

| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
|-----------------|-------------------|---------------------------|---------------------|----------|----------------------|
| Name | UID | Name List | UID List | | |
| Verification | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU/ SCP | None |
| Verification | 1.2.840.10008.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU/ SCP | None |
| Verification | 1.2.840.10008.1.1 | Explicit VR Big Endian | 1.2.840.10008.1.2.2 | SCU/ SCP | None |

5.1.2.2 Real World Activity – Store

Acuson X500 facilitates users to store images and structured reports as they are being created or later in review mode.

Queueing images during acquisition:

“Autostore to DICOM” option in DICOM presets has to be set. One or more of “Print/Store 1”, “Print/Store 2” and “Clip Store” keys on the control panel can be configured for Store (Disk Store, D.Store, Clip capture). When the user presses one of the configured keys, an image or clip is acquired, stored on the hard disk and queued up to be transferred to the storage server. Structured reports, if any, will be stored automatically after the study is closed and each time the report is modified after study close.

Queueing images and structured reports in Review mode:

User can select one or more closed studies and queue them up for Storage. The DICOM Store button is available in Review screen for this operation. All images and structured reports (if any) are stored. The study must be closed to generate a structured report.

Transfer of images to the storage server:

See section 4.1.2.

Associated Real World Activities

When images and/or structured reports are transferred from the hard disk to a DICOM Store SCP, the system establishes an association between the Acuson X500 AE and the configured DICOM device. The association may be used to store multiple images and/or structured reports and is closed when no images or structured reports are available to be stored to the remote device for five seconds.

Proposed Presentation Context

The following Presentation Contexts are presented to the SCP in an A-ASSOCIATE-RQ for DIMSE C-STORE storage services. The storage services utilize C-STORE services, as defined by the DICOM Standard. Table 4 represents all “Store” presentation contexts.

Table 4 Store Presentation Context.

| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
|--|-------------------------------|------------------------------------|------------------------|------|----------------------|
| Name | UID | Name List | UID List | | |
| Ultrasound Multi-Frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | Lossy JPEG 8 Bit Image Compression | 1.2.840.10008.1.2.4.50 | SCU | None |
| Ultrasound Multi-Frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| Ultrasound Multi-Frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Ultrasound Multi-Frame Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.3 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| Ultrasound Multi-Frame Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.3 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Ultrasound Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.6 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| Ultrasound Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.6 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |
| Comprehensive SR | 1.2.840.10008.5.1.4.1.1.88.33 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |

The Acuson X500 system always acts as an SCU for store and is the client in a client-server model.

Stress-Echo multi-frame images are only stored using Ultrasound Multi-Frame Image Storage with Lossy JPEG 8 Bit Image Compression. All other single and multi-frame images can be stored using other SOP classes and transfer syntaxes as described above.

SOP Specific Conformance to Storage Service SOP Classes

The Store Real World Activity provides standard extended conformance as an SCU for the following standard Storage Service Class SOP:

Table 5 Supported SOP Classes.

| Service SOP Class Name | SOP Class UID | Conformance Level |
|--|-------------------------------|-------------------|
| Ultrasound Multi-Frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | Standard Extended |
| Ultrasound Multi-Frame Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.3 | Standard Extended |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | Standard Extended |
| Ultrasound Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.6 | Standard Extended |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | Standard Extended |
| Comprehensive SR | 1.2.840.10008.5.1.4.1.1.88.33 | Standard Extended |

This is accomplished using the DIMSE C-STORE Service. The SCU issues a service request with a SOP instance that meets the requirements of the desired ultrasound, secondary capture, or structured report IOD.

The only Structured Report Template supported by the Acuson X500 is TID 5000 "OB-GYN Ultrasound Procedure Report".

The following table denotes the attributes included in the Ultrasound Image Object as implemented on the Acuson X500. Attributes not listed are not used.

Table 6 Ultrasound Image and Ultrasound Retired Image IOD Attributes

| Module | Attribute | Tag | Notes |
|------------------------|---------------------------------|-------------|---|
| Patient Identification | Patient's Name | (0010,0010) | X500 Patient Data Screen – Last Name, First & Middle fields. Populated from Modality Worklist if used. |
| | Patient ID | (0010,0020) | X500 Patient Data Screen – ID field. Default is today's date & time (e.g., 03_04_2003_17_54_43 = Apr. 3, 2003 at 5:54:43 PM). Populated from Modality Worklist if used. |
| Patient Demographic | Patient's Birth Date | (0010,0030) | X500 Patient Data Screen – DOB field. Default is a zero length attribute. Populated from Modality Worklist if used. |
| | Patient's Sex | (0010,0040) | X500 Patient Data Screen – Gender field. M = male F = female. O = Other Default is a zero length attribute. Populated from Modality Worklist if used. |
| | Patient's Age | (0010,1010) | Calculated from Patient Data Screen DOB field. |
| | Patient's Size | (0010,1020) | X500 Patient Data Screen – Height field. Populated from Modality Worklist if used. |
| | Patient's Weight | (0010,1030) | X500 Patient Data Screen – Weight field. Populated from Modality Worklist if used. |
| Patient Study | Admitting Diagnosis Description | (0008,1080) | X500 Patient Data Screen – Indication field. Populated from Modality Worklist if used. |
| General Study | Study Instance UID | (0020,000D) | Populated from Modality Worklist if used; generated by X500 otherwise |
| | Study Date | (0008,0020) | Date the exam started. |
| | Study Time | (0008,0030) | Time the exam started. |
| | Referring Physician's Name | (0008,0090) | X500 Patient Data Screen – Physician field. Populated from Modality Worklist if used. |
| | Study ID | (0020,0010) | Generated by X500 |
| | Accession Number | (0008,0050) | X500 Patient Data Screen – Accession # field. Populated from Modality Worklist if used. |

| Module | Attribute | Tag | Notes |
|----------------|---|-------------|---|
| | Study Description | (0008,1030) | Populated with the first attribute from Modality Worklist in this list that contains a valid value: Study Description (0008,1030), Scheduled Procedure Step Description (0040,0007), Requested Procedure Description (0032,1060). If Modality Worklist was not used or none of the attributes contains a valid value X500 Patient Data Screen – Indication field is used. |
| General Series | Modality | (0008,0060) | Always set to "US" |
| | Series Instance UID | (0020,000E) | Generated by X500 |
| | Series Number | (0020,0011) | Series Number in study (1-n). |
| | Laterality | (0020,0060) | Always sent as 0 length attribute |
| | ^(b) Series Date | (0008,0021) | Date the series started. |
| | ^(b) Series Time | (0008,0031) | Time the series started. |
| | ^(b) Series Description | (0008,103E) | Populated with Scheduled Procedure Step Description if a value was provided by Modality Worklist. |
| | ^(b) Protocol Name | (0018,1030) | The exam type of the most recent image stored in a particular series. If no images are stored for a series then the value is set to "Ultrasound". |
| | ^(b) Request Attributes Sequence | (0040,0275) | Populated with Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | > ^(b) Scheduled Procedure Step Description | (0040,0007) | Populated with Scheduled Procedure Step Description (0040,0007) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | > ^(b) Scheduled Protocol Code Sequence | (0040,0008) | Populated with Scheduled Protocol Code Sequence (0040,0008) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | > ^(b) Scheduled Procedure Step ID | (0040,0009) | Populated with Scheduled Procedure Step ID (0040,0009) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |

| Module | Attribute | Tag | Notes |
|-------------------|---|-------------|--|
| | > ^(b) Requested Procedure ID | (0040,1001) | Populated with Requested Procedure ID (0040,1001) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | ^(b) Performed Procedure Step Start Date | (0040,0244) | Date the Performed Procedure Step was started. |
| | ^(b) Performed Procedure Step Start Time | (0040,0245) | Time the Performed Procedure Step was started. |
| | ^(b) Performed Procedure Step ID | (0040,0253) | Populated with Scheduled Procedure Step ID (0040,0009) if provided by Modality Worklist. |
| | ^(b) Performed Procedure Step Description | (0040,0254) | Populated with Scheduled Procedure Step Description (0040,0007) if provided by Modality Worklist. |
| | ^(b) Performed Procedure Protocol Code Sequence | (0040,0260) | Populated with Scheduled Protocol Code Sequence (0040,0008) if provided by Modality Worklist. |
| | ^(b) Comments on the Performed Procedure Step | (0040,0280) | Populated with Comments on the Scheduled Procedure Step (0040,0400) if provided by Modality Worklist. |
| General Equipment | Manufacturer | (0008,0070) | Set to "Siemens Ultrasound" |
| | Institution Name | (0008,0080) | X500 System Presets – Organization Name field. |
| | Software Versions | (0018,1020) | Set to the DICOM Software Version |
| | Manufacturer's Model Name | (0008,1090) | Set to "Acuson X500" |
| General Image | Instance Number | (0020,0013) | Image number in study (1 – n) |
| | Patient Orientation | (0020,0020) | Always sent as 0 length attribute |
| Image Pixel | Samples per Pixel | (0028,0002) | Set to 1 for MONOCHROME2 images, 3 for RGB images. |
| | Photometric Interpretation | (0028,0004) | Set to "MONOCHROME2" or "RGB" |
| | Planar Configuration | (0028,0006) | Color-by-pixel. Set to 0 for RGB images, not provided for MONOCHROME2 images. |
| | Rows | (0028,0010) | Set to 480 for NTSC; 547 for PAL. For post-processed images and screen captures, this value may be up to 600. |

| Module | Attribute | Tag | Notes |
|--------------------|---|--------------|--|
| | Columns | (0028,0011) | Set to 640 for NTSC; 692 for PAL. For post-processed images and screen captures, this value may be up to 800. |
| | Bits Allocated | (0028,0100) | Set to 8. |
| | Bits Stored | (0028,0101) | Set to 8. |
| | High Bit | (0028,0102) | Set to 7. |
| | Pixel Representation | (0028,0103) | Set to 0 |
| | Pixel Data | (7FE0, 0010) | |
| US Image | Image Type | (0008,0008) | Always sent as a 0 length attribute. |
| | Heart Rate | (0018,1088) | Only provided if heart rate is > 0 |
| | Lossy Image Compression | (0028,2110) | "00" |
| SOP Common | SOP Class UID | (0008,0016) | 1.2.840.10008.5.1.4.1.1.6.1 or 1.2.840.10008.5.1.4.1.1.6 |
| | SOP Instance UID | (0008,0018) | Generated by X500 |
| | Specific Character Set | (0008,0005) | Set to values as defined in Section 8.4 of this document |
| Image Plane | Pixel Spacing | (0028,0030) | Pixel Spacing information is only provided for single, full screen, 2D image types (2D image types are B-mode, B-mode with color, B-mode with power). |
| Region Calibration | ^(c) Sequence of Ultrasound Regions | (0018,6011) | |
| | > ^(c) Region Spatial Format | (0018,6012) | B-Mode (Tissue or Color) = 0001H M-Mode (Tissue or Color) = 0002H Spectral (CW/PW) Doppler = 0003H |
| | > ^(c) Region Data Type | (0018,6014) | B-Mode, M-Mode = 0001H (Tissue) Spectral Doppler = 0004H (CW Spectral Doppler) Spectral Doppler = 0003H (PW Spectral Doppler) |
| | > ^(c) Region Flags | (0018,6016) | 1st Bit (LSB) = 1 (All images acquired are transparent) 2nd Bit = 1 (All images acquired are automatically scaled) 3rd Bit = 1 for frequency scale 3rd Bit = 0 for velocity scale. The value of the 3rd bit is undefined for any mode other than Doppler. The value for 3rd bit is undefined if both frequency and velocity scales are selected on the Doppler image. 4th Bit is Reserved and value is always 0 |

| Module | Attribute | Tag | Notes |
|--------------------|---|-------------|--|
| | > ^(c) Region Location Min X0 | (0018,6018) | |
| | > ^(c) Region Location Min Y0 | (0018,601A) | |
| | > ^(c) Region Location Max X1 | (0018,601C) | |
| | > ^(c) Region Location Max Y1 | (0018,601E) | |
| | > ^(c) Physical Units X direction | (0018,6024) | B-Mode (Tissue or Color) = 0003H (cm) M-Mode (Tissue or Color) = 0004H (seconds) Spectral (CW/PW) Doppler = 0004H (seconds) |
| | > ^(c) Physical Units Y direction | (0018,6026) | B-Mode (Tissue or Color) = 0003H (cm) M-Mode (Tissue or Color) = 0003H (cm) Spectral (CW/PW) Doppler = 0007H (cm/sec) |
| | > ^(c) Physical Delta X | (0018,602C) | |
| | > ^(c) Physical Delta Y | (0018,602E) | |
| | > ^(c) Reference Pixel X0 | (0018,6020) | Attribute only set for Spectral Doppler Regions |
| | > ^(c) Reference Pixel Y0 | (0018,6022) | Attribute only set for Spectral Doppler Regions |
| | > ^(c) Reference Pixel Physical Value X | (0018,6028) | Attribute only set for Spectral Doppler Regions When provided, value is always 0. |
| | > ^(c) Reference Pixel Physical Value Y | (0018,602A) | Attribute only set for Spectral Doppler Regions When provided, value is always 0. |
| Private Attributes | ^(a) Private Creator | (0011,0010) | Reserves tags 0011,1000 through 0011,10FF for use as private tags. |
| | ^(a) Siemens Medical Solutions Model Name | (0011,1010) | Always set to "Acuson X500". |
| | ^(a) DIMAQ Software Version | (0011,1011) | Set to version of DIMAQ software installed. |
| | ^(a) Private Data | (0011,1020) | For internal X500 use only. |
| | ^(a) Private Data | (0011,1021) | For internal X500 use only. |
| | ^(a) Private Creator | (0013,0010) | Reserves tags 0013,1000 through 0013,10FF for use as private tags. |

| Module | Attribute | Tag | Notes |
|--------|---|-------------|--|
| | ^(a) Siemens Medical Solutions Model Name | (0013,1010) | Always set to "Acuson X500". |
| | ^(a) DIMAQ Software Version | (0013,1011) | Set to version of DIMAQ software installed. |
| | ^(a) Private Data | (0013,1020) | For internal X500 use only. |
| | ^(a) Private Creator | (0015,0010) | This group is populated only if data is available. Reserves tags 0015,1000 through 0015,10FF for use as private tags. |
| | ^(a) Siemens Medical Solutions Model Name | (0015,1010) | Always set to "Acuson X500". |
| | ^(a) DIMAQ Software Version | (0015,1011) | Set to version of DIMAQ software installed. |
| | ^(a) Private Data | (0015,1020) | For internal X500 use only. |
| | ^(a) Private Creator | (0017,0010) | This group is populated only if data is available. Reserves tags 0017,1000 through 0017,10FF for use as private tags. |
| | ^(a) Siemens Medical Solutions Model Name | (0017,1010) | Always set to "Acuson X500". |
| | ^(a) DIMAQ Software Version | (0017,1011) | Set to version of DIMAQ software installed. |
| | ^(a) Private Data | (0017,1020) | For internal X500 use only. |
| | Private Creator | (0019,0010) | Reserves tags 0019,1000 through 0019,10FF for use as private tags. |
| | Import Structured Reports | (0019,1020) | Set to "O" if Obstetric or cardiac SR options was purchased and SR generation was configured. Otherwise set to "No". Instructs SCP that it should attempt to import Obstetric or cardiac measurements from SR. |

^(a)The Attribute is only provided if the image is written to media.

^(b)The Attribute is only provided if the procedure step is queried from the MWL server.

^(c)Region Calibration is provided only for 2D (B-Mode), M-Mode and Spectral Doppler Regions. Region Calibration is not supported on Ultrasound RETIRED images, Screen Captures and post-processed images. Region Calibration is not supported for M-Mode or Spectral Doppler still images taken from Live Imag-

ing.

Table 7 Ultrasound MultiFrame and Ultrasound MultiFrame Retired Image IOD Attributes

| Module | Attribute | Tag | Notes |
|------------------------|---------------------------------|-------------|---|
| Patient Identification | Patient's Name | (0010,0010) | X500 Patient Data Screen – Last Name, First & Middle fields. Populated from Modality Worklist if used. |
| | Patient ID | (0010,0020) | X500 Patient Data Screen – ID field. Default is today's date & time (e.g., 03_04_2003_17_54_43 = Apr. 3, 2003 at 5:54:43 PM). Populated from Modality Worklist if used. |
| Patient Demographic | Patient's Birth Date | (0010,0030) | X500 Patient Data Screen – DOB field. Default is a zero length attribute. Populated from Modality Worklist if used. |
| | Patient's Sex | (0010,0040) | X500 Patient Data Screen – Gender field. M = male F = female. O= Other Default is a zero length attribute. Populated from Modality Worklist if used. |
| | Patient's Age | (0010,1010) | Calculated from Patient Data Screen DOB field. |
| | Patient's Size | (0010,1020) | X500 Patient Data Screen – Height field. Populated from Modality Worklist if used. |
| | Patient's Weight | (0010,1030) | X500 Patient Data Screen – Weight field. Populated from Modality Worklist if used. |
| Patient Study | Admitting Diagnosis Description | (0008,1080) | X500 Patient Data Screen – Indication field. Populated from Modality Worklist if used. |
| General Study | Study Instance UID | (0020,000D) | Populated from Modality Worklist if used; generated by X500 otherwise. |
| | Study Date | (0008,0020) | Date the exam started. |
| | Study Time | (0008,0030) | Time the exam started. |
| | Referring Physician's Name | (0008,0090) | X500 Patient Data Screen – Physician field. Populated from Modality Worklist if used. |
| | Study ID | (0020,0010) | Generated by X500 |
| | Accession Number | (0008,0050) | X500 Patient Data Screen – Accession # field. Populated from Modality Worklist if used. |

| Module | Attribute | Tag | Notes |
|----------------|---|-------------|---|
| | Study Description | (0008,1030) | Populated with the first attribute from Modality Worklist in this list that contains a valid value: Study Description (0008,1030), Scheduled Procedure Step Description (0040,0007), Requested Procedure Description (0032,1060). If Modality Worklist was not used or none of the attributes contains a valid value X500 Patient Data Screen – Indication field is used. |
| General Series | Modality | (0008,0060) | Always set to "US" |
| | Series Instance UID | (0020,000E) | Generated by X500 |
| | Series Number | (0020,0011) | Series Number in study (1-n). |
| | Laterality | (0020,0060) | Always sent as 0 length attribute |
| | ^(b) Series Date | (0008,0021) | Date the series started. |
| | ^(b) Series Time | (0008,0031) | Time the series started. |
| | ^(b) Series Description | (0008,103E) | Populated with Scheduled Procedure Step Description if a value was provided by Modality Worklist. |
| | ^(b) Protocol Name | (0018,1030) | The exam type of the most recent image stored in a particular series. If no images are stored for a series then the value is set to "Ultrasound". |
| | ^(b) Request Attributes Sequence | (0040,0275) | Populated with Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | > ^(b) Scheduled Procedure Step Description | (0040,0007) | Populated with Scheduled Procedure Step Description (0040,0007) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | > ^(b) Scheduled Protocol Code Sequence | (0040,0008) | Populated with Scheduled Protocol Code Sequence (0040,0008) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | > ^(b) Scheduled Procedure Step ID | (0040,0009) | Populated with Scheduled Procedure Step ID (0040,0009) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |

| Module | Attribute | Tag | Notes |
|-------------------|---|-------------|--|
| | > ^(b) Requested Procedure ID | (0040,1001) | Populated with Requested Procedure ID (0040,1001) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | ^(b) Performed Procedure Step Start Date | (0040,0244) | Date the Performed Procedure Step was started. |
| | ^(b) Performed Procedure Step Start Time | (0040,0245) | Time the Performed Procedure Step was started. |
| | ^(b) Performed Procedure Step ID | (0040,0253) | Populated with Scheduled Procedure Step ID (0040,0009) if provided by Modality Worklist. |
| | ^(b) Performed Procedure Step Description | (0040,0254) | Populated with Scheduled Procedure Step Description (0040,0007) if provided by Modality Worklist. |
| | ^(b) Performed Procedure Protocol Code Sequence | (0040,0260) | Populated with Scheduled Protocol Code Sequence (0040,0008) if provided by Modality Worklist. |
| | ^(b) Comments on the Performed Procedure Step | (0040,0280) | Populated with Comments on the Scheduled Procedure Step (0040,0400) if provided by Modality Worklist. |
| General Equipment | Manufacturer | (0008,0070) | Set to "Siemens Ultrasound" |
| | Institution Name | (0008,0080) | X500 System Presets – Organization Name field. |
| | Software Versions | (0018,1020) | Set to the DICOM Software Version |
| | Manufacturer's Model Name | (0008,1090) | Set to "Acuson X500". |
| General Image | Instance Number | (0020,0013) | Image number in study (1 - n) |
| | Patient Orientation | (0020,0020) | Always sent as 0 length attribute |
| Image Pixel | Samples per Pixel | (0028,0002) | Set to 1 for MONOCHROME2 images, 3 for RGB and YBR_FULL_422 images. |
| | Photometric Interpretation | (0028,0004) | "YBR_FULL_422" if sent compressed, "RGB" or "MONOCHROME2" if sent uncompressed. |
| | Planar Configuration | (0028,0006) | Color-by-pixel. Set to 0 for RGB and YBR_FULL_422 images, not provided for MONOCHROME2 images. |

| Module | Attribute | Tag | Notes |
|-------------|---|--------------|--|
| | Rows | (0028,0010) | Set to 480 for NTSC; 547 for PAL. For Stress Echo clips this value may be from 228 to 288. |
| | Columns | (0028,0011) | Set to 640 for NTSC; 692 for PAL. For Stress Echo clips this value may be from 288 to 384. |
| | Bits Allocated | (0028,0100) | Set to 8. |
| | Bits Stored | (0028,0101) | Set to 8. |
| | High Bit | (0028,0102) | Set to 7. |
| | Pixel Representation | (0028,0103) | Set to 0. |
| | Pixel Data | (7FE0, 0010) | |
| US Image | Image Type | (0008,0008) | Sent as a 0 length attribute. |
| | ^(a) View List | (0009,212A) | Private attribute |
| | > ^(a) View Name | (0009,2120) | Private attribute |
| | ^(a) Stage Name | (0008,2120) | |
| | ^(a) Stage Number | (0008,2122) | |
| | ^(a) Number of Stages | (0008,2124) | |
| | ^(a) View Name | (0008,2127) | |
| | ^(a) View Number | (0008,2128) | |
| | ^(a) Number of Views in Stage | (0008,212A) | |
| | ^(a) Trigger Time | (0018,1060) | |
| | ^(a) Nominal Interval | (0018,1062) | |
| | Heart Rate | (0018,1088) | Only provided if heart rate is > 0 |
| | Lossy Image Compression | (0028,2110) | Always set to "01" |
| SOP Common | SOP Class UID | (0008,0016) | 1.2.840.10008.5.1.4.1.1.3.1 or 1.2.840.10008.5.1.4.1.1.3 |
| | SOP Instance UID | (0008,0018) | Generated by X500 |
| | Specific Character Set | (0008,0005) | Set to values as defined in Section 8.4 of this document |
| Image Plane | Pixel Spacing | (0028,0030) | Pixel Spacing information is only provided for single, full screen, 2D image types (2D image types are B-Mode, B-Mode Color, B-Mode with power). |
| Cine | Frame Time | (0018,1063) | |
| Multi-Frame | Number of Frames | (0028,0008) | |
| | Frame Increment Pointer | (0028,0009) | 00181063H |

| Module | Attribute | Tag | Notes |
|----------------------|-----------------------------------|-------------|---|
| Region Calibration | (c)Sequence of Ultrasound Regions | (0018,6011) | |
| | >(c)Region Spatial Format | (0018,6012) | B-Mode (Tissue or Color) = 0001H M-Mode (Tissue or Color) = 0002H Spectral (CW/PW) Doppler = 0003H |
| | >(c)Region Data Type | (0018,6014) | B-Mode, M-Mode = 0001H (Tissue) Spectral Doppler = 0004H (CW Spectral Doppler) Spectral Doppler = 0003H (PW Spectral Doppler) |
| | >(c)Region Flags | (0018,6016) | 1st Bit (LSB) = 1 (All images acquired are transparent) 2nd Bit = 1 (All images acquired are automatically scaled) 3rd Bit = 1 for frequency scale 3rd Bit = 0 for velocity scale. The value of the 3rd bit is undefined for any mode other than Doppler. The value for 3rd bit is undefined if both frequency and velocity scales are selected on the Doppler image. 4th Bit is Reserved and value is always 0 |
| | >(c)Region Location Min X0 | (0018,6018) | |
| | >(c)Region Location Min Y0 | (0018,601A) | |
| | >(c)Region Location Max X1 | (0018,601C) | |
| | >(c)Region Location Max Y1 | (0018,601E) | |
| | >(c)Physical Units X direction | (0018,6024) | B-Mode (Tissue or Color) = 0003H (cm) M-Mode (Tissue or Color) = 0004H (seconds) Spectral (CW/PW) Doppler = 0004H (seconds) |
| | >(c)Physical Units Y direction | (0018,6026) | B-Mode (Tissue or Color) = 0003H (cm) M-Mode (Tissue or Color) = 0003H (cm) Spectral (CW/PW) Doppler = 0007H (cm/sec) |
| | >(c)Physical Delta X | (0018,602C) | |
| >(c)Physical Delta Y | (0018,602E) | | |

| Module | Attribute | Tag | Notes |
|--------------------|---|-------------|---|
| | > ^(c) Reference Pixel X0 | (0018,6020) | Attribute only set for Spectral Doppler Regions |
| | > ^(c) Reference Pixel Y0 | (0018,6022) | Attribute only set for Spectral Doppler Regions |
| | > ^(c) Reference Pixel Physical Value X | (0018,6028) | Attribute only set for Spectral Doppler Regions When provided, value is always 0. |
| | > ^(c) Reference Pixel Physical Value Y | (0018,602A) | Attribute only set for Spectral Doppler Regions When provided, value is always 0. |
| | > ^(a,c) Transducer Frequency | (0018,6030) | |
| | > ^(a,c) Pulse Repetition Frequency | (0018,6032) | |
| | > ^(a,c) Doppler Correction Angle | (0018,6034) | |
| Private Attributes | ^(d) Private Creator | (0011,0010) | Reserves tags 0011,1000 through 0011,10FF for use as private tags. |
| | ^(d) Siemens Medical Solutions Model Name | (0011,1010) | Always set to "Acuson X500". |
| | ^(d) DIMAQ Software Version | (0011,1011) | Set to version of DIMAQ software installed. |
| | ^(d) Private Data | (0011,1020) | For internal X500 use only. |
| | ^(d) Private Data | (0011,1021) | For internal X500 use only. |
| | ^(d) Private Creator | (0013,0010) | Reserves tags 0013,1000 through 0013,10FF for use as private tags. |
| | ^(d) Siemens Medical Solutions Model Name | (0013,1010) | Always set to "Acuson X500". |
| | ^(d) DIMAQ Software Version | (0013,1011) | Set to version of DIMAQ software installed. |
| | ^(d) Private Data | (0013,1020) | For internal X500 use only. |
| | ^(d) Private Creator | (0015,0010) | This group is populated only if data is available. Reserves tags 0015,1000 through 0015,10FF for use as private tags. |
| | ^(d) Siemens Medical Solutions Model Name | (0015,1010) | Always set to "Acuson X500". |
| | ^(d) DIMAQ Software Version | (0015,1011) | Set to version of DIMAQ software installed. |
| | ^(d) Private Data | (0015,1020) | For internal X500 use only. |

| Module | Attribute | Tag | Notes |
|--------|---|-------------|--|
| | ^(d) Private Creator | (0017,0010) | This group is populated only if data is available. Reserves tags 0017,1000 through 0017,10FF for use as private tags. |
| | ^(d) Siemens Medical Solutions Model Name | (0017,1010) | Always set to "Acuson X500". |
| | ^(d) DIMAQ Software Version | (0017,1011) | Set to version of DIMAQ software installed. |
| | ^(d) Private Data | (0017,1020) | For internal X500 use only. |
| | Private Creator | (0019,0010) | Reserves tags 0019,1000 through 0019,10FF for use as private tags. |
| | Import Structured Reports | (0019,1020) | Set to "O" if Obstetric SR options was purchased and SR generation was configured. Otherwise set to "No". Instructs SCP that it should attempt to import Obstetric measurements from SR. |

^(a)The Attribute is only provided for Stress Echo Images.

^(b)The Attribute is only provided if the procedure step is queried from the MWL server.

^(c)Region Calibration is provided only for 2D (B-Mode), M-Mode and Spectral Doppler Regions. Region Calibration is not supported on Ultrasound Retired images, Screen Captures and post-processed images. Region Calibration is not supported for M-Mode or Spectral Doppler still images taken from Live Imaging.

^(d)The Attribute is only provided if the image is written to media.

Table 8 Secondary Capture Image IOD Attributes

| Module | Attribute | Tag | Notes |
|------------------------|---------------------------------|-------------|---|
| Patient Identification | Patient's Name | (0010,0010) | X500 Patient Data Screen – Last Name, First & Middle fields. Populated from Modality Worklist if used. |
| | Patient ID | (0010,0020) | X500 Patient Data Screen – ID field. Default is today's date & time (e.g., 03_04_2003_17_54_43 = Apr. 3, 2003 at 5:54:43 PM). Populated from Modality Worklist if used. |
| Patient Demographic | Patient's Birth Date | (0010,0030) | X500 Patient Data Screen – DOB field. Default is a zero length attribute. Populated from Modality Worklist if used. |
| | Patient's Sex | (0010,0040) | X500 Patient Data Screen – Gender field. M = male F = female. O= Other Default is a zero length attribute. Populated from Modality Worklist if used. |
| | Patient's Age | (0010,1010) | Calculated from Patient Data Screen DOB field. |
| | Patient's Size | (0010,1020) | X500 Patient Data Screen – Height field. Populated from Modality Worklist if used. |
| | Patient's Weight | (0010,1030) | X500 Patient Data Screen – Weight field. Populated from Modality Worklist if used. |
| Patient Study | Admitting Diagnosis Description | (0008,1080) | X500 Patient Data Screen – Indication field. Populated from Modality Worklist if used. |
| General Study | Study Instance UID | (0020,000D) | Populated from Modality Worklist if used; generated by X500 otherwise. |
| | Study Date | (0008,0020) | Date the exam started. |
| | Study Time | (0008,0030) | Time the exam started. |
| | Referring Physician's Name | (0008,0090) | X500 Patient Data Screen – Physician field. Populated from Modality Worklist if used. |
| | Study ID | (0020,0010) | Generated by X500 |
| | Accession Number | (0008,0050) | X500 Patient Data Screen – Accession # field. Populated from Modality Worklist if used. |

| Module | Attribute | Tag | Notes |
|----------------|---|-------------|---|
| | Study Description | (0008,1030) | Populated with the first attribute from Modality Worklist in this list that contains a valid value: Study Description (0008,1030), Scheduled Procedure Step Description (0040,0007), Requested Procedure Description (0032,1060). If Modality Worklist was not used or none of the attributes contains a valid value X500 Patient Data Screen – Indication field is used. |
| General Series | Modality | (0008,0060) | Always set to "US" |
| | Series Instance UID | (0020,000E) | Generated by X500 |
| | Series Number | (0020,0011) | Series Number in study (1-n). |
| | Laterality | (0020,0060) | Always sent as 0 length attribute |
| | ^(b) Series Date | (0008,0021) | Date the series started. |
| | ^(b) Series Time | (0008,0031) | Time the series started. |
| | ^(b) Series Description | (0008,103E) | Populated with Scheduled Procedure Step Description if a value was provided by Modality Worklist. |
| | ^(b) Protocol Name | (0018,1030) | The exam type of the most recent image stored in a particular series. If no images are stored for a series then the value is set to "Ultrasound". |
| | ^(b) Request Attributes Sequence | (0040,0275) | Populated with Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | > ^(b) Scheduled Procedure Step Description | (0040,0007) | Populated with Scheduled Procedure Step Description (0040,0007) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | > ^(b) Scheduled Protocol Code Sequence | (0040,0008) | Populated with Scheduled Protocol Code Sequence (0040,0008) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | > ^(b) Scheduled Procedure Step ID | (0040,0009) | Populated with Scheduled Procedure Step ID (0040,0009) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |

| Module | Attribute | Tag | Notes |
|---------------------|---|-------------|--|
| | > ^(b) Requested Procedure ID | (0040,1001) | Populated with Requested Procedure ID (0040,1001) from Scheduled Procedure Step Sequence (0040,0100) if provided by Modality Worklist. |
| | ^(b) Performed Procedure Step Start Date | (0040,0244) | Date the Performed Procedure Step was started. |
| | ^(b) Performed Procedure Step Start Time | (0040,0245) | Time the Performed Procedure Step was started. |
| | ^(b) Performed Procedure Step ID | (0040,0253) | Populated with Scheduled Procedure Step ID (0040,0009) if provided by Modality Worklist. |
| | ^(b) Performed Procedure Step Description | (0040,0254) | Populated with Scheduled Procedure Step Description (0040,0007) if provided by Modality Worklist. |
| | ^(b) Performed Procedure Protocol Code Sequence | (0040,0260) | Populated with Scheduled Protocol Code Sequence (0040,0008) if provided by Modality Worklist. |
| | ^(b) Comments on the Performed Procedure Step | (0040,0280) | Populated with Comments on the Scheduled Procedure Step (0040,0400) if provided by Modality Worklist. |
| SC Equipment Module | Conversion Type | (0008,0064) | Set to "WSD" |
| General Equipment | Manufacturer | (0008,0070) | Set to "Siemens Ultrasound" |
| | Institution Name | (0008,0080) | X500 System Presets – Organization Name field. |
| | Software Versions | (0018,1020) | Set to the DICOM Software Version |
| | Manufacturer's Model Name | (0008,1090) | Set to "X500". |
| General Image | Instance Number | (0020,0013) | Image number in study (1 – n) |
| | Patient Orientation | (0020,0020) | Always sent as 0 length attribute. |
| Image Pixel | Samples per Pixel | (0028,0002) | Set to 1 for MONOCHROME2 images, 3 for RGB images. |
| | Photometric Interpretation | (0028,0004) | "RGB" or "MONOCHROME2" |
| | Planar Configuration | (0028,0006) | Color-by-pixel. Set to 0 for RGB images, not provided for MONOCHROME2 images. |
| | Rows | (0028,0010) | Set to 480 for NTSC; 547 for PAL. |

| Module | Attribute | Tag | Notes |
|--------------------|---|-------------|---|
| | Columns | (0028,0011) | Set to 640 for NTSC; 692 for PAL. |
| | Bits Allocated | (0028,0100) | Set to 8. |
| | Bits Stored | (0028,0101) | Set to 8. |
| | High Bit | (0028,0102) | Set to 7. |
| | Pixel Representation | (0028,0103) | Set to 0. |
| | Pixel Data | (7FE0,0010) | |
| SOP Common | SOP Class UID | (0008,0016) | 1.2.840.10008.5.1.4.1.1.7 |
| | SOP Instance UID | (0008,0018) | Generated by X500 |
| | Specific Character Set | (0008,0005) | Set to values as defined in Section 8.4 of this document |
| Private Attributes | ^(a) Private Creator | (0011,0010) | Reserves tags 0011,1000 through 0011,10FF for use as private tags. |
| | ^(a) Siemens Medical Solutions Model Name | (0011,1010) | Always set to "Acuson X500". |
| | ^(a) DIMAQ Software Version | (0011,1011) | Set to version of DIMAQ software installed. |
| | ^(a) Private Data | (0011,1020) | For internal X500 use only. |
| | ^(a) Private Data | (0011,1021) | For internal X500 use only. |
| | ^(a) Private Creator | (0013,0010) | Reserves tags 0013,1000 through 0013,10FF for use as private tags. |
| | ^(a) Siemens Medical Solutions Model Name | (0013,1010) | Always set to "Acuson X500". |
| | ^(a) DIMAQ Software Version | (0013,1011) | Set to version of DIMAQ software installed. |
| | ^(a) Private Data | (0013,1020) | For internal X500 use only. |
| | ^(a) Private Creator | (0015,0010) | This group is populated only if data is available. Reserves tags 0015,1000 through 0015,10FF for use as private tags. |
| | ^(a) Siemens Medical Solutions Model Name | (0015,1010) | Always set to "Acuson X500". |
| | ^(a) DIMAQ Software Version | (0015,1011) | Set to version of DIMAQ software installed. |
| | ^(a) Private Data | (0015,1020) | For internal X500 use only. |
| | ^(a) Private Creator | (0017,0010) | This group is populated only if data is available. Reserves tags 0017,1000 through 0017,10FF for use as private tags. |

| Module | Attribute | Tag | Notes |
|--------|---|-------------|--|
| | ^(a) Siemens Medical Solutions Model Name | (0017,1010) | Always set to "Acuson X500". |
| | ^(a) DIMAQ Software Version | (0017,1011) | Set to version of DIMAQ software installed. |
| | ^(a) Private Data | (0017,1020) | For internal X500 use only. |
| | Private Creator | (0019,0010) | Reserves tags 0019,1000 through 0019,10FF for use as private tags. |
| | Import Structured Reports | (0019,1020) | Set to "O" if Obstetric SR options was purchased and SR generation was configured. Otherwise set to "No". Instructs SCP that it should attempt to import Obstetric measurements from SR. |

^(a)The Attribute is only provided if the image is written to media.

^(b)The Attribute is only provided if the procedure step is queried from the MWL server.

The following table denotes the attributes included in the Comprehensive SR Object as implemented on the Acuson X500. Attributes not listed are not used.

Table 9 Comprehensive SR IOD Attributes

| Module | Attribute | Tag | Notes |
|---------|----------------------|-------------|---|
| Patient | Patient's Name | (0010,0010) | X500 Patient Data Screen – Last Name, First & Middle fields. Populated from Modality Worklist if used. |
| | Patient ID | (0010,0020) | X500 Patient Data Screen – ID field. Default is today's date & time (e.g., 03_04_2003_17_54_43 = Apr. 3, 2003 at 5:54:43 PM). Populated from Modality Worklist if used. |
| | Patient's Birth Date | (0010,0030) | X500 Patient Data Screen – DOB field. Default is a zero length attribute. Populated from Modality Worklist if used. |
| | Patient's Sex | (0010,0040) | X500 Patient Data Screen – Gender field. M = male F = female. O = Other Default is a zero length attribute. Populated from Modality Worklist if used. |

| Module | Attribute | Tag | Notes |
|--------------------|--|-------------|---|
| | Patient's Size | (0010,1020) | X500 Patient Data Screen – Height field. Populated from Modality Worklist if used. |
| | Patient's Weight | (0010,1030) | X500 Patient Data Screen – Weight field. Populated from Modality Worklist if used. |
| Patient Study | Admitting Diagnosis Description | (0008,1080) | X500 Patient Data Screen – Indication field. Populated from Modality Worklist if used. |
| General Study | Study Instance UID | (0020,000D) | Populated from Modality Worklist if used; generated by X500 otherwise |
| | Study Date | (0008,0020) | Date the exam started. |
| | Study Time | (0008,0030) | Time the exam started. |
| | Referring Physician's Name | (0008,0090) | X500 Patient Data Screen – Physician field. Populated from Modality Worklist if used. |
| | Study ID | (0020,0010) | Generated by X500 |
| | Accession Number | (0008,0050) | X500 Patient Data Screen – Accession # field. Populated from Modality Worklist if used. |
| | Study Description | (0008,1030) | Populated with the first attribute from Modality Worklist in this list that contains a valid value: Study Description (0008,1030), Scheduled Procedure Step Description (0040,0007), Requested Procedure Description (0032,1060). If Modality Worklist was not used or none of the attributes contains a valid value X500 Patient Data Screen – Indication field is used. |
| SR Document Series | Modality | (0008,0060) | Always set to "SR" |
| | Series Instance UID | (0020,000E) | Generated by X500 |
| | Series Number | (0020,0011) | Series Number in study (2-n). |
| | Series Date | (0008,0021) | Date the series started. |
| | Series Time | (0008,0031) | Time the series started. |
| | Referenced Performed Procedure Step Sequence | (0008,1111) | Populated with MPPS SOP Class UID and MPPS SOP instance UID of MPPS command sent for the procedure step(s) performed. |
| General Equipment | Manufacturer | (0008,0070) | Set to "Siemens Ultrasound" |
| | Institution Name | (0008,0080) | X500 System Presets – Organization Name field. |
| | Software Versions | (0018,1020) | Set to the DICOM Software Version |

| Module | Attribute | Tag | Notes |
|---------------------|---|-------------|---|
| | Manufacturer's Model Name | (0008,1090) | Set to "Acuson X500" |
| SR Document General | Content Date | (0008,0023) | Date the report was created |
| | Content Time | (0008,0033) | Time the report was created |
| | Instance Number | (0020,0013) | Always set to 0. |
| | Completion Flag | (0040,A491) | Always set to "PARTIAL" |
| | Verification Flag | (0040,A493) | Always set to "UNVERIFIED" |
| | Predecessor Documents Sequence | (0040,A360) | Supplied if a previous SR was generated for the study. Populated with SOP Class UID and SOP Instance UID of the previous Obstetric SRs for the study, if any. See table C17-2 in PS 3.3-2004 for sequence definition. |
| | Performed Procedure Code Sequence | (0040,A372) | Populated with contents of Procedure Code Sequence from Modality Worklist if available, empty otherwise. See table C17-2 in PS 3.3-2004 for sequence definition. |
| | Current Requested Procedure Evidence Sequence | (0040,A375) | Lists all images and clips in the study. See table C17-2 in PS 3.3-2004 for sequence definition. |
| SOP Common | SOP Class UID | (0008,0016) | 1.2.840.10008.5.1.4.1.1.88.33 |
| | SOP Instance UID | (0008,0018) | Generated by X500. |
| | Specific Character Set | (0008,0005) | Set to values as defined in Section 8.4 of this document. |
| | Instance Creation Date | (0008,0012) | Date the SOP Instance was created. |
| | Instance Creation Time | (0008,0013) | Time the SOP Instance was created. |
| | Instance Creator UID | (0008,0014) | |
| Private Attributes | Private Creator | (0019,0010) | Reserves tags 0019,1000 through 0019,10FF for use as private tags. |
| | Import Structured Reports | (0019,1020) | Set to "O" if Obstetric SR options was purchased and SR generation was configured. Otherwise set to "No". Instructs SCP that it should attempt to import Obstetric measurements from SR. |

Error Handling

The following table indicates the response status codes that are handled by the Acuson X500 AE, which a SCP may return following the SCU's C-STORE-RSP command.

A successful C-STORE operation will allow the Acuson X500 AE to continue to the next action desired by the user.

Table 10 C-STORE Status Responses.

| Service Status | Further Meaning | Protocol Codes | Related Fields |
|----------------|---|----------------------|----------------|
| Refused | Out of resources. | A7xx | None |
| Error | Data set does not match SOP Class. Cannot understand. | A9xx Cxxx | None |
| Warning | Coercion of data Elements. Data set does not match SOP Class. Elements discarded. | B000 B007 B006 | None |
| Success | | 0000 | None |

If the C-STORE operation is not successful, the image(s) and Structured Report(s), if any, are spooled on the Acuson X500 hard drive. A user-configured number of additional attempts are made to store the image(s) and Structured Report(s). If these attempts fail, the user must select the job and press "Retry Job" on the Store Status page to complete the C-STORE operation.

All image and Structured Report storage on the Acuson X500 system hard drive is temporary in nature. If an attempt is made to store images on a full Acuson X500 system hard drive, the system will attempt to delete studies archived to CD or DICOM. If no deleteable data exists, a "DISK FULL" message is displayed on the Acuson X500 system display. The user must then delete studies not archived in order to store additional images.

5.1.2.3 Real World Activity - Print

Acuson X500 facilitates user to print images as they are being created or later in review mode.

Paging images during acquisition

One or more of "Print/Store 1" and "Print/Store 2" keys on the control panel can be configured for Print (DICOM B/W Print and/or DICOM Color Print). When the user presses one of the configured keys on the control panel, the image is acquired, stored on the hard disk and placed in a page under the respective printer layout (DICOM B/W Printer Layout or DICOM Color Printer Layout).

Paging images in Review mode

User can select either individual images from open or closed studies, or one or more closed studies and queue them up for print. DICOM B/W Printer and DICOM Color Printer buttons are available in Review screen for this operation. When a study is selected for print, all single-frame images belonging to the study will be printed.

Transfer of pages to the Printer

Pages may be immediately transferred or delayed till the end of study using the transfer configuration.

Acuson X500 supports two configurations: “Print At End of Exam” and “Print When Page Is Full”.

If the configuration is set to “Print At End of Exam”, all pages queued to destination devices will be transferred as a batch when the user selects “Close Study” or “New Patient”.

If the configuration is set to “Print When Page Is Full”, a page is transferred to destination devices immediately after it is full.

For both “Print At End of Exam”, and “Print When Page Is Full” settings, image transfer will be delayed if the Acuson X500 is busy performing another DICOM Command (Store/Print/Echo).

Associated Real World Activities

An association is established when the user initiates a “B/W Print” or “Color Print” operation from the Review screen. Individual images or entire exams can be transferred to the selected DICOM Print device. The association is closed no pages are available to be printed for five seconds. An association may also be opened after a network outage or when the system is powered-on if images are queued to be printed.

Proposed Presentation Context to a Grayscale Print Server

Table 11 Grayscale Print Presentation Context.

| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
|---|-----------------------|---------------------------------|---------------------|------|----------------------|
| Name | UID | Name List | UID List | | |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU | None |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

SOP Specific Conformance to Basic Grayscale Print Management Meta SOP Class

The Acuson X500 AE provides standard conformance of the Grayscale Meta SOP classes as an SCU. Specifically, with respect to the Basic Grayscale Print Management Meta SOP Class this means conformance to the underlying SOP classes:

Table 12 Conformance to Grayscale Print Meta SOP Class.

| SOP Class Name | SOP Class UID | Conformance Level |
|-------------------------------------|------------------------|-------------------|
| Basic Film Session SOP Class | 1.2.840.10008.5.1.1.1 | Standard |
| Basic Film Box SOP Class | 1.2.840.10008.5.1.1.2 | Standard |
| Basic Grayscale Image Box SOP Class | 1.2.840.10008.5.1.1.4 | Standard |
| Printer SOP Class | 1.2.840.10008.5.1.1.16 | Standard |

All mandatory elements of these classes are supported.

Specific Conformance to Basic Film Session SOP Class

DICOM specified usage - M = Mandatory; U = User Option

Table 13 Supported DIMSE Services for Basic Film Session SOP Class.

| Name | Usage | Description |
|----------|-------|---------------------------|
| N-Create | M | Creates the Film Session. |
| N-Set | U | Not used. |
| N-Delete | U | Deletes the Film Session. |
| N-Action | U | Not used. |

SOP Specific Conformance to Basic Film Box SOP Class

Table 14 Supported DIMSE Services for Basic Film Box SOP Class.

| Name | Usage | Description |
|----------|-------|--|
| N-Create | M | Creates the Film Box. |
| N-Set | U | Not used. |
| N-Delete | U | Deletes the Film Box. Issued after each film is printed. |
| N-Action | M | PRINT. Sent after each Film Box is filled, and at the end of the exam to force a print of partially filled Film Box. |

Table 15 Attributes set for the Basic Film Box SOP Class.

| Attribute Name | Attribute Tag | Usage | Range | Description |
|---------------------------|---------------|-------|--|--|
| Image Display Format | (2010,0010) | M | STANDARD\ X,Y | Where X, Y can be configured/ selected as 1*1, 1*2, 2*2, 2*3, 3*2, 3*3, 3*5, 4*5, 4*6, 5*6 |
| Film Orientation | (2010,0040) | U | PORTRAIT LANDSCAPE | Range may be limited by print server/printer. |
| Film Size ID | (2010,0050) | U | 8INX10IN 8.5INX11IN 10INX12IN 10INX14IN 11INX14IN 11INX17IN 14INX14IN 14INX17IN 24CMX24CM 24CMX30CM A3 A4 | Range may be limited by print server/printer. |
| Magnification Type | (2010,0060) | U | REPLICATE BILINEAR CUBIC NONE | |
| Min. Density | (2010,0120) | U | 0-99,999,999 | Printer specific |
| Max Density | (2010,0130) | U | 0-99,999,999 | Printer specific |
| Configuration Information | (2010,0150) | U | | Printer specific |
| Smoothing Type | (2010,0080) | U | | Printer specific |
| Border Density | (2010,0100) | U | BLACK WHITE | |
| Empty Image Density | (2010,0110) | U | BLACK WHITE | |
| Trim | (2010,0140) | U | YES NO | |

SOP Specific Conformance to Basic Grayscale Image Box SOP Class

Table 16 Supported DIMSE Services for the Basic Grayscale Image Box SOP.

| Name | Usage | Description |
|-------|-------|--|
| N-Set | M | The SCP for each potential image of the film box creates an image box instance. Only those instances, which actually contain images, will be updated with the N-SET message. |

Table 17 Attributes set for the Basic Grayscale Image Box SOP Class.

| Name | Attribute | Range | Description |
|----------------|-------------|-----------------|---|
| Image Position | (2020,0010) | 1-30 | Value according to Image Display Format |
| Polarity | (2020,0020) | NORMAL, REVERSE | Intensity mapping between display and print |

Table 18 Supported DIMSE Services for the Printer SOP.

| Name | Usage | Description |
|----------------|-------|---|
| N-Event-Report | M | Ignored and not handled. |
| N-Get | U | May be issued by this device at any time to get printer status. |

Table 19 Supported Printer SOP Class Elements.

| Name | Usage | Range | Description |
|----------------------------|-------|--------------------|--|
| Printer Status | U | WARNING FAILURE | During a "Failure" the Print job will be displayed as "Failed" |
| Printer Status Information | U | Vendor specific | Reported to user if printer status = WARNING or FAILURE. |

Proposed Presentation Context to a Color Print Server

Table 20 Color Print Server Presentation Context.

| Name | Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
|---|------------------------|---------------------------------|---------------------|--|------|----------------------|
| | UID | Name List | UID List | | | |
| Basic Color Print Management Meta SOP Class | 1.2.840.10008.5.1.1.18 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 | | SCU | None |
| Basic Color Print Management Meta SOP Class | 1.2.840.10008.5.1.1.18 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | | SCU | None |

SOP Specific Conformance to Basic Color Print Management Meta SOP Class

The Acuson X500 Print AE provides standard conformance to the color printing Meta SOP classes as an SCU. Specifically, with respect to the Basic Color Print Management Meta SOP Class this means conformance to the underlying SOP classes:

Table 21 Conformance to Color Print Meta SOP Class.

| SOP Class Name | SOP Class UID | Conformance Level |
|---------------------------------|-------------------------|-------------------|
| Basic Film Session SOP Class | 1.2.840.10008.5.1.1.1 | Standard |
| Basic Film Box SOP Class | 1.2.840.10008.5.1.1.2 | Standard |
| Basic Color Image Box SOP Class | 1.2.840.10008.5.1.1.4.1 | Standard |
| Printer SOP Class | 1.2.840.10008.5.1.1.16 | Standard |

SOP Specific Conformance to Basic Color Image Box SOP Class

The Basic Color Print Management Meta SOP Class makes identical use of the *Basic Film Session SOP Class*, *Basic Film Box SOP Class* and *Printer SOP Class* elements, which have been previously described for grayscale image printing. Therefore, these will not be described again in this section on color printing. However, it should be noted that certain attributes, such as Medium Type which is defined in the Basic Film Session SOP Class, are highly likely to require printer/print server specific media.

Table 22 Supported DIMSE Services for the Basic Color Image Box SOP Class.

| Name | Usage | Description |
|-------|-------|--|
| N-Set | M | The SCP for each potential image of the film box creates an image box instance. Only those instances, which actually contain images, will be updated with the N-SET message. |

Table 23 Attributes set for the Basic Color Image Box SOP Class.

| Name | Attribute | Range | Description |
|----------------------|-------------|----------------|-------------------------------------|
| Planar Configuration | (0028,0006) | Color-by-plane | Red plane, Green plane, Blue plane. |

The Printer SOP Class behavior is identical to that used for grayscale printing.

Error Handling

The Acuson X500 Print AE supports the following error codes and reports failures to the user.

Table 24 Supported Error Codes for Printer Classes.

| Service Status | Further Meaning | Protocol Codes |
|----------------|---|-------------------|
| Success | Film accepted for Printing | 0000 |
| Warning | Film accepted for Printing, one or more settings ignored. | 107,116,B600,B605 |
| Failure | Printing not successful | C602, C603, C613 |

If the print operation is not successful, the image(s) are spooled on the Acuson X500 hard drive. A user-configured number of additional attempts are made to print the image(s). If these attempts fail, the user must select the job and press "Retry Job" on the Print Status page to complete the print operation.

5.1.2.4 Real World Activity - Worklist

A separate Network association is established by the AE for each Worklist query operation, with only one active query at a time. The association is closed at completion of the query.

Table 25 Worklist Presentation Context Table

| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
|--|------------------------|---------------------------|-------------------|------|----------------------|
| Name | UID | Name List | UID List | | |
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

The Acuson X500 will always act as an SCU and as the client in a client-server model.

SOP Specific Conformance to Modality Worklist Service SOP Classes

The Worklist AE provides conformance to the following DICOM Service SOP Classes as an SCU all at a standard extended level of conformance:

Table 26 Supported SOP Classes

| Supported SOP Class Name | SOP Class UID | Conformance Level |
|--|------------------------|-------------------|
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | Standard Extended |

The following table provides the list of attributes requested in the Modality Worklist Query.

Table 27 Modality Worklist Information Model Attributes

| Attribute Name | Tag |
|--|--------------|
| Specific Character Set | (0008,0005) |
| Accession number | (0008,0050)* |
| Referring Physician's Name | (0008,0090) |
| Study Description | (0008,1030) |
| Admitting Diagnoses Description | (0008,1030) |
| Referenced Study Sequence | (0008,1110) |
| >Referenced SOP Class UID | (0008,1150) |
| >Referenced SOP Instance UID | (0008,1155) |
| Patient's Name | (0010,0010)* |
| Patient ID | (0010,0020)* |
| Patient's Birth Date | (0010,0030) |
| Patient's Sex | (0010,0040) |
| Patient's Size | (0010,1020) |
| Patient's Weight | (0010,1030) |
| Medical Alerts | (0010,2000) |
| Contrast Allergies | (0010,2110) |
| Pregnancy Status | (0010,21C0) |
| Last Menstrual Date | (0010,21D0) |
| Patient Comments | (0010,4000) |
| Study Instance UID | (0020,000D) |
| Requesting Physician | (0032,1032) |
| Requested Procedure Description | (0032,1060) |
| Requested Procedure Code Sequence | (0040,0008) |
| >Code Value | (0008,0100) |
| >Coding Scheme Designator | (0008,0102) |
| >Coding Scheme Version | (0008,0103) |
| >Code Meaning | (0008,0104) |
| Special Needs | (0038,0050) |
| Patient State | (0038,0500) |
| Scheduled Procedure Step Sequence | (0040,0100) |
| >Modality | (0008,0060) |
| >Scheduled Station AE Title | (0040,0001)* |
| >Scheduled Procedure Step Start Date | (0040,0002) |
| >Scheduled Procedure Step Start Time | (0040,0003) |
| >Scheduled Performing Physician's Name | (0040,0006) |
| >Scheduled Procedure Step Description | (0040,0007) |
| >Scheduled Protocol Code Sequence | (0040,0008) |
| >>Code Value | (0008,0100) |

| Attribute Name | Tag |
|--|--------------|
| >>Coding Scheme Designator | (0008,0102) |
| >>Coding Scheme Version | (0008,0103) |
| >>Code Meaning | (0008,0104) |
| >Scheduled Procedure Step ID | (0040,0009) |
| >Comments on the Scheduled Procedure Step | (0040,0400) |
| Requested Procedure ID | (0040,1001)* |
| Reason for the Requested Procedure | (0040,1002) |
| *Indicates parameter may be populated for query. | |

5.1.2.5 Real World Activity - Modality Performed Procedure Step

This operation allows the AE to create an instance of the Modality Performed Procedure Step SOP Class (MPPS) and provide information about a specific real world Performed Procedure Step that is under control of the SCU. This operation is invoked through the DIMSE N-CREATE and N-SET services.

Only the IHE (refer to IHE Rev 5.5) Simple and Abandoned Cases for the relationship between Scheduled Procedure Steps and Performed Procedure Steps is supported. Both cases specify that a 1-to-1 relationship must exist between Scheduled Procedure Step and Performed Procedure Step. In the Simple Case the Performed Procedure Step is completed successfully. In the Abandoned Case the Performed Procedure Step is abandoned before being completed.

A list of scheduled procedures and procedure steps will be accessible from the Worklist and Procedure screens. The Performed Procedure Step User Interface allows the operator to set the status of the performed procedure step. The system shall establish an association for N-CREATE and N-SET, if another N-CREATE or N-SET is available within 5 seconds, it will be sent using the same association.

Starting a Performed Procedure Step

When the user depresses the 'OK' button on the New Patient Screen a performed procedure SOP Class instance will be created using the N-CREATE DIMSE service for the selected scheduled procedure.

Ending a Performed Procedure Step

When the user selects 'Completed' or 'Discontinued' from the MPPS User Interface, the performed procedure step will be closed using the N-SET DIMSE service.

New Patient Request

If the 'New Patient' button is selected and there are opened performed procedure steps, the user shall be prompted for a closure status for the opened procedure step by the MPPS User Interface. Any opened procedure steps must be closed before any 'new patient' data can be entered.

System Shutdown

If the user requests 'System Shutdown' and there is an open performed procedure step, the user will be prompted for a closure status for the open procedure step. All procedure steps should be closed before the system can be shutdown. Failure to close a procedure step will result in the procedure step being set to Discontinued.

Error Handling

If the MPPS operation is not successful, the MPPS command is spooled on the Acuson X500 hard drive. A user-configured number of additional attempts are made to complete the MPPS Commands. If these attempts fail, the user must select this job and press "Retry Job" on the Store Status page to complete the MPPS operation.

Proposed Presentation Context

| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
|-----------------------------------|-------------------------|---------------------------|-------------------|------|----------------------|
| Name | UID | Name List | UID List | | |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

Table 28 MPPS Presentation Context Table

The Acuson X500 system will always act as an SCU and be the client in a client – server model.

SOP Specific Conformance to Modality Performed Procedure Step SOP Classes

The Modality Performed Procedure Step AE provides a conforming implementation of the following DICOM Service SOP Class as an SCU at a standard extended level of conformance.

Table 29 Supported SOP Class

| Supported SOP Class Name | SOP Class UID | Conformance Level |
|-----------------------------------|-------------------------|-------------------|
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Standard Extended |

The following tables provide the list of attributes supported by the AE in the implementation of MPPS SOP Class including N-CREATE, N-SET and Final State attributes. The Acuson X500 sends N-SET only at final state.

Table 30 Modality Performed Procedure Step Attributes in N-CREATE

| Attribute | Tag | Notes |
|---------------------------------------|-------------|---|
| Specific Character Set | (0008,0005) | Set to values as defined in Section 8.4 of this document |
| Scheduled Step Attribute Sequence | (0040,0270) | |
| >Study Instance UID | (0020,000D) | Value obtained from Modality WorkList; generated by X500 in some cases |
| >Referenced Study Sequence | (0008,1110) | Populated with contents of Referenced Study Sequence from Modality Worklist if used, empty otherwise. See table F.7.2-1 in PS 3.4-2004 for sequence definition. |
| >Referenced Patient Sequence | (0008,1120) | Always empty |
| >Accession Number | (0008,0050) | Value obtained from Modality WorkList |
| >Requested Procedure ID | (0040,1001) | Value obtained from Modality WorkList |
| >Requested Procedure Description | (0032,1060) | Value obtained from Modality WorkList |
| >Scheduled Procedure Step ID | (0040,0009) | Value obtained from Modality WorkList |
| >Scheduled Procedure Step Description | (0040,0007) | Value obtained from Modality WorkList |
| >Scheduled Protocol Code Sequence | (0040,0008) | |
| >>Code Value | (0008,0100) | Value obtained from Modality WorkList |
| >>Coding Scheme Designator | (0008,0102) | Value obtained from Modality WorkList |
| >>Coding Scheme Version | (0008,0103) | Value obtained from Modality WorkList |
| >>Code Meaning | (0008,0104) | Value obtained from Modality WorkList |
| Patient's Name | (0010,0010) | Value obtained from Modality WorkList |
| Patient ID | (0010,0020) | Value obtained from Modality WorkList. |
| Patient's Birth Date | (0010,0030) | Value obtained from Modality WorkList |
| Patient's Sex | (0010,0040) | Value obtained from Modality WorkList |
| Performed Procedure Step ID | (0040,0253) | Value obtained from Modality WorkList |
| Performed Station AE Title | (0040,0241) | The AE title of the X500 on which the procedure was performed. |
| Performed Station Name | (0040,0242) | |
| Performed Location | (0040,0243) | |
| Performed Procedure Step Start Date | (0040,0244) | The start date of the performed procedure step. |
| Performed Procedure Step Start Time | (0040,0245) | The start time of the performed procedure step. |
| Performed Procedure Step Status | (0040,0252) | Always set to "In-Progress". |
| Performed Procedure Step Description | (0040,0254) | Value obtained from Modality WorkList |

| Attribute | Tag | Notes |
|--------------------------------------|-------------|---|
| Performed Procedure Type Description | (0040,0255) | Always sent as 0 length attribute |
| Procedure Code Sequence | (0008,1032) | Populated with contents of Requested Procedure Code Sequence from Modality Worklist if used, empty otherwise. See table F.7.2-1 in PS 3.4-2004 for sequence definition. |
| Performed Procedure Step End Date | (0040,0250) | Always sent as 0 length attribute |
| Performed Procedure Step End Time | (0040,0251) | Always sent as 0 length attribute |
| Modality | (0008,0060) | Always set to US |
| Study ID | (0020,0010) | Populated from Requested Procedure ID (0040,1001) if Modality Worklist is used; created by X500 otherwise |
| Performed Protocol Code Sequence | (0040,0260) | Always empty |
| Performed Series Sequence | (0040,0340) | Always empty |

Table 31 Modality Performed Procedure Step Attributes in N-SET

| Attribute | Tag | Notes |
|-----------------------------------|-------------|--|
| Performed Procedure Step Status | (0040,0252) | Set to "Discontinued" or "Completed" based on user selection. |
| Performed Procedure Step End Date | (0040,0250) | Date the procedure step was completed |
| Performed Procedure Step End Time | (0040,0251) | Time the procedure step was completed |
| Performed Series Sequence | (0040,0340) | Shall contain only one series |
| >Performing Physician's Name | (0008,1050) | |
| >Protocol Name | (0008,1030) | Exam type specified by the operator. |
| >Operator's Name | (0008,1070) | |
| >Series Instance UID | (0020,000E) | The Instance UID of the series to which the procedure belongs. |
| >Series Description | (0008,103E) | Always sent as 0 length attribute |
| >Retrieve AE Title | (0008,0054) | Always sent as 0 length attribute |

| Attribute | Tag | Notes |
|---|-------------|---|
| >Referenced Image Sequence | (0008,1140) | List of all the images in the series. |
| >>Referenced SOP Class UID | (0008,1150) | The SOP class UID can be one of: Ultrasound Multi-Frame Image Storage 1.2.840.10008.5.1.4.1.1.3.1 Ultrasound Multi-Frame Image Storage (Retired) 1.2.840.10008.5.1.4.1.1.3 Ultrasound Image Storage 1.2.840.10008.5.1.4.1.1.6.1 Ultrasound Image Storage (Retired) 1.2.840.10008.5.1.4.1.1.6 Secondary Capture Image Storage 1.2.840.10008.5.1.4.1.1.7 |
| >>Referenced SOP Instance UID | (0008,1155) | The SOP instance UID of the image. |
| >Referenced Non-Image Composite SOP Instance Sequence | (0040,0220) | Always empty |

5.1.2.6 Real-World Activity Storage Commitment

This operation allows the AE to create an instance of the Storage Commitment SOP Class and to provide information about a specific Real World Activity that is under the control of the SCU. The AE invokes a request for safekeeping of images by the N-ACTION REQUEST. Referenced in the N-ACTION Request are the SOP class UID(s) and SOP instance UID(s) for all STORE Class objects requesting commitment by the SCU.

Storage Commit

The Storage Commitment (if enabled) command is sent in the following situations:

- a. On series close, when all images and Structured Reports have previously stored successfully.
- b. The series was previously closed, all previous stores have succeeded and the last image or Structured Report stores successfully.
- c. The series was previously closed, at least one store has succeeded, at least one store has failed and the last store with non-zero retry count fails or succeeds.
- d. A series has been partially committed as in c. Later, due to "Retry Job" button press on the Store Status UI screen the store jobs are retried. Another Storage Commit is sent when at least one store has succeeded and the last store with non-zero retry count fails or succeeds.

The Acuson X500 waits for the return of a successful N-ACTION RESPONSE Status Code applicable for the associated request indicating whether the commitment request was successful or a failure. The Acuson X500 waits for the N-EVENT REPORT from the SCP for at most 48 hours.

The Acuson X500 is capable of accepting the N-EVENT REPORT on the association it initiates for the N-ACTION or one initiated by the SCP. Studies with all SOP instances marked as 'successful' in the N-EVENT REPORT will be eligible for deletion from the system hard drive.

The Acuson X500 allows the user to configure a Storage Commitment Server which may be different from the Storage Server. Thus, the Storage Commitment SCP must wait for an appropriate time for the stored images to arrive from the Storage server.

Image-By-Image and Batch Storage Commitment are supported as specified in "Vista DICOM Conformance Requirements for Image Modalities in radiology, Cardiology, Dental, Ophthalmology and other specialties" (Version 2.3).

Storage Commitment of Structured Reports is supported.

Proposed Presentation Context

| Abstract Syntax | | Transfer Syntax | | Role | Extended Negotiation |
|-------------------------------|----------------------|---------------------------|-------------------|------|----------------------|
| Name | UID | Name List | UID List | | |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU | None |

Table 32 Storage Commitment Presentation Context Table

The Acuson X500 system will act as an SCU in the 'Push Model' Storage Commitment SOP Class.

SOP Specific Conformance to Storage Commitment SOP Class

The Storage Commitment AE provides conformance to the following DICOM Service SOP Class as an SCU at a standard level of conformance.

| Supported SOP Class Name | SOP Class UID | Conformance Level |
|-------------------------------|----------------------|-------------------|
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Standard |

Storage Commitment to Storage Media (CD) is not supported.

Table 33 Supported SOP Class

The following table provides the list of attributes supported by the AE in the implementation of Storage Commitment SOP Class:

Table 34 Storage Commitment Request Attributes in N-ACTION REQUEST

| Attribute | Tag | Notes |
|------------------------------|-------------|-------------------|
| Transaction UID | (0008,1195) | Generated by X500 |
| Referenced SOP Sequence | (0008,1199) | |
| >Referenced SOP Class UID | (0008,1150) | |
| >Referenced SOP Instance UID | (0008,1155) | |

5.1.2.7 Error Handling

If the storage commitment operation is not successful, a user-configured number of additional attempts are made. If these attempts fail, the user must select the job and press “Retry Job” on the DICOM Store Queue page to complete the storage commitment operation.

6.0 Removable Media Interchange Specifications

This implementation supports 120mm CD and DVD medium.

6.1 Supported Application Profiles

Acuson X500 provides standard conformance to the following four Ultrasound Application Profiles. A DICOM 3.0 conformant DICOMDIR file is created together with the directory structures and image files.

Table 35 Application Profiles, Real-World Activities, and Roles

| Supported AP | Real-World Activity | Roles | SC Option |
|------------------|---------------------|----------|-------------|
| STD-US-ID-SF-CDR | Create CD-R | FSC, FSR | Interchange |
| STD-US-ID-MF-CDR | Create CD-R | FSC, FSR | Interchange |
| STD-US-SC-SF-CDR | Create CD-R | FSC | Interchange |
| STS-US-SC-MF-CDR | Create CD-R | FSC | Interchange |
| STD-US-ID-SF-DVD | Create DVD | FSC, FSR | Interchange |
| STD-US-ID-MF-DVD | Create DVD | FSC, FSR | Interchange |
| STD-US-SC-SF-DVD | Create DVD | FSC | Interchange |
| STS-US-SC-MF-DVD | Create DVD | FSC | Interchange |

6.2 Supported SOP Classes

6.2.1 Supported SOP Classes and Transfer Syntaxes

This implementation provides standard conformance to the following DICOM 3.0 SOP Classes.

Table 36 Transfer Syntaxes for Media Interchange

| Service SOP Class Name | SOP Class UID | Transfer Syntax Name | Transfer Syntax UID List |
|--------------------------------------|-------------------------------|------------------------------------|--------------------------|
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Ultrasound Multi-Frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | Lossy JPEG 8 Bit Image Compression | 1.2.840.10008.1.2.4.50 |
| Comprehensive SR | 1.2.840.10008.5.1.4.1.1.88.33 | DICOM Explicit VR Little Endian | 1.2.840.10008.1.2.1 |

6.3 Information Object Definition and DICOMDIR Keys

6.3.1 DICOM File Meta Information

The following table denotes the attributes included in the Ultrasound Image Object as implemented on the Acuson X500 in addition to the attributes listed in Table 6.

Table 37 US Image Attributes Used (Refer Table 6 for additional attributes used)

| Attribute Name | Tag | Notes |
|--------------------------------|-------------------------|---|
| File Preamble | No Tag or Length fields | All bytes are set to 00H |
| DICOM Prefix | No Tag or Length fields | Set to DICOM Prefix "DICM" |
| Group length | (0002,0000) | |
| File Meta Information Version | (0002,0001) | Always set to 0001H |
| Media Storage SOP Class UID | (0002,0002) | Always Ultrasound Image 1.2.840.10008.5.1.4.1.1.6.1 |
| Media Storage SOP Instance UID | (0002,0003) | |
| Transfer Syntax UID | (0002,0010) | Always Explicit VR Little Endian 1.2.840.10008.1.2.1 |
| Implementation Class UID | (0002,0012) | Always set to 1.3.12.2.1107.5.5.5 |
| Implementation Version Name | (0002,0013) | Always set to MergeCOM3_351 |

The following table denotes the attributes included in the Ultrasound Multi-Frame Image Object as implemented on the Acuson X500 in addition to the attributes listed in Table 7.

Table 38 USMF Image Attributes Used (Refer Table 7 for additional attributes used)

| Attribute | Tag | Notes |
|--------------------------------|-------------------------|--|
| File Preamble | No Tag or Length fields | All bytes are set to 0 |
| DICOM Prefix | No Tag or Length fields | Set to "DICM" |
| Group length | (0002,0000) | |
| File Meta Information Version | (0002,0001) | Always set to 0001H |
| Media Storage SOP Class UID | (0002,0002) | Always Ultrasound Multi-Frame Image 1.2.840.10008.5.1.4.1.1.3.1 |
| Media Storage SOP Instance UID | (0002,0003) | |
| Transfer Syntax UID | (0002,0010) | Always Lossy JPEG 8 Bit Compression 1.2.840.10008.1.2.4.50 |
| Implementation Class UID | (0002,0012) | Always set to 1.3.12.2.1107.5.5.5 |
| Implementation Version Name | (0002,0013) | Always set to MergeCOM3_351 |

The following table denotes the attributes included in the Comprehensive SR Object as implemented on the Acuson X500 in addition to the attributes listed in Table 8.

Table 39 Comprehensive SR Attributes Used (Refer Table 8 for additional attributes used)

| Attribute | Tag | Notes |
|--------------------------------|-------------------------|--|
| File Preamble | No Tag or Length fields | All bytes are set to 0 |
| DICOM Prefix | No Tag or Length fields | Set to "DICM" |
| Group length | (0002,0000) | |
| File Meta Information Version | (0002,0001) | Always set to 0001H |
| Media Storage SOP Class UID | (0002,0002) | Always Comprehensive SR 1.2.840.10008.5.1.4.1.1.88.33 |
| Media Storage SOP Instance UID | (0002,0003) | |
| Transfer Syntax UID | (0002,0010) | Always Explicit VR Little Endian 1.2.840.10008.1.2.1 |
| Implementation Class UID | (0002,0012) | Always set to 1.3.12.2.1107.5.5.5 |
| Implementation Version Name | (0002,0013) | Always set to MergeCOM3_351 |

6.3.2 Basic Directory Information Object Definitions - File-set Identification Module

| Attribute | Tag | Notes |
|-------------|-------------|--|
| File-Set ID | (0004,1130) | Set to serial number + YYMMDD + 3 digit counter. Volume Label has this same value. |

6.3.3 Basic Directory Information Object Definitions - Directory Identification Module

| Attribute | Tag | Notes |
|--|-------------|-------|
| Offset of the First Directory Record of the Root Directory Entry | (0004,1200) | |
| Offset of the Last Directory Record of the Root Directory Entry | (0004,1202) | |
| File-set Consistency Flag | (0004,1212) | |
| Directory Record Sequence | (0004,1220) | |
| >Offset of the Next Directory Record | (0004,1400) | |
| >Record In-use Flag | (0004,1410) | |
| >Offset of Referenced Lower-Level Directory Entity | (0004,1420) | |
| >Directory Record Type | (0004,1430) | |

| Attribute | Tag | Notes |
|--------------------------------------|-------------|--------------|
| >Referenced File ID | (0004,1500) | |
| >Referenced SOP Class UID in File | (0004,1510) | |
| >Referenced SOP Instance UID in File | (0004,1511) | |

6.3.4 Physical Storage Media and Media Formats

The physical storage media supported are 120mm CD-R, CD-RW, DVD-R, DVD+R, DVD-RW, and DVD+RW medium.

7.0 Communication Profiles

All Acuson X500 system application entities utilize the DICOM 3.0 TCP/IP communication support as defined in PS3.8 (Part 8) of the DICOM 3.0 Standard.

7.1 TCP/IP Stack Supported

Each process inherits its TCP/IP stack from the Acuson X500's operating systems TCP/IP stack. The local AE Port number is always set to 104.

7.1.1 Physical Media Supported

Standard representations of IEEE 802.3 10BaseT/100BaseT ("twisted pair") is supported

7.1.2 Chapter Extensions/Specializations/Privatizations

Pixel Spacing information is only provided for single, full screen, and 2D image types (B-mode, B-mode with color, and B-mode with power).

The private attributes listed in the following table are used by the Acuson X500 AE with Stress Echo images.

Table 40 Private Attributes

| Attribute Name | Tag | VR | Description |
|----------------|-------------|----|---|
| View List | (0009,212A) | SQ | Names of all views represented in the study |
| >View Name | (0009,2120) | SH | Name of a view |

Appendix A lists the DICOM SR mappings used by the Acuson X500 AE in Obstetric Structured Reports. All private concept names use the Coding Scheme Designator "G60S".

8.0 Configuration

Acuson X500 Networking and DICOM parameters can be configured through the Acuson X500 System Presets Menu screens. The following configuration is supported:

- General system
- Network (local and remote)
- DICOM Store
- DICOM Print
- DICOM Modality Worklist
- DICOM Storage Commitment
- DICOM Modality Performed Procedure Step

8.1 General System Configuration

The following system parameter can be configured via the Acuson X500 System Presets Basic Menu screens. This parameter is mapped to a DICOM image attribute:

- Hospital Name

8.1.1 Hospital Name

The user can enter the organization (i.e. hospital, clinic, etc.) as a text string in the Hospital Name field of the System Presets - General menu. The Organization Name field is transferred to DICOM devices as Institution Name - DICOM data element (0008, 0080).

8.2 DICOM Network Configuration

DICOM and networking parameters can be configured for both the local Acuson X500 device and remote DICOM service class providers through the System Presets DICOM Menu.

8.2.1 Local

The Acuson X500 local network parameters are configurable. The following network parameters can be configured for a Acuson X500 device:

- Host Name
- IP address
- Subnet IP mask
- Default Gateway

- DICOM Application Entity Title

8.2.2 Remote

Multiple DICOM service class providers can be configured through the system presets. The following network parameters can be configured for each remote device:

- DICOM Device Application Entity Title
- IP address
- Port Number

8.2.2.1 DICOM Store Configuration

Several configuration settings are provided in addition to those described in Section 8.2.2.

The Image Format setting provides control over the Presentation Contexts proposed during Association negotiation. This is documented in Section 5.1.1.2.

Configuration options are provided to control the Photometric Interpretation of grayscale images stored to a DICOM Store SCP. Grayscale images can be stored as RGB, YBR-FULL-422 or MONOCHROME2.

Many Acuson X500 B-Mode and M-Mode images contain no significant color, the only color is in the Acuson 'a' transducer position marker and the ECG trace, if used. When "Store grayscale images as Monochrome" is selected, single-frame images with no significant color content will be stored as MONOCHROME2. When "No Monochrome conversion" is selected, all multi-frame and single frame images are stored as RGB or YBR-FULL-422. When "Store grayscale images and clips as Monochrome" is selected, B-Mode and M-Mode single- and multi-frame images with no significant color content will be stored as MONOCHROME2. When "Store all images as Monochrome" is selected then all single frame images (with or without significant color) are stored as Monochrome2. When "Store all images and clips as Monochrome" is selected then all single and multi frame images (with or without any significant color) are stored as Monochrome2.

8.2.2.2 DICOM Storage Commitment Configuration

Configuration of DICOM Storage Commitment remote devices must be performed separately from DICOM Store Configuration. The Acuson X500 supports Storage Commitment to the same remote device as Store or to a different device.

8.2.2.3 DICOM Modality Worklist Configuration

Configuration of DICOM Modality Worklist remote devices

8.2.2.4 DICOM Modality Performed Procedure Step Configuration

Configuration of DICOM Modality Performed Procedure Step remote devices must be performed separately from DICOM Modality Worklist Configuration. The Acuson X500 supports MPPS to the same remote device as Modality Worklist or to a different device.

The “Store Image Format” setting controls the Referenced SOP Class UID (0008,1150) in the Referenced Image Sequence (0008,1140) of the MPPS N-SET sent by the Acuson X500. Due to the Acuson X500’s ability to select from multiple Presentation Contexts during Association Negotiation, it is necessary to use this setting.

In the majority of installations the “Store Image Format” should be left at the default setting of “New Ultrasound”. There are two cases when the “Store Image Format” must be set to “Old Ultrasound” or “Secondary Capture”:

1. When the active Storage Server “Image Format” is set to “Old Ultrasound” or “Secondary Capture”.
2. When the active Storage Server “Image Format” is set to “Automatic”, but the Storage Server does not support US Image and US Multi-frame Image.

In both cases the correct setting can be determined by reviewing the DICOM Conformance Statement of the Storage Server and following the instructions below. DICOM Conformance Statements are usually available on the manufacturer’s Web site.

If at least one of US Image and US Multi-frame Image are listed in the DICOM Conformance Statement and the active Storage Server “Image Format” is set to “Automatic” then “New Ultrasound” is the correct setting for “Store Image Format”.

If the above is not true and at least one of US Image (Retired) and US Multi-frame Image (Retired) are listed in the DICOM Conformance Statement then “Old Ultrasound” is the correct setting for “Store Image Format”.

If neither of the above are true then “Secondary Capture” is the correct setting for “Store Image Format”.

8.2.2.5 DICOM Print Configuration

For each DICOM Print server, the following data is configurable by the user using the System Presets DICOM Print User Interface. The effect of changing parameters of the DICOM Print server will be seen at the next created film sheet. The current film sheet is not affected by changing these parameters.

Table 41 User-Configurable Printer Parameters.





| Parameter | Description |
|------------------|--|
| Printer Type: | Color or Black and White - depends on printer |
| Film Size | Select the size of the film - 8x10 inches, 8.5x11 inches, 10x12 inches, 10x14 inches, 11x14 inches, 11x17 inches, 14x14 inches, 14x17 inches, 24x24 centimeters, 24x30 centimeters, A3, or A4. |
| Film Orientation | Select from Portrait: <div style="text-align: center; margin: 10px 0;">  </div> or Landscape: <div style="text-align: center; margin: 10px 0;">  </div> |

Table 41 User-Configurable Printer Parameters. (Continued)

| Parameter | Description |
|----------------------------|--|
| Display Format | <p>You must supply the number of rows and columns of images on the printed sheet.</p> <p>For example, a 6 on 1 print with Landscape mode should have 3 columns and 2 rows:</p>  <p>A 6 on 1 with Portrait mode would have 2 columns and 3 rows:</p>  |
| Print Priority | HIGH, MEDIUM or LOW |
| Medium Type | PAPER, CLEAR FILM, BLUE FILM, TRANSPARENCY or CURRENT (to use the currently loaded media) |
| Film Destination | MAGAZINE, PROCESSOR or CURRENT |
| Max. Density | Used to define the Black value - printer specific |
| Min. Density | Used to define the White value - printer specific |
| Smoothing Type | Printer specific value |
| Border Density | BLACK or WHITE |
| Empty Image Density | BLACK or WHITE |
| Trim | YES/NO to having a border around each image |
| Polarity | Normal/reverse. Normal means black is printed as black. Reverse means the grayscale is inverted so that black comes out as white and white as black. |
| Magnification | Replicate, Bilinear, Cubic, None |
| Configuration Information: | Printer Specific values |

8.3 External Equipment Configuration

The Acuson X500 user can configure “Hard Key” to “Output Device” mapping through the System Presets - Customize Keys. Print images are acquired and sent to the assigned device when the user presses the associated key. The following key assignments are supported:

- **Clip Store** – This key can be assigned to Multi-frame Store Capture, Cine Store or Disk Store
- **Print/Store 1** – This key can be assigned to any configured DICOM Printer, DICOM Store or OEM printer device.
- **Print/Store 2** – This key can be assigned to any configured DICOM Printer, DICOM Store or OEM printer device.

8.4 Support of Extended Character Sets

The “ISO-IR 100”, “GB18030”, “ISO_IR 144” character sets are supported by the Acuson X500 system based on the following language selections:

English, French, Italian, German, Spanish: “ISO_IR 100”

Chinese: “GB18030”

Russian: “ISO_IR 144”

9.0 Security

9.1 Security Profiles

None supported.

9.2 Association Level Security

None supported.

9.3 Application Level Security

None supported.

9.4 Virus Protection

The Acuson X500 computer system's networking has been configured to significantly reduce the possibility of virus and hacking vulnerabilities. On the X500 computer system, all ingress TCP and UDP ports are closed and/or absent of any type of server. The only exception to this is due to the necessity of a DICOM server available at ingress TCP port 104. Additionally, all non-essential computer services and components are disabled to minimize X500 egress network footprint.

Outside of some minimal network exchanges required by the X500's commercial computer operating system, the only network connections initiated by the X500 are for DICOM connectivity and network-share export function.

10.0 Appendix A: Mapping for Obstetric DICOM SR

The following table lists the DICOM Structured Report mappings used by the Acuson X500 1.0 in Obstetric Structured Reports. All private concept names use the Coding Scheme Designator "G60S".

Table 42 DICOM SR Mappings for OB Structured Reports

| Platform Name | Base Measurement Concept Name | Template Section | Inferred From | Modifiers/Properties |
|---------------------------|-----------------------------------|---|---------------|----------------------|
| PATIENT_NAME | | | | |
| PATIENT_ID | | | | |
| SEX | | | | |
| REFERRING_MD | | | | |
| BIRTH_DATE | | | | |
| PATIENT_AGE | | | | |
| HEIGHT | LN 8302-2, "Patient Height" | TID 5001: Patient Characteristics | | |
| WEIGHT | LN 29463-7, "Patient Weight" | TID 5001: Patient Characteristics | | |
| LMP | LN 11955-2, "LMP" | TID 5002: OB-GYN Procedure Summary CID 12003: OB-GYN Dates | | |
| DIAGNOSIS_DESCR IPTION | DCM, 121106, "Comment" | TID 5001: Patient Characteristics | | |
| GRAVIDA | LN 11996-6, "Gravida" | TID 5001: Patient Characteristics | | |
| PARA | LN 11977-6, "Para" | TID 5001: Patient Characteristics | | |
| AB | LN 11612-9, "Aborta" | TID 5001: Patient Characteristics | | |
| ECTOPIC | LN 33065-4, "Ectopic Pregnancies" | TID 5001: Patient Characteristics | | |
| FETAL_AGE | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | | |

| | | | | |
|------------------------|--------------------------------------|--|---|--|
| LMP_ESTIMATE_FETAL_AGE | LN 11885-1, "Gestational Age by LMP" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | | |
| IVF_DATE_TIME | No match in DICOM Spec | TID 5001: Patient Characteristics | Private Code: IVF_Date_Time_G60S | |
| IVF_ESTIMATE_FETAL_AGE | No match in DICOM Spec | TID 5001: Patient Characteristics | Private Code: IVF_Est_Fetal_Age_G60S | |
| BPD | LN 11820-8, "Biparietal Diameter" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_BPD_HADLOCK | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11902-4, BPD, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | |
| MA_BPD_HADLOCK_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11902-4, BPD, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_HADLOCK_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11902-4, BPD, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_MERZ | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Merz_G60S (Equation) | |
| MA_BPD_MERZ_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Merz_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|--------------------|-------------------------------|--|---|--|
| MA_BPD_MERZ_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Merz_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_LASSER | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Lasser_G60S (Equation) | |
| MA_BPD_LASSER_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Lasser_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_LASSER_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Lasser_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_REMPEN | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33083-7, BPD, Rempen 1991 (CID 12013: Gestational Age Equations and Tables) | |
| MA_BPD_REMPEN_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33083-7, BPD, Rempen 1991 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_REMPEN_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33083-7, BPD, Rempen 1991 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_ASUM | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33079-5, BPD, ASUM 1989 (CID 12013: Gestational Age Equations and Tables) | |
| MA_BPD_TOKYO | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33085-2, BPD, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | |

| | | | | |
|-------------------|--|--|--|---|
| MA_BPD_TOKYO_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33085-2, BPD, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_TOKYO_M1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33085-2, BPD, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_OSAKA | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33082-9, BPD, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | |
| MA_BPD_OSAKA_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33082-9, BPD, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_OSAKA_M1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33082-9, BPD, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_BPD_JSUM | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_JSUM_G60S (Equation) | |
| MA_BPD_JSUM_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_JSUM_G60S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| OFD | LN 11851-3, "Occipital-Frontal Diameter" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_OFD_MERZ | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_OFD_Merz_G60S (Equation) | |

| | | | | |
|--------------------|----------------------------------|--|---|--|
| MA_OFD_ASUM | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33119-9, OFD, ASUM 2000 (CID 12013: Gestational Age Equations and Tables) | |
| HC | LN 11984-2, "Head Circumference" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_HC_HADLOCK | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11932-1, HC, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | |
| MA_HC_HADLOCK_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11932-1, HC, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_HC_HADLOCK_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11932-1, HC, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_HC_MERZ | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33115-7, HC, Merz 1988 (CID 12013: Gestational Age Equations and Tables) | |
| MA_HC_MERZ_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33115-7, HC, Merz 1988 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_HC_MERZ_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33115-7, HC, Merz 1988 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_HC_LASSER | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_HC_Lasser_G60 S (Equation) | |

| | | | | |
|--------------------|---------------------------------------|--|---|--|
| MA_HC_LASSER_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_HC_Lasser_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_HC_LASSER_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_HC_Lasser_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| AC | LN 11979-2, "Abdominal Circumference" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_AC_HADLOCK | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11892-7, AC, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | |
| MA_AC_HADLOCK_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11892-7, AC, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_AC_HADLOCK_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11892-7, AC, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_AC_MERZ | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AC_Merz_G60S (Equation) | |
| MA_AC_MERZ_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AC_Merz_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_AC_MERZ_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AC_Merz_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|-------------------|-------------------------------|--|--|--|
| MA_AC_LASSER | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AC_Lasser_G60S (Equation) | |
| MA_AC_LASSER_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AC_Lasser_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_AC_LASSER_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AC_Lasser_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_AC_JSUM | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AC_JSUM_G60S (Equation) | |
| MA_AC_JSUM_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AC_JSUM_G60S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_AC_JSUM_M1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AC_JSUM_G60S (Equation) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| FL | LN 11963-6, "Femur Length" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_FL_JEANTY | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11923-0, FL, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | |
| MA_FL_JEANTY_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11923-0, FL, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|--------------------|-------------------------------|--|---|--|
| MA_FL_JEANTY_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11923-0, FL, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_FL_HADLOCK | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11920-6, FL, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | |
| MA_FL_HADLOCK_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11920-6, FL, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_FL_HADLOCK_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11920-6, FL, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_FL_MERZ | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33542-2, FL, Merz 1988 (CID 12013: Gestational Age Equations and Tables) | |
| MA_FL_MERZ_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33542-2, FL, Merz 1988 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_FL_MERZ_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33542-2, FL, Merz 1988 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_FL_TOKYO | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33103-3, FL, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | |
| MA_FL_TOKYO_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33103-3, FL, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|------------------|---------------------------------|---|--|--|
| MA_FL_TOKYO_M1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33103-3, FL, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_FL_OSAKA | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33101-7, FL, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | |
| MA_FL_OSAKA_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33101-7, FL, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_FL_OSAKA_M1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33101-7, FL, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_FL_JSUM | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_FL_JSUM_G60S (Equation) | |
| MA_FL_JSUM_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_FL_JSUM_G60S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_FL_JSUM_M1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_FL_JSUM_G60S (Equation) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| CRL | LN 11957-8, "Crown Rump Length" | TID 5011: Early Gestation CID 12009: Early Gestation Biometry Measurements | | |
| MA_CRL_HADLOCK | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11910-7, CRL, Hadlock 1992 (CID 12013: Gestational Age Equations and Tables) | |

| | | | | |
|----------------------|-------------------------------|--|---|--|
| MA_CRL_HADLOCK_P2SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11910-7, CRL, Hadlock 1992 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_HADLOCK_M2SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11910-7, CRL, Hadlock 1992 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_ROBINSON | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11914-9, CRL, Robinson 1975 (CID 12013: Gestational Age Equations and Tables) | |
| MA_CRL_ROBINSON_P2SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11914-9, CRL, Robinson 1975 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_ROBINSON_M2SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11914-9, CRL, Robinson 1975 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_HANSMAN | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33540-6, CRL, Hansmann 1986 (CID 12013: Gestational Age Equations and Tables) | |
| MA_CRL_HANSMAN_P2SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33540-6, CRL, Hansmann 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_HANSMAN_M2SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33540-6, CRL, Hansmann 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_LASSER | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | Private Code: MA_CRL_Lasser_G60S (Equation) | |

| | | | | |
|--------------------|-------------------------------|--|--|--|
| MA_CRL_LASSER_P2SD | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | Private Code: MA_CRL_Lasser_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_LASSER_M2SD | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | Private Code: MA_CRL_Lasser_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_ASUM | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33090-2, CRL, ASUM 2000 (CID 12013: Gestational Age Equations and Tables) | |
| MA_CRL_TOKYO | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33096-9, CRL, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | |
| MA_CRL_TOKYO_P1SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33096-9, CRL, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_TOKYO_M1SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33096-9, CRL, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_OSAKA | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33093-6, CRL, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | |
| MA_CRL_OSAKA_P1SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33093-6, CRL, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_OSAKA_M1SD | LN 18185-9, "Gestational Age" | TID 50011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 33093-6, CRL, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|-------------------|---------------------------------|---|---|--|
| MA_CRL_JSUM | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | Private Code: MA_CRL_JSUM_G60 S (Equation) | |
| MA_CRL_JSUM_P1SD | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | Private Code: MA_CRL_JSUM_G60 S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CRL_JSUM_M1SD | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | Private Code: MA_CRL_JSUM_G60 S (Equation) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| HL | LN 11966-9, "Humerus Length" | TID 5006: Long Bones Section CID 12006: Fetal Long Bones Biometry Measurements | | |
| MA_HL_JEANTY | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 11936-2, Humerus, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | |
| MA_HL_JEANTY_P2SD | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 11936-2, Humerus, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_HL_JEANTY_M2SD | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 11936-2, Humerus, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_HL_OSAKA | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 33117-3, Humerus Length, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | |

| | | | | |
|-------------------|-------------------------------|---|--|--|
| MA_HL_OSAKA_P1SD | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 33117-3, Humerus Length, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_HL_OSAKA_M1SD | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 33117-3, Humerus Length, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| UL | LN 11969-3, "Ulna Length" | TID 5006: Long Bones Section CID 12006: Fetal Long Bones Biometry Measurements | | |
| MA_UL_JEANTY | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 11944-6, Ulna, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | |
| MA_UL_JEANTY_P2SD | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 11944-6, Ulna, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_UL_JEANTY_M2SD | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 11944-6, Ulna, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| TL | LN 11968-5, "Tibia Length" | TID 5006: Long Bones Section CID 12006: Fetal Long Bones Biometry Measurements | | |
| MA_TL_JEANTY | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 11941-2, Tibia, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | |
| MA_TL_JEANTY_P2SD | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 11941-2, Tibia, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|--------------------|--|--|---|--|
| MA_TL_JEANTY_M2SD | LN 18185-9, "Gestational Age" | TID 5006: Long Bones Section TID 5008: Fetal Biometry Group | LN 11941-2, Tibia, Jeanty 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| FT | LN 11965-1, "Foot Length" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_FT_MERCER | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11926-3, Foot Length, Mercer 1987 (CID 12013: Gestational Age Equations and Tables) | |
| MA_FT_MERCER_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11926-3, Foot Length, Mercer 1987 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_FT_MERCER_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11926-3, Foot Length, Mercer 1987 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MSD | LN 11850-5, "Gestational Sac Diameter" | TID 5011: Early Gestation CID 12009: Early Gestation Biometry Measurements | | |
| MA_MSD_REMPEN | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11929-7, GS, Rempen 1991 (CID 12013: Gestational Age Equations and Tables) | |
| MA_MSD_REMPEN_P2SD | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11929-7, GS, Rempen 1991 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|---------------------|--------------------------------------|---|---|---|
| MA_MSD_REMPEN_M2SD | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11929-7, GS, Rempen 1991 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_MSD_HELLMAN | LN 18185-9, "Gestational Age" | TID 5011: Early Gestation Section TID 5008: Fetal Biometry Group | LN 11928-9, GS, Hellman 1969 (CID 12013: Gestational Age Equations and Tables) | |
| BN | LN 11629-3, "Outer Orbital Diameter" | TID 5007: Fetal Cranium Section CID 12007: Fetal Cranium | | |
| MA_BN_JEANTY | LN 18185-9, "Gestational Age" | TID 5007: Fetal Cranium Section TID 5008: Fetal Biometry Group | Private Code: MA_BN_Jeanty_G60S (Equation) | |
| MA_BN_JEANTY_P2SD | LN 18185-9, "Gestational Age" | TID 5007: Fetal Cranium Section TID 5008: Fetal Biometry Group | Private Code: MA_BN_Jeanty_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_BN_JEANTY_M2SD | LN 18185-9, "Gestational Age" | TID 5007: Fetal Cranium Section TID 5008: Fetal Biometry Group | Private Code: MA_BN_Jeanty_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_BN_TONGSONG | LN 18185-9, "Gestational Age" | TID 5007: Fetal Cranium Section TID 5008: Fetal Biometry Group | Private Code: MA_BN_Tongsong_G60S (Equation) | |
| MA_BN_TONGSONG_P2SD | LN 18185-9, "Gestational Age" | TID 5007: Fetal Cranium Section TID 5008: Fetal Biometry Group | Private Code: MA_BN_Tongsong_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_BN_TONGSONG_M2SD | LN 18185-9, "Gestational Age" | TID 5007: Fetal Cranium Section TID 5008: Fetal Biometry Group | Private Code: MA_BN_Tongsong_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|------------------|---|--|---|--|
| ASD | LN 11818-2, "Anterior-Posterior Abdominal Diameter" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_ASD_MERZ | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_ASD_Merz_G60S (Equation) | |
| ATD | LN 11862-0, "Transverse Abdominal Diameter" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_ATD_MERZ | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_ATD_Merz_G60S (Equation) | |
| GS | LN 11850-5, "Gestational Sac Diameter" | TID 5011: Early Gestational Section TID 5008: Fetal Biometry Group | | |
| MA_GS_TOKYO | LN 18185-9, "Gestational Age" | TID 5011: Early Gestational Section TID 5008: Fetal Biometry Group | LN 33108-2, GS, Tokyo 1986 | |
| MA_GS_TOKYO_P1SD | LN 18185-9, "Gestational Age" | TID 5011: Early Gestational Section TID 5008: Fetal Biometry Group | LN 33108-2, GS, Tokyo 1986 | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_GS_TOKYO_M1SD | LN 18185-9, "Gestational Age" | TID 5011: Early Gestational Section TID 5008: Fetal Biometry Group | LN 33108-2, GS, Tokyo 1986 | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|----------------|--------------------------------------|--|--|--|
| FTA | LN 33068-8, "Thoracic Area" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_FTA_OSAKA | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_FTA_Osaka_G60S (Equation) | |
| MA_FTA_OSAKA_W | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_FTA_Osaka_W_G60S (Equation) | |
| CL | LN 11962-8, "Clavicle Length" | TID 5006: Long Bones Section CID 12006: Fetal Long Bones Biometry Measurements | | |
| TC | LN 11988-3, "Thoracic Circumference" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| RIGHT_RL | LN 11836-4, "Right Kidney Length" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| LEFT_RL | LN 11834-9, "Left Kidney Length" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |

| | | | | |
|-----------------------------|---|--|--|--|
| RIGHT_RAP | LN 11855-4, "Right Kidney Thickness" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| LEFT_RAP | LN 11853-9, "Left Kidney Thickness" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| AFI | LN 11627-7, "Amniotic Fluid Index" | TID 5010: Amniotic Sac Section | | |
| YOLK_SAC | LN 11816-6, "Yolk Sac Length" | TID 5011: Early Gestation CID 12009: Early Gestation Biometry Measurements | | |
| HW | LN 12170-7, "Width of Hemisphere" | TID 5007: Fetal Cranium Section CID 12007: Fetal Cranium | | |
| TCD | LN 11863-8, "Trans Cerebellar Diameter" | TID 5007: Fetal Cranium Section CID 12007: Fetal Cranium | | |
| LVW | LN 12171-5, "Lateral Ventricular Width" | TID 5007: Fetal Cranium Section CID 12007: Fetal Cranium | | |
| CIST_MAGNA | LN 11860-4, "Cisterna Magna Length" | TID 5007: Fetal Cranium Section CID 12007: Fetal Cranium | | |
| CERVIX_LEN (used for OB) | LN 11961-0, "Cervix Length" | TID 5015: Pelvis and Uterus Section CID 12011: Ultrasound Pelvis and Uterus | | |

| | | | | |
|----------------|---|---|--|--|
| NT | LN 12146-7, "Nuchal Fold Thickness" | TID 5007: Fetal Cranium Section CID 12007: Fetal Cranium | | |
| UMB_VD | SRT-G-0364, "Vessel Lumen Diameter" | TID 5026: Pelvic Vascular Ultrasound Measurement Group CID 12119: Vascular Ultrasound Property CID 12122: Other Vascular Properties | | |
| FHR | LN 11948-7, "Fetal Heart Rate" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | | |
| APTD | LN 11819-0, "Anterior-Posterior Trunk Diameter" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| TTD | LN 11864-6, "Transverse Thoracic Diameter" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| UMB_A_SYSTOLE | LN 11726-7, "Peak Systolic Velocity" | TID 5026: OB-GYN Pelvic Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12140, Pelvic Vasculature Anatomical Location): SRT T-F1810, "Umbilical Artery" | |
| UMB_A_DIASTOLE | LN 11653-3, End Diastolic Velocity" | TID 5026: OB-GYN Pelvic Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12140, Pelvic Vasculature Anatomical Location): SRT T-F1810, "Umbilical Artery" | |

| | | | | |
|--------------|--|---|---|--|
| UMB_A_S_D | LN 12144-2, "Systolic to Diastolic Velocity Ratio" | TID 5026: OB-GYN Pelvic Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12140, Pelvic Vasculature Anatomical Location): SRT T-F1810, "Umbilical Artery" | |
| UMB_A_RI | LN 12023-8, "Resistivity Index" | TID 5026: OB-GYN Pelvic Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12140, Pelvic Vasculature Anatomical Location): SRT T-F1810, "Umbilical Artery" | |
| UMB_A_PI | LN 12008-9, "Pulsatility Index" | TID 5026: OB-GYN Pelvic Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12140, Pelvic Vasculature Anatomical Location): SRT T-F1810, "Umbilical Artery" | |
| MCA_SYSTOLE | LN 11726-7, "Peak Systolic Velocity" | TID 5025: OB-GYN Fetal Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12141, Fetal Vasculature Anatomical Location): SRT T-45600, "Middle Cerebral Artery" | |
| MCA_DIASTOLE | LN 11653-3, End Diastolic Velocity" | TID 5025: OB-GYN Fetal Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12141, Fetal Vasculature Anatomical Location): SRT T-45600, "Middle Cerebral Artery" | |
| MCA_S_D | LN 12144-2, "Systolic to Diastolic Velocity Ratio" | TID 5025: OB-GYN Fetal Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12141, Fetal Vasculature Anatomical Location): SRT T-45600, "Middle Cerebral Artery" | |

| | | | | |
|--------------------|---------------------------------|--|---|--|
| MCA_RI | LN 12023-8, "Resistivity Index" | TID 5025: OB-GYN Fetal Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12141, Fetal Vasculature Anatomical Location): SRT T-45600, "Middle Cerebral Artery" | |
| MCA_PI | LN 12008-9, "Pulsatility Index" | TID 5025: OB-GYN Fetal Vascular Measurement Group CID 12119: Vascular Ultrasound Property CID 12120: Blood Velocity Measurements | Anatomy Group (CID 12141, Fetal Vasculature Anatomical Location): SRT T-45600, "Middle Cerebral Artery" | |
| EFW1_HADLOCK1 | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock1_G6 0S (Equation) | |
| EFW1_HADLOCK1_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock1_G6 0S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_HADLOCK1_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock1_G6 0S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_HADLOCK2 | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock2_G6 0S (Equation) | |
| EFW1_HADLOCK2_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock2_G6 0S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_HADLOCK2_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock2_G6 0S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|--------------------|--------------------------------|--|---|--|
| EFW1_HADLOCK3 | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock3_G6 0S (Equation) | |
| EFW1_HADLOCK3_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock3_G6 0S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_HADLOCK3_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock3_G6 0S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_HADLOCK4 | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock4_G6 0S (Equation) | |
| EFW1_HADLOCK4_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock4_G6 0S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_HADLOCK4_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Hadlock4_G6 0S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_SCMACHER | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Scmacher_G6 0S (Equation) | |
| EFW1_SCMACHER_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Scmacher_G6 0S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_SCMACHER_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Scmacher_G6 0S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|--------------------|--------------------------------|---|---|---|
| EFW1_HANSMANN | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33139-7, EFW by BPD, TTD, Hansmann 1986 (CID 12014: OB Fetal Body Weight Equations and Tables) | |
| EFW1_HANSMANN_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33139-7, EFW by BPD, TTD, Hansmann 1986 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_HANSMANN_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33139-7, EFW by BPD, TTD, Hansmann 1986 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_MERZ | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Merz_G60S (Equation) | |
| EFW1_MERZ_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Merz_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_MERZ_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_Merz_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_SHEPARD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 11739-0, EFW by AC and BPD, Shepard 1982 (CID 12014: OB Fetal Body Weight Equations and Tables) | |

| | | | | |
|--------------------|--------------------------------|---|---|---|
| EFW1_SHEPARD_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 11739-0, EFW by AC and BPD, Shepard 1982 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_SHEPARD_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 11739-0, EFW by AC and BPD, Shepard 1982 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_TOKYO | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33144-7, EFW by BPD, APAD, TAD, FL, Tokyo 1987 (CID 12014: OB Fetal Body Weight Equations and Tables) | |
| EFW1_TOKYO_P1SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33144-7, EFW by BPD, APAD, TAD, FL, Tokyo 1987 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_TOKYO_M1SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33144-7, EFW by BPD, APAD, TAD, FL, Tokyo 1987 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_EFW1_TOKYO | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW1_Tokyo_G60S (Equation) | |
| MA_EFW1_TOKYO_P1SD | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW1_Tokyo_G60S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|--------------------|--------------------------------|---|---|--|
| MA_EFW1_TOKYO_M1SD | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW1_Tokyo_G60S (Equation) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_OSAKA | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33140-5, EFW by BPD, FTA, FL, Osaka 1990 (CID 12014: OB Fetal Body Weight Equations and Tables) | |
| EFW1_OSAKA_P1SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33140-5, EFW by BPD, FTA, FL, Osaka 1990 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_OSAKA_M1SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33140-5, EFW by BPD, FTA, FL, Osaka 1990 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_EFW1_OSAKA | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW1_Osaka_G60S (Equation) | |
| MA_EFW1_OSAKA_P1SD | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW1_Osaka_G60S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_EFW1_OSAKA_M1SD | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW1_Osaka_G60S (Equation) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW1_JSUM | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW1_JSUM_G60S (Equation) | |

| | | | | |
|--------------------|--------------------------------|---|---|---|
| MA_EFW1_JSUM | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW1_JSUM_G60S (Equation) | |
| EFW2_HADLOCK1 | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock1_G60S (Equation) | |
| EFW2_HADLOCK1_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock1_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_HADLOCK1_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock1_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_HADLOCK2 | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock2_G60S (Equation) | |
| EFW2_HADLOCK2_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock2_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_HADLOCK2_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock2_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_HADLOCK3 | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock3_G60S (Equation) | |
| EFW2_HADLOCK3_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock3_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|--------------------|--------------------------------|---|--|--|
| EFW2_HADLOCK3_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock3_G6 0S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_HADLOCK4 | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock4_G6 0S (Equation) | |
| EFW2_HADLOCK4_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock4_G6 0S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_HADLOCK4_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Hadlock4_G6 0S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_SCMACHER | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Scmacher_G6 0S (Equation) | |
| EFW2_SCMACHER_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Scmacher_G6 0S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_SCMACHER_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Scmacher_G6 0S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_HANSMANN | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33139-7, EFW by BPD, TTD, Hansmann 1986 (CID 12014: OB Fetal Body Weight Equations and Tables) | |

| | | | | |
|--------------------|--------------------------------|---|---|---|
| EFW2_HANSMANN_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33139-7, EFW by BPD, TTD, Hansmann 1986 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_HANSMANN_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33139-7, EFW by BPD, TTD, Hansmann 1986 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_MERZ | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Merz_G60S (Equation) | |
| EFW2_MERZ_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Merz_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_MERZ_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_Merz_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_SHEPARD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 11739-0, EFW by AC and BPD, Shepard 1982 (CID 12014: OB Fetal Body Weight Equations and Tables) | |
| EFW2_SHEPARD_P2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 11739-0, EFW by AC and BPD, Shepard 1982 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|--------------------|--------------------------------|---|---|---|
| EFW2_SHEPARD_M2SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 11739-0, EFW by AC and BPD, Shepard 1982 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_TOKYO | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33144-7, EFW by BPD, APAD, TAD, FL, Tokyo 1987 (CID 12014: OB Fetal Body Weight Equations and Tables) | |
| EFW2_TOKYO_P1SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33144-7, EFW by BPD, APAD, TAD, FL, Tokyo 1987 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_TOKYO_M1SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33144-7, EFW by BPD, APAD, TAD, FL, Tokyo 1987 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_EFW2_TOKYO | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW2_Tokyo_G60S (Equation) | |
| MA_EFW2_TOKYO_P1SD | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW2_Tokyo_G60S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_EFW2_TOKYO_M1SD | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW2_Tokyo_G60S (Equation) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |

| | | | | |
|--------------------|--------------------------------|---|---|---|
| EFW2_OSAKA | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33140-5, EFW by BPD, FTA, FL, Osaka 1990 (CID 12014: OB Fetal Body Weight Equations and Tables) | |
| EFW2_OSAKA_P1SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33140-5, EFW by BPD, FTA, FL, Osaka 1990 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_OSAKA_M1SD | LN 11727-5, "Estimated Weight" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary CID 12012: OB Equations and Tables | LN 33140-5, EFW by BPD, FTA, FL, Osaka 1990 (CID 12014: OB Fetal Body Weight Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_EFW2_OSAKA | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW2_Osaka_G60S (Equation) | |
| MA_EFW2_OSAKA_P1SD | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW2_Osaka_G60S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_EFW2_OSAKA_M1SD | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW2_Osaka_G60S (Equation) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| EFW2_JSUM | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: EFW2_JSUM_G60S (Equation) | |
| MA_EFW2_JSUM | LN 18185-9, "Gestational Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: MA_EFW2_JSUM_G60S (Equation) | |

| | | | | |
|--------------------------------|---------------------------|--|--|---|
| RATIO_HC_AC_CAM PBELL | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_HC_AC_Cam pbell_G60S (Equation) | |
| RATIO_HC_AC_CAM PBELL_UPPER | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_HC_AC_Cam pbell_G60S (Equation) | SRT R-0038B, Nor- mal Range Upper Limit |
| RATIO_HC_AC_CAM PBELL_LOWER | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_HC_AC_Cam pbell_G60S (Equation) | SRT R-10041, Normal Range Lower Limit |
| RATIO_FL_AC_HADL OCK | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_FL_AC_Hadlo ck_G60S (Equation) | |
| RATIO_FL_AC_HADL OCK_UPPER | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_FL_AC_Hadlo ck_G60S (Equation) | SRT R-0038B, Nor- mal Range Upper Limit |
| RATIO_FL_AC_HADL OCK_LOWER | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_FL_AC_Hadlo ck_G60S (Equation) | SRT R-10041, Normal Range Lower Limit |
| RATIO_FL_BPD_HO HLER | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_FL_BPD_Hoh ler_G60S (Equation) | |
| RATIO_FL_BPD_HO HLER_UPPER | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_FL_BPD_Hoh ler_G60S (Equation) | SRT R-0038B, Nor- mal Range Upper Limit |
| RATIO_FL_BPD_HO HLER_LOWER | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_FL_BPD_Hoh ler_G60S (Equation) | SRT R-10041, Normal Range Lower Limit |
| RATIO_TCD_AC_ME YER | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_TCD_AC_Me yer_G60S (Equation) | |
| RATIO_TCD_AC_ME YER_UPPER | No match in DICOM Spec | TID 5004: Fetal Biom- etry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_TCD_AC_Me yer_G60S (Equation) | SRT R-0038B, Nor- mal Range Upper Limit |

| | | | | |
|----------------------------|------------------------|--|--|---------------------------------------|
| RATIO_TCD_AC_MEYER_LOWER | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_TCD_AC_Meyer_G60S (Equation) | SRT R-10041, Normal Range Lower Limit |
| RATIO_LVW_HW_JOHNSON | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_LVW_HW_Johnson_G60S (Equation) | |
| RATIO_LVW_HW_JOHNSON_UPPER | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_LVW_HW_Johnson_G60S (Equation) | SRT R-0038B, Normal Range Upper Limit |
| RATIO_LVW_HW_JOHNSON_LOWER | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_LVW_HW_Johnson_G60S (Equation) | SRT R-10041, Normal Range Lower Limit |
| RATIO_CI_HADLOCK | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_CI_Hadlock_G60S (Equation) | |
| RATIO_CI_HADLOCK_UPPER | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_CI_Hadlock_G60S (Equation) | SRT R-0038B, Normal Range Upper Limit |
| RATIO_CI_HADLOCK_LOWER | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_CI_Hadlock_G60S (Equation) | SRT R-10041, Normal Range Lower Limit |
| RATIO_CI_CHITTY | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_CI_CHITTY_G60S (Equation) | |
| RATIO_CI_CHITTY_UPPER | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_CI_CHITTY_G60S (Equation) | SRT R-0038B, Normal Range Upper Limit |
| RATIO_CI_CHITTY_LOWER | No match in DICOM Spec | TID 5004: Fetal Biometry Ratio Section CID 12004: Fetal Biometry Ratios | Private Code: RATIO_CI_CHITTY_G60S (Equation) | SRT R-10041, Normal Range Lower Limit |

| | | | | |
|-------------------|--|--|--|--|
| AXT | LN 33191-8, "APAD * TAD" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |
| MA_AXT_TOKYO | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AXT_Tokyo_G60S (Equation) | |
| MA_AXT_TOKYO_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AXT_Tokyo_G60S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_AXT_TOKYO_M1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_AXT_Tokyo_G60S (Equation) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| USMA_AVARAGE | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | | |
| USMA_HADLOCK1 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock1_G60S (Equation) | |
| USMA_HADLOCK2 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock2_G60S (Equation) | |
| USMA_HADLOCK3 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock3_G60S (Equation) | |
| USMA_HADLOCK4 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock4_G60S (Equation) | |

| | | | | |
|----------------|--|--|---|--|
| USMA_HADLOCK5 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock5_G60S (Equation) | |
| USMA_HADLOCK6 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock6_G60S (Equation) | |
| USMA_HADLOCK7 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock7_G60S (Equation) | |
| USMA_HADLOCK8 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock8_G60S (Equation) | |
| USMA_HADLOCK9 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock9_G60S (Equation) | |
| USMA_HADLOCK10 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock10_G60S (Equation) | |
| USMA_HADLOCK11 | LN 11888-5, "Composite Ultrasound Age" | TID 5003: OB-GYN Procedure Fetus Summary CID 12019: OB-GYN Fetus Summary | Private Code: US_MA_Hadlock11_G60S (Equation) | |
| CORBPD | LN 11824-0, "BPD area corrected" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | | |

| | | | | |
|------------------------|---|--|--|---|
| CORBPD_DOUBLIET | LN 11824-0, "BPD area corrected" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group CID 12005: Fetal Biometry Measurements | Private Code: CORBPD_Doubliet_G60S (Equation) | |
| BIOP_MOVEMENT | LN 11631-0, "Gross Body Movement" | TID 5009: Fetal Biophysical Profile Section | | |
| BIOP_BREATHING | LN 11632-7, "Fetal Breathing" | TID 5009: Fetal Biophysical Profile Section | | |
| BIOP_TONE | LN 11635-0, "Fetal Tone" | TID 5009: Fetal Biophysical Profile Section | | |
| BIOP_FHR | LN 11635-5, "Fetal Heart Reactivity" | TID 5009: Fetal Biophysical Profile Section | | |
| BIOP_AFV | LN 11630-1, "Amniotic Fluid Volume" | TID 5009: Fetal Biophysical Profile Section | | |
| BIOP_TOTAL | LN 11634-3, "Biophysical Profile Sum Score" | TID 5009: Fetal Biophysical Profile Section | | |
| MA_CORBPD_HADLOCK | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11902-4, BPD, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | |
| MA_CORBPD_HADLOCK_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11902-4, BPD, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_HADLOCK_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 11902-4, BPD, Hadlock 1984 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_MERZ | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Merz_G60S (Equation) | |

| | | | | |
|-----------------------|-------------------------------|--|---|--|
| MA_CORBPD_MERZ_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Merz_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_MERZ_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Merz_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_LASSER | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Lasser_G60S (Equation) | |
| MA_CORBPD_LASSER_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Lasser_G60S (Equation) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_LASSER_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_Lasser_G60S (Equation) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_REMPEN | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33083-7, BPD, Rempen 1991 (CID 12013: Gestational Age Equations and Tables) | |
| MA_CORBPD_REMPEN_P2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33083-7, BPD, Rempen 1991 (CID 12013: Gestational Age Equations and Tables) | SRT R-00387, 2 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_REMPEN_M2SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33083-7, BPD, Rempen 1991 (CID 12013: Gestational Age Equations and Tables) | SRT R-00388, 2 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_ASUM | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33079-5, BPD, ASUM 1989 (CID 12013: Gestational Age Equations and Tables) | |

| | | | | |
|----------------------|-------------------------------|--|--|---|
| MA_CORBPD_TOKYO | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33085-2, BPD, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | |
| MA_CORBPD_TOKYO_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33085-2, BPD, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_TOKYO_M1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33085-2, BPD, Tokyo 1986 (CID 12013: Gestational Age Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_OSAKA | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33082-9, BPD, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | |
| MA_CORBPD_OSAKA_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33082-9, BPD, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_OSAKA_M1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | LN 33082-9, BPD, Osaka 1989 (CID 12013: Gestational Age Equations and Tables) | SRT R-00347, 1 Sigma Lower Value of population (CID 226, Population Statistical Descriptors) |
| MA_CORBPD_JSUM | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_JSUM_G60S (Equation) | |
| MA_CORBPD_JSUM_P1SD | LN 18185-9, "Gestational Age" | TID 5005: Fetal Biometry Section TID 5008: Fetal Biometry Group | Private Code: MA_BPD_JSUM_G60S (Equation) | SRT R-00346, 1 Sigma Upper Value of population (CID 226, Population Statistical Descriptors) |

11.0 Appendix B: Mapping for Cardiac DICOM SR

The following table lists the DICOM Structured Report mappings used by the Acuson X500 1.0 in Adult Echocardiography Structured Reports. All private concept names use the Coding Scheme Designator "99SIEMENSUS".

Table 43 DICOM SR Mappings for Adult Echocardiography Structured Reports

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|---------------|---|--|-----------------------------------|---|
| AV/LA_AO | Aortic Root Diameter | (18015-8, LN, "Aortic Root Diameter") | (T-42000, SRT, "Aorta") | (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole") |
| AV/LA_ACS | Aortic Valve Cusp Separation | (17996-0, LN, "Aortic Valve Cusp Separation") | (T-35400, SRT, "Aortic Valve") | |
| AV/LA_LA_diam | Left Atrial Antero-posterior Systolic Dimension by M-Mode | (29469-4, LN, "Left Atrium Antero-posterior Systolic Dimension") | (T-32300, SRT, "Left Atrium") | |
| AV/LA_LA_diam | Left Atrial Antero-posterior Systolic Dimension by 2-D | (29469-4, LN, "Left Atrium Antero-posterior Systolic Dimension") | (T-35400, SRT, "Aortic Valve") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| AV/LA_AO/LA | | (CA_LA_AO2LA, 99SIEMENSUS, "Aortic Root Diameter to Left Atrium Diameter Ratio") | (T-35400, SRT, "Aortic Valve") | |
| AV/LA_RV_diam | | (CA_RV_INTDIA_AV2LA, 99SIEMENSUS, "Right Ventricular Internal Diastolic Dimension at AV/LA") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| AV/LA_LVET | Aortic Valve | (18041-4, LN, "Aortic Valve Ejection Time") | (T-35400, SRT, "Aortic Valve") | |
| AV/LA_LVPEP | Pre-Ejection Period | (18068-7, LN, "Pre-Ejection Period") | (T-35400, SRT, "Aortic Valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| AV/LA_LVSTI | | (CA_AV_PEP2ET, 99SIEMENSUS, "Aortic Valve Pre-ejection Period to Ejection Time") | (T-35400, SRT, "Aortic Valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| AR_Decel_Rate | Aortic Valve Regurgitant Diastolic Deceleration Slope | (20216-8, LN, "Deceleration Slope") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| AR_AR_Vmax | | (20351-3, LN, "Maximum Velocity") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| AR_AR_Pgmax | | (20247-3, LN, "Peak Gradient") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-------------------|--|---|----------------------------------|---|
| AR_AR_V_ed | | (11653-3, LN, "End Diastolic Velocity") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| AR_AR_Time | | (17998-6, LN, "Aortic Valve Regurgitant Diastolic Deceleration Time") | (T-35400, SRT, "Aortic Valve") | |
| AR_AI_Vmax | | (20351-3, LN, "Maximum Velocity") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| AR_AI_Pgmax | | (20247-3, LN, "Peak Gradient") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| AR_AI_PHT | Aortic Valve Regurgitant Diastolic Pressure Half-time | (20280-4, LN, "Pressure Half-Time") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| AR_AI_Dec_Time | | (20217-6, LN, "Deceleration Time") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| AR_AI_Dec_Slope | | (20216-8, LN, "Deceleration Slope") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| Cubed(2D)_LVIDd | Left Ventricular Internal End Diastolic Dimension by 2-D | (29436-3, LN, "Left Ventricle Internal End Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Cubed(2D)_LVIDs | Left Ventricular Internal Systolic Dimension by 2-D | (29438-9, LN, "Left Ventricle Internal Systolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Cubed(2D)_FS | Left Ventricular Fractional Shortening by 2-D | (18051-3, LN, "Left Ventricular Fractional Shortening") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Cubed(2D)_IVSd | Interventricular Septum Diastolic Thickness by 2-D | (18154-5, LN, "Interventricular Septum Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Cubed(2D)_IVSs | Interventricular Septum Systolic Thickness by 2-D | (18158-6, LN, "Interventricular Septum Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Cubed(2D)_LVPWd | Left Ventricular Posterior Wall Diastolic Thickness by 2-D | (18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-------------------|---|--|-----------------------------------|--|
| Cubed(2D)_LVPWs | Left Ventricle Posterior Wall Systolic Thickness by 2-D | (18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Cubed(2D)_EDV | Left Ventricular End Diastolic Volume | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| Cubed(2D)_ESV | Left Ventricular End Systolic Volume | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| Cubed(2D)_EF | Left Ventricular EF | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| Cubed(2D)_SV | Left Ventricular Stroke Volume | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| Cubed(2D)_CO | Left Ventricular Cardiac Output | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| Cubed(2D)_CI | Left Ventricular Cardiac Index | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| Cubed(2D)_RVDd | Right Ventricular Internal Diastolic Dimension by 2-D | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| Cubed(2D)_RVAWd | Right Ventricle Anterior Wall Diastolic Thickness | (18153-7, LN, "Right Ventricle Anterior Wall Diastolic Thickness") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| Cubed(2D)_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-------------------|---|--|----------------------------------|---|
| Cubed(2D)_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | |
| LV/Cubed(M)_LVIDd | Left Ventricle Internal End Diastolic Dimension | (29436-3, LN "Left Ventricle Internal End Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Cubed(M)_LVIDs | Left Ventricle Internal Systolic Dimension | (29438-9, LN, "Left Ventricle Internal Systolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Cubed(M)_FS | Left Ventricular Fractional Shortening | (18051-3, LN, "Left Ventricular Fractional Shortening") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Cubed(M)_IVSd | Interventricular Septum Diastolic Thickness | (18154-5, LN, "Interventricular Septum Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Cubed(M)_IVSs | Interventricular Septum Systolic Thickness | (18158-6, LN, "Interventricular Septum Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Cubed(M)_LVPWd | Left Ventricle Posterior Wall Diastolic Thickness | (18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Cubed(M)_LVPWs | Left Ventricle Posterior Wall Systolic Thickness | (18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Cubed(M)_EDV | Left Ventricular End Diastolic Volume | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_ESV | Left Ventricular End Systolic Volume | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_EF | Left Ventricular EF | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_SV | Left Ventricular Stroke Volume | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_CO | Left Ventricular Cardiac Output | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-----------------------|--|--|-----------------------------------|---|
| LV/Cubed(M)_CI | Left Ventricular Cardiac Index | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_LV_Mass | Left Ventricular Mass by M-Mode | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_LV_Mass-c | | (CA_LV_MASSc, 99SIEMENSUS, "Left Ventricle Mass corrected") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_LV_Mass-l | | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area"); (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_RVDd | Right Ventricular Internal Diastolic Dimension by M-Mode | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_LVET | | (18041-4, LN, "Aortic Valve Ejection Time") | (T-32600, SRT, "Left Ventricle") | |
| LV/Cubed(M)_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Cubed(M)_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| LV/Cubed(M)_mVcf | | (CA_LF_VCF, 99SIEMENSUS, "systolic index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125206, DCM, "Cube Method") |
| Teichholz(2D)_LVIDd | Left Ventricular Internal End Diastolic Dimension by 2-D | (29436-3, LN, "Left Ventricle Internal End Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Teichholz(2D)_LVIDs | Left Ventricular Internal Systolic Dimension by 2-D | (29438-9, LN, "Left Ventricle Internal Systolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|---------------------|--|--|----------------------------------|---|
| Teichholz(2D)_FS | Left Ventricular Fractional Shortening by 2-D | (18051-3, LN, "Left Ventricular Fractional Shortening") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Teichholz(2D)_IVSd | Interventricular Septum Diastolic Thickness by 2-D | (18154-5, LN, "Interventricular Septum Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Teichholz(2D)_IVSs | Interventricular Septum Systolic Thickness by 2-D | (18158-6, LN, "Interventricular Septum Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Teichholz(2D)_LVPWd | Left Ventricular Posterior Wall Diastolic Thickness by 2-D | (18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Teichholz(2D)_LVPWs | Left Ventricle Posterior Wall Systolic Thickness by 2-D | (18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Teichholz(2D)_EDV | Left Ventricular End Diastolic Volume by Teichholz Method | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz"); (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Teichholz(2D)_ESV | Left Ventricular End Systolic Volume by Teichholz Method | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Teichholz(2D)_SV | Left Ventricular Stroke Volume by Teichholz Method | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Teichholz(2D)_EF | Left Ventricular EF by Teichholz Method | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Teichholz(2D)_CO | Left Ventricular Cardiac Output by Teichholz Method | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-----------------------|---|--|-----------------------------------|--|
| Teichholz(2D)_CI | Left Ventricular Cardiac Index by Teichholz Method | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Teichholz(2D)_RVDd | Right Ventricular Internal Diastolic Dimension by 2-D | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Teichholz(2D)_RVAWd | Right Ventricle Anterior Wall Diastolic Thickness | (18153-7, LN, "Right Ventricle Anterior Wall Diastolic Thickness") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Teichholz(2D)_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Teichholz(2D)_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| LV/Teichholz(M)_LVIDd | Left Ventricle Internal End Diastolic Dimension | (29436-3, LN "Left Ventricle Internal End Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Teichholz(M)_LVIDs | Left Ventricle Internal Systolic Dimension | (29438-9, LN, "Left Ventricle Internal Systolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Teichholz(M)_FS | Left Ventricular Fractional Shortening | (18051-3, LN, "Left Ventricular Fractional Shortening") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Teichholz(M)_IVSd | Interventricular Septum Diastolic Thickness | (18154-5, LN, "Interventricular Septum Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Teichholz(M)_IVSs | Interventricular Septum Systolic Thickness | (18158-6, LN, "Interventricular Septum Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Teichholz(M)_LVPWd | Left Ventricle Posterior Wall Diastolic Thickness | (18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Teichholz(M)_LVPWs | Left Ventricle Posterior Wall Systolic Thickness | (18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Teichholz(M)_EDV | Left Ventricular End Diastolic Volume by Teichholz Method | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|---------------------------|--|---|----------------------------------|---|
| LV/Teichholz(M)_ESV | Left Ventricular End Systolic Volume by Teichholz Method | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| LV/Teichholz(M)_EF | Left Ventricular EF by Teichholz Method | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| LV/Teichholz(M)_CO | Left Ventricular Cardiac Output by Teichholz Method | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| LV/Teichholz(M)_CI | Left Ventricular Cardiac Index by Teichholz Method | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| LV/Teichholz(M)_LV_Mass | Left Ventricular Mass by M-Mode | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| LV/Teichholz(M)_LVET | | (18041-4, LN, "Aortic Valve Ejection Time") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| LV/Teichholz(M)_LV_Mass-c | | (CA_LV_MASSc, 99SIEMENSUS, "Left Ventricle Mass corrected") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| LV/Teichholz(M)_LV_Mass-l | | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area"); (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| LV/Teichholz(M)_SV | Left Ventricular Stroke Volume | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|----------------------|---|---|-----------------------------------|---|
| Teichholz(M)_RVDd | Right Ventricular Internal Diastolic Dimension by M-Mode | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Teichholz(M)_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Teichholz(M)_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Teichholz(M)_mVcf | | (CA_LF_VCF, 99SIEMENSUS, "systolic index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (125209, DCM, "Teichholz") |
| Simpson_BP_EDV | Left Ventricular End Diastolic Volume by 2-D Biplane by Method of Disks | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane") |
| Simpson_BP_EDV2 | | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Simpson_BP_EDV4 | | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Simpson_BP_ESV | Left Ventricular End Systolic Volume by 2-D Biplane by Method of Disks | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane") |
| Simpson_BP_ESV2 | | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Simpson_BP_ESV4 | | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-------------------|---|---|----------------------------------|---|
| Simpson_BP_EF | Left Ventricular EF by 2-D Biplane by Method of Disks | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (1252087, DCM, "Method of Disks, Biplane") |
| Simpson_BP_EF2 | Left Ventricular EF | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Simpson_BP_EF4 | Left Ventricular EF | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Simpson_BP_SV | Left Ventricular Stroke Volume by 2-D Biplane by Method of Disks | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane") |
| Simpson_BP_SV2 | | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Simpson_BP_SV4 | | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Simpson_BP_CO | Left Ventricular Cardiac Output by 2-D Biplane by Method of Disks | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane") |
| Simpson_BP_CO2 | | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Simpson_BP_CO4 | | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|----------------------|--|---|----------------------------------|---|
| Simpson_BP_CI | Left Ventricular Cardiac Index by 2-D Biplane by Method of Disks | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane") |
| Simpson_BP_CI2 | | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Simpson_BP_CI4 | | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Simpson_BP_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane") |
| Simpson_BP_SI2 | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Simpson_BP_SI4 | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Simpson_BP_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | |
| Simpson_BP_difD(2ch) | | (CA_LV_DIF, 99SIEMENSUS, LV difference) | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (R-4089A, SRT, "Cardiac Cycle Point") = (F-32011, SRT, "End Diastole") |
| Simpson_BP_difD(4ch) | | (CA_LV_DIF, 99SIEMENSUS, LV difference) | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (R-4089A, SRT, "Cardiac Cycle Point") = (F-32011, SRT, "End Diastole") |
| Simpson_BP_difS(2ch) | | (CA_LV_DIF, 99SIEMENSUS, LV difference) | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (109070, DCM, "Cardiac Cycle Point") = (F-32011, SRT, "End Systole") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|----------------------|--|--|----------------------------------|--|
| Simpson_BP_difS(4ch) | | (CA_LV_DIF, 99SIEMENSUS, LV difference) | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (109070, DCM, "Cardiac Cycle Point") = (F-32011, SRT, "End Systole") |
| Simpson_SP_EDV | Left Ventricular End Diastolic Volume by 2-D Single Plane by Method of Disks (4-Chamber) | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber"); (G-C036, SRT, "Measurement Method") = (125208, DCM, "Method of Disks, Single Plane") |
| Simpson_SP_ESV | Left Ventricular End Systolic Volume by 2D Single Plane by Method of Disks (4-Chamber) | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber"); (G-C036, SRT, "Measurement Method") = (125208, DCM, "Method of Disks, Single Plane") |
| Simpson_SP_EF | Left Ventricular EF by 2D Single Plane by Method of Disks (4-Chamber) | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") (G-C036, SRT, "Measurement Method") = (125208, DCM, "Method Of Disks, Single Plane") |
| Simpson_SP_SV | Left Ventricular Stroke Volume by 2-D Single Plane by Method of Disks (4-Chamber) | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (1110321 DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber"); (G-C036, SRT, "Measurement Method") = (125208, DCM, "Method of Disks, Single Plane") |
| Simpson_SP_CO | Left Ventricular Cardiac Output by 2-D Single Plane by Method of Disks (4-Chamber) | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber"); (G-C036, SRT, "Measurement Method") = (125208, DCM, "Method of Disks, Single Plane") |
| Simpson_SP_CI | Left Ventricular Cardiac Index by 2-D Single Plane by Method of Disks (4-Chamber) | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber"); (G-C036, SRT, "Measurement Method") = (125208, DCM, "Method Of Disks, Single Plane") |
| Simpson_SP_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber"); (G-C036, SRT, "Measurement Method") = (125208, DCM, "Method Of Disks, Single Plane") |
| Simpson_SP_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|------------------------|---|---|----------------------------------|--|
| Mitral_Valve_SV | Left Ventricular Stroke Volume by Doppler Volume Flow | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left Ventricle Outflow Tract") |
| Mitral_Valve_CO | Left Ventricular Cardiac Output by Doppler Volume Outflow | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left Ventricle Outflow Tract") |
| Mitral_Valve_CI | Left Ventricular Cardiac Index by Doppler Volume Flow | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow") |
| Mitral_Valve_MV_Vmean | | (20352-1, LN, "Mean Velocity") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Mitral_Valve_MV_PGmean | | (20256-4, LN, "Mean Gradient") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Mitral_Valve_MV_VTI | | (20354-7, LN, "Velocity Time Integral") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Mitral_Valve_IVRT | Left Ventricular Isovolumic Relaxation Time | (18071-1, LN, "Left Ventricular Isovolumic Relaxation Time") | (T-32600, SRT, "Left Ventricle") | |
| Mitral_Valve_MV_E_pt | | (18037-2, LN, "Mitral Valve E-Wave Peak Velocity") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_MV_A_pt | | (17978-8, LN, "Mitral Valve A-Wave Peak Velocity") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_Dec_Time | | (G-0384, SRT, "Mitral Valve E-Wave Deceleration Time") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_Dec_Slope | | (20216-8, LN, "Deceleration Slope") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_MV_Vmax | | (20351-3, LN, "Maximum Velocity") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Mitral_Valve_MV_PGmax | | (18057-0, LN, "Mitral Valve Diastolic Peak Instantaneous Gradient") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_MV_PHT | | (20280-4, LN, "Pressure Half-Time") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-------------------------|----------------------------|---|----------------------------------|--|
| Mitral_Valve_MV_diam | | (G-038F, SRT, "Cardiovascular Orifice Diameter") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_LVIMP | | (G-037F, SRT, "Left Ventricular Index of Myocardial Performance") | (T-32600, SRT, "Left Ventricle") | |
| Mitral_Valve_LVET | | (18041-4, LN, "Aortic Valve Ejection Time") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left ventricle outflow tract") |
| Mitral_Valve_MV_C_Odur | | (G-0387, SRT, "Mitral Valve Closure to Opening Time") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_A_duration | | (G-0385, SRT, "Mitral Valve A-Wave Duration") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_E_duration | | (CA_MV_DURe, 99SIEMENSUS, "Mitral Valve E-Wave Duration") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left ventricle outflow tract") |
| Mitral_Valve_E/A | | (18038-0, LN, "Mitral Valve E to A Ratio") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_A/E | | (CA_MV_A2E, 99SIEMENSUS, "Mitral Valve A to E Ratio") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_MVA(PHT) | | (G-038E, SRT, "Cardiovascular Orifice Area") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow"); (G-C036, SRT, "Measurement Method") = (125210, DCM, "Area by PHT") |
| Mitral_Valve_LVOT_VTI | | (20354-7, LN, "Velocity Time Integral") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site")= (T-32650, SRT, "Left Ventricle Outflow Tract"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Mitral_Valve_LVOT_diam | | (G-038F, SRT, "Cardiovascular Orifice Diameter") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left Ventricular Outflow Tract") |
| Mitral_Valve_MVA(VTI) | | (G-038E, SRT, "Cardiovascular Orifice Area") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow"); (G-C036, SRT, "Measurement Method") = (125212, DCM, "Continuity Equation") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|---------------------------|---|--|----------------------------------|--|
| Mitral_Valve_MVA(Trace) | | (G-038E, SRT, "Cardiovascular Orifice Area") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125220, DCM, "Planimetry"); (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Mitral_Valve_HR | | (8867-4, LN, "Heart Rate") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_CE_amp | | (CA_MV_EWaveAmpl, 99SIEMENSUS, "Amplitude E Wave MMode") | (T-35300, SRT, "Mitral Valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| Mitral_Valve_CA_amp | | (CA_MV_AWaveAmpl, 99SIEMENSUS, "Amplitude A Wave MMode") | (T-35300, SRT, "Mitral Valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| Mitral_Valve_DE_excursion | | (17997-8, LN, "Mitral Valve D-E Excursion") | (T-35300, SRT, "Mitral Valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| Mitral_Valve_DE_amp | | (CA_MV_DEWaveAmpl, 99SIEMENSUS, "Amplitude D-E Wave MMode") | (T-35300, SRT, "Mitral Valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| Mitral_Valve_EPSS | | (18036-4, LN, "Mitral Valve EPSS, E wave ") | (T-35300, SRT, "Mitral Valve") | |
| Mitral_Valve_EF_slope | | (18040-6, LN, "Mitral Valve E-F Slope by M-Mode") | (T-35300, SRT, "Mitral Valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| Mitral_Valve_CA/CE | | (CA_MV_A2E, 99SIEMENSUS, "Mitral Valve A to E Ratio") | (T-35300, SRT, "Mitral Valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| Aortic_Valve_SV | Left Ventricular Stroke Volume by Doppler Volume Flow | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site") = (T-32650, SNM3, "Left Ventricle Outflow Tract") |
| Aortic_Valve_CO | Left Ventricular Cardiac Output by Doppler Volume Outflow | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left Ventricle Outflow Tract") |
| Aortic_Valve_CI | Left Ventricular Cardiac Index by Doppler Volume Flow | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left Ventricle Outflow Tract") |
| Aortic_Valve_IVRT | Left Ventricular Isovolumic Relaxation Time | (18071-1, LN, "Left Ventricular Isovolumic Relaxation Time") | (T-32600, SRT, "Left Ventricle") | |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-------------------------|----------------------------|--|----------------------------------|---|
| Aortic_Valve_LVPEP | Pre-Ejection Period | (18068-7, LN, "Pre-Ejection Period") | (T-35400, SRT, "Aortic Valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| Aortic_Valve_LVET | | (18041-4, LN, "Aortic Valve Ejection Time") | (T-32600, SRT, "Left Ventricle") | |
| Aortic_Valve_LVSTI | | (CA_AV_PEP2ET, 99SIEMENSUS, "Aortic Valve Pre-ejection Period to Ejection Time") | (T-35400, SRT, "Aortic Valve") | (G-0373, SRT, "Image Mode") = (DOPPLER, 99SIEMENSUS, "Doppler Mode") |
| Aortic_Valve_AV_VTI | | (20354-7, LN, "Velocity Time Integral") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_AV_Vmax | | (11726-7, LN, "Peak Systolic Velocity") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_AV_Vmean | | (20352-1, LN, "Mean Velocity") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_AV_PGmax | | (20247-3, LN, "Peak Gradient") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_AV_PGmean | | (20256-4, LN, "Mean Gradient") | (T-35400, SRT, "Aortic Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_LVOT_VTI | | (20354-7, LN, "Velocity Time Integral") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left ventricle outflow tract"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_LVOT_Vmax | | (11726-7, LN, "Peak Systolic Velocity") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left ventricle outflow tract"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_LVOT_Vmean | | (20352-1, LN, "Mean Velocity") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left ventricle outflow tract"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_LVOT_PGmax | | (20247-3, LN, "Peak Gradient") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left ventricle outflow tract"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|--------------------------|---|---|--|---|
| Aortic_Valve_LVOT_PGmean | | (20256-4, LN, "Mean Gradient") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left ventricle outflow tract"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_LVOT_diam | | (G-038F, SRT, "Cardiovascular Orifice Diameter") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left Ventricular Outflow Tract") |
| Aortic_Valve_AVA(VTI) | | (G-038E, SRT, "Cardiovascular Orifice Area") | (T-35400, SRT, "Aortic Valve"); | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow"); (G-C036, SRT, "Measurement Method") = (125215, DCM, "Continuity Equation by Velocity Time Integral") |
| Aortic_Valve_AVA(Vmax) | | (G-038E, SRT, "Cardiovascular Orifice Area") | (T-35400, SRT, "Aortic Valve"); | (G-C036, SRT, "Measurement Method") = (125214, DCM, "Continuity Equation by Peak Velocity"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Aortic_Valve_AVA(Trace) | | (G-038E, SRT, "Cardiovascular Orifice Area") | (T-35400, SRT, "Aortic Valve") | (G-C036, SRT, "Measurement Method") = (125220, DCM, "Planimetry") |
| Aortic_Valve_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site") = (T-32650, SRT, "Left ventricle outflow tract"); |
| Aortic_Valve_HR | | (8867-4, LN, "Heart Rate") | (T-35400, SRT, "Aortic Valve") | |
| Aortic_Valve_VSD_Vmax | | (11726-7, LN, "Peak Systolic Velocity") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (D4-31150, SRT, "Ventricular Septal Defect") |
| Aortic_Valve_VSD_PGmax | | (20247-3, LN, "Peak Gradient") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (D4-31150, SRT, "Ventricular Septal Defect") |
| Aortic_Valve_Qp/Qs | | (29462-9, LN, "Pulmonary-to-Systemic Shunt Flow Ratio") | (P5-30031, SRT, "Cardiac Shunt Study") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow") |
| Aortic_Valve_Qp-Qs | | (CA_CS_PF-SF, 99SIEMENSUS, "Pulmonary minus Systemic Shunt Flow") | (P5-30031, SRT, "Cardiac Shunt Study") | |
| TR_TR_Vmax | Tricuspid Valve Regurgitant Peak Velocity | (11726-7, LN, "Peak Systolic Velocity") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-------------------|--|---|------------------------------------|---|
| TR_TR_Vmean | | (20352-1, LN, "Mean Velocity") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| TR_TR_PGmax | | (20247-3, LN, "Peak Gradient") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| TR_TR_PGmean | | (20256-4, LN, "Mean Gradient") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| TR_TR_VTI | Tricuspid Regurgitation Velocity Time Integral | (20354-7, LN, "Velocity Time Integral") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| TR_TR_RAP | Right Atrium Systolic Pressure | (18070-3, LN, "Right Atrium Systolic Pressure") | (T-32200, SRT, "Right Atrium") | |
| TR_RVSP | | (G-0380, SRT, "Right Ventricular Peak Systolic Pressure") | (T-32500, SRT, "Right Ventricle") | |
| PR_PR_Vmax | Pulmonary Valve Regurgitant Peak Velocity | (20351-3, LN, "Maximum Velocity") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PR_PR_Vmean | | (20352-1, LN, "Mean Velocity") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PR_PR_PGmax | | (20247-3, LN, "Peak Gradient") | (T-35200, SRT, "Pulmonary valve"); | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PR_PR_PGmean | | (20256-4, LN, "Mean Gradient") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PR_PR_Ved | | (11653-3, LN, "End Diastolic Velocity") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PR_PR_VTI | Pulmonary Regurgitation Velocity Time Integral | (20354-7, LN, "Velocity Time Integral") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PR_PAEDP | | (F-0212C, SRT, "Pulmonary Artery Pressure") | (T-44000, SRT, "Pulmonary artery") | (R-4089A, SRT, "Cardiac Cycle Point") = (F-32011, SRT, "End Diastole") |
| PR_RAP | Right Atrium Systolic Pressure | (18070-3, LN, "Right Atrium Systolic Pressure") | (T-32200, SRT, "Right Atrium") | |
| MR_MR_Vmax | | (11726-7, LN, "Peak Systolic Velocity") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| MR_MR_PGmax | | (20247-3, LN, "Peak Gradient") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| MR_dt | | (20217-6, LN, "Deceleration Time") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-----------------------------|-------------------------------------|--|------------------------------------|---|
| MR_dP/dt | | (18035-6, LN, "Mitral Regurgitation dP/dt derived from Mitral Regurgitation velocity") | (T-35300, SRT, "Mitral Valve") | |
| Pulmonary_Valve_PV_Vmax | | (11726-7, LN, "Peak Systolic Velocity") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Pulmonary_Valve_PV_PGmax | | (20247-3, LN, "Peak Gradient") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Pulmonary_Valve_PV_Vmean | | (20352-1, LN, "Mean Velocity") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Pulmonary_Valve_PV_PGmean | | (20256-4, LN, "Mean Gradient") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Pulmonary_Valve_PV_VTI | | (20354-7, LN, "Velocity Time Integral") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Pulmonary_Valve_RVET | | (18042-2, LN, "Pulmonic Valve Ejection Time") | (T-35200, SRT, "Pulmonary valve") | |
| Pulmonary_Valve_RV_Acc_Time | | (20168-1, LN, "Acceleration Time") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Pulmonary_Valve_PA_Acc_Time | | (20168-1, LN, "Acceleration Time") | (T-44000, SRT, "Pulmonary artery") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Pulmonary_Valve_RVPEP | Right Ventricle Pre-Ejection Period | (20301-8, LN, "Right Ventricle Pre-Ejection Period") | (T-35200, SRT, "Pulmonary valve") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); |
| Pulmonary_Valve_PV_diam | | (G-038F, SRT, "Cardiovascular Orifice Diameter") | (T-32500, SRT, "Right Ventricle") | (G-C0E3, SRT, "Finding Site") = (T-32550, SRT, "Right Ventricular Outflow Tract"); (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Pulmonary_Valve_SV | | (F-32120, SRT, "Stroke Volume") | (T-32500, SRT, "Right Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site") = (T-32550, SRT, "Right Ventricular Outflow Tract") |
| Pulmonary_Valve_CO | | (F-32100, SRT, "Cardiac Output") | (T-32500, SRT, "Right Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site") = (T-32550, SRT, "Right Ventricular Outflow Tract") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-----------------------------|---|--|--|--|
| Pulmonary_Valve_SI | | (F-00078, SRT, "Stroke Index") | (T-32500, SRT, "Right Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site")= (T-32550, SRT, "Right Ventricle Outflow Tract") |
| Pulmonary_Valve_CI | | (F-32110, SRT, "Cardiac Index") | (T-32500, SRT, "Right Ventricle") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow"); (G-C0E3, SRT, "Finding Site")= (T-32550, SRT, "Right Ventricle Outflow Tract") |
| Pulmonary_Valve_RVSTI | | (CA_PV_PEP2ET, 99SIEMENSUS, "Pulmonic Valve Pre-ejection Period to Ejection Time") | (T-35200, SRT, "Pulmonary valve") | (G-0373, SRT, "Image Mode") = (DOPPLER, 99SIEMENSUS, "Doppler Mode") |
| Pulmonary_Valve_RV_Acc_T/ET | | (G-0388, SRT, "Ratio of Pulmonic Valve Acceleration Time to Ejection Time") | (T-35200, SRT, "Pulmonary valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Pulmonary_Valve_HR | | (8867-4, LN, "Heart Rate") | (T-35200, SRT, "Pulmonary valve") | |
| Pulmonary_Valve_Qp/Qs | Pulmonary-to-Systemic Shunt Flow Ratio by Doppler Volume Flow | (29462-9, LN, "Pulmonary-to-Systemic Shunt Flow Ratio") | (P5-30031, SRT, "Cardiac Shunt Study") | (G-C036, SRT, "Measurement Method") = (125219, DCM, "Doppler Volume Flow") |
| Pulmonary_Valve_Qp-Qs | | (CA_CS_PF-SF, 99SIEMENSUS, "Pulmonary minus Systemic Shunt Flow") | (P5-30031, SRT, "Cardiac Shunt Study") | |
| Pulmonary_Vein_PVs1_Vel | | (29450-4, LN, "Pulmonary Vein Systolic Peak Velocity") | (T-48581, SRT, "Pulmonary Venous Structure") | |
| Pulmonary_Vein_PVs2_Vel | | (29450-4, LN, "Pulmonary Vein Systolic Peak Velocity") | (T-48581, SRT, "Pulmonary Venous Structure") | |
| Pulmonary_Vein_PVd_Vel | | (29451-2, LN, "Pulmonary Vein Diastolic Peak Velocity") | (T-48581, SRT, "Pulmonary Venous Structure") | |
| Pulmonary_Vein_PVa_Vel | | (29453-8, LN, "Pulmonary Vein Atrial Contraction Reversal Peak Velocity") | (T-48581, SRT, "Pulmonary Venous Structure") | (R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-----------------------------|---|---|--|--|
| Pulmonary_Vein_PVa_Dur | | (G-038B, SRT, "Pulmonary Vein A-Wave Duration") | (T-48581, SRT, "Pulmonary Venous Structure") | |
| Pulmonary_Vein_PVs_VTI | | (G-038C, SRT, "Pulmonary Vein S-Wave Velocity Time Integral") | (T-48581, SRT, "Pulmonary Venous Structure") | |
| Pulmonary_Vein_PVd_VTI | | (G-038D, SRT, "Pulmonary Vein D-Wave Velocity Time Integral") | (T-48581, SRT, "Pulmonary Venous Structure") | |
| Pulmonary_Vein_PVs2/PVd | | (29452-0, LN, "Pulmonary Vein Systolic to Diastolic Ratio") | (T-48581, SRT, "Pulmonary Venous Structure") | |
| Pulmonary_Vein_PVd_Decel_T | | (20217-6, LN, "Deceleration Time") | (T-48581, SRT, "Pulmonary Venous Structure") | (R-4089A, SRT, "Cardiac Cycle Point") = (F-32010, SRT, "Diastole") |
| Pulmonary_Vein_Sys_Fraction | | (CA_PVE_SF, 99SIEMENSUS, "Pulmonary Vein Systolic Fraction") | (T-48581, SRT, "Pulmonary Venous Structure") | |
| RV_RV_diam | | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension") | (T-32500, SRT, "Right Ventricle") | |
| LV_Mass_T-E_A_Sax_Epi | Left Ventricle Epicardial Diastolic Area, psax pap view | (G-0379, SRT, "Left Ventricle Epicardial Diastolic Area, psax pap view") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125222, DCM, "Left Ventricle Mass by Truncated Ellipse") |
| LV_Mass_T-E_A_Sax_Endo | | (G-0375, SRT, "Left Ventricular Diastolic Area") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125222, DCM, "Left Ventricle Mass by Truncated Ellipse"); (111031, DCM, "Image View") = (G-039B, SRT, "Parasternal short axis at the Papillary Muscle level") |
| LV_Mass_T-E_a | | (G-0377, SRT, "Left Ventricle Semi-major Axis Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125222, DCM, "Left Ventricle Mass by Truncated Ellipse") |
| LV_Mass_T-E_d | | (G-0378, SRT, "Left Ventricle Truncated Semi-major Axis Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125222, DCM, "Left Ventricle Mass by Truncated Ellipse") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|---------------------------|----------------------------|---|----------------------------------|--|
| LV_Mass_T-E_t | | (CA_LV_MyoTh, 99SIEMENSUS, "Myocardial thickness calculated from short axis epicardial and cavity areas") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125222, DCM, "Left Ventricle Mass by Truncated Ellipse") |
| LV_Mass_T-E_b | | (CA_LV_RadSAX, 99SIEMENSUS, "Short axis radius calculated from short axis cavity area") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125222, DCM, "Left Ventricle Mass by Truncated Ellipse") |
| LV_Mass_T-E_ LV_Mass | | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125222, DCM, "Left Ventricle Mass by Truncated Ellipse") |
| LV_Mass_T-E_ LV_Mass-I | | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area"); (G-C036, SRT, "Measurement Method") = (125222, DCM, "Left Ventricle Mass by Truncated Ellipse") |
| LV_Mass_A-L_A_Sax_Epi | | (G-0379, SRT, "Left Ventricle Epicardial Diastolic Area, psax pap view") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (LVMassAL, 99SIEMENSUS, "Left Ventricle Mass by Area Length") |
| LV_Mass_A-L_A_Sax_Endo | | (G-0375, SRT, "Left Ventricular Diastolic Area") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (111031, DCM, "Image View") = (G-039B, SRT, "Parasternal short axis at the Papillary Muscle level"); (G-C036, SRT, "Measurement Method") = (LVMassAL, 99SIEMENSUS, "Left Ventricle Mass by Area Length") |
| LV_Mass_A-L_LVL | | (18077-8, LN, "Left Ventricle diastolic major axis") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (LVMassAL, 99SIEMENSUS, "Left Ventricle Mass by Area Length") |
| LV_Mass_A-L_t | | (CA_LV_MyoTh, 99SIEMENSUS, "Myocardial thickness calculated from short axis epicardial and cavity areas") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (LVMassAL, 99SIEMENSUS, "Left Ventricle Mass by Area Length") |
| LV_Mass_A-L_LV_Mass | | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (LVMassAL, 99SIEMENSUS, "Left Ventricle Mass by Area Length") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-----------------------|----------------------------|--|----------------------------------|---|
| LV_Mass_A-L_LV_Mass-I | | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area"); (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (G-C036, SRT, "Measurement Method") = (LVMassAL, 99SIEMENSUS, "Left Ventricle Mass by Area Length") |
| DTI_Ea(m) | | (G-037A, SRT, "Left Ventricular Peak Early Diastolic Tissue Velocity") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0391, SRT, "Medial Mitral Annulus") |
| DTI_Aa(m) | | (G-037C, SRT, "LV Peak Diastolic Tissue Velocity During Atrial Systole") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0391, SRT, "Medial Mitral Annulus") |
| DTI_Ea/Aa(m) | | (CA_MV_Ea2Aa, 99SIEMENSUS, "Mitral Valve DTI Ea to Aa Ratio") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0391, SRT, "Medial Mitral Annulus") |
| DTI_E/Ea(m) | | (G-037B, SRT, "Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0391, SRT, "Medial Mitral Annulus") |
| DTI_ATa(m) | | (CA_MV_ATa, 99SIEMENSUS, "Mitral Valve DTI acceleration time of Ea") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0391, SRT, "Medial Mitral Annulus") |
| DTI_ARa(m) | | (CA_MV_ARa, 99SIEMENSUS, "Mitral Valve DTI acceleration rate of Ea") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0391, SRT, "Medial Mitral Annulus") |
| DTI_DTa(m) | | (CA_MV_DTa, 99SIEMENSUS, "Mitral Valve DTI deceleration time of Ea") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0391, SRT, "Medial Mitral Annulus") |
| DTI_DRa(m) | | (CA_MV_DRa, 99SIEMENSUS, "Mitral Valve DTI deceleration rate of Ea") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0391, SRT, "Medial Mitral Annulus") |
| DTI_Sa(m) | | (G-037D, SRT, "Left Ventricular Peak Systolic Tissue Velocity") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0391, SRT, "Medial Mitral Annulus") |
| DTI_Ea(l) | | (G-037A, SRT, "Left Ventricular Peak Early Diastolic Tissue Velocity") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0392, SRT, "Lateral Mitral Annulus") |
| DTI_Aa(l) | | (G-037C, SRT, "LV Peak Diastolic Tissue Velocity During Atrial Systole") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0392, SRT, "Lateral Mitral Annulus") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|---------------------------|----------------------------|--|-----------------------------------|---|
| DTI_Ea/Aa(l) | | (CA_MV_Ea2Aa, 99SIEMENSUS, "Mitral Valve DTI Ea to Aa Ratio") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0392, SRT, "Lateral Mitral Annulus") |
| DTI_E/Ea(l) | | (G-037B, SRT, "Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0392, SRT, "Lateral Mitral Annulus") |
| DTI_ATa(l) | | (CA_MV_ATa, 99SIEMENSUS, "Mitral Valve DTI acceleration time of Ea") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0392, SRT, "Lateral Mitral Annulus") |
| DTI_ARa(l) | | (CA_MV_ARa, 99SIEMENSUS, "Mitral Valve DTI acceleration rate of Ea") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0392, SRT, "Lateral Mitral Annulus") |
| DTI_DTa(l) | | (CA_MV_DTa, 99SIEMENSUS, "Mitral Valve DTI deceleration time of Ea") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0392, SRT, "Lateral Mitral Annulus") |
| DTI_DRa(l) | | (CA_MV_DRa, 99SIEMENSUS, "Mitral Valve DTI deceleration rate of Ea") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0392, SRT, "Lateral Mitral Annulus") |
| DTI_Sa(l) | | (G-037D, SRT, "Left Ventricular Peak Systolic Tissue Velocity") | (T-32600, SRT, "Left Ventricle") | (G-C0E3, SRT, "Finding Site") = (G-0392, SRT, "Lateral Mitral Annulus") |
| Tricuspid_Valve_TV_VTI | | (20354-7, LN, "Velocity Time Integral") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Tricuspid_Valve_TV_Vmax | | (20351-3, LN, "Maximum Velocity") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Tricuspid_Valve_TV_PGmax | | (20247-3, LN, "Peak Gradient") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Tricuspid_Valve_TV_Vmean | | (20352-1, LN, "Mean Velocity") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Tricuspid_Valve_TV_PGmean | | (20256-4, LN, "Mean Gradient") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Tricuspid_Valve_TV_E_pt | | (18031-5, LN, "Tricuspid Valve E Wave Peak Velocity") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Tricuspid_Valve_TV_A_pt | | (18030-7, LN, "Tricuspid Valve A Wave Peak Velocity") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| Tricuspid_Valve_E/A | | (18039-8, LN, "Tricuspid Valve E to A Ratio") | (T-35100, SRT, "Tricuspid valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|---------------------------|----------------------------|--|-----------------------------------|--|
| Tricuspid_Valve_RVIMP | | (G-0381, SRT, "Right Ventricular Index of Myocardial Performance") | (T-32500, SRT, "Right Ventricle") | |
| Tricuspid_Valve_RVET | | (18042-2, LN, "Pulmonic Valve Ejection Time") | (T-32500, SRT, "Right Ventricle") | |
| Tricuspid_Valve_TV_C-Odur | | (G-0389, SRT, "Tricuspid Valve Closure to Opening Time") | (T-35100, SRT, "Tricuspid valve") | |
| PISA(MR)_Radius | | (CA_DOME, 99SIEMENSUS, "Dome Radius") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PISA(MR)_Aliasing_Vel | | (20252-1, LN, "Mean Velocity") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PISA(MR)_MR_VTI | | (20354-7, LN, "Velocity Time Integral") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PISA(MR)_MR_Vmax | | (11726-7, LN, "Peak Systolic Velocity") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PISA(MR)_Flow_Rate | | (34141-2, LN, "Peak Instantaneous Flow Rate") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|-----------------------|----------------------------|--|--------------------------------|---|
| PISA(MR)_EO_Area | | (G-038E, SRT, "Cardiovascular Orifice Area") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow"); (G-0373, SRT, "Image Mode") = (DOPPLER, 99SIEMENSUS, "Doppler Mode") |
| PISA(MR)_Flow_Vol | | (33878-0, LN, "Volume Flow") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42E61, SRT, "Regurgitant Flow") |
| PISA(MS)_Radius | | (CA_DOME, 99SIEMENSUS, "Dome Radius") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| PISA(MS)_Aliasing_Vel | | (20252-1, LN, "Mean Velocity") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| PISA(MS)_MS_VTI | | (20354-7, LN, "Velocity Time Integral") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| PISA(MS)_MS_Vmax | | (20351-3, LN, "Maximum Velocity") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| PISA(MS)_MS_PGmax | | (20247-3, LN, "Peak Gradient") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|--------------------|----------------------------|---|--------------------------------|---|
| PISA(MS)_Angle | | (CA_MS_Angle, 99SIEMENSUS, "Angle measured at Mitral Valve Stenosis") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| PISA(MS)_Flow_Rate | | (34141-2, LN, "Peak Instantaneous Flow Rate") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| PISA(MS)_Flow_Vol | | (33878-0, LN, "Volume Flow") | (T-35300, SRT, "Mitral Valve") | (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow") |
| PISA(MS)_MVA | | (G-038E, SRT, "Cardiovascular Orifice Area") | (T-35300, SRT, "Mitral Valve") | (G-C048, SRT, "Direction of Flow") = (R-42047, SRT, "Antegrade Flow"); (G-C036, SRT, "Measurement Method") = (125216, DCM, "Proximal Isovelocity Surface Area"); (G-0373, SRT, "Image Mode") = (DOPPLER, 99SIEMENSUS, "Doppler Mode") |
| LA_Vol(A-L)_A1-A4C | | (17977-0, LN, "Left Atrium Systolic Area") | (T-32300, SRT, "Left Atrium") | (G-C036, SRT, "Measurement Method") = (125220, DCM, "Planimetry"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| LA_Vol(A-L)_A2-A2C | | (17977-0, LN, "Left Atrium Systolic Area") | (T-32300, SRT, "Left Atrium") | (G-C036, SRT, "Measurement Method") = (125220, DCM, "Planimetry"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| LA_Vol(A-L)_L | | (29469-4, LN, "Left Atrium Antero-posterior Systolic Dimension") | (T-32300, SRT, "Left Atrium") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|---------------------------|----------------------------|---|--------------------------------|---|
| LA_Vol(A-L)_LA_Vol(A-L) | | (122408, DCM, "Left Atrial ES Volume") | (T-32300, SRT, "Left Atrium") | (G-C036, SRT, "Measurement Method") = (125204, DCM, "Area-Length Biplane"); (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| LA_Vol(A-L)_LA_Vol-I | | (122408, DCM, "Left Atrial ES Volume") | (T-32300, SRT, "Left Atrium") | (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area"); (G-C036, SRT, "Measurement Method") = (125204, DCM, "Area-Length Biplane"); (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| LA_Vol(Simp)_4CH | | (122408, DCM, "Left Atrial ES Volume") | (T-32300, SRT, "Left Atrium") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Single Plane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| LA_Vol(Simp)_2CH | | (122408, DCM, "Left Atrial ES Volume") | (T-32300, SRT, "Left Atrium") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Single Plane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| LA_Vol(Simp)_LA_Vol(Simp) | | (122408, DCM, "Left Atrial ES Volume") | (T-32300, SRT, "Left Atrium") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); |
| LA_Vol(Simp)_LA_Vol-I | | (122408, DCM, "Left Atrial ES Volume") | (T-32300, SRT, "Left Atrium") | (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area"); (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| RA_Vol(A-L)_A1-A4C | | (17988-7, LN, "Right Atrium Systolic Area") | (T-32200, SRT, "Right Atrium") | (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| RA_Vol(A-L)_A2-A2C | | (17988-7, LN, "Right Atrium Systolic Area") | (T-32200, SRT, "Right Atrium") | (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") ; (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole") |
| RA_Vol(A-L)_L | | (CA_RA_DIA, 99SIEMENSUS, "Right Atrium Systolic Dimension") | (T-32200, SRT, "Right Atrium") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|---------------------------|----------------------------|--|----------------------------------|---|
| RA_Vol(A-L)_RA_Vol(A-L) | | (CA_RA_VOL, 99SIEMENSUS, "Right Atrium Systolic Volume") | (T-32200, SRT, "Right Atrium") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (R-4089A, SRT, "Cardiac Cycle Point") = (F-32020, SRT, "Systole") |
| RA_Vol(A-L)_RA_Vol-I | | (CA_RA_VOL, 99SIEMENSUS, "Right Atrium Systolic Volume") | (T-32200, SRT, "Right Atrium") | (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area"); (G-C036, SRT, "Measurement Method") = (125204, DCM, "Area-Length Biplane") |
| RA_Vol(Simp)_4CH | | (CA_RA_VOL, 99SIEMENSUS, "Right Atrium Systolic Volume") | (T-32200, SRT, "Right Atrium") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| RA_Vol(Simp)_2CH | | (CA_RA_VOL, 99SIEMENSUS, "Right Atrium Systolic Volume") | (T-32200, SRT, "Right Atrium") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| RA_Vol(Simp)_RA_Vol(Simp) | | (CA_RA_VOL, 99SIEMENSUS, "Right Atrium Systolic Volume") | (T-32200, SRT, "Right Atrium") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane") |
| RA_Vol(Simp)_RA_Vol-I | | (CA_RA_VOL, 99SIEMENSUS, "Right Atrium Systolic Volume") | (T-32200, SRT, "Right Atrium") | (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area"); (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane") |
| Axius-EF_EF2 | | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (1252087, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Axius-EF_SV2 | | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Axius-EF_CO2 | | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Axius-EF_HR2 | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|--------------------------|-------------------------------------|--|----------------------------------|--|
| Axius-EF_EDV2 | | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Axius-EF_ESV2 | | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19B, SRT, "Apical Two Chamber") |
| Axius-EF_EF4 | | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (1252087, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Axius-EF_SV4 | | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Axius-EF_CO4 | | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Axius-EF_HR4 | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Axius-EF_EDV4 | | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Axius-EF_ESV4 | | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-C036, SRT, "Measurement Method") = (125207, DCM, "Method of Disks, Biplane"); (111031, DCM, "Image View") = (G-A19C, SRT, "Apical Four Chamber") |
| Mod._Simpson_LVLd_apical | Left Ventricle Diastolic Major Axis | (18077-8, LN, "Left Ventricle Diastolic Major Axis") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|--------------------------|---------------------------------------|--|----------------------------------|---|
| Mod._Simpson_LVAd_sax_MV | Left Ventricular Diastolic Area | (G-0375, SRT, "Left Ventricular Diastolic Area") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (111031, DCM, "Image View") = (G-039A, SRT, "Parasternal short axis at the Mitral Valve level") |
| Mod._Simpson_LVAd_sax_PM | Left Ventricular Diastolic Area | (G-0375, SRT, "Left Ventricular Diastolic Area") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (111031, DCM, "Image View") = (G-039B, SRT, "Parasternal short axis at the Papillary Muscle level") |
| Mod._Simpson_LVLs_apical | Left Ventricle systolic major axis | (18076-0, LN, "Left Ventricle systolic major axis") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Mod._Simpson_LVAs_sax_MV | Left Ventricular Systolic Area | (G-0374, SRT, "Left Ventricular Systolic Area") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (111031, DCM, "Image View") = (G-039A, SRT, "Parasternal short axis at the Mitral Valve level") |
| Mod._Simpson_LVAs_sax_PM | Left Ventricular Systolic Area | (G-0374, SRT, "Left Ventricular Systolic Area") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (111031, DCM, "Image View") = (G-039B, SRT, "Parasternal short axis at the Papillary Muscle level") |
| Mod._Simpson_EDV | Left Ventricular End Diastolic Volume | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Mod._Simpson_ESV | Left Ventricular End Systolic Volume | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Mod._Simpson_SV | Left Ventricular Stroke Volume | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Mod._Simpson_CO | Left Ventricular Cardiac Output | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Mod._Simpson_EF | Left Ventricular EF | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Mod._Simpson_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Mod._Simpson_CI | Left Ventricular Cardiac Index | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Mod._Simpson_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | |
| Bullet_LVLd_apical | Left Ventricle Diastolic Major Axis | (18077-8, LN, "Left Ventricle Diastolic Major Axis") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Bullet_LVAd_sax_MV | Left Ventricular Diastolic Area | (G-0375, SRT, "Left Ventricular Diastolic Area") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (111031, DCM, "Image View") = (G-039A, SRT, "Parasternal short axis at the Mitral Valve level") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|--------------------|--|--|-----------------------------------|---|
| Bullet_LVLs_apical | Left Ventricle systolic major axis | (18076-0, LN, "Left Ventricle systolic major axis") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Bullet_LVAs_sax_MV | Left Ventricular Systolic Area | (G-0374, SRT, "Left Ventricular Systolic Area") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode"); (111031, DCM, "Image View") = (G-039A, SRT, "Parasternal short axis at the Mitral Valve level") |
| Bullet_EDV | Left Ventricular End Diastolic Volume | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Bullet_ESV | Left Ventricular End Systolic Volume | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Bullet_SV | Left Ventricular Stroke Volume | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Bullet_CO | Left Ventricular Cardiac Output | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Bullet_EF | Left Ventricular EF | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Bullet_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Bullet_CI | Left Ventricular Cardiac Index | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Bullet_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | |
| Gibson(2D)_RVAWd | Right Ventricle Anterior Wall Diastolic Thickness | (18153-7, LN, "Right Ventricle Anterior Wall Diastolic Thickness") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_RVDd | Right Ventricular Internal Diastolic Dimension by 2-D | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_IVSd | Interventricular Septum Diastolic Thickness by 2-D | (18154-5, LN, "Interventricular Septum Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_LVIDd | Left Ventricular Internal End Diastolic Dimension by 2-D | (29436-3, LN, "Left Ventricle Internal End Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_LVPWd | Left Ventricular Posterior Wall Diastolic Thickness by 2-D | (18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|--------------------|--|--|-----------------------------------|--|
| Gibson(2D)_IVSs | Interventricular Septum Systolic Thickness by 2-D | (18158-6, LN, "Interventricular Septum Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_LVIDs | Left Ventricular Internal Systolic Dimension by 2-D | (29438-9, LN, "Left Ventricle Internal Systolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_LVPWs | Left Ventricle Posterior Wall Systolic Thickness by 2-D | (18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_EDV | Left Ventricular End Diastolic Volume | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_ESV | Left Ventricular End Systolic Volume | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_SV | Left Ventricular Stroke Volume | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_CO | Left Ventricular Cardiac Output | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_EF | Left Ventricular EF | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_FS | Left Ventricular Fractional Shortening by 2-D | (18051-3, LN, "Left Ventricular Fractional Shortening") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_CI | Left Ventricular Cardiac Index | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-03A2, SRT, "2D mode") |
| Gibson(2D)_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | |
| LV/Gibson(M)_RVDd | Right Ventricular Internal Diastolic Dimension by M-Mode | (20304-2, LN, "Right Ventricular Internal Diastolic Dimension") | (T-32500, SRT, "Right Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_IVSd | Interventricular Septum Diastolic Thickness | (18154-5, LN, "Interventricular Septum Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_LVIDd | Left Ventricle Internal End Diastolic Dimension | (29436-3, LN "Left Ventricle Internal End Diastolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_LVPWd | Left Ventricle Posterior Wall Diastolic Thickness | (18152-9, LN, "Left Ventricle Posterior Wall Diastolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|----------------------|--|---|----------------------------------|---|
| LV/Gibson(M)_IVSs | Interventricular Septum Systolic Thickness | (18158-6, LN, "Interventricular Septum Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_LVIDs | Left Ventricle Internal Systolic Dimension | (29438-9, LN, "Left Ventricle Internal Systolic Dimension") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_LVPWs | Left Ventricle Posterior Wall Systolic Thickness | (18156-0, LN, "Left Ventricle Posterior Wall Systolic Thickness") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_EDV | Left Ventricular End Diastolic Volume | (18026-5, LN, "Left Ventricular End Diastolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_ESV | Left Ventricular End Systolic Volume | (18148-7, LN, "Left Ventricular End Systolic Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_SV | Left Ventricular Stroke Volume | (F-32120, SRT, "Stroke Volume") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_CO | Left Ventricular Cardiac Output | (F-32100, SRT, "Cardiac Output") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_EF | Left Ventricular EF | (18043-0, LN, "Left Ventricular Ejection Fraction") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_FS | Left Ventricular Fractional Shortening by 2-D | (18051-3, LN, "Left Ventricular Fractional Shortening") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_SI | | (F-00078, SRT, "Stroke Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_CI | Left Ventricular Cardiac Index | (F-32110, SRT, "Cardiac Index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_HR | | (8867-4, LN, "Heart Rate") | (T-32600, SRT, "Left Ventricle") | |
| LV/Gibson(M)_mVcf | | (CA_LF_VCF, 99SIEMENSUS, "systolic index") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (MeasPackage, 99SIEMENSUS, Measurement Package Membership) = (MPMMode, 99SIEMENSUS, "MMode Gibson Package) |
| LV/Gibson(M)_LVET | | (18041-4, LN, "Aortic Valve Ejection Time") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |
| LV/Gibson(M)_LV_Mass | Left Ventricular Mass by M-Mode | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode") |

| X500 Label | Name of ASE Concept | Base Measurement Concept Name | Echo Section | Concept or Acquisition Context Modifiers |
|------------------------|---------------------|---|----------------------------------|---|
| LV/Gibson(M)_LV_Mass-c | | (CA_LV_MASSc, 99SIEMENSUS, "Left Ventricle Mass corrected") | (T-32600, SRT, "Left Ventricle") | (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (G-C036, SRT, "Measurement Method") = (MethodMassASEcor, 99SIEMENSUS, "LeftVentricle corrected by M-Mode"); (MeasPackage, 99SIEMENSUS, Measurement Package Membership) = (MPMMode, 99SIEMENSUS, "MMode Gibson Package) |
| LV/Gibson(M)_LV_Mass-l | | (18087-7, LN, "Left Ventricle Mass") | (T-32600, SRT, "Left Ventricle") | (121425, DCM, "Index") = (8277-6, LN, "Body Surface Area"); (G-0373, SRT, "Image Mode") = (G-0394, SRT, "M mode"); (MeasPackage, 99SIEMENSUS, Measurement Package Membership) = (MPMMode, 99SIEMENSUS, "MMode Gibson Package) |