

# Device Errors

## Biliblanket phototherapy light

**FOLLOW THESE TIPS TO USE IT SAFELY.**

BY EILEEN K. WOO, RN, BSN  
Nurse-Consultant

Center for Devices and Radiological Health  
Food and Drug Administration • Rockville, Md.

A 6-DAY-OLD, 34-WEEK-GESTATION INFANT WAS BEING treated with a biliblanket phototherapy light for **hyperbilirubinemia**. **When** the nurse removed the biliblanket from the infant, she noticed a scratch and some blood on the infant's knee.

The biliblanket, used to treat jaundice, consists of a fiber-optic pad and a portable illuminator from which therapeutic light is delivered to the infant. The pad should be covered completely with a disposable cover and placed underneath or wrapped around the infant. Any alteration of the cover, such as cutting it in half, may contribute to an adverse event.

### What went wrong?

Investigation of the infant's scratch revealed that the disposable cover had been cut crosswise in half, exposing the infant to a **small** amount of dried adhesive where the cable was connected to the pad. Fortunately, the injury was minor.

### What precautions can you take?

Follow these tips to ensure that you use the biliblanket safely:

- Inspect the biliblanket phototherapy system to make sure the fiber-optic pad is smooth and free from flaws. Make sure the cable and illuminator work properly.
- Check that the disposable cover completely protects the fiber-optic pad, including the area where the cable is connected to the pad. Don't cut or alter the cover.
- Secure the disposable cover around the pad cable with self-adhesive tabs.
- Expose as much of the infant's **skin** as possible to the illuminating side of the fiber-optic pad.
- Remove everything except the disposable cover from between the infant's **skin** and the light pad. The infant may wear a diaper, clothing, or a blanket over the pad.
- Shield the infant's eyes from direct exposure to the light with an opaque eye mask. ▀

Although you need to support the adverse event-reporting policy of your health care facility, you may voluntarily report a medical device that doesn't perform as intended by calling MedWatch at 1-800-FDA-1088 (fax: 1-800-FDA-0178). The opinions and statements contained in this report are those of the author and may not reflect the views of the Department of Health and Human Services. *Device Errors* is coordinated by Chris Parmentier, RN.