

DNA Preparation from Paraffin Tissue

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Chloroform

Mallinckrodt, Cat. 4440

EDTA, 0.5 M

Ethanol, absolute

Isoamyl alcohol

Sigma, Cat. I-3643

Phenol

Proteinase K

RNase A

Boehringer, Cat. 109 169

Sodium acetate, pH 5.2

Sodium chloride, 5 M

Sodium thiocyanate (NaSCN), 1 M

Sigma, Cat. S 7757, 250 g

TE Buffer (Tris-EDTA), pH 7.4

Tween 20

Xylene

Preparation

Chloroform/Isoamyl alcohol 24:1

Chloroform 24 ml

Isoamyl alcohol 1 ml

DNA extraction buffer

1.5 ml 5M NaCl

5.0 ml 0.5M EDTA

0.5 ml Tween 20

Fill up to 100 ml with sterile water.

Proteinase K (10 mg/ml)

Dissolve 100 mg Proteinase K in 10 ml TE for 30 min at RT.

Aliquot and store at -20° .

RNase A (20 mg/ml)

Dissolve 200 mg RNase A in 5 ml sterile water.

Boil for 15 min.

Cool to room temperature

Aliquot and store at -20° .

Procedure

1. Cut 50 μm slices of formalin-fixed and paraffin-embedded tumor samples.
(Note: Cut 4 μm slices before and after each 50 μm slice for Haematoxylin-Eosin staining to insure that the tissue is still representative.)
2. Incubate in xylene at 45°C for 15 min.
3. Centrifuge 10 min at 14,000 rpm.
4. Pipet off supernatant.
5. Repeat steps 2-4.
6. Add 1 ml 100% ethanol to tissue pellet, vortex, centrifuge for 10 min at 14,000 rpm, and pipet off supernatant.
7. Add 1 ml 90% ethanol to tissue pellet, vortex, centrifuge for 10 min at 14,000 rpm, and pipet off supernatant.
8. Add 1 ml 70% ethanol to tissue pellet, vortex, centrifuge for 10 min at 14,000 rpm, pipet off supernatant, and dry pellet in speed vac.
9. Resuspend pellet in 1 ml NaSCN (1 M) and incubate at 37°C overnight.
10. Centrifuge for 10 min at 14,000 rpm; pipet off supernatant.
11. Resuspend pellet in 1 ml of DNA extraction buffer.
12. Centrifuge for 10 min at 14,000 rpm, pipett off supernatent.
13. Repeat steps 11-12 twice.
14. Add 5 μl RNase (20 mg/ml) and incubate for 1 hr at 37°C (RNase treatment is optional, paraffin material often does not contain large amounts of RNA).

15. Add 40 μ l of Proteinase K (10 mg/ml), vortex briefly, and incubate at 55°C overnight (if tissue is not completely dissolved, add additional proteinase K and continue incubating; tissue should be dissolved).
16. Add 440 μ l of phenol, shake vigorously by hand for 5 min, and centrifuge for 5 min at 8000 rpm.
17. Pipet supernatant into a new tube, add a solution of 220 μ l phenol plus 220 μ l chloroform/isoamyl alcohol (24:1), shake vigorously by hand for 5 min, and centrifuge for 5 min at 8000 rpm.
18. Pipet supernatant into a new tube, add 440 μ l chloroform/isoamyl alcohol (24:1), shake vigorously by hand for 5 min, and centrifuge for 5 min at 8000 rpm.
19. Pipet supernatant into a new tube (2 ml eppendorf tube), add 1/10 volume of sodium acetate (pH 5.2), add 3 volumes of ice cold 100 % ethanol, and keep tube for 1 hr at -80°C or overnight at -20°C.
20. Centrifuge for 30 min at 14,000 rpm at 4°C.
21. Remove and save supernatant (optional, in case you are not sure of the precipitation).
22. Wash pellet in 70% ethanol, spin for 15 min at 14,000 rpm at 4°C, remove supernatant.
23. Dry pellet in speed vac.
24. Add 20-50 μ l sterile H₂O (depending on the amount of DNA you expect, which is subjective to experience).
25. Shake gently in thermomixer at 37°C for 2 hr (DNA should be dissolved, but if you have doubts put on a rotating shaker in the cold room overnight).
26. Measure DNA concentration with a spectrophotometer and run around 200 ng on a 1% agarose gel.