

Protocol for Amplification of RNA –Modified from Eberwine

Version 1.4 Nov. 2001

Reagents

- Superscript II RT, Gibco Life technologies # 18064-014
- *E. coli* DNA Polymerase I, Gibco Life technologies # 18010-025
- *E. coli* DNA Ligase, Gibco Life technologies # 18052-019
- *E. coli* RNase H, Gibco Life technologies # 18021-014
- T4 DNA Polymerase, Gibco Life technologies # 18005-025
- 5X second strand buffer, Gibco Life technologies # 10812-014
- dNTP Set, 10mM Gibco Life
- T7-(dT)24 Primer, GENSET Corp), HPLC purified DNA-5'-GGC-CAG-TGA-ATT-GTA-ATA-CGA-CTC-ACT-ATA-GGG-AGG-CGG-(dT)24-3'
- DEPC treated water, Ambion #9902 or Water, Molecular Biology Grade, BioWhittaker, #16-001Y
- RNAsin, Promega #N2515
- EDTA 0.5 M
- NaOH 1 M
- Phase Lock gel, 5 prime to 3 prime, Inc #pl-188233 for 200 or #pl-175850 for 50
- Phenol/Chloroform/isoamyl alcohol, Ambion #9732
- 7.5M Ammonium Acetate, Sigma #A2706
- 95% Ethanol
- 80% Ethanol
- Linear acrylamide, Ambion #9520
- T7 Megascript Kit, Ambion #1334
- Random hexamers, Gibco Life technologies #48190-011
- Qiagen Rneasy Mini Kit (50 samples) #74104

Check Before Starting

- Heating Block at 70°C
- 95% and 80% ethanol in -20°C
- Cold bath at 16°C

First Strand Synthesis

1. Primer hybridization:
 - 1-3 μg Total RNA
 - (Variable) Nuclease free H_2O (to complete 1st strand final reaction volume of 20 μl)
 - 1 μl T7-(dT)₂₄ primer (2 $\mu\text{g}/\mu\text{l}$)
2. Incubate at 70°C for 10 minutes
3. Quick spin and put on ice
4. Add:
 - 4 μl 5X First strand cDNA buffer
 - 2 μl 0.1M DTT
 - 2 μl 10mM dNTP mix
 - 1 μl RNAsin
5. Add 2 μl Superscript II
6. Mix well and incubate at 42°C for 1 hour
7. Centrifuge briefly, place on ice

Second Strand Synthesis

8. Add:
 - 91 μl DEPC treated H_2O
 - 30 μl Second strand buffer
 - 3 μl 10 mM dNTP mix
 - 4 μl DNA Polymerase I (10U/ μl)
 - 1 μl DNA Ligase (10U/ μl)
 - 1 μl RNase H (2U/ μl)Final Volume: 150 μl
9. Gently tap tube to mix, briefly centrifuge
10. Incubate at 16°C for 2 hours
11. Add 2 μl (10U) T4 DNA Polymerase
12. Cool for 5 minutes at 16°C
13. Stop reaction with 10 μl of 0.5 M EDTA, then add 10 μl of 1M NaOH
14. Incubate at 65°C for 10 minutes, then neutralize solution with 25 μl Tris-HCl (pH=7.5)

Clean Up of Double Stranded cDNA

15. Pellet the Phase Lock Gel in a microcentrifuge at maximum speed for 30 seconds
16. Add 198 μl (equal volume) of (25:24:1) Phenol : chloroform : isoamyl alcohol (saturated with 10 mM Tris-HCL pH 8.0/1mM EDTA) to the final DNA synthesis preparation (198 μl) to a final volume of 396 μl . Mix well by pipetting up & down vigorously.

17. Transfer the entire cDNA-phenol/chloroform mixture to the PLG tube
18. Do not vortex. Microcentrifuge at maximum speed for 2 minutes
19. Transfer the aqueous supernatant to a new 1.5 ml tube
20. Add 1 μ l linear acrylamide
21. Add 0.5 volumes of 7.5M Ammonium Acetate + 2.5 volumes (include the added Ammonium Acetate) of 95% ethanol stored at -20 to the sample and vortex. (may stop here, keep sample at -20°C and continue next day)
22. Centrifuge at maximum speed in a microcentrifuge at room temperature for 20 minutes
23. Remove supernatant. Wash pellet (may be invisible) with 0.5 ml of 80% ethanol
24. Centrifuge at maximum speed for 5 minutes at room temperature
25. Carefully pour-off the 80% ethanol.
26. Repeat the 80% ethanol wash once again
27. Air dry the pellet (~ 15 min.)
28. Resuspend the pellet in 16 μ l of Nuclease-free water

In Vitro Transcription

29. Ambion T7 Megascript kit- Follow manufacturer's instructions for a total 40 μ l reaction (double the 20 μ l standard reaction, 37°C for 4-5 hours)

SET UP REACTION AT ROOM TEMPERATURE

If done on ice, spermidine in the buffer can lead to template precipitation

- 16 μ l of template DS DNA (from step 28).
 - 4 μ l of 10x Reaction Buffer
 - 4 μ l of ATP solution (75mM T7)
 - 4 μ l of CTP solution (75mM T7)
 - 4 μ l of GTP solution (75mM T7)
 - 4 μ l of UTP solution (75mM T7)
 - 4 μ l of Enzyme Mix
 - Incubate at 37°C for 5-6 hours.
30. Add 60 μ l Nuclease free water to bring total volume up to 100 μ l.
 31. Follow "RNA clean-up" protocol in the Qiagen RNeasy mini handbook
 32. Elute with 30 μ l of Nuclease free water.
 33. Check O.D. and ratio

Expected yield: ~10 X starting amount of total RNA

34. Proceed using the Liu lab protocol for generating probe for microarrays using total RNA. Use 3 to 5 μ g amplified RNA and **random hexamer primers (2 μ l of 3 μ g/ μ l instead of oligo-dT primer)**

For 2nd Round Amplifications:

- Resuspend 0.5-1.0µg of amplified RNA in 11 µl ultrapure water

First Strand Synthesis

1. Add 1 µl Random hexamer (1 mg/ml)
2. Incubate 70°C for 10 minutes, then chill on ice
3. Equilibrate at room temperature for 10 minutes
4. Add:
 - 4 µl 5X First strand cDNA buffer
 - 2 µl 0.1M DTT
 - 2µl 10mM dNTP mix
 - 1µl RNAsin
5. Mix and incubate at 42°C for 2 minutes
6. Add 2 µl Superscript II
7. Mix well and incubate at 42°C for 1 hour
8. Add 1 µl RNase H and incubate at 37°C for 20 minutes
9. Heat to 95°C for 2 minutes and chill on ice

Second Strand Synthesis

10. Add 1 µl T7-oligodT primer (0.5 mg/ml), incubate 70°C for 5 minutes and at 42°C for 10 minutes
11. Then add:
 - 91 µl DEPC treated H₂O
 - 30 µl Second strand buffer
 - 3 µl 10 mM dNTP mix
 - 4 µl DNA Polymerase I (10U/ µl)
 - 1 µl DNA Ligase (10U/µl)
 - 1 µl RNase H (2U/µl)Final Volume: 150 µl
12. Gently tap tube to mix, briefly centrifuge
13. Incubate at 16°C for 2 hours
14. Add 2 µl (10U) T4 DNA Polymerase
15. Cool for 5 minutes at 16°C
16. Stop reaction with 10 µl of 0.5 M EDTA, then add 10 µl of 1M NaOH
17. Incubate at 65°C for 10 minutes, then neutralize solution with 25 µl Tris-HCl (pH=7.5)

Clean Up of Double Stranded cDNA and In Vitro Transcription: *Proceed as above*