RT-PCR REACTION

Section of Cancer Genomics, Genetics Branch, NCI National Institutes of Health

Reagents

cDNA (see cDNA Synthesis for RT-PCR Reaction Protocol)
Molecular Biology Grade Water
5μM Forward Primer
5μM Reverse Primer
Power SYBR GREEN PCR Master Mix [2x]
Applied Biosystems, Cat. 4367659

Materials and Equipment

RT-PCR instrument Applied Biosystems Sequence Detection System 7000 (ABI Prism 7000) **MicroAmp Optical Tubes** Applied Biosystems, Cat. N801-0533

OR

Optical Tubes (8 tubes/strip) Applied Biosystems, Cat. 4316567

OR

96-well Optical Reaction Plate with Barcode (code 128) Applied Biosystems, Cat. 4306737 **Optical Caps** Applied Biosystems, Cat. 4323032

Procedure

- 1. Program Sequence Detection System 7000 using Well Inspector.
- 2. Calculate volume of cDNA for 300 ng (\sim 3 µL) for each sample.
- 3. Calculate the number of tubes/wells planned for each gene. On ice, make Master master mix for each gene by combining the following reagents. For each tube/well, the volumes required are (includes 15% extra for pipetting error):
 - 14.3 µl Power SYBR GREEN PCR Master Mix [2x]
 - 1.44 µl Forward Primer
 - 1.44 µL Reverse Primer
 - 8.05 µl water

- 4. Aliquot 22 µl Master master mix into tubes/wells on ice.
- 5. Add 300 ng cDNA (~3 $\mu L)$ to each tube/well as appropriate and pipet up and down to mix.
- 6. Centrifuge tubes/plate for 1 min at 1200 rpm.
- 7. Load tubes/plate into Sequence Detection System 7000 and cap each tube/well with Optical Caps.
- 8. Set appropriate reaction conditions (number of cycles, temperatures, times, etc) and start reaction.