

# ProteinChip® Technology: A new tool for proteome analysis and biomarker discovery

**Kelly A. Conrads, Ph.D.**

Ciphergen Biosystems, Inc.  
kconrads@ciphergen.com




**CIPHERGEN®**

# Outline


---

- Background on technology – What is SELDI?
- Applications
- Example biomarker discover study – Ovarian Cancer


# ProteinChip® Platform: Components



ProteinChip Arrays  
» Separation




ProteinChip Reader  
» Detection




ProteinChip Software  
» Analysis


# ProteinChip® Arrays: A Variety of Surfaces



## Chemical Surfaces – Protein Expression Profiling:



**Hydrophobic**




**Anionic**



**Cationic**




**Metal Chelate**




**Normal Phase**


## Biological Surfaces – Protein Interaction Assays:



**PS-10 or PS-20**



**Antibody - Antigen**




**Receptor - Ligand**




**DNA - Protein**

# ProteinChip® Arrays: Retentate Chromatography




# ProteinChip® Proteomics: Complexity of Samples



## Tumor Tissue & Cultures

**Biopsy**  
**LCM**  
**Nuclear Extracts**  
**Cell lysates**  
**Plant Extracts**



## Tumor Fluid

**Seminal Plasma**  
**Nipple Aspirates**  
**Fine Needle**  
**Aspirates**



## Fluids

**Serum**  
**Plasma**  
**Urine**  
**CSF**  
**Blood eluates**

# ProteinChip® Array Preparation

**Step 1:** Complex protein sample is placed on a ProteinChip Array


- Proteins bind to chemical or biological sites on the ProteinChip surface

**Step 2:** Remove unbound proteins


- Wash the ProteinChip with appropriate stringency buffer
- Bound proteins are retained

**Step 3:** Add Energy Absorbing Molecules or "Matrix"



- EAM is applied to each spot to facilitate desorption and ionization in the TOF-MS Chip Reader



# Detection Using SELDI



- Retained proteins are “eluted” from the ProteinChip Array by laser desorption/ionization
- Ionized proteins are detected and their mass accurately determined by Time-of-Flight Mass Spectrometry
- ProteinChip® software produces map of proteins showing their accurate mass and relative ion intensity





# ProteinChip® Array Preparation: Sample Workflow

---

- **Sample deposit (2 min)**
- **Allow proteins to bind to surface (45-60 min)**
- **Wash (5-15 min)**
- **Matrix deposition (10 min)**
- **Chip introduction (1 min)**
- **Chip reading (15 min)**
- **Data analysis (30 min to ...)**

# Protein Discovery and Characterization

## ProteinChip<sup>®</sup> Advantages

---

- **Capture, detect and analyze proteins directly from crude biological samples without tagging or labeling**
- **Wash step – analysis of protein in MS unfriendly mixtures**
- **Lower limit of detection: 1-50 fmole of protein**
- **Small sample size: 0.5  $\mu$ l - 400  $\mu$ l**
- **Medium to high throughput: 8, 24, 96, or 192 samples**
- **Rapid results: 10 minutes**
- **Significant result: Accurate molecular weight  
(500 Da- 300 kDa)**

# ProteinChip® Applications

---

## Research Proteomics Basic Research

- Immunoassay development and screening
- Receptor-ligand assays
- Protein-protein interactions
- DNA-protein interactions

## Process Proteomics Protein Purification

- Purification development and monitoring
- Peptide mapping and protein identification
- Analysis of post-translational modifications
- Epitope mapping


## Clinical Proteomics Profiling

- Target discovery and validation
- Disease monitoring (drug efficacy studies)
- Toxicology
- Diagnostics


# Interaction Discovery Mapping™

## Protein-Protein Interactions


# ProteinChip® Arrays: A Variety of Surfaces




## Biological Surfaces – Protein Interaction Assays:



PS-10 or PS-20



Antibody - Antigen




Receptor - Ligand




DNA - Protein


# Protein-Protein Interaction Analysis




1. Preactivated ProteinChip




2. Bind 'Capture' Molecule




3. Block unused sites




4. Analyte Capture




5. Wash



6. Add EAM




7. Laser Ionization / Desorption



8. Direct Mass Detection


# Antibody Capture of IL8



Two different isoforms detected using antibody capture

**\*\*Advantage over ELISA\*\***

# Protein-Protein Interaction Analysis





# Protein Interactions: Beads & Arrays



## Arrays

- **Quick**
  - <4 hrs for typical capture experiment
- **High-throughput**
  - Up to 12 Arrays in parallel
- **Easily automated**
  - Simple to transfer protocols to robotics
- **Ability to work with small volumes**
  - 1-2 ul sample can be applied directly to the spot



## Beads

- **High capacity**
  - Large Dynamic Range
  - Ability to work with larger volumes
- **Ability to multiplex**
  - More than one antibody can be immobilized
- **Amenable to Orthogonal sample processing**
  - Coupling of affinity chromatography with subsequent purification steps
- **Semi-preparative Method**
  - Ability to ID captured proteins

# Process Proteomics Applications

## Protein Purification

# Expression Optimization of a His-tagged Protein




**Objective:** Determine the best conditions of expression of a recombinant 23.8 kDa His-tagged protein in a bacterial host strain

**Conditions tested:** 1. Strain, 2. Expression vector construct, 3. Temperature.

**ProteinChip Array:** Selective retention of the His-tagged target is found on IMAC-Cu+2

**Results:** C-Terminal vector construct in Strain 2 at 37°C attained the highest expression levels. Time spent for monitoring expression in 15 fermentation samples: Half a day.

# IMAC-Ni ProteinChip® Array



- Various amounts of imidazole were added to the binding step to determine what amount would improve the specific capture
- Protein binds to the IMAC-Ni surface, seen at 48 kDa
- The presence of low amounts of imidazole (10 mM) doesn't interfere with binding
- There are a significant number of endogenous impurity proteins binding to the IMAC surface
- The number of impurity proteins decreases as the amount of imidazole increases
- 8.5-85 kDa MW region is shown

# Process Proteomics: Key Benefits

---




- **A versatile technology**
  - Analytical capabilities for crude fermentation/CC samples through purification and final formulation
- **Significantly reduce analysis/purification time**
  - Rapid, 'on-chip' analytical and purification method development
  - Small sample requirement allows development to begin earlier in the process
  - Reduce number of SDS-PAGE, Westerns, ELISAs, HPLC runs
- **Get more information from the sample**
  - Accurate MW tracking of target and impurities

# Expression Difference Mapping™ Applications

## Biomarker Discovery

# Pattern Track™ Process



# Pattern Track™ Process

## Discovery

**Discovery Pilot Study: biomarker scouting study to find biomarker candidates**

**Multiple Biomarker Candidates**



**Validation Study: biomarker study using optimal conditions from Pilot Study to validate biomarker candidates**

**Validated Biomarkers**



- Maximize number of conditions to discover biomarkers

- Univariate analysis to identify biomarker candidates

- Minimize number of conditions to validate biomarkers

- Multivariate analysis to classify patient groups based on multiple biomarker patterns

- Biomarker Patterns™ Software selects the best biomarker candidates

## Purification

**Purified Biomarkers**



- SELDI Assisted Purification

## Identification

**Identified Biomarkers**



- SEND ID on-chip peptide mapping

- SELDI Tandem Mass Spectrometry


## Assay

**Quantitative Predictive Multi Biomarker Pattern Assays**

- ProteinChip System, Series 4000




# Fractionation



# Benefits of Fractionation

- Fractionation expands the investigated part of the proteome:
- Neat = 48 peaks
- 6 Fractions = 180 unique peaks, or 132 new peaks.



# Discovery: Surface Specificity

H4


IMAC-Cu

WCX, pH 4

WCX, pH 7


SAX, pH 9

SAX, pH 5




# Discovery: Visualizing Biomarkers


## Spectral View




## Gel View



## Overlay View





# 7.8 kDa Peak is Differentially expressed







# 7.8 kDa Peak is Differentially expressed

**7.8 kDa peak discovered on IMAC30-Cu<sup>2+</sup> array, Low Mass**

**p-value < 0.001  
AUROC = 0.947**





# Principle Component Analysis







# Multiple Biomarker Diagnostic Assays




ID the biomarkers,  
Link to biology of  
disease

	Clinical Dx Sensitivity "True Positives"	Clinical Dx Specificity "True Negatives"
<b>Single Biomarker</b>	<b>65%</b>	<b>35%</b>
<b>Multi-Biomarker Pattern</b>	<b>&gt;90%</b>	<b>&gt;90%</b>

Improve Predictive Accuracy and Diagnostic Value

# Host Response Protein Amplification Cascade: Common Proteins Yield Uncommon Fragments

- **Apolipoprotein A1**
- **Transthyretin**
- **Inter alpha-trypsin inhibitor 4**
- **Haptoglobin a**
- **Serum amyloid A**
- **Vitamin D Binding Protein**
- **C3 anaphylotoxin**



Specific Disease Processes Cleave and Modify Common Proteins into Uncommon Fragments with Diagnostic Utility

# Common Proteins



# Disease Processes




# Specific Diagnostic Fragments

Inter Alpha Trypsin Inhibitor 4 (ITIH4)



## Pancreatic Cancer

nvhsgstffkyyylqgakipkpeasfspr



## Ovarian Cancer

mnfrpgvlssrqlglpgppdvpdhaayhpfr



## Diabetes

srqlglgppdvpdhaayhpfr

# Biomarkers found using SELDI technology

- **Cancer**
  - Breast
  - Prostate
  - Bladder
  - Leukemia
  - Lung
  - Brain
- **Pharmaceutical testing**
  - Toxicity markers
  - Non-responders
- **Other**
  - Acute renal failure
  - Acute heart failure
  - Exposure to airborne toxins
- **Neuropsychiatric**
  - Depression
  - Schizophrenia
  - Alzheimer's disease
  - Parkinson's
  - Huntington's
- **Infectious diseases**
  - Yersinia pestis
  - Mycobacterium
  - Caulobacter
  - Streptococcus
  - Botulism
  - Prions

# A Novel Multi-Marker Assay for Ovarian Cancer

Zhang, Z. *et al.* (2004) *Cancer Research* **64**, 5882-5890.



CIPHERGEN®

# Study participants

**Johns Hopkins University:** Zhen Zhang, Ph.D., Jinong Li, Ph.D., Lori J. Sokoll, Ph.D., Alex J. Rai, Ph.D., Jason M. Rosenzweig, Bonnie Cameron, Daniel W. Chan, Ph.D.

**MD Anderson:** Robert C. Bast Jr., M.D., Yinhua Yu, M.D.

**Duke:** Andrew Berchuck, M.D.

**Royal Hospital for Women (Sydney):** Carolien van Haaften-Day, Ph.D., Neville F. Hacker, M.D.

**Groningen University Hospital:** Henk W. A. de Bruijn, Ph.D., Ate G. J. van der Zee, M.D.

**Bart's and The London, Queen Mary School of Medicine, London University:** Ian J. Jacobs, M.D.

**Ciphergen:** Xiao-Ying Meng, M.Sc., Eric T. Fung, M.D., Ph.D.



**CIPHERGEN**<sup>®</sup>

# Basic statistics of ovarian cancer

---

- **Prevalence 40/100,000 (1 in 2500)**
- **23,000 new cases diagnosed annually**
- **14,000 deaths annually**
- **Overall 5 year survival 20-30%**
- **75% of cases are diagnosed in late stage (stage III/IV)**
- **90% cure rate in stage I/IIa**
- **Therefore, detection in earlier stages critical in improving overall survival**


# Diagnosing early stage ovarian cancer

---

- **Low prevalence of disease places extremely high requirements on screening test**
- **CA-125 is the only generally accepted tumor marker for ovarian cancer**
  - Elevated in ~80% of all patients with ovarian cancer
  - Elevated in <50% of patients with early stage ovarian cancer
  - Not useful for screening
  - Used typically for:
    - *Initial diagnosis in symptomatic women in conjunction with other tests (e.g., transvaginal ultrasound, computed tomography)*
    - *Monitoring of response to therapy*



# Discovery: Expression Differences




# Methodology Outline

---

- **Protein expression profiling for all samples was performed simultaneously**
- **Protocol: anion exchange fractionation of serum followed by binding of fractions to ProteinChip arrays**
- **Each sample was bound to each array in triplicate**
- **Instrument performance was monitored using**
  - Insulin (test for instrument resolution)
  - Immunoglobulin (test for instrument sensitivity)
- **Assay performance was monitored by processing a standard serum sample**

# Fractionation



# Data mining

---


- **Data pre-processing**

- Mass calibration, baseline subtraction, total ion current normalization
- Peak selection – automated followed by manual confirmation


- **Data analysis**

- Unified Maximum Separability Analysis (UMSA): Modified support vector machine algorithm
- Bootstrap re-sampling technique on site A and site B data separately to determine best set of peaks from independent analysis
- Biomarkers for inclusion are those that
  - *Are deemed significant repeatedly during the bootstrap process*
  - *Change in the same direction between both sites*
- Validate biomarkers on independent data set


# Study design for biomarker discovery



# Reproducibility



Mass Range = 2-100 kDa




Mass Range = 35-62 kDa

Average CV's range from 10-25%

# Performance of individual biomarkers in discovery set & validation set

<b>Comparison</b>	<b>All epithelial ovarian cancer</b>	<b>Stage I/II invasive ovarian cancer</b>
<b>CA125, Discovery set</b>	<b>&lt;.00001</b>	<b>&lt;.00001</b>
<b>CA125, Validation set</b>	<b>.00001</b>	<b>&lt;.00001</b>
<b>Marker 1, Discovery set</b>	<b>&lt;.00001</b>	<b>&lt;.00001</b>
<b>Marker 1, Validation set</b>	<b>&lt;.00001</b>	<b>&lt;.00001</b>
<b>Marker 2, Discovery set</b>	<b>.00002</b>	<b>.00004</b>
<b>Marker 2, Validation set</b>	<b>&lt;.00001</b>	<b>&lt;.00001</b>
<b>Marker 3, Discovery set</b>	<b>&lt;.00001</b>	<b>.059178</b>
<b>Marker 3, Validation set</b>	<b>&lt;.00001</b>	<b>.079999</b>

P value estimated for two-group T test between healthy controls and ovarian cancer patient population



**CA125 ROC curve, area=0.76961, std = 0.054821,**

**Model w/ CA125 ROC curve, area=0.88526, std = 0.037679,**

**Model w/o CA125 ROC curve, area=0.9195, std = 0.032696,**




# Purification and identification



# ProteinChip® Assisted purification and ID

- **Spin column chromatography, with monitoring of fractions using the ProteinChip® Biomarker System**
- **Protease digestion**
- **Matching of peptide fingerprint with database**
- **Confirmation using the CIPHERGEN ProteinChip Interface installed on the ABI/Sciex QSTAR™ tandem MS**



# Apolipoprotein A1

---

- **Major lipoprotein in HDL's**
- **Negative acute phase reactant protein**
- **Total apolipoprotein A has previously been reported to be decreased in patients with ovarian cancer (Kuesel et al, 1992)**

# Transthyretin

---

- **Marker identified by profiling is an N-terminal truncation of transthyretin (pre-albumin)**
- **Transthyretin has been reported to be decreased in ovarian cancer patients previously (Mahlck et al, 1994)**
- **Biological function**
  - Homotetramer
  - Transports thyroid hormone, retinol binding protein (RBP)
  - Decreases in pre-albumin lead to decreases in RBP
  - Decreased levels of RBP leads to malignant transformation of ovarian epithelial cells (in culture)
- **No reports of truncated forms**

# Peptide fragment of inter alpha trypsin inhibitor IV (ITIH4)

---


- **ITIH4 Known to be cleaved by Kallikrein enzymes**
- **Kallikrein enzyme levels known to change in Ovarian Cancer**
- **Hypothesis: kallikreins, or other cancer-related proteases, may be responsible for generating fragments found in this study**
- **Test hypothesis using IDM platform**

# Specificity of marker panel

	Ovarian cancer	Breast cancer	Prostate cancer	Colon cancer
CA125	0.000000	0.39	0.07	0.86
Apolipoprotein A1 (intact)	0.0004	0.77	0.02	0.69
Transthyretin (intact)	0.00005	0.51	0.22	0.01


Additional specificity conferred by:  
Use of combination of markers  
Nature of cleavage products

# Interacting Proteins




# Secondary interaction discovery


*Confirm existing interactors, discover new ones*



## Analysis on ProteinChip® Array



Low mass range



High mass range




# Ongoing research

---

- **Expansion of sample sets (>1500 total)**
  - University of Leuven (prognosis and treatment response)
  - Groningen (prognosis and treatment response)
  - Bart's College (screening study)
  - Mayo Clinic/National Cancer Institute (initial diagnosis)
- **Validation of existing marker set**
- **Combination of markers, novel and previously characterized**
- **Discovery of additional markers by**
  - New generation protein expression profiling protocol
  - Interaction mapping
  - Kallikrein substrate experiments
- **Development of multiplexed ProteinChip assays for markers**

# Pattern Track™ Process



# ProteinChip® System, Series 4000

*Personal edition*



*Enterprise edition*



# ProteinChip System, Series 4000

<b>Superior Quantitation</b>	<ul style="list-style-type: none"><li>■ Raster laser design for maximum spot coverage</li><li>■ Auto laser energy setting</li><li>■ Extended Linear Dynamic Range</li><li>■ Improved ProteinChip Arrays and protocols</li></ul>
<b>Highest Sensitivity</b>	<ul style="list-style-type: none"><li>■ New detector</li><li>■ New Ion Source increases ion efficiency</li><li>■ Detector blanking reduces noise</li><li>■ Innovative flight tube design</li></ul>
<b>Highest Resolution/Resolving power</b>	<ul style="list-style-type: none"><li>■ Improved fractionation tools increase resolution up to 3000 proteins</li></ul>
<b>Highest Throughput</b>	<ul style="list-style-type: none"><li>■ Unattended runs of up to 168 ProteinChip Arrays</li></ul>
<b>Improved <i>low abundance</i> biomarker discovery</b>	<ul style="list-style-type: none"><li>■ New Deep Proteome™ Tools to explore the last 1% of the serum proteome</li></ul>

# Improved Resolution/Resolving power

*MultiSelect™ Fractionation*

**Load Sample**



**Depletion Chemistries**

**High Specificity Chemistries**



**Low Specificity,  
High Capacity Chemistries**

**Disassembly Chemistries**



# Improved sensitivity with multiple chemistry approach


## Multiple Chemistry Pre-fractionation

100 fmol/ul


10 fmol/ul

1 fmol/ul

Control




## Single Chemistry Pre-fractionation



# Protein Equalizer™ Beads

*Affinity Bioseparation of Proteins via Combinatorial Ligands*

## Define Equalizer Bead Library



- Each bead binds a unique protein
- Capacity of every bead is the same

Bind

Wash away Excess

Elute Equalized Sample

4 Yellow

4 Blue


3 Green

1 Red

1 Pink

# Effect of Protein Equalizer™ Beads on serum proteins

Unfractionated Serum profiled Before and After processing on Protein Equalizer Beads











# SELDI-MS Platform Protein Identification Using ProteinChip Interface

1. Enrich the proteins  
on BioSeptra resins and  
digest with trypsin



Protein Fragments  
obtained after trypsin  
digestion

2. Capture the peptides  
on NP-20 ProteinChip Array




3. MS-MS analysis on a  
Micromass Q-TOF II  
outfitted with CIPHERGEN  
ProteinChip Interface



MS-MS data with  
Ion assignments

# Protein Identification (Mascot Search)

**Search title** : 25 kDa  
**Database** : SwissProt 41.0 (206209 sequences; 113194615 residues)  
**Taxonomy** : Rattus (6684 sequences)  
**Timestamp** : 21 Mar 2003 at 13:53:43 GMT  
**Top Score** : 56 for **P08009**, Glutathione S-transferase Yb3 (EC 2.5.1.18) (Chain 4) (GST class-mu 3)



1.	<a href="#">P08009</a>	Mass: 25763	Total score: <b>56</b>	Peptides matched: 10
	Glutathione S-transferase Yb3 (EC 2.5.1.18) (Chain 4) (GST class-mu 3)			
2.	<a href="#">P08010</a>	Mass: 25785	Total score: <b>51</b>	Peptides matched: 10
	Glutathione S-transferase Yb2 (EC 2.5.1.18) (Chain 4) (GST class-mu 2) ST			
	<a href="#">P08010-00-00-01</a>	Mass: 25686	Total score: 45	Peptides matched: 9
	Glutathione S-transferase Yb2 (EC 2.5.1.18) (Chain 4) (GST class-mu 2) ST			
7.	<a href="#">P04905</a>	Mass: 25996	Total score: 34	Peptides matched: 8
	Glutathione S-transferase Yb1 (EC 2.5.1.18) (Chain 3) (GST M1-1) (GST class-mu 1)			
	<a href="#">P04905-00-00-01</a>	Mass: 25997	Total score: 34	Peptides matched: 8
	Glutathione S-transferase Yb1 (EC 2.5.1.18) (Chain 3) (GST M1-1) (GST class-mu 1)			
	<a href="#">P04905-00-00-02</a>	Mass: 26069	Total score: 27	Peptides matched: 7
	Glutathione S-transferase Yb1 (EC 2.5.1.18) (Chain 3) (GST M1-1) (GST class-mu 1)			

Match to: **P08010**; Score: 51  
 Glutathione S-transferase Yb2 (EC 2.5.1.18) (Chain 4) (GST class-mu 2)

Nominal mass ( $M_r$ ): 25785; Calculated pI value: 7.30

NCBI BLAST search of [P08010](#) against nr  
 Unformatted [sequence string](#) for pasting into other applications

Taxonomy: [Rattus norvegicus](#)

Fixed modifications: Propionamide (C)  
 Cleavage by Trypsin: cuts C-term side of KR unless next residue is P  
 Number of mass values searched: 46  
 Number of mass values matched: 10  
 Sequence Coverage: **43%**

Matched peptides shown in **Bold Red**

```

1  PMLTGYWDIR GLAHAIRLFL EYTDTSYEDK KYSMGDAPDY DRSQWLSEKF
51  KLGLDFFNLP YLIDGSHKIT QSNAILRYLG RKHNLGGETE EERIRVDVLE
101 NQAMDTRLQL AMVCYSPDFE RKKPEYLEGL PEKMKLYSEF LGKQPFWAGN
151 KITYVDFLVY DVLDQHRIFE PKCLDAFPNL KDFVARFEGF KKISDYMKSG
201 RFLSKPIFAK MAFWNPK
    
```

Match to: **P08009**; Score: 56  
 Glutathione S-transferase Yb3 (EC 2.5.1.18) (Chain 4) (GST class-mu 3)

Nominal mass ( $M_r$ ): 25763; Calculated pI value: 7.27

NCBI BLAST search of [P08009](#) against nr  
 Unformatted [sequence string](#) for pasting into other applications

Taxonomy: [Rattus norvegicus](#)

Fixed modifications: Propionamide (C)  
 Cleavage by Trypsin: cuts C-term side of KR unless next residue is P  
 Number of mass values searched: 46  
 Number of mass values matched: 10  
 Sequence Coverage: **48%**


Matched peptides shown in **Bold Red**

```

1  PMLTGYWDIR GLAHAIRLLL EYTDSSYEK RYTMGDAPDF DRSQWLNEKF
51  KLGLDFFNLP YLIDGSHKIT QSNAILRYLG RKHNLGGETE EERIRVDVLE
101 NQLMDNRMLV ARLCYNPDFE KLPKGYLEQL PGMRLYSEF LGKRPWFAGD
151 KITFVDFIAY DVLERNQVFE ATCLDAFPNL KDFIARFEGF KKISDYMKSS
201 RFLPRPLFTK MAIWGSK
    
```


# Fractionation - scheme

## ■ Employed fractionation scheme:




# Results – Fractions

- **Fractionation expands the investigated part of the proteome:**
- **Neat = 48 peaks**
- **6 Fractions = 180 unique peaks, or 132 new peaks.**



# The ProteinChip® Bioprocessor



# SELDI Citations

Over 200 SELDI Papers:  
Nature, Science, The Lancet.....

<http://www.ciphergen.com/pub/searchPub.asp>

CIPHERGEN®			
HOME	PRINTER VIEW	SEARCH: <input type="text"/>	GO
+1-888-864-3770		+1-510-505-2100	
About Us	Products & Services	Sales & Support	Technology & Apps
Publications	Investors		
Publications: <a href="#">Search database</a>			
Search Again			
190 total documents returned: <a href="#">ProteinChip Citations</a> (✓ denotes recommended readings)			
ProteinChip Citations			
<a href="#">view</a>	190 <b>A data-analytic strategy for protein biomarker discovery: profiling of high-dimensional proteomic data for cancer detection</b>	Yasui, Y., Pepe, M., Thompson, M.L., Adam, B.-L., Whight, G.L., Qu, Y., Potter, J.D., Winget, M., Thornquist, M., Feng, Z.	2003
<a href="#">view</a>	189 <b>Haptoglobin-alpha Subunit As Potential Serum Biomarker in Ovarian Cancer: Identification and Characterization Using Proteomic Profiling and Mass Spectrometry</b>	Ye, B., Cramer, D.W., Skates, S.J., Gygi, S.P., Pratomo, V., Fu, L., Horick, N.K., Licklider, L.J., Schorge, J.O., Berkowitz, R.S., Mok, S.C.	2003
<a href="#">view</a>	188 <b>A panel of cerebrospinal fluid potential biomarkers for the diagnosis of Alzheimer's disease</b>	Carrette, O., Demalle, I., Scherl, A., Valkinoglu, O., Cortinals, G., Burkhard, P., Hochstrasser, D.F., Sanchez, J.-C.	2003
<a href="#">view</a>	187 <b>Stereoselective Inactivation of <i>Torpedo californica</i> Acetylcholinesterase by Isomaltation: Inhibitory Reactions with (1R)- and (1S)-Isomers Proceed by Different Mechanisms</b>	Doorn, J.A., Thompson, C.M., Christner, R.B., Richardson, R.J.	2003
<a href="#">view</a>	186 <b>Estradiol Control of Expression and Levels of Estradiol-Binding Proteins in the Medial Preoptic Area, Medial Hypothalamus and Pituitary</b>	Gao, G., Herbert, Z., Kong, J., Gabrielson, N., Mautz, A., Wu, D., Jankowski, G.F., Caldwell, J.D.	2003
<a href="#">view</a>	185 <b>Modern Tumor Marker Discovery in Urology: Surface Enhanced Laser Desorption and Ionization (SELDI)</b>	Gretzer, M.B., Parin, A.W., Chan, D.W., Veltri, R.V.	2003
<a href="#">view</a>	184 <b>Putative protein markers in the sera of men with prostatic neoplasms</b>	Lehrer, S., Roboz, J., Ding, H., Zhao, S., Diamond, E.J., Holland, J.F., Stone, N.N., Drollier, M.J., Stock, R.G.	2003
<a href="#">view</a>	183 <b>Extraction of oxytocin and arginine-vasopressin from serum and plasma for radioimmunoassay and surface-enhanced laser desorption/ionization time-of-flight mass spectrometry</b>	Cool, D.R., DeRosier, D.	2003
<a href="#">view</a>	182 <b>Presenilin-1, Nicastrin, Amyloid Precursor Protein, and gamma-Secretase Activity Are Co-localized in the Lysosomal Membrane</b>	Pasternak, S.H., Bagshaw, R.D., Gural, M., Zhang, S., Ackersley, C.A., Pak, B.J., Callahan, J.W., Mahuran, D.J.	2003
<a href="#">view</a>	181 <b>beta-Secretase Cleavage at Amino Acid Residue 31 in the Amyloid beta Peptide Is Dependent upon gamma-Secretase Activity</b>	Shi, X.-P., Tugusheva, K., Bruce, J.E., Lucka, A., Wu, G.-X., Chen-Dodson, E., Price, E., Li, Y., Xu, M., Huang, Q., Sardana, M.K., Hazuda, D.J.	2003
<a href="#">view</a>	180 <b>Use of ProteinChip™ Array Surface Enhanced Laser Desorption/Ionization Time-of-Flight Mass Spectrometry (SELDI-TOF MS) to Identify Thymosin Beta-4, A Differentially Secreted Protein from Lymphoblastoid Cell Lines</b>	Diamond, D.L., Zhang, Y., Gaiger, A., Smithgall, M., Vedvick, T.S., Carter, D.	2003
<a href="#">view</a>	179 <b>Analysis of complex autoantibody repertoires by surface-enhanced laser desorption/ionization-time of flight mass spectrometry</b>	Grus, F.H., Joachim, S.C., Pfeiffer, N.	2003
<a href="#">view</a>	178 <b>Applications of high-throughput methods to cancer metastases</b>	Keller, E.T., Yao, Z.	2002
<a href="#">view</a>	177 <b>Macrophage proteomic fingerprinting predicts HIV-1-associated cognitive impairment</b>	Luo, X., Carlson, K.A., Wojna, V., Mayo, R., Biskup, T.M., Stoner, J., Anderson, J., Gendelman, H.E., Meléndez, L.M.	2003



**AACR** American Association  
for Cancer Research

## 95th Annual Meeting

March 27 - 31, 2004  
Orange County Convention Center  
Orlando, Florida

Over 51 papers featuring SELDI at the 2004 AACR meeting (NCI, JHU, MD Anderson, Stanford etc. covering biomarker applications for all majors cancers



# 20 Million Possible Diagnostic Patterns from 12 Common Proteins

- **Apolipoprotein A1**
- **Transthyretin**
- **Inter alpha-trypsin inhibitor 4**
- **Haptoglobin**
- **Serum amyloid A**
- **Vitamin D Binding Protein**
- **C3 anaphylotoxin**
- **Apolipoprotein A4**
- **C4 Anaphylatoxin Des-Arg**
- **Alpha-1 antichymotrypsin**
- **Apolipoprotein C1**
- **Human Serum Albumin**

12 Most Common Proteins

X 5 Fragments (estimate)

---

- **Fragments !**


60x59x58x57...

20 Million Possible  
Diagnostic Patterns


# Individual markers have low correlation to each other

	CA125	Marker 1	Marker 2
CA125	1		
Marker 1	-0.26	1	
Marker 2	-.015	.32	1
Marker 3	.05	-0.42	-.07

# Flowchart of data analysis




**Preprocess**  
 Baseline subtract  
 Normalize  
 Calibrate



## Univariate analysis

Biomarker Wizard or  
 CIPHERGENExpress™  
 Cluster Analysis,  
 Detect peaks,  
 Calculate P-values,  
 ROC/AUC



## Multivariate analysis

CIPHERGENExpress  
 Hierarchical  
 Clustering, PCA;  
 Biomarker Patterns  
 Software

