



WIENER LABORATORIOS S.A.I.C. - Riobamba 2944 - 2000 Rosario - Argentina Phone +54 (341) 432-9191/6 - Fax +54 (341) 432-5454/5555 Internet: http://www.wiener-lab.com.ar

# **Section 6 – Summary**

# 510(k) Summary "This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21CFR 807.92"

"The assigned 510(k) number is:  $\frac{102/334}{}$ Introduction According to the requirements of 21 CFR 862.1145, the following information provides sufficient details to understand the basis of a determination of substantial equivalence. 6-1 Submitter Wiener Lab Group Name, Address, Riobamba 2944 2000 - Rosario - Argentina Contact Contact person: Viviana Cétola Date Prepared: July, 2002 6-2 Device Name Proprietary name: WIENER LAB. CA-COLOR ARSENAZO III AA Common name: Calcium test system. Classification name: Azo Dye, Calcium Device Class II 6-3 Predicate We claim substantial equivalence to the currently marketed WIENER LAB. CA-COLOR AA (Cat. Nº 1152002) kit. **Device** 

### 6-4 Device Description

Calcium reacts with Arsenazo III, yielding a blue colored complex, which is photocolorimetrically measured at 650 nm. 8-hvdroxvauinoline is added to remove magnesium interference.

#### 6-5 Intended Use

The WIENER LAB CA-COLOR ARSENAZO III AA test system is a quantitative in vitro diagnostic device intended to be used in the quantitative determination of calcium in human sera, heparinized plasmas and urine on both manual and automated systems. Measurements of calcium are used in the diagnosis and treatment of parathyroid diseases, a variety of bone diseases, chronic renal diseases and tetany (intermittent muscular contractions or spasms).

# and Differences

6-6 Equivalencies The WIENER LAB. CA-COLOR ARSENAZO III AA test system is substantially equivalent to other products in commercial distribution intended for similar use. Most notably it is substantially equivalent to the currently marketed WIENER LAB. CA-COLOR AA test system.

> The following table illustrates the similarities and differences between the WIENER LAB. CA-COLOR ARSENAZO III AA test system and the currently marketed WIENER LAB. CA-COLOR AA test system.

	CA-COLOR AA	CA-COLOR ARSENAZO III AA		
Intended use	Quantitative determination of calcium in human serum, heparinized plasma and urine.			
		Continued on next page		

	CA-COLOR AA	CA-COLOR ARSENAZO III AA			
Test principle	Calcium reacts with o- Cresolphtalein complexone (o-CPC) at pH 10.8, yielding a purple colored complex, which is photocolorimetrically measured at 570 nm. 8-hydroxyquinoline is added to remove magnesium interference	Calcium reacts with Arsenazo III, yielding a blue colored complex, which is photocolorimetrically measured at 650 nm. 8-hydroxyquinoline is added to remove magnesium interference			
Essential Components	o-CPC 8-hydroxyquinoline	Arsenazo III 8-hydroxyquinoline			
Reagent Storage	Room temperature	2 – 10°C			
Reagent Deterioration	Reagent Blank > 0.400 O.D.	Reagent turbidity Reagent Blank > 0.800 O.D.			
Preparation of Working Reagent	Mixture of R1 and R2 (1:1) or they can be used separately.	None			
Working Reagent Stability	Stable 4 days at 2-10°C	Same as kit stability			
Precautions	All glassware should be cleaned with diluted hydrochloric acid and rinsed with distilled water.				
Working Temperatures	Room temperature - 37°C				
Wavelength of reading.	560 – 590 nm	620 – 650 nm			
Continued on next page					

	CA-COLOR AA	CA-COLOR ARSENAZO III AA		
Linearity	20 mg/dl			
Expected values	Serum: 8.5-10.5 mg/dl Urine: 60-200 mg/24hr			
Within-run precision	Normal Level Serum: CV = 1.28%	Normal Level Serum: CV = 1.93%		
	High Level Serum: CV = 1.30%	High Level Serum: CV = 0.88%		
	Normal Level Urine CV = 1.06%	Normal Level Urine CV = 2.57%		
	High Level Urine CV = 0.68% High Level Urine CV = 1.79%			
Run-to-run precision	Normal Level Serum: CV = 1.74%	Normal Level Serum: CV = 1.74%		
	High Level Serum: CV = 1.70%	High Level Serum: CV = 1.29%		
	Normal Level Urine CV = 2.50%	Normal Level Urine CV = 2.44%		
	High Level Urine CV = 1.34%	High Level Urine CV = 2.62%		

## 6-7 Conclusion

Based on the data above mentioned, we believe that the extended claims continue to support substantial equivalence to other products in commercial distribution intended for similar use.



Food and Drug Administration 2098 Gaither Road Rockville MD 20850

SEP 1 0 2002

Dr. Viviana Cetola QC/QA Manager Weiner Laboratorios S.A.I.C. Riobamba 2944 Rosario, Santa Fe Argentina

Re: k021334

Trade/Device Name: Ca-Color Arsenazo III AA

Regulation Number: 21 CFR 862.1145 Regulation Name: Calcium test system

Regulatory Class: Class II

Product Code: CJY
Dated: August 1, 2002
Received: August 13, 2002

Dear Dr. Cetola:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsma/dsmamain.html".

Sincerely yours,

Steven I. Gutman, M.D., M.B.A.

Director

Division of Clinical Laboratory-Devices

Steven Butman

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known):	1334			•
Device Name: Wiener lab.				
Ca-color, Arsen	iazo III A	4		
Indications For Use:				FOA)
		·		FRA/CORH/OUT/
The "Wiener lab. Ca-Color quantitative in vitro diagnos quantitative determination of plasmas and urine on be Measurements of calcium are parathyroid diseases, a variable diseases and tetany (intermitted).	stic device of calcium oth manua e used in the ariety of be tent muscula	intended to in human seal and autor he diagnosis a one diseases ar contractions	be used/ era, hepa nated sy and treati , chronic s or spasi	in the arinized ystems. ment of c renal ms).
Concurrence of CDR	UL, Office of D	evice Evaluation (	ODE)	-
		•		
Prescription Use V (Per 21 CFR 801.109)	OR	Over-The-Co	unter Use_,	
(Division Sign-Off) Division of Clinical Laboratory 510(k) Number	Devices	Opao	na i otinat i	,

54.34