

National Health and Nutrition Examination Survey 2003-2004

Documentation, Codebook, and Frequencies

MEC Laboratory Component: Biochemistry Profile

Survey Years:
2003 to 2004

SAS Export File:
L40_C.XPT



First Published: January 2006
Last Revised: October 2007

NHANES 2003–2004 Data Documentation

Laboratory Assessment: Lab 40 - Standard Biochemistry Profile

Years of Coverage: 2003–2004

First Published: January 2006

Last Revised: October 2007

The file was updated to correct the **Blood Urea Nitrogen (BUN)** method description.

Component Description

This battery of measurements are used in the diagnosis and treatment of certain liver, heart, and kidney diseases, acid-base imbalance in the respiratory and metabolic systems, other diseases involving lipid metabolism and various endocrine disorders as well as other metabolic or nutritional disorders.

1. Alanine Aminotransferase (ALT)

Alanine aminotransferase measurements are used in the diagnosis and treatment of certain liver diseases (e.g., viral hepatitis and cirrhosis) and heart diseases. Elevated levels of the transaminases can indicate myocardial infarction, hepatic disease, muscular dystrophy, or organ damage. Serum elevations of ALT activity are rarely observed except in parenchymal liver disease, since ALT is a more liver-specific enzyme than aspartate aminotransferase (AST).

2. Albumin

Albumin measurements are used in the diagnosis and treatment of numerous diseases primarily involving the liver or kidneys.

3. Alkaline Phosphatase (ALP)

Increased ALP activity is associated with two groups of diseases: those affecting liver function and those involving osteoblastic activity in the bones. In hepatic disease, an increase in ALP activity is generally accepted as an indication of biliary obstruction. An increase in serum phosphatase activity is associated with primary hyperparathyroidism, secondary hyperparathyroidism owing to chronic renal disease, rickets, and osteitis deformans juvenilia due to vitamin D deficiency and malabsorption or renal tubular dystrophies. Increased levels of ALP are also associated with Von Recklinghausen's disease with bone involvement and malignant infiltrations of bone. Low levels are associated with hyperthyroidism, and with the rare condition of idiopathic hypophosphatasia associated with rickets and the excretion of excess phosphatidyl ethanolamine in the urine.

4. Aspartate Aminotransferase (AST)

AST measurements are used in the diagnosis and treatment of certain types of liver and heart disease. Elevated levels of the transaminases can signal myocardial infarction, hepatic disease, muscular dystrophy, or organ damage.

5. Bicarbonate (HCO₃)

Together with pH determination, bicarbonate measurements are used in the diagnosis and treatment of numerous potentially serious disorders associated with acid-base imbalance in the respiratory and metabolic systems.

6. Blood Urea Nitrogen (BUN)

BUN measurements are used in the diagnosis of certain renal and metabolic diseases. The determination of serum urea nitrogen is the most widely used test for the evaluation of kidney function. The test is frequently requested in conjunction with the serum creatinine test for the differential diagnosis of prerenal, renal, and postrenal uremia. High BUN levels are associated with impaired renal function, increased protein catabolism, nephritis, intestinal obstruction, urinary obstruction, metallic poisoning, cardiac failure, peritonitis, dehydration, malignancy, pneumonia, surgical shock, Addison's disease, and uremia. Low BUN levels are associated with amyloidosis, acute liver disease, pregnancy, and nephrosis. Normal variations are observed according to a person's age and sex, the time of day, and diet, particularly protein intake.

7. Calcium

Elevated total serum calcium levels are associated with idiopathic hypercalcemia, vitamin D intoxication, hyperparathyroidism, sarcoidosis, pneumocystic carinii pneumonia, and blue diaper syndrome. Low calcium levels are associated with hypoparathyroidism, pseudo-hypoparathyroidism, chronic renal failure, rickets, infantile tetany, and steroid therapy.

8. Cholesterol

An elevated cholesterol level is associated with diabetes, nephrosis, hypothyroidism, biliary obstruction, and those rare cases of idiopathic hypercholesterolemia and hyperlipidemia; low levels are associated with hyperthyroidism, hepatitis, and sometimes severe anemia or infection.

9. Creatinine

Creatinine measurement serves as a test for normal glomerular filtration. Elevated levels are associated with acute and chronic renal insufficiency and urinary tract obstruction. Levels below 0.6 mg/dL are of no significance.

10. Gamma Glutamyl Transaminase (GGT)

GT measurement is principally used to diagnose and monitor hepatobiliary disease. It is currently the most sensitive enzymatic indicator of liver disease, with normal values rarely found in the

presence of hepatic disease. It is also used as a sensitive screening test for occult alcoholism. Elevated levels are found in patients who chronically take drugs such as phenobarbital and phenytoin.

11. Glucose

Glucose measurements are used in the diagnosis and treatment of pancreatic islet cell carcinoma and of carbohydrate metabolism disorders, including diabetes mellitus, neonatal hypoglycemia, and idiopathic hypoglycemia.

12. Iron

Iron (non-heme) measurements are used in the diagnosis and treatment of diseases such as iron deficiency anemia, chronic renal disease, and hemochromatosis (a disease associated with widespread deposit in the tissues of two iron-containing pigments, hemosiderin and hemofuscin, and characterized by pigmentation of the skin).

13. Lactate Dehydrogenase (LDH)

LDH measurements are used in the diagnosis and treatment of liver diseases such as acute viral hepatitis, cirrhosis, and metastatic carcinoma of the liver; cardiac diseases such as myocardial infarction; and tumors of the lungs or kidneys.

14. Phosphorus

There is a reciprocal relationship between serum calcium and inorganic phosphorus. Any increase in the level of inorganic phosphorus causes a decrease in the calcium level by a mechanism not clearly understood. Hyperphosphatemia is associated with vitamin D hypervitaminosis, hypoparathyroidism, and renal failure. Hypophosphatemia is associated with rickets, hyperparathyroidism, and Fanconi syndrome. Measurements of inorganic phosphorus are used in the diagnosis and treatment of various disorders, including parathyroid gland, kidney diseases, and vitamin D imbalance.

15. Sodium, Potassium, and Chloride

Hyponatremia (low serum sodium level) is associated with a variety of conditions, including severe polyuria, metabolic acidosis, Addison's disease, diarrhea, and renal tubular disease. Hypernatremia (increased serum sodium level) is associated with Cushing's syndrome, severe dehydration due to primary water loss, certain types of brain injury, diabetic coma after therapy with insulin, and excess treatment with sodium salts.

Hypokalemia (low serum potassium level) is associated with body potassium deficiency, excessive potassium loss caused by prolonged diarrhea or prolonged periods of vomiting and increased secretion of mineralocorticosteroids. Hyperkalemia (increased serum potassium

level) is associated with oliguria, anuria, and urinary obstruction.

Low serum chloride values are associated with salt-losing nephritis, Addisonian crisis, prolonged vomiting, and metabolic acidosis caused by excessive production or diminished excretion of acids. High serum chloride values are associated with dehydration and conditions causing decreased renal blood flow, such as congestive heart failure.

16. Total Bilirubin

Elevated levels are associated with hemolytic jaundice, paroxysmal hemoglobinuria, pernicious anemia, polycythemia, icterus neonatorum, internal hemorrhage, acute hemolytic anemia, malaria, and septicemia. Low bilirubin levels are associated with aplastic anemia, and certain types of secondary anemia resulting from toxic therapy for carcinoma and chronic nephritis.

17. Total Protein

Total protein measurements are used in the diagnosis and treatment of a variety of diseases involving the liver, kidney, or bone marrow, as well as other metabolic or nutritional disorders.

18. Triglycerides

Triglyceride measurements are used in the diagnosis of diabetes mellitus, nephrosis, liver obstruction, and other diseases involving lipid metabolism and various endocrine disorders and in the treatment of patients with these diseases.

19. Uric Acid

Uric acid measurements are used in the diagnosis and treatment of numerous renal and metabolic disorders, including renal failure, gout, leukemia, psoriasis, starvation, or other wasting conditions and in the treatment of patients receiving cytotoxic drugs.

Eligible Sample

Participants aged 12 year and older are tested.

Description of Laboratory Methodology

The 19 analytes described in this method constitute the routine biochemistry profile. The analyses are performed with a Beckman Synchron LX20. Each analyte is described separately within each pertinent section of this document. NOTE: Glucose, cholesterol, and triglycerides were analyzed as part of this profile, but the results do not replace the formalized reference methods data from NHANES 2003-

2004 samples analyzed at other institutions.

1. Alanine Aminotransferase (ALT)

The LX20 uses an enzymatic rate method to measure ALT activity in serum or plasma. In the reaction, ALT catalyzes the reversible transamination of L-alanine and α -ketoglutarate to pyruvate and L-glutamate. The pyruvate is then reduced to lactate in the presence of lactate dehydrogenase (LDH) with the concurrent oxidation of NADH to NAD. The system monitors the rate of change in absorbance at 340 nm over a fixed-time interval. The rate of change in absorbance is directly proportional to the ALT activity in the sample.

2. Albumin

The method used to measure the albumin concentration on the LX20 is a bichromatic digital endpoint method. In the reaction, the albumin combines with Bromcresol Purple (BCP) reagent to form a complex. The system monitors the change in absorbance at 600 nm. The change in absorbance is directly proportional to the concentration of albumin in the sample.

3. Alkaline Phosphatase (ALP)

The LX system uses an enzymatic rate using a 2-Amino-2-Methyl-1-Propanol (AMP) buffer to measure ALP activity in serum or plasma. In the reaction, the ALP catalyzes the hydrolysis of the colorless organic phosphate ester substrate, p-Nitrophenylphosphate, to the yellow colored product p-Nitrophenol and phosphate. This reaction occurs at an alkaline pH of 10.3. The system monitors the rate of change in absorbance at 410 nm over a fixed-time interval. This rate of change in absorbance is directly proportional to the ALP activity in the serum.

4. Aspartate Aminotransferase (AST)

The LX20 uses an enzymatic rate method to measure the AST activity in serum or plasma. In the reaction, the AST catalyzes the reversible transamination of L-aspartate and α -ketoglutarate to oxaloacetate and L-glutamate. The oxaloacetate is then reduced to malate in the presence of malate dehydrogenase with the concurrent oxidation of NADH to NAD. The system monitors the rate of change in absorbance at 340 nm over a fixed-time interval. The rate of change in absorbance is directly proportional to the AST activity in the sample.

5. Bicarbonate (HCO₃)

The LX20 system uses indirect (or diluted) ISE methodology to measure the total CO₂ level in serum, plasma or urine. The system

measures the rate of pH change as CO₂ ions diffuse across a membrane. The electrode used for CO₂ determination is actually a pH electrode with the tip covered by a silicone rubber membrane and lowers the pH of a bicarbonate solution between the tip of the membrane and the tip of the pH electrode. The rate of pH change is directly proportional to the carbon dioxide (CO₂) in the sample.

6. Blood Urea Nitrogen (BUN)

The LX20 modular chemistry (BUNm) is used to quantitatively determine the concentration of blood urea nitrogen in serum or plasma by means of the enzymatic conductivity rate method. A precise volume of sample is injected into the urease reagent in a reaction cup containing an electrode that responds to changes in solution conductivity. Electronic circuits determine the rate of increase in conductivity, which is directly proportional to the concentration of urea in the sample.

7. Calcium

The LX20 system uses indirect (or diluted) ISE methodology to measure calcium concentration in serum, plasma, or urine. The system determines calcium concentration by measuring calcium ion activity in solution. When the sample buffer mixture contacts the electrode, calcium ions complex with the ionophore at the electrode surface. Changes in potential develop at the electrode surface as the reaction occurs. These changes in potential are referenced to a sodium reference electrode. The reference signal is used in calculating the analyte concentrations based on the Nernst equation.

8. Cholesterol

The LX20 uses the timed-endpoint method to measure the cholesterol concentration in serum or plasma. In the reaction, the cholesterol esterase hydrolyzes cholesterol esters to free cholesterol and fatty acids. The free cholesterol is oxidized to cholesten-3-one and hydrogen peroxide by cholesterol oxidase. Peroxidase catalyzes the reaction of hydrogen peroxide with 4-aminoantipyrine and phenol to produce a colored quinoneimine product. The system monitors the change in absorbance at 520 nm at a fixed-time interval. The change in absorbance is directly proportional to the concentration of cholesterol in the sample.

9. Creatinine

The LX20 modular chemistry side uses the Jaffe rate method (kinetic alkaline picrate) to determine the concentration of creatinine in serum, plasma, or urine. A precise volume of sample is introduced into a

reaction cup containing an alkaline picrate solution. Absorbance readings are taken at both 520 nm and 560 nm. Creatinine from the sample combines with the reagent to produce a red color complex. The observed rate measurement at 25.6 seconds after sample introduction has been shown to be a direct measure of the concentration of the creatinine in the sample.

10. Gamma Glutamyltransaminase (-GT)

The LX uses an enzymatic rate method to determine the GGT activity in serum or plasma. In the reaction, the GGT catalyzes the transfer of a gamma-glutamyl group from the colorless substrate, gamma-glutamyl-p-nitroaniline, to the acceptor, glycylglycine with production of the colored product, p-nitroaniline. The system monitors the rate of change in absorbance at 410 nm over a fixed-time interval. The rate of change in absorbance is directly proportional to the activity of GGT in the sample.

11. Glucose

On the Modular Chemistry side of the LX20, glucose concentration in biologic fluids is determined by the oxygen rate method employing a Beckman Oxygen electrode. A precise volume of sample is introduced in a reaction cup containing an electrode that responds to oxygen concentration. Electronic circuits determine the rate of oxygen consumption, which is directly proportional to the concentration of glucose in the sample.

12. Iron

The method used to measure the iron concentration is a timed-endpoint method. In the reaction, iron is released from transferrin by acetic acid and is reduced to the ferrous state by hydroxylamine and thioglycolate. The ferrous ion is immediately complexed with the FerroZine Iron Reagent. The system monitors the change in absorbance at 560 nm at a fixed-time interval. This change in absorbance is directly proportional to the concentration of iron in the sample

13. Lactate Dehydrogenase (LDH)

The LX20 with LD reagent (using lactate as substrate) utilizes an enzymatic rate method to measure LD activity in biological fluids. In the reaction, the LD catalyzes the reversible oxidation of L-Lactate to Pyruvate with the concurrent reduction of β -Nicotinamide Adenine Dinucleotide (NAD) to β -Nicotinamide Adenine Dinucleotide (reduced form) (NADH). The system monitors the rate of change in absorbance at 340 nm over a fixed-time interval. The rate of change in absorbance is directly proportional to the activity of LD in the sample.

14. Phosphorus

The LX system uses a timed-rate method to determine the concentration of phosphorus in serum, plasma and urine. In the reaction, inorganic phosphorus reacts with ammonium molybdate in an acidic solution to form a colored phosphomolybdate

15. Sodium, Potassium, and Chloride

The LX system utilizes indirect (or diluted) I.S.E. methodology to determine the concentration of sodium in biological fluids. The LX determines sodium ion concentration by measuring electrolyte activity in solution. When the sample/buffer mixture contacts the electrode, sodium ions undergo an ion exchange in the hydrated outer layer of the glass electrode. As the ion exchange takes place, a change in voltage (potential) is developed at the face of the electrode. The potential follows the Nernst equation and allows the calculation of sodium concentration in a solution.

The LX system uses indirect (or diluted) I.S.E. methodology to measure potassium in biological fluids. The system determines potassium ion concentration by measuring electrolyte activity in solution. The potassium electrode consists of valinomycin membrane. The voltage (potential) change that takes place within the membrane follows the Nernst equation and allows the calculation of potassium concentration in solution.

The LX system uses indirect (or diluted) I.S.E. methodology to determine chloride concentration in biological fluids. Chloride is measured using an Ag/AgCl electrode. At the face of the electrode, solid AgCl dissolves to the extent as to saturate the solution around the tip with silver (Ag^+) and Chloride (Cl^-) ions until equilibrium is established. The product of the ion concentrations in solution, at equilibrium, with an excess of the slightly soluble AgCl is defined as the solubility product constant (K_{sp}). When chloride sample is added, the K_{sp} of the solution at the tip is disrupted as AgCl precipitates out of solution. To reestablish the equilibrium, Ag^+ ions are generated from the tip causing a change in the potential. According to the Nernst equation, this change is proportional to the concentration of chloride in the sample.

16. Total Bilirubin

The LX20 uses a timed-endpoint Diazo method to measure the concentration of total bilirubin in serum or plasma. In the reaction, bilirubin reacts with diazo reagent in the presence of caffeine, benzoate, and acetate as accelerators to form azobilirubin. The system monitors the change in absorbance at 520 nm at a fixed-time interval. This

change in absorbance is directly proportional to the concentration of total bilirubin in the sample.

17. Total Protein

The LX20 uses a timed rate biuret method to measure the concentration of total protein in serum or plasma. Proteins in the sample combine with the reagent producing alkaline copper-protein chelate. The rate change in absorbance is monitored by a detector at 545 nm. The observed rate of chelate formation is directly proportional to the total protein concentration in the sample.

18. Triglycerides

The LX uses a timed-endpoint method to determine the concentration of triglycerides in serum or plasma. Triglycerides in the sample are hydrolyzed to glycerol and free fatty acids by the action of lipase. A sequence of three coupled enzymatic steps using glycerol kinase (GK), glycerophosphate oxidase (GPO), and horseradish peroxidase (HPO) causes the oxidative coupling of 3,5-dichloro-2-hydroxybenzenesulfonic acid (DHBS) with 4-aminoantipyrine to form a red quinoneimine dye. The system monitors the change in absorbance at 520 nm for a fixed-time interval. The change in absorbance is directly proportional to the concentration of triglycerides in the sample.

19. Uric acid

The LX20 uses a timed endpoint method to measure the concentration of uric acid in serum, plasma or urine. Uric acid is oxidized by uricase to produce allantoin and hydrogen peroxide. The hydrogen peroxide reacts with 4-aminoantipyrine (4-AAP) and 3,5-dichloro-2-hydroxybenzene sulfonate (DCHBS) in a reaction catalyzed by peroxidase to produce a colored product. The system monitors the change in absorbance at 520 nm at a fixed time interval. The change in absorbance is directly proportional to the concentration of uric acid in the sample.

There were changes to the lab site from the 1999 to 2004. Coulston Foundation at Alamogordo, New Mexico performed testing from 1999 to 2001 and Collaborative Laboratory Services at Ottumwa, Iowa performed testing from 2002 to 2004.

A detailed description of the laboratory method used can be found at NHANES website.

Laboratory Quality Control and Monitoring

The NHANES quality control and quality assurance protocols (QA/QC) meet the 1988 Clinical Laboratory Improvement Act mandates. Detailed quality control and quality assurance instructions are discussed in the NHANES Laboratory/Medical Technologists Procedures Manual (LPM).

Read the LABDOC file for detailed QA/QC protocols

A detailed description of the quality assurance and quality control procedures can be found on the NHANES website.

Data Processing and Editing

Specimens were processed, stored and shipped to Collaborative Laboratory Services in Ottumwa, Iowa. Detailed specimen collection and processing instructions are discussed in the NHANES Laboratory/Medical Technologists Procedures Manual (LPM). Read the LABDOC file for detailed data processing and editing protocols. The analytical methods are described in the Description of the Laboratory Methodology section.

Many derived variables were created in this data file. The formula for their derivation is as follows:

LBDSALSI:

The albumin in g/dL (LBXSAL) was converted to g/L (LBDSALSI) by multiplying by 10.

LBDSBUSI:

The blood urea nitrogen (BUN) in mg/dL (LBXSBU) was converted to mmol/L (LBDSBUSI) by multiplying by 0.357

LBDESCASI

The calcium in mg/dL (LBXSCA) was converted to mmol/L (LBDESCASI) by multiplying by 0.250

LBDSCHSI

The cholesterol in mg/dL (LBXSCH) was converted to mmol/L (LBDSCHSI) by multiplying by 0.02586.

LBDSGLSI

The glucose in mg/dL (LBXSGL) was converted to mmol/L (LBDSGLSI) by multiplying by 0.05551.

LBDSIRSI

The iron in $\mu\text{g/dL}$ (LBXSIR) was converted to $\mu\text{mol/L}$ (LBDSIRSI) by multiplying by 0.1791.

LBDSPHSI

The phosphorus in mg/dL (LBXSPH) was converted to mmol/L (LBDSPHSI) by multiplying by 0.3229.

LBDSTBSI

The total bilirubin in mg/dL (LBXSTB) was converted to $\mu\text{mol/L}$ (LBDSTBSI) by multiplying by 17.1.

LBDSTPSI

The total protein in g/dL (LBXSTP) was converted to g/L (LBDSTPSI) by multiplying by 10.

LBDSTRSI

The triglycerides in mg/dL (LBXSTR) were converted to mmol/L (LBDSTRSI) by multiplying by 0.01129.

LBDSUASI

The uric acid in mg/dL (LBXSUA) was converted to $\mu\text{mol/L}$ (LBDSUASI) by multiplying by 59.48.

LBDSCRSI

The creatine in mg/dL (LBXSCR) was converted to $\mu\text{mol/L}$ (LBDSCRSI) by multiplying by 88.4.

LBDSGBSI

The globulin in g/dL (LBXSGB) was converted to g/L (LBDSGBSI) by multiplying by 10.

Detailed instructions on specimen collection and processing can be found at the NHANES website.

Analytic Notes

Analytical Note for Serum Creatinine (Lab 40) for NHANES 2003-2004:

Correction for Serum Creatinine for NHANES 2003-2004 is not necessary:

Serum creatinine is not standardized in many laboratories. The National Kidney Disease Education Program is attempting to have all laboratories standardize serum creatinine to reference methods (Myers, GL, et al. Recommendations for Improving Serum Creatinine Measurement: A Report from the Laboratory Working Group of the National Kidney Disease Education Program. Clin. Chem. 2006; 5-18). Equations for estimating glomerular filtration rate (GFR) from standardized creatinine have been published (Stevens LA, et al. N Engl J Med. 2006 Jun 8;354(23):2473-83). Serum creatinine assays on 190 stored specimens from NHANES 2003-2004 were used to determine if serum creatinine needed to be adjusted when compared to a method traceable to a "gold" standard reference method. The Cleveland Clinic Foundation (CCF) laboratory analyzed the serum creatinine specimens using a Roche coupled enzymatic assay (creatininase, creatinase, sarcosine oxidase, kits # 1775677 and 1775766) performed on a Roche P Module instrument. The Roche method calibrators were traceable to an isotope dilution mass spectrometric method for serum creatinine using standard references methods (NIST SRM 967) and confirmed by analysis of CAP LN-24 linearity set based on NIST assigned values. Serum creatinine by the Roche method was then compared to the original

NHANES 2003-2004 measurements which used the Jaffe kinetic alkaline picrate method performed on a Beckman LX-20 analyzer. There were no significant differences in results between these two measurements. The comparison of values revealed the mean (SD) serum creatinine at NHANES, CCF, and their difference were 0.977 (0.341), 0.970 (0.337), and 0.007 (0.0583) mg/dL, respectively (paired t-test, $p=0.09$). The regression of CCF (Y) on NHANES (X) had an intercept and slope of 0.02 (0.01) and 0.98 (0.01) mg/dL and a correlation of 0.99. The difference between the methods was within the analytical error of the NHANES method (CV of 4.3% at 0.67 mg/dL or 0.029 mg/dL). Thus, no correction is necessary for serum creatinine values in NHANES 2003-2004. The analysis of NHANES 2003-2004 laboratory data must be conducted with the key survey design and basic demographic variables. The NHANES 2003-2004 Household Questionnaire Data Files contain demographic data, health indicators, and other related information collected during household interviews. They also contain all survey design variables and sample weights for these age groups. The phlebotomy file includes auxiliary information such as the conditions precluding venipuncture. The household questionnaire and phlebotomy files may be linked to the laboratory data file using the unique survey participant identifier SEQN.

LBXSTR:

This value was obtained from the standard battery of biochemical assessments. Use of the laboratory test result from the reference method (LBXTR), rather than the (LBXSTR) value, is generally recommended. For most triglyceride analyses, the appropriate variable to use is (LBXTR). The value from the biochemistry profile (LBXSTR) should not be used routinely.

LBXSCH:

This value was obtained from the standard battery of biochemical assessments. Use of the laboratory test result from the reference method (LBXTC), rather than the (LBXSCH) value, is generally recommended. For most analyses of serum cholesterol, the appropriate variable to use will be (LBXTC). The (LBXSCH) value from the biochemistry profile should not be used routinely.

LBXSGL

This value was obtained from the standard battery of biochemical assessments. Use of the laboratory test result from the reference method (LBXGLU), rather than the (LBXSGL) value, is generally recommended. These serum glucose values (LBXSGL) reported in this release should not be used to determine undiagnosed diabetes or prediabetes. Instead, plasma glucose values (LBXGLU) should be used based on the reference analytic method of this analyte. Use the special weights included in this data file when analyzing data.

References

1. N/A

Locator Fields

Title: Biochemistry Profile

Contact Number: 1-866-441-NCHS

Years of Content: 2003–2004

First Published: January 2006

Revised: October 2007

Access Constraints: None

Use Constraints: None

Geographic Coverage: National

Subject: Biochemistry Profile

Record Source: NHANES 2003–2004

Survey Methodology: NHANES 2003–2004 is a stratified multistage probability sample of the civilian non-institutionalized population of the U.S.

Medium: NHANES Web site; SAS transport files

**National Health and Nutrition Examination Survey
Codebook for Data Production (2003-2004)**

**Biochemistry Profile (L40_C)
Person Level Data**

First Published: January 2006

Last Revised: August 2006



SEQN	Target
	B(12 Yrs. to 150 Yrs.)
Hard Edits	SAS Label
	Respondent sequence number
English Text: Respondent sequence number.	
English Instructions:	

LBXSAL	Target
	B(12 Yrs. to 150 Yrs.)
Hard Edits	SAS Label
	Albumin (g/dL)
English Text: Albumin (g/dL)	
English Instructions:	

Code or Value	Description	Count	Cumulative	Skip to Item
1.9 to 5.5	Range of Values	6492	6492	
.	Missing	498	6990	

LBDSALSI	Target
	B(12 Yrs. to 150 Yrs.)
Hard Edits	SAS Label
	Albumin (g/L)
English Text: Albumin (g/L)	
English Instructions:	

Code or Value	Description	Count	Cumulative	Skip to Item
19 to 55	Range of Values	6492	6492	
.	Missing	498	6990	

LBXSATSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Alanine aminotransferase ALT (U/L)			
English Text: Alanine aminotransferase(U/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
5 to 1997	Range of Values	6489	6489	
.	Missing	501	6990	

LBXSASSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Aspartate aminotransferase AST (U/L)			
English Text: Aspartate aminotransferase(U/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
7 to 1672	Range of Values	6488	6488	
.	Missing	502	6990	

LBXSAPSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Alkaline phosphatase (U/L)			
English Text: Alkaline phosphatase(U/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
16 to 1378	Range of Values	6492	6492	
.	Missing	498	6990	

LBXSBU	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Blood urea nitrogen (mg/dL)			
English Text: Blood urea nitrogen (mg/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
2 to 91	Range of Values	6492	6492	
.	Missing	498	6990	

LBDSBUSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Blood urea nitrogen (mmol/L)			
English Text: Blood urea nitrogen (mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
0.71 to 32.49	Range of Values	6492	6492	
.	Missing	498	6990	

LBXSCA	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Total calcium (mg/dL)			
English Text: Total calcium (mg/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
6.7 to 12.5	Range of Values	6492	6492	
.	Missing	498	6990	

LBDSCASI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Total calcium (mmol/L)			
English Text: Total calcium (mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
1.675 to 3.125	Range of Values	6492	6492	
.	Missing	498	6990	

LBXSCH	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Cholesterol (mg/dL)			
English Text: Cholesterol (mg/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
72 to 712	Range of Values	6492	6492	
.	Missing	498	6990	

LBDSCHSI		Target		
		B(12 Yrs. to 150 Yrs.)		
Hard Edits		SAS Label		
		Cholesterol (mmol/L)		
English Text: Cholesterol (mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
1.862 to 18.412	Range of Values	6492	6492	
.	Missing	498	6990	

LBXSC3SI		Target		
		B(12 Yrs. to 150 Yrs.)		
Hard Edits		SAS Label		
		Bicarbonate (mmol/L)		
English Text: Bicarbonate(mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
12 to 35	Range of Values	6492	6492	
.	Missing	498	6990	

LBXSGTSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Gamma glutamyl transferase (U/L)			
English Text: Gamma glutamyl transferase : SI (U/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
3 to 2274	Range of Values	6491	6491	
.	Missing	499	6990	

LBXSGL	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Glucose, serum (mg/dL)			
English Text: Glucose, serum (mg/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
46 to 512	Range of Values	6492	6492	
.	Missing	498	6990	

LBDSGLSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Glucose, serum (mmol/L)			
English Text: Glucose, serum (mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
2.55 to 28.42	Range of Values	6492	6492	
.	Missing	498	6990	

LBXSIR	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Iron, refrigerated (ug/dL)			
English Text: Iron, refrigerated (ug/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
5 to 278	Range of Values	6486	6486	
.	Missing	504	6990	

LBDSIRSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Iron, refrigerated (umol/L)			
English Text: Iron, refrigerated (umol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
0.9 to 49.8	Range of Values	6486	6486	
.	Missing	504	6990	

LBXSLDSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Lactate dehydrogenase LDH (U/L)			
English Text: Lactate dehydrogenase LDH (U/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
39 to 1292	Range of Values	6487	6487	
.	Missing	503	6990	

LBXSPH	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Phosphorus (mg/dL)			
English Text: Phosphorus (mg/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
2.2 to 6.8	Range of Values	6491	6491	
.	Missing	499	6990	

LBDSPHSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Phosphorus (mmol/L)			
English Text: Phosphorus (mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
0.71 to 2.196	Range of Values	6491	6491	
.	Missing	499	6990	

LBXSTB	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Total bilirubin (mg/dL)			
English Text: Total bilirubin (mg/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
0.1 to 12.9	Range of Values	6487	6487	
.	Missing	503	6990	

LBDSTBSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Bilirubin, total (umol/L)			
English Text: Bilirubin, total (umol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
1.71 to 220.59	Range of Values	6487	6487	
.	Missing	503	6990	

LBXSTP	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Total protein (g/dL)			
English Text: Total protein (g/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
5.4 to 10	Range of Values	6489	6489	
.	Missing	501	6990	

LBDSTPSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Total protein (g/L)			
English Text: Total protein (g/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
54 to 100	Range of Values	6489	6489	
.	Missing	501	6990	

LBXSTR	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Triglycerides (mg/dL)			
English Text: Triglycerides (mg/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
20 to 5210	Range of Values	6488	6488	
.	Missing	502	6990	

LBDSTRSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Triglycerides (mmol/L)			
English Text: Triglycerides (mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
0.226 to 58.821	Range of Values	6488	6488	
.	Missing	502	6990	

LBXSUA	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Uric acid (mg/dL)			
English Text: Uric acid (mg/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
1.4 to 14.9	Range of Values	6490	6490	
.	Missing	500	6990	

LBDSUASI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Uric acid (umol/L)			
English Text: Uric acid (umol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
83.3 to 886.3	Range of Values	6490	6490	
.	Missing	500	6990	

LBXSCR	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
0 to 99.9	Creatinine (mg/dL)			
English Text: Creatinine (mg/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
0.3 to 13.7	Range of Values	6492	6492	
.	Missing	498	6990	

LBDSRSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Creatinine (umol/L)			
English Text:				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
26.52 to 1211.08	Range of Values	6492	6492	
.	Missing	498	6990	

LBXSNASI		Target		
		B(12 Yrs. to 150 Yrs.)		
Hard Edits		SAS Label		
		Sodium (mmol/L)		
English Text: Sodium(mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
126 to 153	Range of Values	6491	6491	
.	Missing	499	6990	

LBXSKSI		Target		
		B(12 Yrs. to 150 Yrs.)		
Hard Edits		SAS Label		
		Potassium (mmol/L)		
English Text: Potassium (mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
2.5 to 5.9	Range of Values	6491	6491	
.	Missing	499	6990	

LBXSCLSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Chloride (mmol/L)			
English Text: Chloride(mmol/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
86 to 116	Range of Values	6492	6492	
.	Missing	498	6990	

LBXSOSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Osmolality (mmol/Kg)			
English Text: Osmolality(mmol/Kg)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
250 to 306	Range of Values	6491	6491	
.	Missing	499	6990	

LBXSGB	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Globulin (g/dL)			
English Text: Globulin (g/dL)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
1.7 to 7.4	Range of Values	6489	6489	
.	Missing	501	6990	

LBDSGBSI	Target			
	B(12 Yrs. to 150 Yrs.)			
Hard Edits	SAS Label			
	Globulin (g/L)			
English Text: Globulin (g/L)				
English Instructions:				
Code or Value	Description	Count	Cumulative	Skip to Item
17 to 74	Range of Values	6489	6489	
.	Missing	501	6990	