

**APPENDIX V**

**FSIS LABORATORY RESIDUE ANALYTICAL  
CAPABILITY**

# **APPENDIX V. FSIS LABORATORY RESIDUE ANALYTICAL CAPABILITY**

## **INTRODUCTION**

The Food Safety and Inspection Service (FSIS) requires practical analytical methods for detecting, quantifying, and identify residues that may be present in meat, poultry, and their processed products. These methods can be used by the Agency for monitoring and surveillance activities to determine whether product is adulterated and for human risk assessment (exposure) purposes. The Agency uses available methodology to take appropriate regulatory action against adulterated products, consistent with the reliability of the analytical data. However, because of the large number of potential residues that may occur in the food chain, practical methods are not available for many compounds of interest. This section describes the types of methods used by FSIS to conduct analyses and their suitability for regulatory use. A list of key terms precedes the method descriptions. Note that the chemistry method descriptions with few exceptions, are referenced to the latest edition of the FSIS Analytical Chemistry Laboratory Guidebook.

## **CRITERIA FOR PRACTICAL METHODS**

The following criteria have been identified as primary concerns for methods suitable for regulatory use.

1. The method requires no more than 2-4 hours of analytical time per sample.
2. A quality assurance plan has been developed for the method.

## **KEY TO ABBREVIATIONS**

**AAS** -- Atomic Absorption Spectrometry

**AOAC** -- Association of Official Analytical Chemists International

**CELIA, CA** -- Competitive Enzyme Labeled Immunoassay for Chloramphenicol: a laboratory test that detects and identifies chloramphenicol residues in cattle and pork muscle

**ECD** -- Electron capture detector

**EI** -- Electron impact

**ELISA** -- Enzyme-Linked Immuno Sorbent Assay

**E-Z SCREEN** -- A proprietary immunoassay system for rapid detecting and identifying various antibiotics and other residues in tissue extracts

**GC** -- Gas chromatography

**GLC** -- Gas liquid chromatography (same as GC)

**GPC** -- Gel Permeation Chromatography

**HFAA** -- Heptafluorobutyric acid anhydride

**HPLC** -- High pressure liquid chromatography

**ICP**-- Inductively Coupled Plasma Spectrophotometry

**JAOAC** -- Journal of the Association of Official Analytical Chemists

**LDL** -- Lowest detectable limit. The smallest amount of individual residue or sample component that can be reliably observed or found in the sample matrix by the current appropriate methodology. Qualitative number. Not applicable to compounds with established tolerances.

**METHOD STATUS** -- See discussion above

**MIC** -- Minimum inhibitory concentration: the minimum amount of antimicrobial compound present in a buffer extract of tissue that will inhibit bacterial growth in a cell culture media

**MPL** -- Minimum proficiency level: the minimum amount of analyte expected to be identified and quantified by a laboratory and upon which ongoing capability will be evaluated. It is the smallest concentration for which the predicted coefficient of variation for reproducibility (CV) is less than or equal to 20 percent and the upper 90 percent confidence level for the predicted CV is less than 30 percent

**MS** -- Mass spectrometry

**NADA** -- New Animal Drug Application, issued by the Center for Veterinary Medicine Food and Drug Administration (CVM-FDA)

**NE** -- Level not established

**NICI** -- Negative ion chemical ionization

**PICI** -- Positive ion chemical ionization

**PP** -- Processed product

**QUANTIFICATION** -- The determination of the amount of residue present in a sample

**ppb** -- Parts per billion

**ppm** -- Parts per million

**QAP**—Quality Assurance Plan

**REFERENCE METHODS** -- Analytical procedures by which other methods may be evaluated and for which standards are established. These methods are considered suitable for regulatory use in the National Residue Program

**RESIDUE** -- Any compound present in edible or target tissues of the animal that results from that compound's use or inadvertent introduction into the animal. "Residue" includes the compound itself, its metabolite, and other substances formed in or on food because of the compound's use or inadvertent

introduction

**SOS** -- Sulfa-on-Site: a rapid in-plant chemical screening test for detecting sulfonamide residues in food animal urine or serum that provides same-day results

**STOP** -- Swab Test on Premises: an overnight in-plant laboratory microbiological screen test for detecting antibiotic residues in edible tissues

**SWAB** -- STOP precursor: an overnight laboratory microbiological screen test for detecting antibiotic residues in edibles tissues

**TLC** -- Thin layer chromatography

**UV** -- Ultraviolet spectrophotometric techniques for detection and quantification

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Albendazole (aminosulfone metabolite)</b>	Marker residue detected and quantified by HPLC-fluorescence detection	20 ppb	50 ppb	Cattle/liver	1
	Extraction with organic solvents followed by HPLC with UV detection; confirmed by GC/EI/MS	0.05 ppm	NE	Red meat, liver, muscle	1
<b>Aldicarb and Metabolites</b>	GPC plus HPLC with post-column fluorescence detection; extracts verified by oxidation to the sulfone	5 ppb	10 ppb	All/liver	1
<b>Aldrin</b>	GPC plus GLC	0.02 ppm	0.1 ppm	All/fat	1
	Extracts from GPC or Mills confirmed by GC/MS.	0.03 ppm	NE	All/fat, PP	1
<b>Amoxicillin Trihydrate</b>	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.02 ppm	0.02 ppm	Cattle, swine/kidney, liver, muscle	2
	Tissue extracts quantified by HPLC using fluorometer	0.01 ppm	0.01 ppm	Cattle, swine/kidney,liver, muscle	1
<b>Ampicillin Trihydrate</b>	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.01 ppm	0.01 ppm	Cattle, swine/all	3
<b>Apramycin</b>	Sample extraction TLC; bioautographed using <i>Bacillus subtilis</i> as a test organism	0.05 ppm	0.1 ppm	Swine/kidney, muscle	1
<b>Atrazine</b>	Fat extracted using C18 columns and quantified by capillary GC with nitrogen/ phosphorous detector	5 ppb	NE	All/fat	1
	Extracts confirmed by GC/MS	5 ppb	NE	All/fat	1

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Bacitracin methylene disalicylate</b> <b>Bacitracin, zinc</b>	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.05 ppm	NE	All/kidney, liver, muscle	4
<b>Bambermycins</b>	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	25 ppb	NE	All/kidney, liver, muscle	5
<b>Bendiocarb</b>	GPC plus HPLC with post-column fluorescence detection.	5 ppb	10 ppb	All/liver	1
<b>Benomyl (Benzimidazole)</b>	pH extraction with organic solvents; followed by HPLC ppm with UV detection; extracts derivatized and confirmed by GC/EI/MS	0.05 ppm	50 ppb	Poultry/liver, muscle	1
	Micro alumina assay: column chromatography plus GLC	0.01 ppm	NE	All/fat, PP	1
<b>BHC</b>	GPC plus GLC	0.01 ppm	0.1 ppm	All/fat	1
	Beta and delta isomers: GPC plus GLC	Identification only	Identification only	All/fat	1
	Extracts from GPC or Mills confirmed by GC/MS	0.02 ppm	NE	All/fat, PP	1
<b>Bufencarb</b>	GPC plus HPLC with post- column fluorescence detection. Extracts are subjected to reverse phase chromatography, derivatized and confirmed by GC/MS	5 ppb	10 ppb	All/liver	1
<b>Cacodylic acid</b>	Dry ashed tissue is dissolved and reacted to produce arsine gas, which reacts to form a blue complex for colorimetric quantification	0.05 ppm	0.20 ppm	All/kidney, liver, muscle	6

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Cadmium</b>	Dry ashed tissue is dissolved and quantified by ICP		0.002 :g/ml	All/kidney, liver, muscle	1
<b>Calcium</b>	Tissue is wet ashed and titrated with specific indicator	0.03%	0.03%	All/muscle	8
	Wet ashed tissue is quantified by AAS	NE	NE	All	
<b>Cambendazole (Benzimidazol)</b>	Extraction with organic solvents followed by HPLC with UV detection; extracts confirmed by GC/EI/MS	0.05 ppm	50 ppb	Red meat/liver, muscle, PP	1
<b>Captan</b>	GPC plus GLC	0.02 ppm	0.05 ppm	Red meat/fat	1
<b>Carbadox</b>	Tissue extract is hydrolyzed and a derivative is prepared and separated by ion exclusion chromatography, then quantified by GC-ECD	7.5 ppb	30 ppb	Swine/liver	1
	Extracts confirmed by GC/EI/MS	7.5 ppb	NE	Swine/liver, muscle	1
<b>Carbarsone</b>	Dry ashed tissue is dissolved and reacted to produce arsine gas, which reacts to form blue complex for colorimetric quantification	0.05 ppm	0.20 ppm	All/kidney, liver, muscle	6
<b>Carbaryl</b>	GPC plus HPLC with post-column fluorescence detection	5 ppb	10 ppb	All/liver	1
	Extracts are subjected to reverse phase chromatography, derivatized and confirmed by GC/MS	10 ppb	NE		

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Carbofuran and metabolite</b>	GPC plus HPLC with post-column fluorescence detection.	5 ppb	10 ppb	All/liver	1
	Extracts are subjected to reverse phase chromatography, derivatized and confirmed by GC/MS	10 ppb	NE		
<b>Carbophenothion</b>	Tissue extracts are quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	0.10 ppm	NE	All/liver, muscle	
	GPC plus GLC	0.03 ppm	0.20 ppm	All/fat	1
	GPC extracts are confirmed by CG/EI/MS	0.01 ppm	NE	Red meat/fat	1
<b>Chloramphenicol</b>	Tissue extracts are screened by E-Z screen	25 ppb	NE	Calf/muscle, kidney	9
	Tissue extract screened for chloramphenicol by CELIA CA	5 ppb	NE	Calf/muscle	10
	Extraction of parent and glucuronide using C18 columns with GC capillary quantification as the trimethylsilyl derivative	0.5 ppb	1.0 ppb	Calf/muscle	1
	Extracts are confirmed using NICI/MS	0.15 ppb	NE	Calf/muscle	1
	C18 cleanup of the hydrolyzed extract with GC capillary quantification as the trimethyl derivative	2.5 ppb	5.0 ppb	Calf/urine	1
<b>Chloramphenicol (continued)</b>	Extracts are confirmed using NICI/MS	5 ppb	NE	Calf/urine	1
<b>Chlordane</b>	GPC followed by GC-ECD	0.10 ppm	0.30 ppm	All/fat	1



Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
	Extracts from GPC or Mills are confirmed by GC/MS	NE	NE	All/fat, PP	1
<b>Chlordecone (Kepone)</b>	GPC plus GLC	0.03 ppm	0.20 ppm	All/fat	1
	GPC extracts are confirmed by GC/EI/MS	0.05 ppm (poultry) 0.20 ppm (cattle)	NE	Poultry, red meat/fat	1
<b>2-Chloro-1-(2,4,-di chlorophenyl) vinyl diethyl phosphate (chlorfenvinphos)</b>	GPC plus GLC	0.03 ppm	0.10 ppm	All/fat	1
	GPC extracts are confirmed by CG/EI/MS	0.01 ppm (poultry) 0.10 ppm (red meat)	NE NE	Poultry, red meat/fat	1
<b>2-Chloro-1-(2,4,5-trivinyl dimethyl phosphate (stirophos)</b>	GPC plus GLC	0.05 ppm	0.30 ppm	All/fat	1
	GPC extracts are confirmed by GC/EI/MS	0.05 ppm (poultry) 0.10 ppm (cattle)	NE	Poultry, red meat/fat	1
<b>Chlorpyrifos</b>	GPC followed by GC-ECD	0.05 ppm	0.20 ppm	All/fat	1
	GPC extracts are confirmed by GC/EI/MS	0.05 ppm (poultry) 0.50 ppm (swine)	NE	Poultry, swine/fat	1
<b>Chlortetracycline</b>	Antibiotic screen test (Swab): ability of tissue fluids containing anti-microbial activity to inhibit microbial growth	0.01 ppm	NE	All/kidney	11
<b>Chlortetracycline (contiued)</b>	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.01 ppm	NE	All/kidney, liver, muscle	12
	Extraction using C18 columns followed by HPLC with UV detection	0.05 ppm	0.10 ppm	All/kidney, liver, muscle	1 (Guide-book draft)

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Chromium</b>	Dry ashed tissue is extracted with organic reagent and quantified using ICP	0.01 mg/l	0.037 mg/l	All/kidney, liver, muscle	1
<b>Clenbuterol</b>	Liquid-liquid extraction with internal standard, detected by GC/MS of the oxazolidin-3-one derivatives.		1.0 ppb	Cattle, swine, sheep/liver, muscle.	1 (9/95)
<b>Clopidol</b>	Organic solvent extraction with HPLC-UV detection	0.1 ppm	NE	Poultry/liver	1
	Organic solvent extraction with GC-EC detection	0.1 ppm	NE	Poultry/liver	1
<b>Clorsulon</b>	Tissue extracts are quantified by HPLC-UV detection	0.25 ppm	0.50 ppm	Red meat/kidney, muscle, liver, PP	1
	Tissue extracts for HPLC are derivatized and confirmed by GC/MS	0.5 ppm	NE	Red meat/kidney, muscle, liver, PP	1
<b>Cloxacillin</b>	Antibiotic screen test(Swab): ability of tissue fluids containing anti- microbial activity to inhibit microbial growth	0.16 ppm	NE	All/kidney	14
	Microbiological assay combined with HPLC separation and quantified by microbial inhibition	0.02 ppm	NE	Dairy cows/ kidney, liver, muscle	15
<b>Cobalt</b>	Dry ashed tissue is dissolved and quantified by ICP	0.009 mg/l	0.03 mg/l	All/kidney, liver, muscle	
<b>Copper</b>	Dry ashed tissue is dissolved and quantified by ICP	0.006 mg/l	0.02 mg/l	All/kidney, liver, muscle	1
<b>Coumaphos and oxygen analog</b>	GPC followed by GC-ECD	0.15 ppm	0.30 ppm	All/fat	1
	GPC extracts are confirmed by CG/EI/MS	0.20 ppm	0.30 ppm	Red meat/fat	1

<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Cresylic acid</b>	Tissue extracts are derivatized and determined by GC-ECD	NE	NE	Poultry/fat	1
<b>Crufomate (Ruelene)</b>	Tissue extracts are quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	0.10 ppm	NE	All/liver, muscle	1
<b>Cyanide salts</b>	Aqueous extraction followed by a colorimetric determination	0.5 ppm	NE	All/all	1
	For confirmation, cyanogen chloride is produced and determined by GC/ECD	0.5 ppm	NE	All/all	1
<b>Cyano (3-phenoxy chlorophenyl) methyl- 4-a-(methylethyl) - benzeneacetate (Fenvalerate)</b>	Organic solvent extracts are quantified as the sum of both isomers by GC/EC; extracts are confirmed by GC/EI/MS	0.03 ppb	NE	All/fat	1
<b>DDE (metabolites of DDT collectively reported as DDT)</b>	GPC followed by GC-ECD	0.02 ppm	0.05 ppm	All/fat	1
	Extracts from GPC are confirmed by GC/MS (LRC)	0.02 ppm	NE	All/fat, PP	1
<b>DDT (isomers of DDT collectively reported as DDT)</b>	GPC followed by GC-ECD	0.04 ppm	0.15 ppm	All/fat	1
<b>Decoquate</b>	Zymark Pytech System; organic extraction followed by HPLC with fluorescence detection	0.20 ppm	0.50 ppm	Cattle, poultry/liver, muscle	1
<b>Dexamethasone</b>	Extract is partitioned on an SPE column with quantitation by HPLC/UV		5.0 ppb	Cattle, swine, veal, sheep/liver, muscle	1 (guide-book draft)

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Deltamethrin</b>	Organic solvent extracts are quantified by GC/ECD;		20 ppb	Bovine, poultry/fat	1
	Solvent extraction followed by a competitive ELISA determination	0.5 ppm	NE	All/fat	1
<b>Diazinon</b>	Tissue extracts are quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	0.1 ppm	0.15 ppm	All/liver, muscle	1
<b>Dibutyltin dilaurate</b>	Tissue extraction acid hydrolysis-morin derivatization--HPLC-UV	0.25 ppm	NE	Turkey/liver	1
<b>Dieldrin</b>	GPC plus GLC	0.01 ppm	0.05 ppm	All/fat	1
	Extracts from GPC are confirmed by GC/MS	0.02 ppm	NE	All/fat, PP	1
<b>Diethylstilbestrol (DES)</b>	Solid phase extraction technique using an internal standard followed by methylsilation for GC/MS quantification and confirmation		0.25 ppb	Cattle, sheep/liver, muscle	1
<b>Dihydro-streptomycin</b>	Antibiotic screen test (Swab): ability of tissue fluids containing anti-microbial activity to inhibit microbial growth	0.25 ppm	NE	All/kidney	11
	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.25 ppm	NE	All/kidney, liver, muscle	16

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>3, 5-Dimethyl-4-(methylthio)phenyl ethylcarbamate and metabolite</b>	GPC plus HPLC with post- column fluorescence detection. Extracts are subjected to reverse phase chromatography, derivatized and confirmed by GC/MS	5 ppb	10 ppb	All/liver	1
<b>Dimetridazole and hydroxy metabolite</b>	Extracts are quantified by HPLC/UV	1.0 ppb	NE	Turkey, swine/ muscle	1
	Tissue extracts from HPLC are confirmed by GC/NICI/MS	1.0 ppb	NE	Turkey, swine/ muscle	1
<b>Dioxacarb</b>	GPC plus HPLC with post- column fluorescence detection	5 ppb	10 ppb	All/liver	1
	Extracts are subjected to reverse phase chromatography, derivatized and confirmed by GC/MS	10 ppb	NE		
<b>Dioxathion</b>	Tissue extracts are quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	0.10 ppm	NE	All/liver, muscle	
<b>Dodecachloro-octahydro-1,3,4-metheno-2H-cyclobuta(cd)-pentalene[Mirex]</b>	GPC plus GLC	0.04 ppm	0.10 ppm	All/fat	1
	Extracts from GPC or Mills are confirmed by GC/MS	0.05 ppm	NE	All/fat, PP	1
<b>Doramectin</b>	Tissue extracts are quantified by HPLC fluorescence detection	2 ppb	7.5 ppb	Red meat/liver, muscle	1

<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Endosulfan I</b>	GPC plus GLC	0.01 ppm	0.10 ppm	All/fat	1
	GPC extracts are confirmed by GC/EI/MS	0.02 ppm	NE	Red meat/fat	
<b>Endosulfan II</b>	GPC plus GLC	0.02 ppm	0.20 ppm	All/fat	1
<b>Endrin</b>	GPC plus GLC	0.03 ppm	0.05 ppm	All/fat	1
	Extracts from GPC or Mills are confirmed by GC/MS	0.05 ppm	NE	All/fat, PP	1
<b>Erythromycin</b>	Antibiotic screen test (Swab): ability of tissue fluids containing antimicrobial activity to inhibit microbial growth	25 ppb	NE	All/kidney	11
	Microbiological assay procedure: ability of tissue fluids containing antimicrobial activity to inhibit microbial growth	25 ppb	NE	All/kidney, liver, muscle	17
<b>Ethion and oxygen analog</b>	Tissue extracts are quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	0.10 ppm	0.15 ppm	All/liver, muscle	1
	Tissue extracts are quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	NE	NE	All/liver muscle	1
<b>Ethylene dibromide</b>	Residue is co-distilled from aqueous suspension and quantified by GLC	0.5 ppb	1.0 ppb	All/fat	1
	MS by NICI to determine bromine	1 ppb	NE	All/fat	1
<b>Fenbendazole</b>	Extraction with organic solvents followed by HPLC with UV detection; extracts derivatized and confirmed by GC/EI/MS		0.5 ppb	All/ liver, muscle	1

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Fenbendazole (continued)</b>	Tissue extracts are quantified by HPLC	200 ppb	400 ppb	Cattle, calf/ liver	18
	Quantification extract purified by TLC, derivatized and identified by HPLC fluorescence	200 ppb	NE	Cattle, calf/ liver	18
<b>Fenitrothion</b>	Tissue extracts are quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	0.10 ppm	0.15 ppm	All/liver, muscle	1
<b>Flucythrinate</b>	Organic solvent extracts are quantified as the sum of both isomers by GC/EC; extracts are confirmed by GC/EI/MS			All/fat	1
	Solvent extraction followed by a competitive ELISA determination		NE	All/fat	1
<b>Gasoline</b>	Fat from product is heated in a sealed vial and gasoline components are identified by pattern recognition using GC/flame ionization detection	0.1 ppm	1.0 ppm	Canned meat	1
<b>Gentamicin sulfate</b>	Tissue extracts are screened by E-Z Screen	50 ppb	NE	All/muscle, liver, kidney	19
	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	NE	NE	Swine/kidney	20
	Extraction followed by detection by HPLC with fluorescence detector	0.2 ppm	0.4 ppm	Swine/kidney	1

<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Halofuginone</b>	Tissue extracts are quantified by HPLC-UV	0.05 ppm	0.05 ppm	Chicken/liver	1
	Tissue extracts are confirmed by GC/MS/MS	0.05 ppm	0.05 ppm	Chicken/liver	21
<b>HCB</b>	GPC plus GLC	0.01 ppm	0.01 ppm	All/fat	1
	Extracts from GPC are confirmed by GC/ MS	0.01 ppm	NE	All/fat, PP	1
<b>Heptachlor and heptachlor epoxide</b>	GPC plus GLC	0.01 ppm	0.05 ppm	All/fat	1
	Extracts from GPC are confirmed by GC/MS		0.01 ppm, Heptachlor; 0.1 ppm, Heptachlor Epoxide	All/fat, PP	1
<b>Hetacillin, Potassium</b>	Antibiotic screen test (Swab): ability of tissue fluids containing anti-microbial activity to inhibit microbial growth	NE	NE	All/kidney	14
<b>5-Hydroxy-thiabendazole</b>	Extraction with organic solvents followed by HPLC with UV detection; extracts confirmed by derivatized GC/EI/MS	0.05 ppm	50 ppb	Red meat/liver, muscle	1
<b>Hygromycin B</b>	Antibiotic screen test (Swab): ability of tissue fluids containing antimicrobial activity to inhibit microbial growth	5.00 ppm	NE	All/kidney	14
<b>Ipronidazole and hydroxy metabolite</b>	Tissue extracts are quantified by HPLC/UV	1.0 ppb	NE	Turkey, swine/ muscle	1
	Tissue extracts from HPLC are confirmed by GC/NICI/MS	1.0 ppb	NE	Turkey, swine/ muscle	1
<b>Iron</b>	Dry ashed tissue is dissolved and quantified by ICP	0.009 mg/l	0.14 µg/ml	All/kidney, liver, muscle	1



Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Ivermectin</b>	Tissue extracts are quantified by HPLC fluorescence	2 ppb	7.5 ppb	Red meat/liver, muscle	1
	APCI/MS to confirm Ivermectin and Doramectin	7.5 ppb	15 ppb	Red meat/liver, muscle	1
<b>Lasalocid</b>	Tissue extracts are quantified by HPLC fluorescence detector	0.025 ppm	0.35 ppm	Cattle/liver Poultry/fat, skin	1 22
	Tissue extraction followed by bioautography	0.005 ppm	0.01 ppm	Poultry/fat, skin	22
	GC pyrolysis of the HPLC extract with MS identification of the fragments	0.2 ppm	NE	Cattle/liver Poultry/fat, skin	22
<b>Lead</b>	Dry ashed tissue is dissolved and quantified by ICP		0.05 µg/ml	All/kidney, liver, muscle	1
	Dry ashed tissue is quantified by anodic stripping voltammetry	1.0 ppb	NE	Poultry/kidney, liver	7
<b>Levamisole</b>	Tissue extracts are quantified by GLC flame photometric detection	0.05 ppm	0.1 ppm	Red meat/liver, muscle	1
	Tissue extracts are subjected to GC/MS	0.05 ppm	NE	Red meat/liver, muscle	1
<b>Lincomycin hydrochloride</b>	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.10 ppm	0.10 ppm	Poultry, swine/all	23
<b>Lindane</b>	GPC plus GLC	0.01 ppm	0.01 ppm	All/fat	1
	Extracts from GPC or Mills are confirmed by GC/MS		0.01 ppm	All/fat, PP	1

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Linuron</b>	GPC plus GLC	0.25 ppm	0.50 ppm	All/fat	1
	Extracts are confirmed by GC/EI/MS	0.05 ppm	0.50 ppm	Red meat/fat	1
<b>Lysergic acid diethylamide</b>	Tissue extracts are spotted for TLC and detected with specific chromagenic reagent	NE	NE	All/kidney, liver, muscle	1
<b>Macrolide antibiotics</b>	Qualitative identification by APCI-MS/MS full Scan			All/kidney, liver, muscle	1
<b>Malathion</b>	Tissue extracts were quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	0.10 ppm	0.2 ppm	All/liver, muscle	1
<b>Manganese</b>	Dry ashed tissue is dissolved and quantified by ICP		0.002 µg/ml	All/kidney, liver, muscle	1
<b>Mebendazole</b>	Extraction with organic solvents followed by HPLC with UV detection; extracts are confirmed by GC/EI/MS	0.05 ppm	50 ppb	Red meat/liver, muscle, processed products	1
<b>Melengestrol acetate (MGA)</b>	Tissue extract is column chromatographed on Florisil and quantified by GLC	5.0 ppb	10.0 ppb	Cattle/fat	1
	Extracts are derivatized with HFB and confirmed by GC/EI/MS	5.0 ppb	10.0 ppb	Cattle/fat	1
<b>Mercury</b>	Tissue is digested in acid. Mercury is reduced to its vapor and quantified by flameless AAS	0.01 ppm	0.02 ppm	All/kidney, liver, muscle	1
<b>Methanearsonic acid</b>	Dry ashed tissue is dissolved and reacted to produce arsine gas, which is quantified by AAS	0.05 ppm	NE	All/kidney, liver, muscle	1
	The same as above, but arsine gas reacts to form blue complex for colorimetric quantification	0.05 ppm	0.20 ppm	All/kidney, liver, muscle	1

<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Methomyl</b>	GPC plus HPLC with post-column fluorescence detection	5 ppb	10 ppb	All/liver	1
<b>Methoxychlor</b>	GPC plus GLC	0.15 ppm	0.50 ppm	All/fat	1
	Extracts from GPC are confirmed by GC/MS	0.15 ppm	0.15 ppm	All/fat,	1
<b>Methyl parathion</b>	Tissue extracts are quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	0.10 ppm	0.20 ppm	All/liver, muscle	1
<b>Monensin</b>	Tissue extract is partitioned by TLC and semi-quantified by inhibition of micro-organism growth		NE	Poultry/fat	1
<b>Morantel tartrate</b>	Tissue extract is hydrolyzed and a derivative is quantified by GLC	0.25 ppm 0.50 ppm	0.50 ppm 0.50 ppm	Cattle/liver Cattle/muscle	1 1
	Identification of a structurally significant hydrolyzed fragment by GC/MS	0.25 ppm	NE	Cattle/liver, muscle	24
<b>Narasin</b>	Tissue extracts are spotted on TLC and quantified with a bio-autographic overlay		NE	Cattle, poultry/liver, kidney, fat	1
<b>Neomycin</b>	Antibiotic screen test (Swab): ability of tissue fluids containing anti-microbial activity to inhibit microbial growth	0.25 ppm	NE	All/kidney	14
	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.25	NE	All/kidney, liver, muscle	25

<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Nequinat</b>	Zymark Pytechnology System: Tissue extracts are screened by HPLC/UV	NE	NE	Cattle/liver, muscle	13
<b>Nicarbazin</b>	Tissues are extracted with ethyl acetate; the dinitro-carbanilide moiety is quantified by HPLC-UV. Extracts verified by photodiode array detection	0.1 ppm	0.4 ppm	Chicken/liver muscle	1
<b>Nickel</b>	Dry ashed tissue is dissolved and quantified by ICP	0.015 mg/l	0.014 µg/ml	All/kidney, liver, muscle	1
<b>Nonachlor</b>	GPC plus GLC	0.03 ppm	0.15 ppm	All/fat	1
<b>Novobiocin</b>	Extracts from GPC are confirmed by GC/MS	0.05 ppm	NE	All/fat, PP	
	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.125 ppm	NE	All/kidney, liver, muscle	4
	Zymark Pytechnology System: organic solvent extraction followed by HPLC/UV detection.	0.50 ppm	NE	All/kidney, liver, muscle	1
	Manual system organic solvent extraction followed by HPLC/UV detection	0.50 ppm	1.0 ppm	All/kidney, liver, muscle	1
<b>Oleandomycin</b>	Antibiotic screen test (Swab): ability of tissue fluids containing antimicrobial activity to inhibit microbial growth	0.25 ppm	NE	All/kidney	14

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Oxfendazole</b>	Extraction with organic solvents followed by HPLC with UV detection; extracts are derivatized and confirmed by GC/EI/MS	0.05 ppm	50 ppb	Red meat/liver, muscle, PP	1
<b>Oxytetracycline hydrochloride</b>	Antibiotic screen test (Swab): ability of tissue fluids containing anti-microbial activity to inhibit microbial growth	0.08 ppm	NE	All/kidney	14
	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.08 ppm	0.08 ppm	All/kidney, liver muscle	14
	Extraction using C18 columns followed by HPLC with UV detection	0.05 ppm	0.10 ppm	All/kidney, liver muscle	1 (Guide-book draft)
<b>Parathion</b>	Tissue extracts are quantified by GLC with flame photometric or nitrogen-phosphorous flame ionization detector	NE	NE	All/liver, muscle	
<b>PBB's</b>	GPC cleanup followed by GLC/ECD and MS/EI	0.05 ppm	NE	All/fat	26
<b>PCB's (reported as Aroclor 1242, 1248, 1254, 1260, etc.)</b>	Column chromatography plus GLC.	0.30 ppm	0.50 ppm	All/fat, PP	1
	GPC plus GLC	0.30 ppm	0.50 ppm	All/fat	
<b>Penicillin, procaine and procaine G</b>	Antibiotic screen test (Swab): ability of tissue containing anti-microbial activity to inhibit microbial growth.	12.5 ppb	NE	All/kidney	14

<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Pentachloroanisole</b>	GPC plus GC	NE	NE	Poultry/fat	
	Extracts confirmed by GC/EI/MS	NE	NE	Poultry/fat	
<b>Pentachlorophenol (PCP)</b>	Tissue extracts for GLC are confirmed by GC/MS.	0.03 ppm	50 ppb	All/liver, muscle	1
	Extracts confirmed by GC/EI/MS	NE	NE	Poultry/fat	1
<b>Permethrin (cis &amp; trans)</b>	Solvent extraction followed by a competitive ELISA determination.	1.0 ppb	50 ppb	All/fat	1
	Organic solvent extracts are quantified as the sum of both isomers; extracts are confirmed by GC/EI/MS	1.0 ppb	50 ppb	All/fat	1
<b>Phencyclidine</b>	Tissue extracts are spotted for TLC with specific chromogenic agent	NE	NE	All/kidney, liver, muscle	1
<b>Phenothrin</b>	Solvent extraction followed by a competitive ELISA determination		NE	All/fat	1
<b>Phosalone</b>	GPC plus GLC	0.01 ppm	0.05 ppm	All/fat	1
	GPC extracts are confirmed by GC/EI/MS	0.02 ppm	0.04 ppm	Red meat/fat	1
<b>Propazine</b>	Fat extracted using C18 columns and quantified by capillary GC detector with nitrogen-phosphorous flame ionization detector.	5 ppb	10 ppb	All/fat	1
	Extracts confirmed by GC/MS	5 ppb	10 ppb	All/fat	1
<b>Promecarb</b>	GPC plus HPLC with post- column fluorescence detection. Extracts are subjected to reverse phase chromatography, derivatized and confirmed by GC/MS	5 ppb	10 ppb	All/liver	1

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Propoxur</b>	GPC plus HPLC with post- column fluorescence detection. Extracts are subjected to reverse phase chromatography, derivatized and confirmed by GC/MS	5 ppb	10 ppb	All/liver	1
<b>Pyrantel tartrate</b>	Tissue extract is hydrolyzed and a derivative is quantified by GLC	0.25 ppm	0.50 ppm	Swine/liver, muscle	1
	Identification of a structurally significant hydrolyzed fragment by GC/MS	0.25 ppm	NE	Swine/liver, muscle	27
<b>Ronnel</b>	GPC plus GLC	0.02 ppm	0.05 ppm	All/fat	1
	Extracts are confirmed by GC/EI/MS	0.01 ppm (poultry) 0.10 ppm (red meat)	0.05 ppm 0.05 ppm	Poultry, red meat/fat	1
<b>Roxarsone</b>	Dry ashed tissue is dissolved and reacted to produce arsine gas, which is quantified by AAS.	0.05 ppm	NE	All/kidney, liver muscle	1
	Dry ashed tissue is dissolved and reacted to produce arsine gas, which reacts to form blue complex for colorimetric quantification	0.05 ppm	0.20 ppm	All/kidney, liver, muscle	1
<b>Simazine</b>	Fat extracted using C18 columns and quantified by capillary GC with nitrogen-phosphorous detector.	5 ppb	10 ppb	All/fat	1
	Extracts confirmed by GC/MS	5 ppb	NE	All/fat	1

<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Simazine</b>	Fat extracted using C18 columns and quantified by capillary GC with nitrogen-phosphorous detector.	5 ppb	10 ppb	All/fat	1
	Extracts confirmed by GC/MS	5 ppb	NE	All/fat	1
<b>Spectinomycin hydrochloride</b>	Microbiological assay: tissue extracts are quantified using a turbidimetric assay.	2.8 ppm	NE	All/kidney, liver, muscle	28
<b>Streptomycin</b>	Antibiotic screen test (Swab): ability of tissue fluids containing antimicrobial activity to inhibit microbial growth	0.25 ppm	NE	All/kidney	14
	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.25 ppm	NE	All/kidney, liver, muscle	30
<b>Styrene</b>	Tissues are subjected to GC/MS head space analysis.		NE	All/kidney, liver, muscle fat, PP	1
<b>Sulfachloro-pyridazine</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry	0.02 ppm	0.05 ppm	Red meat/liver, muscle	1
<b>Sulfadiazine</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry	0.02 ppm	0.05 ppm	Red meat/liver, muscle	1
	Extraction followed by GC/CI and EI/MS	0.05 ppm	0.05 ppm	Red meat/liver, muscle	1
<b>Sulfadimethoxine</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry	0.02 ppm	0.05 ppm	All/liver, muscle	1



<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Sulfadoxine</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry.	0.05 ppm	0.05 ppm	Red meat/liver, muscle	1
	Tissue extracts are confirmed by GC/EI/MS	NE	0.05 ppm	Red meat/liver, muscle	1
<b>Sulfaethoxy-pyridazine</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry	0.02 ppm	0.05 ppm	Red meat/liver, muscle	1
<b>Sulfamethazine</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry.	0.02 ppm	0.05 ppm	All/liver, muscle	1
	Tissue extracts are confirmed by GC/EI/MS	0.05 ppm	0.05 ppm	All/liver, muscle	1
	Tissue extracts are screened by E-Z Screen	50 ppb	NE	All/liver	1
<b>Sulfamethoxy-pyridazine</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry	0.02 ppm	0.05 ppm	Red meat/liver, muscle	1
<b>Sulfaphenazole</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry	0.05 ppm	0.05 ppm	Red meat/liver, muscle	1
	Tissue extracts are confirmed by GC/EI/MS	NE	0.05 ppm	Red meat/liver, muscle	1
<b>Sulfapyridine</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry.	0.02 ppm	0.05 ppm	All/liver, muscle	1
	Tissue extracts are confirmed by GC/EI/MS	0.05	0.05 ppm	All/liver, muscle, PP	1

<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Sulfaquinoxaline</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry.	12.5 ppb	0.05 ppm	Poultry/liver, muscle	1
	Tissue extracts are confirmed by GC/EI/MS	25 ppb	0.05 ppm	Poultry/liver, muscle	1
<b>Sulfathiazole</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry	0.02 ppm	0.05 ppm	Red meat/liver, muscle	1
	Tissue extracts are confirmed by GC/EI/MS	0.05 ppm	0.05 ppm	Red meat/liver, muscle, PP	1
<b>Sulfatroxazole</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry.	0.05 ppm	0.05 ppm	Red meat/liver, muscle	1
	Tissue extracts are confirmed by GC/EI/MS	NE	0.05 ppm	Red meat/liver, muscle	1
<b>Sulfisoxazole</b>	TLC fluorescence: tissue extracts are partitioned by TLC and quantified by densitometry.	0.05 ppm	0.05 ppm	Red meat/liver, muscle	1
	Tissue extracts are confirmed by GC/EI/MS	NE	0.05 ppm	Red meat/liver, muscle	1
<b>TDE (metabolite of DDT reported as DDT)</b>	GPC plus GLC	0.03 ppm	0.15 ppm	All/fat	1
	Extracts from GPC are confirmed by GC/MS	0.02 ppm	0.04 ppm	All/fat, PP	1
<b>Terbuthylazine</b>	Fat extracted using C18 columns and quantified by capillary GC with nitrogen-phosphorous detector Extracts confirmed by GC/MS	5 ppb	10 ppb	All/fat	1

<b>Compound</b>	<b>Method Description</b>	<b>LDL/MIC</b>	<b>MPL</b>	<b>Species/Tissues</b>	<b>Reference</b>
<b>Terpene polychlorinates (Strobane)</b>	Micro alumina assay: column chromatography plus GLC	0.50 ppm	1.0 ppm	All/fat	1
	Mills method: Florisil column chromatography plus GLC	0.50 ppm	1.0 ppm	All/fat	1
<b>Tetracycline hydrochloride</b>	Antibiotic screen test (Swab): ability of tissue chloride fluids containing antimicrobial activity to inhibit microbial growth	0.08 ppm	NE	All/kidney	14
	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth.	0.08 ppm	NE	All/kidney, liver, muscle	29
<b>Tetracycline hydrochloride</b>	Extraction using C18 columns followed by HPLC with UV detection.	0.05 ppm	0.10 ppm	All/kidney, liver, muscle	1 (Guide-book draft)
<b>Thiabendazole</b>	pH extraction with organic solvents followed by HPLC with UV detection; extracts derivatized and confirmed by GC/EI/MS	0.05 ppm	50 ppb	Red meat/PP, liver, muscle	1
<b>Tiamulin</b>	Organic solvent extraction followed by GC of the 8-hydroxymutinin metabolite.	0.2 ppm	NE	Swine/liver	30
	Extracts confirmed by GC/MS	NE	NE	Swine/liver	30

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
Tilmicosin	Antibiotic screen test (SWAB): ability of tissue chloride fluids containing antimicrobial activity to inhibit microbial growth.	NE	NA	All/kidney	28
	Microbial assay procedure: ability of tissue extract containing antimicrobial activity to inhibit microbial growth.	NE	NA	All/kidney	29
	Liquid-liquid extraction with quantitation by ion pairing HPLC/UV		0.3 ppm M 0.6 ppm L	Beef/liver, kidney, muscle	1 (Guidebook draft)
Tin	Tissue is dry ashed and dissolved and quantified by AAS (used to screen for organotin compounds)	0.02 ppm	0.1 ppm	All/kidney, liver, muscle	1
Toxaphene	GPC plus GLC	0.50 ppm	1.0 ppm	All/fat	1
	Tissue extracts are screened by E-Z screen.	50 ppm	NE	All/muscle, liver, kidney	9
Tylosin	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.2 :g/mL	NA	Cattle/muscle	31
Virginiamycin	Microbiological assay procedure: ability of tissue extracts containing antimicrobial activity to inhibit microbial growth	0.64 ppm	NA	Swine/kidney, liver, muscle	31
	Organic solvent extraction followed by HPLC/UV quantification	0.1 ppm	0.2 ppm	All/kidney, liver, muscle	1

Compound	Method Description	LDL/MIC	MPL	Species/Tissues	Reference
<b>Zeranol and metabolite taleranol</b>	Extraction followed by radioimmunoassay		NE	Cattle/liver, muscle	1
	Solid phase extraction using an internal standard followed by polymethyl-silation for GC/MS quantification and confirmation		NE	Cattle, sheep/liver, muscle	1
<b>Zinc</b>	Dry ashed tissue is dissolved and quantified by ICP	0.006 mg/l	0.006 µg/ml	All/kidney, liver, muscle	1

### **Residue Analytical Capability References**

1. FSIS Chemistry Laboratory Guidebook, Residue Chemistry.
2. NADA 55-580 & 550-089, Beecham.
3. NADA 55-030, Squibb.
4. Kramer et al., FDA (1974).
5. NADA 44-759, Hoechst.
6. AOAC Book of Methods, 14<sup>th</sup> Edit., 25.050.
7. AOAC, 60,826-832 (1977).
8. AOAC Book of Methods, 15<sup>th</sup> Edit., Pg. 941; Analytical Chemistry Laboratory Guidebook, Food Chemistry, (1993).
9. Environmental Diagnostics, Card Test
10. Antimic. Ag. Chemo., 25.2, 205-211 (1984)
11. J. Food Prot., 1981, 44, 828-831
12. Sec. 6.312 FSIS Microbiology Laboratory Guidebook
13. Copy available upon request from FSIS Eastern Laboratory
14. J. Food Prot., 1981, 44, 828-831
15. NADA 55-069, Beecham-Masengill
16. Sec. 6.315 FSIS Microbiology Laboratory Guidebook

17. Sec. 6.316 FSIS Microbiology Laboratory Guidebook
18. NADA 128-620, American Hoechst
19. Environmental Diagnostics, Gentamicin Card
20. NADA 103-037, 91-191 Schering
21. NADA 130-951, American Hoechst
22. NADA 96-298V, Hoffman-LaRoche
23. NADA 97-505, Upjohn
24. NADA 92-444, 93-903 Pfizer
25. Sec. 6.317 FSIS Microbiology Laboratory Guidebook
26. Guidbook CHC<sub>3</sub> method as screen only. Confirmed by full scan Mass Spectroscopy.
27. NADA 47-244, Upjohn
28. Sec. 6.315 FSIS Microbiology Laboratory Guidebook
29. Sec. 6.312 FSIS Microbiology Laboratory Guidebook
30. Investigational New Animal Drug 1776, Diamond-Shamrock Corp.
31. NADA 91-467, 91-513 Smith Kline & French