

UNITED STATES OF AMERICA
FOOD AND DRUG ADMINISTRATION

+ + + + +

DENTAL PRODUCTS ADVISORY PANEL

+ + + + +

TUESDAY,
JULY 13, 2004

+ + + + +

The above-entitled Meeting was conducted at 8:30 a.m., at the Hilton Washington DC North/Gaithersburg, Salons A and B, 620 Perry Parkway, Gaithersburg, Maryland, Dr. Jon B. Suzuki, Chairman, presiding.

PANEL MEMBERS PRESENT:

JON B. SUZUKI, DDS, PhD, MBA, Chairman, Professor at the University of Pittsburgh School of Dental Medicine

MICHAEL E. ADJODHA, MChE, Executive Secretary, Department of Health and Human Services, FDA, Center for Devices and Radiological Health, Office of Device Evaluation, Division of Anesthesiology, General Hospital, Infection Control, and Dental Devices

SALOMON AMAR, DDS, PhD, Voting Member, Professor of Periodontology at Boston University School of Dental Medicine

DAVID L. COCHRAN, DDS, PhD, Voting Member (Non-Voting for this Meeting), Professor and Chairman of Periodontology at the University of Texas Health Science Center, San Antonio

ELIZABETH S. HOWE, Consumer Representative, Outreach

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

Coordinator, National Foundation for
Ectodermal Dysplasias, Auburn, Washington

ALLISON F. LAWTON, MBA, Drug Industry
Representative, Senior Vice President of
the Genzyme Corporation, Cambridge,
Massachusetts

WILLIAM J. O'BRIEN, MS, PhD, Voting Member (Non-
Voting for this Meeting), Professor of
Materials Science at the University of Michigan
School of Dentistry, Ann Arbor

DANIEL R. SCHECHTER, JD, Device Industry
Representative, General Counsel for
Parkell, Incorporated, Farmingdale, New
York

INDER SHARMA, PhD, Consultant, Deputized to Vote,
Associate Professor of Biostatistics at
the Morehouse School of Medicine,
Department of Community Health and
Preventative Medicine, Atlanta, Georgia

DOMENICK T. ZERO, DDS, MS, Voting Member, Professor
and Chairman of Preventative Dentistry at
Indiana University School of Dentistry,
Indianapolis

JOHN R. ZUNIGA, PhD, DMD, Voting Member, Professor
and Graduate Program Director of Oral
Surgery at the University of North
Carolina School of Dentistry, Chapel Hill

SPONSOR PRESENTERS:

MARK CITRON, Vice President, Regulatory Affairs,
BioMimetic Pharmaceuticals, Inc,
Franklin, TN

ROBERT GENCO, DDS, PhD, Vice Provost, State
University of New York at Buffalo

WILLIAM V. GIANNOBLE, DDS, DMSc, Associate Professor
at University of Michigan

SAMUEL E. LYNCH, DMD, DMSc, President and CEO,
BioMimetic Pharmaceuticals

MYRON NEVINS, DDS, Associate Professor, Harvard
University

FDA PRESENTERS:

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

ANGELA E. BLACKWELL, MS, Biomedical Engineer, Dental
Devices Branch, DHHS/FDA/CDRH/ODE

JUDY S. CHEN, MS, Mathematical Statistician
(Biomedical), Division of Biostatistics,
DHHS/FDA/CDRH/OSB

M. SUSAN RUNNER, DDS, MS, Captain, USPHS, Deputy
Division Director, DAGID and Chief,
Dental Devices Branch, DHHS/FDA/CDRH/ODE

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

A-G-E-N-D-A

Introductions and opening.....5

Presentation on the product GEM 21S
 Mark Citron.....13

Brief Overview
 Dr. Samuel Lynch.....16

Mode of Action
 Dr. William Giannoble.....24

Animal and Human Histology Data
 Dr. Myron Nevins.....37

Results of Randomized Control Clinical Trial
 Dr. Robert Genco.....48

Concluding Remarks
 Dr. Samuel Lynch.....

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22

P-R-O-C-E-E-D-I-N-G-S

8:30 a.m.

CHAIRMAN SUZUKI: The Dental Products Panel of the CDRH Medical Devices Advisory Committee. My name is Jon Suzuki. I'm serving as the Chairman of the Dental Panel. And I would like to call this meeting to order.

The Executive Secretary, Michael Adjodha, will make some introductory remarks.

Mr. Adjodha?

EXECUTIVE SECRETARY ADJODHA: Thank you, Chairman Suzuki.

My name is Michael Adjodha, Executive Secretary of the Dental Products Panel.

Allow me to introduce the members of our panel. Please raise your hand as I call your name.

The Chairman of the panel is Dr. Jon B. Suzuki. Chairman Suzuki is a periodontist and immunologist, and is the Associate Dean of the School of Dental Medicine at Temple University in Philadelphia, Pennsylvania. Note that change from the agenda. This change is recent.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Joining him are the following panel
2 members:

3 Dr. Salomon Amar is a periodontist and is
4 Professor at the Department of Periodontology and
5 Oral Biology of Boston University, Boston,
6 Massachusetts.

7 Dr. David L. Cochran i s a periodontist
8 and is Chair of the Department of Periodontics at the
9 Health Science Center at the University of Texas, San
10 Antonio, Texas.

11 Ms. Elizabeth Howe is a consumer
12 representative and is the Outreach Coordinator for
13 the National Foundation for Ectodermal Dysplasias in
14 Auburn, Washington.

15 Ms. Allison F. Lawton is our drug
16 industry representative and is Senior Vice President
17 for Genzyme Corporation, Cambridge, Massachusetts.

18 Dr. William J. O'Brien is a materials
19 engineer and is Professor at the School of Dentistry
20 at the University of Michigan, Ann Arbor, Michigan.

21 Mr. Daniel R. Schechter is the Device
22 Industry Representative and is General Counsel for

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Parkell, Inc., Farmingdale, New York.

2 Dr. Domenick T. Zero is a cariologist and
3 is Chairman of the Department of Preventative and
4 Community Dentistry at Indiana University,
5 Indianapolis, Indiana.

6 And Dr. John R. Zuniga is an oral surgeon
7 and is Professor at the School of Dentistry of the
8 University of North Carolina at Chapel Hill, Chapel
9 Hill, North Carolina. Dr. Zungia is recovering from
10 an automobile accident and we're pleased he could be
11 with us today.

12 Joining the Panel members if the
13 following consultant: Dr. Inder J. Sharma is a
14 biostatistics consultant and is an Associate
15 Professor at the Department of Community Health and
16 Preventative Medicine of Morehouse School of
17 Medicine, Atlanta, Georgia.

18 Joining us at the table is Dr. Susan
19 Runner, Deputy Director of FDA's Division of
20 Anesthesiology, Infection Control, General Hospital,
21 and Dental Devices.

22 I will now read into the record a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 memorandum from the Center Director regarding voting
2 status of our Panel Consultant.

3 Pursuant to the authority granted under
4 the Medical Devices Advisory Committee charter, dated
5 October 27, 1990 and as amended on April 20, 1995, I
6 appoint the following consultant as a voting members
7 of the Dental Products Panel for the meeting to be
8 held on Tuesday, July 13, 2004. Inder J. Sharma, PhD

9 For the record, this individual is a
10 special government employee and is a consultant to
11 this Panel under the Medical Advisory Committee. He
12 has undergone customary conflict of interest review
13 and he has reviewed the material to be considered for
14 the meeting. Signed Daniel G. Schultz, MD, Acting
15 Director Center for Devises and Radiological Health,
16 July 8, 2004.

17 Next I'll read into the record a conflict
18 of interest statement for the this meeting.

19 The following announcement addresses
20 conflict of interest issues associated with this
21 meeting and is to be a part of the record to preclude
22 even the appearance of impropriety.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 To determine if any conflict existed, the
2 agency reviewed the submitted agenda for this meeting
3 and all financial interests reported by the Committee
4 participants. The conflict of interest statutes
5 prohibits special government employees from
6 participating in matters that could affect their or
7 their employer's financial interests. The Agency has
8 determined, however, that participation of certain
9 members and consultants, the need for whose services
10 that waives the potential of conflict of interest
11 involved is in the best interest of the government.

12 Therefore, waivers have been granted for
13 Drs. Cochran, O'Brien and Sharma for their interests
14 in firms that could potentially effect the panel's
15 recommendations.

16 Dr. Cochran's waiver involves a grant to
17 his institution for the sponsor study for which he
18 had no knowledge of the funding and had no
19 involvement in the data generation or analysis.

20 Dr. Cochran's waiver is limited in that
21 it allows him to participate in the panel discussion
22 but excludes him from voting.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Dr. O'Brien's waiver involves a grant to
2 his institution for the sponsor;s study for which he
3 had no knowledge of the funding and no involvement in
4 the data generation or analysis. Dr. O'Brien's
5 waiver is limited in that it allows him to
6 participate in the panel discussion but excludes him
7 from voting.

8 Dr. Inder Sharma's waiver involves a
9 philanthropic contribution from the firm at issue at
10 his institution for which he has no involvement and
11 is uncompensated.

12 Dr. Sharma's waiver allows him to
13 participate fully in today's deliberation. Copies of
14 these waivers may be obtained from the Agency's
15 Freedom of Information Office, Room 12A-15 of the
16 Parklawn Building.

17 We would like to note for the record, the
18 Agency took into consideration on other matters
19 regarding Dr. Domenick Zero. This panelist reported
20 past and current interest involving firms at issue,
21 but are matters that are not related to today's
22 agenda. The Agency has determined, therefore, that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 this panelist may participate fully in all
2 discussions.

3 In event that the discussions involve any
4 other products or firms not already on the agenda for
5 which a FDA participant has had financial interests,
6 the participant should excuse him or herself from
7 such involvement and exclusion should be noted for
8 the record.

9 With respect to all participants we ask
10 in the interest of fairness that all persons making
11 statements or presentations disclose any current or
12 previous financial involvement of any firm whose
13 products they may wish to comment on.

14 I'd like to request that everyone in
15 attendance at this meeting take the time to sign the
16 attendance sheet available at the front door.

17 Now transmitting you back to Chairman
18 Suzuki.

19 CHAIRMAN SUZUKI: Okay. Thank you. I
20 note for the record that voting members resent
21 constitute a quorum as required by 21 CFR Part 14.

22 We will now proceed the first of two open

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 public hearing sessions for this meeting. The
2 second open public session will follow the panel
3 discussion this afternoon. At these times public
4 attendees are given an opportunity to address the
5 panel to present data or views relevant to the
6 panel's activities. No individual has given advance
7 notice of wishing to address this panel. If there's
8 anyone now wishing to address the panel, because
9 identify yourselves at this time. Okay. Thank you.

10 I'd like to remind public observers at
11 this meeting that while a portion of this meeting is
12 open to the public observation, public attendees may
13 not participate except at the specific request of the
14 Chair. You will be given no more than 10 minutes for
15 your presentation.

16 I would like to ask at this time that
17 persons addressing the panel come forward to the
18 microphone and speak clearly, as the transcriptist
19 is dependent on this as a means for providing an
20 articulate transcription of the proceedings of this
21 meeting.

22 If you have a hard copy of your talk

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 available, please provide it to the Executive
2 Secretary for use by the transcriptist to help
3 provide an accurate recording of these proceedings.

4 We're also requesting that all persons
5 making statements during the open public hearings
6 disclose if they have financial interests with the
7 sponsor of the products under consideration.

8 Before making your presentation to the
9 panel, in addition to stating your name and
10 affiliation, please state the nature of your
11 financial interest in the product under
12 consideration, including who is paying for your
13 attendance at this meeting.

14 Okay. At this time we'll follow the
15 agenda and we will present with the sponsor
16 presentation on the product GEM 21S. Mr. Mark
17 Citron.

18 MR. CITRON: Good morning. My name is
19 Mark Citron. I'm Vice President of Regulatory
20 Affairs at BioMimetic Pharmaceuticals.

21 On behalf of BioMinetics we would like to
22 thank the panel and the FDA for the time and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 attention that the FDA and the panel have spent in
2 reviewing our PMA and meeting today to provide your
3 recommendation regarding approval of our device.

4 We have the privilege today to present to
5 you the results of decades of what began as
6 scientific research, progressed to product
7 development and clinical trials leading to today's
8 presentation of the GEM 21S control comparison
9 randomized study results. For the next 60 minutes
10 we will present these preclinical and clinical
11 results and respond to any questions you may have.

12 I will begin by introducing today's
13 speakers and our agenda.

14 First, Dr. Samuel Lynch, President and
15 CEO of BioMinetic will provide the brief overview of
16 the GEM 21S device and the development of the device.

17 Dr. Lynch is a periodontist and has conducted
18 extensive scientific research on PDGF as well as
19 other growth factors involved in tissue repair
20 covering many years.

21 Next, Dr. William Giannoble of the
22 University of Michigan will speak on the mode of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 action of GEM 21S with particular emphasis on the
2 growth factor component recombinant human platelet-
3 derived growth factor.

4 Dr. Ron Nevins, a clinical professor at
5 the Harvard School of Dental Medicine, who is also in
6 private practice, will present the animal and human
7 histology data.

8 Dr. Bob Genco, who is currently the
9 Distinguished Professor of Oral Biology and
10 Microbiology at the State University of New York at
11 Buffalo and recently appointed the Vice President of
12 Research at the State University of New York at
13 Buffalo will present the results of the randomized
14 control clinical trial.

15 Finally, Dr. Lynch will provide
16 concluding remarks to the formal presentations.

17 We welcome the panel's questions, and we
18 have available today several of the key scientific
19 researchers who have been involved in the GEM 21S
20 program, and they are prepared to respond to your
21 questions. These include the study statistician Dr.
22 Phil Lavin. He's an Associate Professor of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Biostatistics at Harvard Medical School and President
2 of Averion, which is a biostatistics consulting firm.

3 We have Dr. Charles Hart, our Vice
4 President and Chief Scientific Officer.

5 Dr. Jeffrey Hollinger, the Director of
6 the Carnegie Mellon University's Bone and Tissue
7 Engineering Center.

8 Dr. Michael Reddy, a clinical professor
9 at the University of Alabama, Birmingham.

10 And finally Dr. Mark Reynolds from the
11 University of Maryland Dental School.

12 Dr. Lynch will now begin.

13 DR. LYNCH: Thank you, Mark. And good
14 morning to the panelists, members of the audience and
15 the FDA.

16 I would also like to thank the panel for
17 your time and consideration today as well as the FDA
18 for their support and recommendations during the
19 development of GEM 21S.

20 I believe it is important to note that
21 our meeting today is the culmination of over 15 years
22 of scientific research by multiple investigational

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 groups working independently and sometimes
2 collaborative. We are fortunate to have many of
3 these research groups represented here today.

4 Two persons who are not here this morning
5 but who deserve substantial credit for the
6 development of the GEM 21 product, and who I would
7 like to take this opportunity to acknowledge and
8 thank, are Dr. Ray Williams, Chairman of
9 Periodontology at the University of North Carolina
10 and formally Chairman of Periodontics at Harvard. My
11 mentor, counselor and friend.

12 And posthumously, Professor Harry
13 Antaniales, whose lab conducted much of the early
14 research on PDGF, who inspired much more of the
15 scientific work in this field and in whose lab I
16 trained.

17 Finally, I would wish to acknowledge my
18 appreciation to the Biomedics Clinical and Regulatory
19 team for their hours of preparing the PMA submission
20 before you today as well as the entire GEM 21 group
21 of clinical investigators who rigorously conducted
22 the pivotal clinical study from both academic

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 research centers and private clinical practices
2 thereby providing us robust data from both important
3 clinical environments.

4 We are fortunate today to have three
5 individuals who are involved in the pivotal clinical
6 trial and who are widely recognized for their
7 expertise in clinical and basic scientific research
8 to speak in favor of the approval of GEM 21.

9 Our first speaker today is Dr. William
10 Giannoble of the University of Michigan and Director
11 of the Michigan Center for Oral Health Research. Dr.
12 Giannoble was a clinical investigator in the GEM 21
13 pivotal clinical trial, and is a recognized expert on
14 the biology of growth factors including platelet-
15 derived growth factor or PDGF.

16 As Mark mentioned, Dr. Giannoble will
17 discuss the mode of action of GEM 21S with particular
18 emphasis on the protein growth factor component.

19 Next Dr. Ron Nevins, a former President
20 of the American Academy of Periodontology and
21 currently the editor and chief of the *International*
22 *Journal of Periodontics and Restorative Dentistry* who

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 also finds time for a busy private practice will
2 present will present the animal and human histology
3 data demonstrating the effectiveness of GEM 21 in
4 promoting periodontal regeneration including new
5 cementum and periodontal ligament coronal to the
6 original apical extent of calculus.

7 Dr. Nevins is uniquely qualified for this
8 presentation, having participated in a GEM 21S
9 pivotal trial also as well as having been the lead
10 investigator for many studies evaluating the human
11 histological response to a number of different
12 drafting materials including PDGF and periodontal
13 bone defects.

14 And finally, Dr. Bob Genco, past
15 President of the International Association of Dental
16 Research and editor and chief of the *Journal of*
17 *Periodontology* will present the results of our
18 randomized control double blinded prospective multi-
19 center pivotal clinical trial.

20 Dr. Genco was the independent medical
21 director for the overall GEM 21S clinical program and
22 has many years of experience in designing, conducting

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 and evaluating the scientific integrity of clinical
2 trials related to periodontology including having
3 served as the formal of this august FDA Advisory
4 Panel.

5 Let me now set the stage for these
6 speakers by briefly describing the GEM 21 product,
7 it's development history and the unmet clinical need
8 that it is designed to satisfy.

9 Next.

10 GEM 21S, as we have alluded to,
11 principally consists of two main components. One
12 component is a particulate beta-tricalcium phosphate
13 or Beta-TCP, which is filled into a cup and
14 terminally sterilized.

15 The other principal component is a
16 physiologic solution containing recombinant platelet-
17 derived growth factor, which is aseptically-filled
18 into a syringe just to facilitate handling of the
19 material.

20 At the time of the surgical procedure,
21 the surgeon or surgical assistant simply peels back
22 the lid of the cup, adds the growth factor solution

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 to fully wet the graft particles. After a few
2 minutes sitting on the surgical tray, specifically
3 we're recommending approximately 10 minutes, the
4 material then forms a cohesive mass of particles
5 which are then packed into the alveolar bone defect.

6 Next, please.

7 One of the main and principle attributes
8 that we would like to stress today is the extensive
9 scientific research known about both principal
10 components of this product, both the PDGF and the
11 Beta-TCP.

12 There are well over 200 publications on
13 PDGF that deal specifically with its beneficial
14 effect on wound healing. These studies have been
15 conducted in a variety of models and systems
16 including in vitro self-culture systems using primary
17 cultures of osteoblast or well qualified osteoblast
18 like cellnoids, primary cultures of periodontal
19 ligament cells and gingival fibroblast cells and
20 many, many other cell types.

21 All of these studies in vitro have
22 clearly demonstrated the receptor binding of the BDGT

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 to the receptors, as you will hear from Dr.
2 Giannoble.

3 In addition, there are multiple
4 publications showing the beneficial effect of PDGF on
5 a wound healing in vivo in mice, rats, rabbits,
6 canines, swine, nonhuman primates and human clinical
7 trials. As you can see, it's a very well studied
8 molecule.

9 In addition, PDGF was the first
10 recombinant human growth factor to be FDA approved as
11 a wound healing agent. It is currently marketed under
12 the trade name Regranex by Johnson & Johnson. Has
13 been on the market for over 5 years and is absolutely
14 well documented safety record with no elicitation of
15 antibodies or any adverse responses in commercial
16 use.

17 In addition, the beta-tricalcium
18 phosphate has is an FDA cleared bone augmentation
19 device. It is the Beta-TCP that we incorporate into
20 GEM 21S. Is on the market in a larger particle form
21 under the trade name Vitoss by Orthovita for
22 orthopedic bone regeneration procedures.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Next please.

2 Thus, as you will hear this morning the
3 benefits of GEM 21S are it is a fully synthetic bone
4 regeneration system supported by hears of research
5 that have elucidated its mechanism of action and
6 demonstration a strong safety profile. And again,
7 rigorously conducted clinical trials and commercial
8 use.

9 The PDGT component has specifically been
10 shown to enhance periodontal regeneration in both
11 animals and humans. Our pivotal clinical trial has
12 demonstrated that the product accelerates the
13 attachment level gain and enhances or improves
14 significantly radiographic evidence for bone
15 regeneration.

16 Finally, we hope to show today that this
17 product demonstrates minimal risk and has the
18 potential for strong benefits in clinical practice.

19 Thank you very much. And I would now
20 like to turn the presentation over to Dr. William
21 Giannobbe to discuss the biology mechanism of action
22 and highlights of some preclinical data on the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 product. Thank you.

2 DR. GIANNOBLE: Thank you, Dr. Lynch.

3 And I'd also like to thank the FDA and
4 the FDA panel members for the opportunity for me to
5 present to you this morning some of the basic biology
6 in the extent of preclinical data that have
7 demonstrated some of the safety and effectiveness of
8 platelet-derived growth factor the GEM 21S system for
9 the promotion of periodontal regeneration.

10 So as we look at periodontal disease,
11 which typically it's a disease that results from a
12 microbial infection that leads to the resorption of
13 alveolar bone through to its cementum and periodontal
14 ligament. There are a variety of different factors
15 that appear to be critically important to the
16 reconstruction of periodontal wounds; those being the
17 appropriate cells within the lesion, they can
18 repopulate the wounds such as osteoblast, cemental
19 blasts, periodontal ligament fibroblasts within the
20 presence of the appropriate scaffold that will then
21 allow cell ingrowth and vascular invasion into the
22 lesion.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 And then the usage of signaling molecules
2 or growth factors that can direct the migration of
3 cells into the wounds, promote proliferation of the
4 cell types within the defect and stimulate matrix
5 biosyntheses.

6 In addition, given that the structure is
7 a vascular, it is critical to provide an angiogenic
8 environment to promote new blood vessel formation to
9 reconstruct these periodontal wounds. And so in my
10 presentation this morning I will focus on platelet-
11 derived growth factor and the scaffold, the
12 osteoconductive scaffold beta-tricalcium phosphate
13 for use in promoting periodontal regeneration.

14 So as we look at the two key components
15 of the GEM 21S system, the first being the Beta-TCP,
16 as this is an osteoconductive scaffold that promotes
17 cell attachment ingrowth, it also has been
18 demonstrated to prevent soft tissue collapse into the
19 soft tissue defects, and also facilitates blood clot
20 stabilization during the initial wound repair
21 process.

22 Recombinant human platelet-derived growth

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 factor BB there is a very extensive profile in terms
2 of its demonstrated ability to promote chemotaxis of
3 the key cell types involved in tissue repair. It is
4 also mitogenetic or promotes proliferation of these
5 various cell types such as periodontal ligament
6 fiberblast and osteoblast.

7 And PDGF has also been demonstrated to be
8 to be an angiogenic molecule by recruiting smooth
9 muscle cells that are important in the formation of
10 new blood vessels.

11 Next slide.

12 So to go into a bit more depth on beta-
13 tricalcium phosphate, this is a synthetic purified
14 calcium phosphate ceramic that has a very extensive
15 history in the FDA as well as a device used in
16 dentistry and in orthopedic applications as a bone
17 void filler. And in this long history of usage there
18 have been no demonstrated adverse events utilizing
19 beta-tricalcium phosphate as a bone void filler in
20 these varieties of applications. And recently the
21 FDA Advisory Panel recommended a reclassification of
22 Beta-TCP from a high risk device to a lower risk

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 device for use in dental applications.

2 This slide demonstrates scanning electron
3 microscopic views of beta-tricalcium phosphate at
4 lower power magnification here and a higher power
5 magnification. The low power view demonstrates the
6 beta-tricalcium phosphate granules which in the
7 formulation for the GEM 21S system range in particle
8 size from 250 to 1,000 microns in diameter. This
9 higher magnification view demonstrates the very open
10 pore structure of the Beta-TCP used in the GEM 21S.
11 It has a 90 percent open pore structure which then
12 this porosity, this ranging from 1 to 1,000 microns
13 in diameter thus allows cellular ingrowth and
14 vascular invasion. This lower panel demonstrates at
15 a different microscopic view the growth of osteoblast
16 like cells on top of the beta-tricalcium phosphate
17 demonstrating that it does promote cell attachment
18 and proliferation on the device.

19 Recombinant human platelet-derived growth
20 factor has been an extensively studied molecule in
21 the area of wound healing. So it's a natural wound
22 healing hormone released from platelets during normal

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 wound repair.

2 The scientific established mode of action
3 is that PDGF has been demonstrated to promote
4 connective tissue formation, also osteogeneses and
5 angiogeneses by the induction of vascular endothelia
6 growth factor in the recruitment of smooth muscle
7 cells.

8 This diagram depicts the binding of
9 platelet-derived growth factor, which is a dyncric
10 protein which binds to cell surface associated
11 tyrosine kinase receptors. These receptors dimerize
12 and then elicit autophosphorylation of the receptor.

13 This autophosphorylation event then leads to a
14 variety of different signal transduction pathways
15 which will then led to the elicitation of the variety
16 of different biological effects such a mitogenesis or
17 cellular proliferation, directed cell mitigation or
18 chemotaxis, and also the blocking of program cell
19 death or promoting cell survival.

20 So PDGF more specifically as we examine
21 its ability to promote periodontal regeneration
22 within the periodontia, platelet-derived growth

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 factor and its associated receptors are naturally
2 induced during normal tissue repair, both soft tissue
3 repair and during the fracture healing procedure.

4 PDGF has been demonstrated to be
5 chemotactic for a variety of cells derived from the
6 periodontia as well as promoting cellular
7 proliferation and matrix biosynthesis. And there is
8 a large body of work supporting the variety of
9 effects as shown here.

10 PDGF also promotes cell survival since a
11 PDGF alpha receptor encodes for a growth arrest
12 specific gene. So PDGF will promote or prevent
13 apoptosis or programmed cell death.

14 PDGF also enhances angiogenesis
15 specifically by promoting the proliferation of smooth
16 muscle cells or parasites around the newly formed
17 blood vessels and it compliments the actions of VEGF
18 or vascular endothelia growth factor that's
19 critically important for blood vessel formation and
20 maturation.

21 This slide published by the San Antonio
22 group demonstrates the effects of recombinant human

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 PDGF in an artificial wound model on promoting cell
2 repopulation. And so what we can see in this slide
3 looking at percent wound fill or cell repopulation of
4 periodontal ligament fibroblast versus a low serum
5 control, this graphic demonstrates that over a period
6 of ten days the significant increase in cellular
7 repopulation into artificial wound defects by the
8 application of recombinant human platelet-derived
9 growth factor.

10 Next slide.

11 This slide demonstrates the ability of
12 platelet-derived growth factor applied onto the beta-
13 tricalcium phosphate osteoconductive device for its
14 release and then subsequent biological activity of
15 the release PDGF. And so this slide shows treated
16 Thymidine incorporation as a method to determine DNA
17 synthesis over time when PDGF has been applied to the
18 beta-tricalcium phosphate device. And so what one
19 can note is that there is a rapid release over the
20 first 24 hours and the PDGF that is released is
21 indeed biologically active as measured of the
22 promotion of DNA synthesis.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 This slide demonstrates results from an
2 in vivo animal study done in Beagle dogs where
3 fenestration bony defects were created on two
4 surfaces and then autoradiography was performed to
5 look at cells that were demonstrating active
6 proliferation within the periodontal wound
7 compartment.

8 So the variety of different cell types
9 examined were those important in periodontal repair
10 such as fibroblasts, cementoblasts, osteoblasts,
11 perivascular and endothelial cells. And what was
12 noted that it was compared to control or surgery
13 alone defects, PDGF promoted at least a three to five
14 full increase in cellular DNA synthesis as noted by
15 the autoradiography. And you can see this in a
16 multitude of different cell types that were found
17 within the lesions, thus demonstrating that the PDGF
18 has pleiotropic effects on promoting a variety of
19 parameters associated with periodontal regeneration.

20 This slide published by Bob Genco's group
21 several years ago in a canine model of surgically
22 created critical size defects in dogs. These are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 class 3 furcation defects that do not typically heal.

2 The defects were treated with guided
3 tissue regeneration, a standard treatment modality
4 for periodontal regeneration versus recombinant human
5 platelet-derived growth factor applied to the tooth
6 root surface combined with GTR. And looking at
7 histomorphen metric analysis to determine the amount
8 of regeneration that occurred within the defects,
9 what was noted was that PDGT strongly augmented the
10 degree of newly formed bone and periodontal ligament,
11 while at the same time blocking really the production
12 of the granulation tissue or scar formation that
13 resulted after this healing period.

14 This slide demonstrates the potent
15 effects of platelet-derived growth factor on
16 promoting osteogenesis. This is a study published
17 several years ago that examined an osteoporosis model
18 where female rats were ovariectomized which induced a
19 rapid bone loss. And the slides on the left
20 demonstrate the metathesis of the tibia in these
21 animals either in an osteoporosis saline control or
22 animals that were delivered a three times per week

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 infusion of 2 milligrams of recombinant human
2 platelet-derived growth factor.

3 What one can note from these long bones
4 was that there was a significant increase in the
5 boning trabecular in both the primary and secondary
6 spongiosa in these bones that were treated with --
7 these animals that were treated with recombinant
8 human PDGF.

9 Using histomorphologic metric analysis of the
10 vertebral body and then tibial metatarsals one could
11 also note a statistically significant improvement in
12 bone density measures nearly two-fold in both of
13 these different bony sites.

14 This slide now demonstrates the platelet-
15 derived growth factor's ability to promote
16 periodontal regeneration. This is a natural disease
17 model in the Beagle dog that will result in loss of
18 connective tissue and alveolar bone. So this slide
19 demonstrates a through and through class 3 furcation
20 defect that typically will not heal on its own.
21 These animals were delivered a single application of
22 recombinant human platelet-derived growth factor in a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 beta-tricalcium phosphate carrier, and this slide
2 shows six weeks after this single application of PDGF
3 plus the Beta-TCP. The promotion of new alveolar
4 bone, a periodontal ligament and tissue consistent
5 with cementum.

6 These various preclinical animal studies
7 performed in dogs were also followed by in nonhuman
8 primates in the monkey model *Macaca Fascicularis*.
9 And what this side is demonstrating is the
10 consistency of effects in the animal model *Macaca*
11 *Fascicularis* versus humans when platelet-derived
12 growth factor was combined with insulin like growth
13 factor one. So this study looked at animals that
14 were treated with a single application. If you look
15 at the parameter of ostis defect fill, there is a
16 striking similarity between the monkey model and this
17 is -- the human data here is derived from a multi-
18 center trial done, it was a phase 1 phase 2 trial
19 done at the Harvard School of Dental Medicine and at
20 the University of North Carolina. And essentially
21 the bottom line of this study was demonstrating that
22 similarity between the animal model and the human.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Next slide.

2 The next few slides will now demonstrate
3 the extensive track record for the various components
4 used for the GEM 21S product. And with the Regranex
5 product that has been FDA approved for, it's been
6 over five years now, it has a very extensive safety
7 record. And so the results shown here are actually a
8 compilation of six randomized controlled trials where
9 the Regranex product demonstrated extensive safety.
10 There was no neutralizing antibodies that were
11 developed. And these patients received the treatment
12 of the Regranex every other day for up to 140 days of
13 a concentration of 100 milligrams per mil of the
14 PDGF. And so to date there have been at least 17
15 million doses applied of Regranex, demonstrating its
16 safety.

17 Also you have provided to you very
18 extensive confirmatory biocompatibility tests. As
19 you can see on the list here, in terms of
20 cytotoxicity, sensitization, acute systemic toxicity,
21 genotoxicity and muscle implantation for the GEM 21S
22 product. And so all of these tests have demonstrated

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 that GEM 21S is both biocompatible and safe.

2 So what I would like to summarize for you
3 this morning is that we have demonstrated that the
4 mechanism of action of platelet-derived growth factor
5 is well established as shown in vitro studies as well
6 as in vivo applications demonstrating its potent
7 ability to promote periodontal regeneration, i.e.,
8 tooth group cementum, periodontal ligament and
9 alveolar bone.

10 This is also a very long history of
11 safety for both of the components, the beta-
12 tricalcium phosphate in both dental and orthopedic
13 applications and the platelet-derived growth factor
14 component, i.e, in the Regranex product for the
15 treatment of neuropathic diabetic ulcers.

16 The results have also been demonstrated
17 to be quite consistent amongst the large body of
18 research done with a variety of different clinical
19 investigator reclinically that bridge and demonstrate
20 consistency to some of the human clinical studies
21 that have been performed.

22 I would like to thank you for your

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 attention and I look forward to the discussion this
2 afternoon.

3 I will now introduce Dr. Myron Nevins who
4 will present the proof of principle data on the
5 ability of platelet-derived growth factor to promote
6 periodontal regeneration in humans.

7 DR. NEVINS: Good morning.

8 I'd like to take this opportunity to
9 thank the FDA and the panel by allowing us to
10 demonstrate the evidence of regeneration, periodontal
11 regeneration that we've been able to achieve with GEM
12 21S.

13 The definition of periodontal
14 regeneration is histologic. It has evolved from
15 proceedings of two world workshops in clinical
16 periodontics and it is inclusive of information of
17 new bone, new cementum connected by a functional
18 periodontal ligament on a root surface that has
19 previously been pathologically exposed.

20 The hierarchy of evidence in periodontal
21 regeneration has taken years to evolve, but because
22 of the histologic definition, it's clear that the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 most compelling evidence are human studies that have
2 histologic evaluation. In lieu of the obvious
3 difficulties in obtaining this information, other
4 means have surfaced to measure success, including
5 RCTs with measuring clinical invaded ethic
6 parameters. Perhaps the more contemporary benchmark
7 has been the use of, be it surgical reopening, which
8 would more closely mimic the radiographs in terms of
9 interpretation.

10 Proven principle assessment is
11 established to demonstrate the safety and the
12 effectiveness of a product. Safety would determine
13 histologic tissue reactions, healing response and
14 provide a clinical assessment for safety.

15 Effectiveness provides human histologic evidence of
16 regeneration, in this case for vertical intrabony
17 defects and also for Class II furcation invasion
18 problems.

19 This design include 11 intra-osseous
20 defects around teeth scheduled for extraction, six
21 intrabony and five Class I furcation defects were
22 treated. They were treated with a combination of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 recombinant BDGF plus a carrier. At nine month post-
2 operative follow-up recordings of the CAL, the pocket
3 depth gingival recession and linear bone -- were
4 recorded. At that time the teeth were abstracted
5 with a small amount of surrounding tissues and
6 submitted to blind histologic analysis to assess
7 regeneration.

8 I should mention at this moment that
9 informed consent was obtained from the patients. The
10 patients were rehabilitated from the site with bone
11 crafting, dental implants and a prothesis to
12 reestablish or in all senses to provide them with a
13 dental solution that they would not have been able to
14 have otherwise.

15 The intrabony defect results demonstrate
16 a pocket depth, a mean pocket depth of 9.7
17 millimeters and at nine months at time of the
18 harvesting block, 3.3 millimeters. This is a change
19 from baseline of 6.42.

20 The importance of this is related to
21 length of roots in a human model. The smallest, the
22 shortest roots are the incisors, the central incisors

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 with 11 millimeters and of course the longest would
2 be the cuspids which approximate 18 millimeters.

3 If we accomplish a 6 millimeter
4 correction, this definitely changes the prognoses of
5 the tooth.

6 The CAL gain started -- level started at
7 11.1 baseline and was 4.9 at nine months. Once
8 again, for a change from baseline of 6.7, which is
9 consistent with the pocket depth reduction.

10 Bone height change shows radiographically
11 a 2.14 improvement.

12 This will become in a few minutes when we
13 look at the histologic measurements, because there
14 will be a correlation between what was here
15 radiographically and histologically.

16 The furcation defects from a pocket depth
17 began at 6.2, in nine months were 2.8 for a change
18 from baseline of 3.4. And the clinical attachment
19 level changed with a change from a baseline of 4.

20 Since these were horizontal as well as
21 vertical probing depths, these are very significant
22 in reversing the invasion of the furcation by

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 inflammatory periodontal disease.

2 When the teeth are described as being
3 candidates for extraction, it's most effective when
4 we look at a clinical photograph. And here we see a
5 maxillary cuspid with bone defects both in vertical
6 dimension and bone morphology that would be serious
7 candidates for extraction.

8 We determine the level of the root that
9 has been exposed to disease by the presence of
10 calculus. So at the base of the calculus a notch is
11 made with a small burr to designate that we actually
12 have evidence that disease occurred at that point.

13 The calculus, of course, is removed
14 before we continue on to place the crafting material.

15 Next.

16 After nine months when the block was
17 removed or harvested, we now have an opportunity to
18 witness the histology in evidence.

19 Here we see a lower power and we're going
20 to observe three different situations. One, the area
21 of the notch where we can see new cementum, new bone
22 and a new mature well vascularized periodontal

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 ligament. If you look closely, you can see sparky
2 fiber attachments on both the bone and the cementum
3 side. So the area of the notch which was placed at
4 the base of the calculus has responded appropriately.

5 Now, in the next observation we'll look
6 at mid-root and then we'll look at the mouth of the
7 defect.

8 As we move occlusally we again witness
9 new bone, new periodontal ligament and a functional
10 vascular periodontal ligament with clear evidence of
11 sparky fiber attachments on both sides indicating its
12 function.

13 Coming to the mouth of the defect, the
14 new cementum has come all the way to the beginning of
15 the bone defect and we have new bone and, once again,
16 the functional periodontal ligament with supercrestal
17 fibers that show very little evidence of any
18 inflammatory infiltrate.

19 This completes the picture of that cuspid
20 that we witnessed.

21 Looking at a second vertical defect just
22 to demonstrate quickly that this occurred more than

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 one time, we again observe a notch. We see in the
2 notch new cementum, new bone connected by a
3 functional vascular periodontal ligament, and in fact
4 we have new cementum and new bone all the way to the
5 top of the defect.

6 We have studied several different
7 materials, but this astounding to see complete
8 regeneration of the defect.

9 Next.

10 However, the most exciting observation
11 that we encountered was the response in Class II
12 furcations where there has been evidence to suggest
13 that we fulfil the definition of periodontal
14 regeneration with any of the materials that are
15 presently available.

16 The notch designated the extent of the
17 calculus and if we take the excerpt from the box, we
18 see new cementum and new bone connected by a new
19 functional periodontal ligament, again with evidence
20 of sparky fiber attachments and no evidence of
21 epithelium.

22 The outstanding observation in my

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 estimation is that even at the fungus of the
2 furcation there is no evidence of epithelium and
3 we've completely resolved the definition of
4 periodontal regeneration with new cementum, new bone
5 and new periodontal ligament. This offers us the
6 opportunity clinically to provide resolution for a
7 problem that escaped clinicians indefinitely.

8 Next.

9 The results and conclusions of this human
10 histologic evidence demonstrate safety; there's
11 normal bone and ligament remodeling. The clinical
12 measurements were demonstrated to be significantly
13 improved. Radiographs were consistent with bone
14 fill. We have no evidence of root resorption or
15 ankylosis.

16 Actually, the histo micromophy that was
17 performed very closely related in size or dimension
18 to the radiographic analysis that was performed.

19 There is no evidence of root resorption
20 or ankylosis whatsoever, so there is nothing to
21 discuss along those lines. And the histologic
22 evaluation we just saw revealed regeneration in both

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 the intrabony and Class II furcation defects.

2 It became time to design a pivotal study.

3 And in doing so, the transition was made to GEM 21S.

4 There were two issues: One, select a carrier and
5 the other to give some consideration to dosing.

6 Allograft was used for the histologic
7 study. Since it's not formally approved by the FDA,
8 a lot of questions -- a lot of producability remained
9 and consideration was given to trying alternatives.
10 Since beta-tricalcium phosphate and allograft were
11 shown to provide comparable delivery properties of
12 recombinant PDGF. The kinetics are similar and the
13 BDGF release from both matrices simulated bone cell
14 proliferation.

15 The study objectives were to compare the
16 in vivo performance of PDGF with the two carriers,
17 beta-tricalcium phosphate and allograft.

18 It was also to access the dose response
19 of the recombinant PDGF.

20 The study designed is a randomized
21 control blind trial and in canine this critical size
22 periodontal defects. Six defects were made in each

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 group and there is an eight week follow-up.

2 Looking at the results we can see that
3 the beta-tricalcium phosphate by itself demonstrated
4 new bone formation, but particles of the carrier
5 remained and obviously it has left a significant
6 portion of the furcation without periodontal
7 regeneration. However, when the product GEM 21S is
8 used the combination of the recombinant PDGF with the
9 beta-tricalcium phosphate received a notch and
10 received complete regeneration with cementum and
11 periodontal ligament indicating a much more favorable
12 response in the type of clinical end point goal that
13 we would hope to achieve for our patients.

14 Next.

15 Evaluating the results of the canine
16 study we see results with both TCP and allograft.
17 And it's clear that the dosage of .3 mg/ml with the
18 TCP and PDGF outperformed the other possibilities.

19 Next.

20 This led to the overall conclusions that
21 GEM 21S, a truly synthetic system is safe and
22 biocompatible with no risk of disease transmission.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 The BDGF when used with beta-tricalcium phosphate or
2 allograft significantly improved the periodontal
3 condition. This is measured in the formation of new
4 bone, new cementum connected by a functional
5 periodontal ligament.

6 There was sufficient evidence to now
7 initiate a pivotal clinical trial and the decision
8 was made to use the .3 mg/ml of BDGF because of the
9 greater effectiveness that was shown.

10 And now I have the pleasure of
11 introducing Dr. Robert Genco, the Director of the
12 Periodontal Disease Clinical Research Center at the
13 State University of New York in Buffalo. Bob has
14 carried out five phase three and pivotal trials of
15 periodontal products that have previously been
16 accepted the FDA, so he's an old hand at it.

17 DR. GENCO: Thank you, Dr. Nivens.

18 And I, too, would like to thank the panel
19 for your special efforts in reading that mass of
20 material that was submitted to you. I have the file;
21 files and files of those submissions and I know what
22 a tremendous effort it is.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 I would like to also thank the FDA, as
2 Sam did, for their help during the design and
3 analyses of the pivotal trial.

4 Now, I've worked in this area for about
5 15 years. One of the groups that Sam mentioned was
6 the Buffalo group that looked at BDGF and other
7 growth factors, and I have a tremendous interest in
8 seeing this come to the benefit of society. And I'm
9 very pleased to present this material today.

10 I have an official role with the company.
11 I'm the Chairman of their Scientific Advisory Board.

12 And longstanding interaction with Dr. Lynch.

13 I'd like to talk about the pivotal trial
14 and share some of the results with you, the
15 highlights the results. The next slide shows the
16 nature of the trial. It was a double -- prospective
17 randomized control trial with 180 patients randomized
18 to three treatment groups.

19 Group one was Beta-TCP plus 0.3
20 milligrams per mil of recombinant PDGF beta subunit.

21 Group two TCP plus 1 milligram per mil of
22 recombinant PDGF.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 And then group three is interesting in
2 that it's an active control. It's TCP with buffer,
3 no recombinant PDGF. And it's an active control in
4 that it's a product that's already on the market for
5 bone regeneration used in orthopedics extensively.
6 And we used a super fine fraction of that Vitoss.
7 And it's a newly formulated form with increase
8 porosity and increased surface area. And really it
9 hadn't been systematically tested in periodontal
10 disease. So, the design is really -- puts a high
11 hurdle to show an adjunctive or additional affect of
12 recombinant BDGF, and you'll see some of the results
13 that bare that out.

14 It's a six month follow-up study. And we
15 looked both at clinical and radiographic pinpoints.

16 The study was carried out in 11 centers,
17 four university centers and 7 private clinical
18 offices.

19 Next slide, please.

20 The investigators, patients, sponsor and
21 monitors and radiographic assessment was all masked.

22 The patients were randomized to one of the three

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 groups by a variable block design, and all of the
2 investigators; that is the examiners who were
3 separate from the operators, separate from the
4 surgeons, the examiners were calibrated, both at
5 baseline and at six months to ensure inter and intra
6 examiner standardization.

7 Next slide, please.

8 The study was independently monitored for
9 quality and safety performed by Target Health, and it
10 was independently analyzed by both Target Health and
11 Averion, Dr. Phil Lavin's company, and he's here.

12 Next slide please.

13 The key inclusion criteria were:

14 Age, 25 to 75 years; the pocket depth of the
15 treatment site had to be at least 7 millimeters deep
16 and had to have an intrabony defect at the time of
17 surgery of at least 4 millimeters.

18 Any configuration of pocket was allowed.
19 It could be 1, 2, 3 combination, combination with
20 circumferential and combination with Class I or II
21 furcation. So these are real live complicated complex
22 intrabony lesions.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 We allowed smokers who smoked up to one
2 pack per day. The rationale was that many, many
3 patients who suffer periodontal disease are smokers.

4 So we wanted to make sure that this was a treatment
5 that would work in smokers who are known to heal less
6 well than nonsmokers.

7 Next slide.

8 The key exclusion criteria included
9 pregnant women or women intending to become pregnant
10 during the study. This was not excluding women of
11 childbearing age. Only those that were pregnant,
12 lactating or intending to become pregnant.

13 History of oral cancer or HIV. Signs of
14 acute infection or abscess at the site, the test
15 site, Class III furcations, surgery under study too
16 from the last year; all of these were exclusion
17 criteria.

18 Next slide, please.

19 The outcome measures are very important
20 to comment to. Clinical attachment level was
21 assessed both at three months and at six months.
22 Linear bone growth and percent of bone fill are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 quantitative measurements of radiographs. These were
2 assessed as companion outcomes. And as we have heard
3 from the previous presentations, in this complex
4 disease, periodontal disease, the pathology involves
5 both soft tissue and hard tissue, so it's reasonable
6 from the clinical pathologic standpoint to assess
7 both tissues, hard and soft for clinical outcome.

8 Then we also used a composite outcome
9 where we blended or we merged, melded both the
10 clinical and the radiographic technique. And the
11 rationale for that is to get at this question of
12 clinical significance. To try to address the issue
13 of what percent of the target population benefitted
14 from therapy. It wasn't meant to look at statistical
15 significance to prove the efficacy. It was to get at
16 this very difficult question. I know I was on the
17 panel for a number of years. We always wrestled with
18 the question of is significant. Did it benefit a
19 significant portion of the target population? And
20 that's why we used this composite outcome.

21 We also looked the pocket depth
22 reduction, gingival recession and wound healing. And

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 at the direction or suggestion of the FDA we compared
2 it to currently approved FDA therapies that were sort
3 of comparable.

4 Next slide.

5 This is the study time table. I draw you
6 attention to day zero to the day the surgery was
7 carried out. At least or less than two weeks prior
8 to that baseline examination, examiner calibration
9 and radiographs had to be made. At least two months
10 prior to that the patients had to be screened,
11 informed consent obtained and an initial preparation
12 carried out.

13 After surgery the patients were followed
14 at three months and at six months, radiographs were
15 taken and all of the clinical measurements made both
16 at three and six months.

17 The next is a videoclip of the actual
18 preparation of the material. And this shows the dry
19 material in a dappen dish and the PDGF solution added
20 to it from the sterile syringe. And then the
21 material is next. And this is done approximately ten
22 minutes, at least ten minutes before the material is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 applied. You can see how the particles adhere to
2 each other, and it's actually a very easily managed
3 material to place in the mouth when you're having to
4 place the material in upper lesions or mandibular
5 lesions. It's actually quite easy to work with.

6 The next is a videoclip of the actual
7 surgical procedure. And this is from one of the
8 clinical sites. You can see the initial probing was
9 carried out. And the -- we'll get some -- there we
10 go. The videoclip shows the depth of the pocket. You
11 see the tissue is quite firm after the initial
12 preparation.

13 Then the root is thoroughly debrided.
14 The issue is removed. All the granulation tissue is
15 removed from the lesion. The lesion, you can see the
16 dimensions here. It's a 3 wall intrabony defect. The
17 root is cleaned absolutely clean. And then the root
18 is treated with tetracycline to condition it. And
19 then material is placed in the lesion to fill the
20 lesion to the brim. And it's packed gently into the
21 region.

22 You can see how easy it is to handle.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 The operators, the surgeons were all
2 standardized. They were standardized to a standard
3 way of making an incision, incision design, to a
4 standard debridement of the root, to a standard use
5 of tetracycline concentration, duration.

6 And as you can see here, see the incision
7 is a scalloped incision and we standardized the
8 suture technique so that the buckle flaps could be
9 opposed to get primary tension healing. Very, very
10 important in these regenerative techniques to make
11 sure that we get full coverage inasmuch as possible
12 of the lesion with the soft tissue.

13 Now the examiners, as I mentioned, they
14 were different than the surgeons. Different set of
15 people. They were calibrated. They were calibrated
16 to look for reproduceability of their own
17 measurements, and that's inter-examiner calibration.

18 And they were calibrated against a gold standard.
19 Someone on the research team who had an intrinsic low
20 error, all of the other examiners were calibrated
21 against that person to ensure consistency across
22 sites so we had more confidence to prove the data.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 The next slide shows the actual results.
2 The Kappa for the intra-examiner reproducibility was
3 0.94 and for the inter-examiner consistency was .89.
4 Both very, very high levels of reproducibility and
5 consistency exhibited by those Kappas.

6 Now the radiographic analysis. Care was
7 taken in that also. For example, the films were
8 taken under a uniform height quality field conditions
9 using the renperil system, and every investigator's
10 team was standardized to take these x-rays at a high
11 quality uniform way.

12 Then the films were sent to a central
13 site, University of Alabama, and Dr. Reddy and his
14 team used standardized techniques and validated
15 measurements which they and Dr. Genco for a decade
16 had developed over the years to measure both linea
17 bone growth as well as percent bone fill. And this
18 is really percent linear bone fill. It's not a
19 volume. It's a two dimensional measure.

20 And I'll show you those measures on the
21 next slide. This is a graphic diagram of the
22 radiograph and the landmarks that were measured at a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 synitho enamel junction, or if there is a filling, it
2 would be the apical portion of the filling,
3 restoration, root apex, crust of bone, base of defect
4 at baseline, and similar measurements at six months.

5 And then the next slide shows how these calculations
6 were made.

7 First, linear bone growth is very simply
8 the measurement from the CEJ to the base of the
9 pocket at baseline, and subtracted from that is the
10 measurement from CEJ at the base of the pocket six
11 months later. Now in this instance, it turned out to
12 be the original value of 6 millimeters, and at six
13 months 3 millimeters, so we have three millimeters of
14 linear bone growth. It's that simple, but very
15 precisely measured.

16 Now percent bone fill-in is simply the
17 linear bone growth divided by the initial depth of
18 the lesion. In this instance, it would be 50
19 percent, so the original depth of the lesion is 6
20 millimeters and the linear bone growth was 3, so you
21 had a 50 percent bone fill. It's a linear bone fill.

22 Now as a matter, because these are field

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 x-rays, they're not taken with stents or any other
2 precaution except good technique. There was a
3 control on elongation or foreshortening, and that was
4 the measurement of the CEJ to the apex, and that was
5 measured on all x-rays pre and post. If they varied
6 by 15 percent either way, then the x-rays were
7 adjusted. They were normalized, a very standard
8 technique used in radiographic analysis. We've used
9 it for years and it works quite well. Now in fact
10 that happened in less than 5 percent of the cases
11 which test to the quality of the x-rays site by site.
12 Next slide, please.

13 Now once the x-rays are sent to the site,
14 then there's a whole other set of calibrations and
15 measurement variability assessed; that is, the actual
16 measurement of the x-rays at the site. The
17 technician made repeated measurements on randomly
18 selected cases, and there was a less than 3 percent
19 variability between measurements, which is very good.

20 Following assurance, all radiographs were
21 then looked at, reviewed by an independent
22 periodontist. And, of course, all the radiographs

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 are blinded anyway so nobody in Alabama knew which
2 group they were from, but the independent
3 periodontist who was not connected to the study
4 looked at the x-rays and looked at the measurements
5 to see if they made sense; were there any really odd-
6 ball measurements. And they were occasionally some
7 measurements that didn't -- so those were remeasured,
8 so there was another level of control placed on the
9 measurements. Next slide.

10 Now the results. I'll first summarize
11 all the results in the next slide, and then go into
12 them individually. First of all, there were no
13 device-related serious adverse effects, an expected
14 result, but it had to be proven. You're using two
15 FDA approved products, put them together, both are
16 safe, together they're safe, but it had to be proven.

17 GEM-21S significantly improved, that is
18 statistically significantly improved CAL at three
19 months. It significantly improved CAL gained between
20 zero and six months. And the area under the curve
21 assessment showed that the three month gain was
22 maintained, it was a really accelerated healing which

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 was maintained at six months. The LPG, that's the
2 linear bone growth, was significantly improved at six
3 months, as was the percent bone fill at six months,
4 and these were highly significant in the .001 range.

5 And the GEM-21S exceeded the benchmarks of
6 effectiveness as compared to Emdogain, an FDA
7 approved product, PepGen P-15, an FDA approved
8 product, Allograft which is FDA allowed, not
9 necessarily approved but it's allowed, and open flat
10 debridement. Next slide, please.

11 DR. SHARMA: Excuse me. I want to just
12 clarify one thing there. These results you're
13 talking about, they are baseline to certain time
14 point.

15 DR. GENCO: That's right.

16 DR. SHARMA: Not to compare it with
17 different groups. Right?

18 DR. GENCO: I'll get into which group,
19 yes. It's the .3 milligram group. That dose group
20 showed these differences, and not the one. Right.
21 But I'll get into that in some detail. I just wanted
22 to give an overview, the result of my judgment and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 education; that Aristotle technique of tell them what
2 you're going to say, say it, and tell them what you
3 said. So I just told you what I'm going to say. Now
4 I'm going to say it.

5 DR. SHARMA: All right.

6 DR. GENCO: The number of subjects were
7 180, 178 finished with a 1 percent drop-off rate,
8 which is amazing for such a study. Forty-three
9 smokers, mean age 51, gender slightly more males than
10 females, approximately 60 percent Caucasian, the rest
11 distributed among Asian, African American and
12 Hispanic. Next slide.

13 Now baseline defect characteristics, the
14 message here - there were no significant differences
15 among treatment groups, and you can see this in the
16 data, this inspect pocket clinical attachment level,
17 bone defect, percent one wall, percent two wall,
18 percent three wall, circumferential. They're all
19 approximately the same, which you'd expect that
20 random variation you get by randomizing. No
21 statistically significant differences. And this is
22 extremely important as we'll see later, because the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 deeper the pocket, the more healing you're going to
2 get, so you really must have all the pocket depths at
3 the baseline comparable. Next slide.

4 Now let's look at total adverse events.
5 No significant differences among the treatment groups
6 with respect to any adverse events, serious,
7 potentially related, unrelated. For example,
8 subjects with at least one adverse event ranged
9 around 70 percent. Well, they all had surgery, and
10 what was that adverse event; pain after surgery,
11 which is to be expected. Not different between the
12 surgical control and the other treatment groups, so
13 there's no effect here of increasing adverse events
14 by adding PDGF to the TCP. Let's look at the serious
15 adverse events. They were present. They were not
16 different among the groups, but they were present.
17 Let's look at them. Next slide.

18 There were four serious adverse events,
19 none related to the study device; including
20 bronchitis, basal cell carcinoma, spinal fusion
21 surgery, and diabetic complications. These are
22 things as we all know in a six month study with 180

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 patients, you're going to get these adverse events -
2 not related to the device.

3 Now let's look at some of the measurement
4 data. One assessment is clinical attachment level
5 gained over three months and over six months
6 comparing the .3 milligram, the 1 milligram and the
7 TCP, and you can see that the .3 milligram was
8 statistically significantly different than the
9 control at three months. That 3.8 millimeter gain
10 was more or less maintained at six months. However,
11 what happened, I think, is that the TCP control
12 gained - and we saw this with the dog study too. The
13 control actually catches up to the treatment over
14 time, so now the difference between .3 milligram and
15 TCP is not statistically significant. And we'll look
16 at this another way looking at the area under the
17 curve analysis. Next slide, please.

18 Now if we compare the gain of GEM-21 with
19 the CAL gain of Emdogain and PepGen using the three
20 studies for Emdogain that were submitted to the FDA
21 and the two studies for PepGen, using those studies
22 as a baseline, you can see that 3.7 CAL gain versus

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 2.7 versus 1.7, at least it's comparable, highly
2 unlikely that GEM-21 under-performs, but at least
3 they're comparable. Next slide.

4 Now one of the problems with such
5 analysis, you have to really be careful as you all
6 know, is that the studies were not done head-to-head.
7 They are separate studies. We're talking about
8 three, two, in our study five different, six
9 different studies compared. And the possibilities
10 for making misinterpretations are great.

11 For example, if you look at the baseline
12 pocket depth, they're pretty comparable between our
13 study and the Emdogain studies, but look at the
14 PepGen study. They started out with shallower
15 pockets, so the comparison with PepGen is fraught
16 with difficulties, because they started with lower
17 pocket depth so they're going to get less healing.
18 And, in fact, that's what we saw. So we really have
19 to be very careful about the interpretations compared
20 to products on the market, and we are.

21 So what we say is they're comparable,
22 very unlikely that GEM-21 under-performs relative to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 the others, so I think that's a conservative way of
2 stating those comparisons. Next slide, please.

3 Now the area under the curve is commonly
4 used in wound healing studies, and its purpose is to
5 detect differences in CAL gain among the treatment
6 groups between baseline and six months. We're using
7 data from zero, three, and six months. Next slide,
8 please.

9 And here are the curves. The green line
10 is the 0.3 milligrams, and you can see it out-
11 performs the other two groups at 3 and at 6. Our
12 interpretation is that this is an early gain over the
13 other two groups, statistically significantly
14 different at .3 milligrams and the control, and then
15 it's maintained. And if you look at the area under
16 the curve, there is a difference between 0.3
17 milligrams and the other two, which is statistically
18 significant at the 0.54 level. There is one subject
19 here who suffered an abscess during the healing phase
20 who lost four millimeters of attachment. If you
21 remove that attachment that subject just becomes
22 0.033. However, we didn't remove that subject

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 because it's an intent to treat analysis. Next
2 slide.

3 Now let's look at radiographic linear
4 bone growth as the companion. Here the .03 milligram
5 per mil out-performs the 1 milligram per mil, and
6 both out-perform the TCP alone in terms of
7 radiographic bone analysis. And these P-values are
8 very, very strong, even taking into consideration
9 multiple variable comparisons using Yates Correction
10 or other techniques. These P-values are extremely
11 powerful. Next slide, please.

12 Now if we compare against current
13 therapies again with the caveats I mentioned before,
14 clearly GEM-21 is comparable to Emdogain - 2.5 to 1.1
15 and probably better than Allograft, certainly better
16 than surgery alone. Next slide, please.

17 Now let's look at this radiographic
18 percent bone fill. That's this derived ratio of
19 linear bone growth as related to the original pocket
20 depth. And here the mean percent bone fill is on the
21 X-axis, and the three groups are depicted by the
22 bars. Again, the green bar in the 0.3 milligram per

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 mil out-performs both the 1 milligram per mil, and
2 both out-perform the TCP alone, and the P-values
3 again are quite low, showing high levels of
4 statistical significance. Next slide.

5 And now comparing radiographic bone fill
6 with the predicate products, you can see that GEM-21
7 is comparable to and probably performs better than
8 the other products with respect to percent bone fill.

9 But certainly it's comparable too. Next slide,
10 please.

11 Now if we look now, drill down into the
12 data and look at the various types of lesions we're
13 treating, you know the one and two walls are very
14 difficult to treat as to the more contained three-
15 wall and circumferential. And you can see that in
16 the data. You look at all of these bars for the one
17 and two-wall are lower than all of these bars for the
18 three-wall and circumferential, so in general, these
19 lesions heal better than these; yet, the 0.3
20 milligram per mil GEM-21 gave 50 percent bone fill in
21 over half the subjects in those very difficult to
22 treat one and two-wall lesions. Again, out-

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 performing the 1 milligram per mil and out-performing
2 the TCP. And if you look at the three-wall and
3 circumferential defects, 65 percent of the lesions
4 were filled, or the lesions were filled 65 percent of
5 the bone with the 0.3 milligram, as compared to 34
6 and 21 for the other one control. Again, the .3
7 milligram out-performs the TCP and even these defects
8 that heal on their own.

9 DR. SHARMA: Is this all this
10 radiographic data at three months or six?

11 DR. GENCO: Well, it's at six months.
12 All the radiographic data is at six months. The CAL
13 data is at three and six. Now the reason the
14 radiographic data at three months is not used is
15 because the material is in the lesion at three
16 months. You can see it on the radiograph. And from
17 the histology, which we've done extensive histology
18 both in man and animals, it's usually gone by three
19 months histologically, you can't see it any more. So
20 we felt safe in looking at the six month x-rays.
21 Next slide, please.

22 Now this is a distribution of the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 cumulative proportion of bone fill and the curve to
2 the right, right is better, left is worse. You can
3 see the curve to the left is the control, and if you
4 look at the proportion, let's say 50 percent of the
5 subjects with control, 20 percent bone fill. In
6 other words, in 50 percent of the subjects given the
7 control, you got 20 percent bone fill, not very good.

8 And 33 percent and 50 percent of the subjects given
9 the 1 milligram per mil, they got 33 percent bone
10 fill, but in 50 percent of the subjects with .3
11 milligram, we got 50 percent fill, so half the
12 subjects, half the pockets were filled with .3
13 milligram, 20 percent of the subjects or half the
14 subjects, 20 percent fill with the control. Again,
15 statistically significantly different. This starts
16 to address the issue of clinical significance. What
17 proportion of the population actually benefits from
18 this? It's not meant to be the statistical proof for
19 efficacy, but some indication of how many people in
20 the population are actually benefitting from this
21 treatment which addresses clinical significance.

22 Now the rationale for the composite

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 outcomes is just as I said, to get some indication of
2 clinical significance. So we use the two primary
3 end-points, CAL and bone, and then made a composite.

4 And this is done in rheumatology and other areas of
5 wound healing. As a matter of fact, we're carrying
6 out a study of cardiovascular disease and we're using
7 a composite of six different cardiovascular
8 variables, so composite variables I think are gaining
9 in attention in the clinical trial methodology area,
10 and are extremely useful when done properly. Next
11 slide.

12 So how do we define a successful outcome?

13 Well, what we did is we took the PMAs for Emdogain
14 and PepGen and took their best results, and we said
15 all right, that's the attachment gain achieved by
16 either Emdogain or PepGen, and the best bone fill or
17 linear bone growth will accept as the cut-off point.

18 Okay. So we put that together. The best attachment
19 gain for Emdogain or PepGen was 2.7 millimeters, and
20 the best bone fill for either was 14.1, so we put
21 those together. That's our composite. If you reach
22 that, we define that as "success". Same for CAL and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 LPG, 2.7 millimeters and 1.1. Let's see what the
2 results are.

3 Now look at this overall. Success was
4 achieved with .3 milligram dose using one of those
5 composites in 70 percent of the subjects, and with
6 the other composite, it's 60 percent of the subjects.

7 Now did the percent success out-perform the control?

8 Yes, but that's not the intent here. The intent
9 here is not the statistical significance, although it
10 was highly significant, but what percent of the
11 population achieved this definition of success?

12 We're quite pleased that 60 to 70 percent of the
13 population benefitted from this product. Next slide,
14 please.

15 Now let's look at one of these lesions.
16 This is a lesion that after breaking the code, we
17 know got the .03 milligram per mil, a 44-year old
18 female with a lesion above the two-wall defect which
19 was circled around the lingual with a class two
20 furcation, and you could see the defect on the
21 distal. You could see it here. And if you look now
22 at the six month x-ray, this is new bone. This not

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 TCP, this is new bone. The TCP is little particles
2 that you could see on the x-ray, very
3 characteristically different than the new bone. So
4 when you make the measurements from the CEJ to the
5 base of the defect, in this instance it's something
6 like 6 millimeters, and in this instance it's
7 something like 3, so we had about a 50 percent fill,
8 just rounded off with a linear bone gain of about 3
9 millimeters, so this is a typical result of 50
10 percent fill and 3 millimeter bone gain.

11 Now as Dr. Nevins mentions, this markedly
12 changes the prognosis of that tooth. The lingual
13 furcation is filled and the distal lesion is pretty
14 much 100 percent healed, and the mesial lesion is
15 about 50 percent healed. This is a very good result
16 clinically. Chances are we've gone from a 6
17 millimeter, 7 millimeter pocket to a 3 millimeter, 4
18 millimeter site. That's maintainable. Next slide,
19 please.

20 So in summary then we had 180 patients
21 who were fully masked in this perspective multi-
22 center trial. Quality was assured by multiple

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 mechanisms, including CRO, blinded investigators,
2 blinded design, arm's length statistical analysis by
3 CRO and another firm. There were no device-related
4 serious adverse events. The GEM-21 has statistically
5 improved CAL at three months, and the CAL under the
6 curve between zero and six months, and the
7 interpretation is it's rapid healing which persisted
8 at six months that was induced by the .3 milligram
9 per mil growth factor.

10 The linear bone growth was improved at
11 six months as was the percent bone fill at six months
12 in a highly statistically significant manner. We
13 feel that the GEM-21 exceeded benchmarks of
14 effectiveness but for caution we'll say it was
15 comparable to the benchmarks, very unlikely to be
16 less effective than already existing products on the
17 market. Next slide, please.

18 Now let's look again at the comparability
19 to the approved products, the GEM-21S, CAL gain 3.7,
20 radiographic fill 2.5 are in the ballpark if not
21 better than the other approved products. Next slide.

22 So overall then, we feel that GEM-21S, a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 fully synthetic and safe product. Take it off the
2 shelf, don't have to worry about contamination with
3 Bovine contaminants, or with Allograft problems,
4 although that's not a major problem, but it is an
5 issue in patient's minds and some clinicians, so this
6 is fully synthetic and safe as a known mechanism of
7 action demonstrated by over a decade and a half of
8 very intense high quality research a mechanism of
9 action of PDGF. And on the characteristics of the
10 Vitoss, and in fact the Vitoss is a new product
11 developed in 1999.

12 The recommitant PDGF-BB component
13 enhances periodontic regeneration in animals and
14 humans, and this is very reproducible result seen in
15 many species. I personally have been involved in
16 three dog studies, and they all show the same thing;
17 complete fill of Class 3 furcations. Accelerates
18 attachment level gain and radiographic evidence of
19 bone regeneration, quite well documented in human
20 study I've just mentioned, and demonstrates a
21 favorable risk to benefit relations, so I would say
22 in general, to sum up, my view is that this product

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 gives a favorable clinical result in 60 to 70 percent
2 of patients when used as indicated with very few side
3 effects that we're not expecting. Thank you very
4 much for your attention.

5 DR. LYNCH: Thank you very much, Dr.
6 Genco, Dr. Nevins and Dr. Giannoble for those
7 presentations. We had asked the FDA, Mr. Adjodja,
8 Dr. Runner for just a few minutes at the conclusion
9 of the presentation just to be available to answer
10 any burning questions on the methodology. We don't
11 want to pre-empt the discussion this afternoon. We
12 understand there may be some more global questions.
13 I think those might be more relevant for this
14 afternoon, and I think that's where that discussion
15 is planned, but we didn't want to leave any -- if
16 there were any lingering burning questions on
17 methodology or specific aspects of the presentation
18 to let those sort of fester in your mind. So we'd be
19 happy to again entertain any specific burning
20 questions relative to methodology that you might have
21 now, or proceed forward and we can address further
22 questions this afternoon.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 CHAIRMAN SUZUKI: Before we begin with
2 the questions, I'd like to remind the audience that
3 I'd like you to reserve any questions regarding the
4 particular hearings until after the presentations.
5 And then secondly, the FDA panel members will have
6 the prerogative of asking questions first, including
7 procedural questions. And I wanted to thank the
8 presenters and the sponsor for presenting such a
9 precise presentation and keeping us on time.

10 As Chair, I'd like to take the
11 prerogative of asking perhaps the first two
12 procedural questions, and that is with respect to the
13 radiographic benefits - and I know Dr. Genco
14 mentioned looking at the composite outcomes in total,
15 but focusing in on the radiographic interpretations,
16 I notice that there is a mean improvement of about
17 2.1 millimeters. I'd like an explanation as to why
18 you think that this is clinically significant.

19 DR. LYNCH: Okay, certainly. And I'm
20 going to moderate the session and probably defer that
21 to many of our panel members since your question
22 specifically refers to the clinical relevance of the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 radiographic bone gain. I think it's appropriate
2 that the clinicians and the panel answer that, and I
3 might refer to Dr. Genco and to Dr. Nevins.

4 DR. GENCO: The radiograph especially at
5 six months under-estimates the healing. You saw the
6 attachment gain was more like 3.8 with a bone gain
7 for the radiograph of 2.1. I think that reflects the
8 under-estimate. 3.8 millimeters in a 7 millimeter
9 pocket, it's a 4 millimeters gain, 7 millimeters to
10 begin with, you're down to 3 millimeters. And I
11 think we and others have done studies showing that if
12 a pocket is 5 millimeters or greater, it has a 6 to 7
13 fold greater chance of losing attachment in the next
14 two years, so if you can get it below 5 or 6, that
15 bodes well for the future. And I think this is what
16 this study has shown, that the pockets are reduced
17 from 7 to approximately 3, 3 and a half. The bone
18 doesn't quite reflect that because it under-estimates
19 the healing, but the pocket reduction and the
20 attachment gain I think really are telling from the
21 clinical standpoint.

22 CHAIRMAN SUZUKI: Okay. Thank you. My

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 second question is a procedural one regarding the
2 surgery itself. In the video clip that we saw, there
3 were just a couple of procedure questions that I
4 had. The first is that in your presentation of
5 materials, you indicated that you use Tetracycline of
6 the preparation of the root surface, but the video
7 clip did not show that. Is there a reason why it was
8 omitted or was that standardized?

9 DR. LYNCH: Yes, and Dr. Giannoble, who
10 was an investigator, could comment, but I could
11 certainly comment on that one, as well. It was
12 omitted from the video clip simply to make that a
13 very concise video clip for no other reasons. It was
14 standardized as to the amount of Tetracycline that
15 was used and the duration of the root conditioning so
16 that was all pre-specified and the examiners or the
17 surgeons were trained on that.

18 CHAIRMAN SUZUKI: Okay. The last
19 question I had about the procedure is that frequently
20 surgeons fenestrate the osseous lesion, and I noticed
21 in the video clip that you did not. Was there a
22 reason for selectively not using that particular

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 step?

2 DR. LYNCH: Again, Dr. Giannoble, if you
3 want to come up and comment, you could certainly feel
4 free to. We did, I believe, allow fenestration at
5 the discretion of the investigators, a very hard
6 cerotic bone, the bony walls there. It wasn't
7 necessary in the particular case that you saw, but if
8 in the judgment of the investigator that the bone is
9 very cerotic and sort of avascular, it was permitted
10 for them to do perforations.

11 CHAIRMAN SUZUKI: So that portion of the
12 surgical procedure was not standardized in terms of
13 the fenestration.

14 DR. LYNCH: Yes, we felt that clinically
15 certainly not all cases would require it, again as
16 the case that you saw, but we felt like certainly
17 some would. So again, in the investigator meeting
18 prior to study initiation, we discussed certainly
19 this very point, and the investigators that we chose,
20 of course, are all very, very highly regarded, very
21 seasoned clinicians, and they felt like you couldn't
22 predetermine that all the lesions should have

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 perforation because it was necessary in all the
2 lesions, so we discussed the parameters or
3 characteristics of the lesion that would require
4 perforation and left that then to the surgeon.

5 CHAIRMAN SUZUKI: Okay. AT this time I'd
6 like to ask the panel if they have any points of
7 clarification that they would like answered. And
8 before you do so, I'd like to ask that you identify
9 yourself in the microphone for the transcriptionist,
10 as well as presenters identifying themselves into the
11 microphone, as well, before you respond. Okay. Dr.
12 Cochran.

13 DR. COCHRAN: David Cochran. Dr. Lynch,
14 I'd like to ask a couple of questions. First of all,
15 in the documentation you provided for us, you used a
16 couple of papers to reference for linear bone growth
17 and percent bone fill. Dr. Genco talked about using
18 the PMAs that were submitted prior for these other
19 products. How did you go about choosing those,
20 particularly there was a study from Greece, and then
21 there was Rutger Persant was another one that you
22 used.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 DR. LYNCH: I'll turn this also over to
2 the rest of the panel, but wherever there was data
3 available from the PMA submissions, specifically the
4 summary of safety and effectiveness for previously
5 approved products, we utilized those, so to translate
6 that means Emdogain and PepGen P-15. As has been
7 mentioned, Allograft has never been formally
8 "approved" or cleared for dental uses at any rate by
9 the agency and so, of course, there are no formal FDA
10 submissions that were available, so we did go back
11 through the literature and did a very extensive
12 literature search on specifically looking again at
13 radiographic assessment of bone fill following
14 Allograft treatment. And we used what data was
15 available in the literature. As Dr. Genco mentioned,
16 we're certainly not by any means claiming superiority
17 to any of those materials that were used. We were
18 just trying to, and at the agency's request, get some
19 comparison of the effectiveness of the results seen
20 in this trial benchmarked against other materials
21 that the clinicians are using.

22 We also utilized or carefully reviewed

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 the paper out of the San Antonio group that clearly
2 shows that the radiograph assessments do under-
3 estimate bone fill as compared to re-entry
4 assessments.

5 DR. COCHRAN: That was a good choice.
6 The second question is what are your thoughts on
7 comparing some of your results to the enamel matrix
8 proteins being that that's a protein-only therapy,
9 and you're using protein plus graft material. Would
10 you comment on that?

11 DR. LYNCH: Again just to stress that the
12 comparisons that we did were simply to compare the
13 results to other benchmarks of effectiveness that
14 were available. The TCP that's used as a carrier, as
15 has been mentioned, is fully resorbed within about
16 three to four months and, therefore, we did not think
17 that that would affect, for example, the radiographic
18 assessment at six months, as was mentioned. That's
19 one reason we did not do the radiographic assessment
20 at three months.

21 Given that the PDGF is clearly gone by
22 six months, given that the matrix is, as far as we

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 can determine, certainly mostly, if not totally,
2 resorbed at six months, as well; that would make that
3 site then somewhat comparable, and I don't want to
4 overstate this to the Emdogain where there was no
5 matrix observed.

6 We know in clinical practice that many
7 people do mix Emdogain, as again has been
8 investigated in your university, with other bone
9 substitute materials to try to contain it in the bone
10 defect. And I think that that was one of the
11 rationales for using a matrix in our product, was to
12 provide the clinician a standard matrix that they
13 could use to mix with the recombinant protein, as
14 compared to Emdogain where the clinicians are often
15 just kind of taking whatever they have on the shelf,
16 so to speak, or whatever grafting material they like
17 and mixing it with the Emdogain, so we feel this
18 provides a more standardized product.

19 DR. COCHRAN: The last question would be
20 in the documentation there's a product mentioned
21 called Vitoss Plus, which is a similar product or an
22 approved product. What is that?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 DR. LYNCH: It's just a different name
2 that -- sorry for the confusion there. Actually, the
3 names of this product that we're reviewing today have
4 sort of transitioned from Beta TCP Plus, at some
5 point I think in the documentation it was called
6 Vitoss Plus, now it's called GEM-21S. It's the same
7 product.

8 CHAIRMAN SUZUKI: Dr. Amar.

9 DR. AMAR: Salomo Amar. I'm going to ask
10 a more general question. Dr. Genco mentioned that at
11 a certain point the control catches up with the rest
12 of the experimental. And my general question, if at
13 six months the control or the experimental -- the
14 control catches up with the experimental, what would
15 be the added benefit of using this molecule as
16 compared to TCP? Is it just for the early reading
17 improvement parameters or are we talking about long-
18 term maintenance?

19 DR. LYNCH: Well, Dr. Genco, why don't
20 you come up here and I'll provide my interpretation,
21 but I would like -- I think Dr. Amar would like to
22 hear your's, as well. I think what Dr. Genco was

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 referring there was for the clinical attachment level
2 of the soft tissue. Certainly, the bone fill as
3 measured radiographically never catches up in the
4 control versus the treatment group, so there's always
5 a strongly significant improvement or benefit in the
6 radiographic bone fill.

7 DR. GENCO: That's one point. We didn't
8 see that phenomena in the bone. Of course, we only
9 looked at 1.2. With respect to early healing it's,
10 of course, benefit to get that healing pretty much
11 underway in the first three months, and you can get
12 on with the rest of the therapy. We think that if
13 you're involved in a complex case that requires
14 implants and other treatment, that to have this early
15 result at three months is of great benefit, so it's
16 an accelerated treatment that fits in with the
17 treatment of advanced case, and it's a definite
18 benefit.

19 DR. AMAR: If the clinical attachment
20 level is the primary outcome and it catches up, aside
21 from the bone failing, it looks pretty similar.

22 DR. GENCO: Well the point is if you --

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 let's say if you have to put a crown on the tooth,
2 you could start putting the crown on in maybe two to
3 three month rather than waiting six.

4 DR. AMAR: The other question that I had
5 is regarding the resorption. I saw sections by Dr.
6 Nevins, at nine months I believe on the furcations
7 showing probably some deposit of Beta Tricalcium
8 Phosphate.

9 DR. GENCO: Myron, do you want to --

10 DR. NEVINS: If you're referring to the -
11 -

12 CHAIRMAN SUZUKI: Can you identify
13 yourself, please.

14 DR. NEVINS: I'm sorry. Myron Nevins.
15 If you're referring to the human histology, we didn't
16 use Beta Tricalcium Phosphate. That was an Allograft
17 study. The Allograft was the carrier for that. The
18 only thing I showed with Beta Tricalcium Phosphate
19 was the one slide at the end on a K-9 study, and that
20 was at eight weeks. And on the GEM-21 there was no
21 evidence of Tricalcium Phosphate at all in the
22 control which was the Tricalcium Phosphate by itself.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 There were pieces of Tricalcium Phosphate.

2 DR. COCHRAN: To follow-up on that, I
3 think in that result didn't you get 70 percent
4 regeneration with the TCP alone in that K-9 study?

5 DR. NEVINS: I would have to --

6 DR. COCHRAN: It's 37 percent, not 70.

7 DR. NEVINS: It's the blue column over
8 TCP.

9 CHAIRMAN SUZUKI: Okay. Any other
10 questions, Dr. Cochran? Okay. Any other questions
11 from the panel? Dr. Sharma.

12 DR. SHARMA: According to the protocol,
13 this is Inder Sharma. According to the protocol, the
14 primary comparison was to be between the high dose
15 and the control. Only if it was significant, then
16 0.3 which is low dose, was to compare with control.
17 But I see the focus of presentation have been mostly
18 on the low dose, so I'm wondering what happened that
19 we are now focusing on what we said in the protocol
20 that this will be about primary comparison, because
21 primary comparison is not significant whether you
22 look at three months or six months.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 DR. LYNCH: I think you're looking at an
2 older version of the protocol. There were formal
3 amendments that were submitted to the agency and
4 approved by the agency throughout the study for
5 various things that we had under discussion with
6 them, so the primary comparator was the .3 mg/ml dose
7 level versus the TCP control.

8 DR. SHARMA: The second question I have
9 is about the composite end-point. Was this a pre-
10 planned comparison or was it decided to do that after
11 the fact?

12 DR. LYNCH: I'm sorry. I'm not sure I
13 understand your question.

14 DR. SHARMA: The composite end-point.

15 DR. LYNCH: Oh, the composite. Okay.

16 DR. SHARMA: Was it a pre-planned
17 comparison using composite end-point, or was it later
18 decided to be --

19 DR. LYNCH: I think Dr. Phil Laven will
20 address that. He's our biostatistician.

21 DR. LAVEN: Hi, Philip Laven,
22 biostatistical consultant to Biomonetics. My company

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 is Avarion. When we came up with the idea for the
2 composite end-point, it was done before the database
3 lock, and it was planned at the same time that we
4 planned looking at the AUC. That end-point was
5 reflective of the fact that we knew that the disease
6 was more extensive than just looking at the delta CAL
7 measurement, and that the composite end-point had to
8 address both the radiographic end-points, as well as
9 the clinical end-points, so this was all done
10 prospectively before the database lock, but was not
11 in the original protocol.

12 CHAIRMAN SUZUKI: Dr. Zero.

13 DR. ZERO: Domenick Zero. I have a
14 question about how the statistical analysis was done,
15 although that's not my main expertise. On looking at
16 the distribution of females, smokers, African
17 Americans, and the CAL values, there are although not
18 statistically significant differences, there are some
19 numerical differences that are noticeable just in
20 looking at the different groups. The data reported,
21 was that an adjusted statistics, or was that just the
22 raw statistic?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 DR. LYNCH: Dr. Laven.

2 DR. LAVEN: Yes, Philip Laven. Those
3 gain that you are seeing there in the report are
4 unadjusted statistics. At the request of the FDA
5 over the last month, we did prepare additional
6 analyses where we did look at controlling for those
7 factors and looking for treatment interactions with
8 factors like smoking, the location of the tooth,
9 whether it was a molar or not, and we did assess
10 those analyses. And those analyses, just to give a
11 sense for where they turned out, there was no
12 treatment interaction with any baseline co-variates,
13 so the treatment advantages that you're materially
14 seeing there, although they're uncorrected, do
15 represent the state-of-the-art for what happened in
16 those groups.

17 DR. ZERO: Thank you.

18 CHAIRMAN SUZUKI: Dr. Zuniga.

19 DR. ZUNIGA: John Zuniga. A couple of
20 very simple, hopefully, questions on the study
21 protocol, just more for clarification. I notice in
22 your management of your post-operative patients, you

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 included the use of NSAIDs for any analgesia. Was
2 that a -- why did you do that, and do you have
3 concern about NSAIDs in this product?

4 DR. LYNCH: No. Sam Lynch. There are no
5 specific concerns about NSAIDs related to this
6 product. We just knew that NSAIDs as a general class
7 of drugs had ability to affect wound healing, and we
8 didn't want some patients to be on NSAIDs by some
9 investigators, and other patients not to be, because
10 we thought that that might affect the -- especially
11 the immediate post-op healing, and we did have an
12 end-point that looked at wound healing over that
13 first three week period, so in order to standardize
14 that regimen we just elected not to use NSAIDs.

15 DR. ZUNIGA: And has that been explored
16 using NSAIDs?

17 DR. LYNCH: It, again, was not in the
18 pivotal clinical trial. We would have to look at the
19 patients in the human histologic study to see if they
20 were given NSAIDs or not.

21 DR. ZUNIGA: And then the second
22 question, again relative to the study protocol, is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 the antibiotic use and some of the post-operative
2 instructions for the patients. Were there any
3 patients that did not comply with the antibiotics?
4 And if so, were there any difference in effects? And
5 then finally, you have pretty strict protocol for
6 soft foods and diet eating on the other side. How
7 compliant was that regarding effects on the patient
8 outcomes?

9 DR. LYNCH: Right. I think as is
10 relatively customary, post-op instructions were given
11 to these patients. I don't believe that there were
12 anything unusual about our post-op instructions
13 compared to many that we give our periodontal
14 patients. In terms of any specific, like protocol
15 violations that were reported where the patient
16 reported chewing on the site of the surgery or that
17 kind of thing, and I don't know - Mark, do you care
18 to comment on that? I don't think there was any --
19 there wasn't certainly any significant violations
20 along that line. There may have been isolated cases.

21
22 MR. CITRON: No protocol violations.

1 DR. LYNCH: Okay. So there were no
2 protocol violations.

3 CHAIRMAN SUZUKI: I have a couple of
4 final questions. With respect to the dose, .3 versus
5 1.0 milligrams per mil, were dose response curves
6 completed prior to your selection of these doses, or
7 were these doses taken from the literature? And
8 secondly, why isn't more better?

9 DR. LYNCH: I think the second half of
10 your question, Dr. Suzuki, we might want to table to
11 this afternoon. It's certainly a very excellent
12 question, and I don't mean to put it off. We could
13 address it here, but for broader questions in terms
14 of that, we might defer those to this afternoon at
15 your discretion.

16 In terms of how we selected the .3 and
17 the 1 mg/ml, as Dr. Giannoble showed one slide from
18 the study by how co-workers at Harvard and UNC
19 several years ago, that study utilized .05 mg/ml, and
20 a .15 mg/ml of PDGF. And then that study used also
21 accommodation with the insulin like growth factor,
22 and it showed that there was no effect, no beneficial

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 effect at the .05 mg/ml. There was a significant
2 beneficial effect on bone similar to what we frankly
3 saw in this study at the .15 mg/ml, so we used that
4 as information. We also then conducted a canine
5 study that looked at the .3 mg/ml, and you say well
6 how did you get from .15 mg/ml in that initial
7 clinical study a few years ago to .3. And the
8 rationale, right or wrong there was that because that
9 initial study utilized a combination of two growth
10 factors, we felt like we might need to utilize the
11 total growth factor dose, if you will, and so that
12 would be .3 mg/ml. So that was the justification for
13 the low dose in our pivotal trial. And the
14 justification for the high dose was taken just as a
15 XXX multiple from that. And as was reported, I
16 believe by Dr. Nevins, we did see absolutely
17 consistent results in the canine study that the .3
18 mg/ml provided the most beneficial response, so that
19 was again the reason for determining that that was
20 our primary comparator, was based upon the canine
21 study.

22 CHAIRMAN SUZUKI: Thank you. John Suzuki

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 again. My last question is with respect to the T-
2 inclusion criteria, and the age range of your
3 periodontitis patients were from age 25 to 75, I
4 believe. In the submitted materials, you indicated
5 that the aggress of periodontitis patients were
6 excluded from this patient group; yet in the oral
7 presentation that was not in the particular slide.
8 Is there a reason for that?

9 DR. LYNCH: Sam Lynch. That was just
10 omitted off the slide just for sake of brevity, and
11 we couldn't include all of the criteria on the slide.

12 But certainly, patients that were considered to have
13 aggressive periodontitis, what we used to call
14 juvenile periodontitis, were excluded from the study.

15 CHAIRMAN SUZUKI: Dr. Amar.

16 DR. AMAR: I have just one more question.
17 You're going to do x-ray analysis and Dr. Genco
18 mentioned that there was no stent, am I correct?
19 There as no stent. And what was the percentage of
20 elongation accepted, and you mentioned 15 percent.
21 Am I correct?

22 DR. GENCO: I can start the answer. Bob

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Genco. Maybe Dr. Reddy can continue. But what I
2 presented was that the elongation -- or for
3 shorthand, 15 percent, it was seen in less than 5
4 percent of the x-rays. Is that the question?

5 DR. AMAR: I guess my question is if it's
6 really 15 percent, the cut-off value on a tooth
7 that's -- a root that is say 10 millimeters, 15
8 percent is 1.5 millimeters, that the effect size -
9 it's about the effect size that we would see on bone
10 fill. And I would have some kind of concerns about
11 that.

12 DR. GENCO: Well, that was used to then
13 adjust the x-rays to normalize.

14 DR. AMAR: So there was no more than 5
15 percent elongation.

16 DR. GENCO: Well, let's let -- well, 5
17 percent of the cases exceeded the criteria of 15
18 percent elongation or foreshortening. Therefore,
19 required normalization, so the extent to which there
20 were over 15 percent - I think that Mike could answer
21 that. Dr. Reddy.

22 DR. REDDY: I think I understand what

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 Salomo is asking. Hi, Michael Reddy. I'm an
2 investigator from the University of Alabama. I did
3 radiograph analysis. You want to know what percent -
4 - what was the highest range of elongation and
5 foreshortening. And I have to look at the database
6 to tell you exactly, but some of them were up to
7 about 25 percent, but there were very few x-rays.
8 Remember, this was an intent to treat analysis, so we
9 simply couldn't say that that didn't make our
10 radiographic criteria, so we included them, and then
11 retrospectively corrected it mathematically. You
12 have to remember that a 15 percent increase in the
13 overall root length, which may be 15 millimeters in
14 length, may just vary the measurement of the bone
15 growth by about 10 percent even if you didn't correct
16 for it, even though we did correct for it. So if you
17 have 2 millimeters of bone growth, you're really only
18 going to change it to about 2.2. But the case that
19 did happen, we did mathematically apply a formula and
20 run an algorithm to correct those. Those are very
21 few. We had a great fear, the same fear you had at
22 the start of the study. That's why we incorporated

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 that into the analysis, this 15 percent cut-off
2 because we were afraid that these are field x-rays,
3 and we may have 50 percent of them with elongation
4 and foreshortening, and it turned out that they were
5 actually very clinical radiographs. Of course,
6 again, we had to get a good x-ray of one tooth --

7 DR. AMAR: So if I understand you
8 correctly, there was an area of elongation about 15
9 percent of max that was corrected.

10 DR. REDDY: No, only if the area was over
11 15 percent was it corrected. There were some sites
12 that were 25 percent foreshortened or elongated.

13 DR. AMAR: That could translate into if
14 the root is 10 millimeters into 1.5 millimeters
15 change?

16 DR. REDDY: It could if it wasn't
17 corrected for, and that's the reason why we put the
18 correction in, exactly.

19 DR. AMAR: Can you elaborate on the
20 correction?

21 DR. REDDY: Yes, the correction simply to
22 go ahead and correct everything back to the baseline.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 We consider whatever length we measured at the
2 baseline from CEJ to apex as the gold standard. If
3 it differed by more than 15 percent, we used a ratio
4 of the original CEJ to apex measurement to new CEJ to
5 apex measurement to mathematically correct all
6 measurements that were subsequently taken, so we
7 wouldn't lose the data.

8 DR. AMAR: Thank you.

9 DR. SHARMA: I have one.

10 CHAIRMAN SUZUKI: Okay. Dr. Sharma.

11 DR. SHARMA: Inder Sharma. I have one
12 final question. Interim analysis were planned for
13 the study and I was wondering where they conducted?
14 And if they were conducted, was there a DSMB or who
15 had access to those results?

16 DR. LYNCH: Would you mind repeating the
17 question, please?

18 DR. SHARMA: The interim analysis for the
19 study --

20 DR. LYNCH: Interim analysis --

21 DR. SHARMA: They were planned, and my
22 question is if those interim analysis were conducted,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 and who had access to the results, was it a DSME
2 independent diversity of monitoring the worker?

3 DR. LYNCH: Yes, I understand. Sam
4 Lynch. There was an interim analysis conducted per
5 the protocol. The analysis was conducted on the
6 first 90 patients to complete three-month follow-up.

7 This was, again, an analysis that was agreed upon
8 with the agency. It was done in a fully blinded
9 fashion by the independent clinical research
10 organization, the CRO that was responsible for
11 monitoring the study, so there certainly was no
12 breaking of the blind or anything.

13 The only data that we got back was that
14 and the reason the FDA had asked us to do that
15 interim analysis was to one of sample size, should we
16 adjust the sample size at that point? Do we increase
17 the number of patients, because we had agreed not to
18 decrease the number of patients because we did not
19 want to take a statistical penalty for the interim
20 analysis, so we had all along said we're going to do
21 180 patients, even if that result was just fantastic.

22 But the question was would we need to add patients

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701