

Media Recipes for *Humulus*

Humulus growth medium I (HGI) (solid) - 1000 ml

- ✓ To a small volume of double distilled water (ddH₂O) add:

potassium nitrate (KNO ₃)	1.9 g
ammonium nitrate (NH ₄ KNO ₃)	1.65 g
magnesium sulfate (MgSO ₄)	0.1807g
calcium chloride, <i>dihydrate</i> (CaCl ₂)	0.332 g
potassium phosphate, <i>monobasic</i> (KH ₂ PO ₄)	0.17 g
Sequestrene 138 ^{TM1} , iron chelate (EDDHA)	0.1 g
Iron stock ²	20.0 ml
MS micronutrients ^{2,3}	10.0 ml
MS vitamins ^{2,3}	10.0 ml
glucose	20.0 g
BA (6-benzylaminopurine)	0.1 mg

- ✓ Stir until well blended
- ✓ Bring to final volume (1000 ml) with ddH₂O
- ✓ Adjust pH to 5.7
- ✓ Add:

gellan gum (Phytigel ^{TM4})	1.75 g
agar (Sigma ^{®4} A7002)	3.5 g

- ✓ Heat and stir until boiling
- ✓ Dispense into stacked Magenta^{®5} GA7 culture vessels (40 ml/vessel)
- ✓ Autoclave

¹ Becker Underwood Inc., Ames, IA

² Recipe follows

³ Murashige & Skoog, 1962

⁴ Sigma-Aldrich, St. Louis, MO

⁵ Magenta Corp., Chicago, IL.

Iron stock solution (100x) (liquid) – 500 ml

- ✓ To a small volume of double distilled water (ddH₂O) add:
 NaEDTA, disodium salt, *dihydrate* (NA₂EDTA*2H₂O) 1.86 g
- ✓ Stir until NaEDTA is completely dissolved
- ✓ In a separate vessel containing a small volume of ddH₂O add:
 ferric sulfate (FeSO₄*7H₂O) 1.39 g
- ✓ Heat and stir until the ferric sulfate is completely dissolved. Allow solution to cool completely
- ✓ Combine the NaEDTA solution with the ferric sulfate solution
- ✓ Bring to volume (500 ml) and stir until the solution turns yellow
- ✓ Dispense into an amber vessel to prevent photodegradation. Store at 2-4 °C

MS³ micronutrient stock solution (100x) (liquid) – 500 ml

- ✓ To a small volume of double distilled water (ddH₂O) add:

boric acid (H ₃ BO ₃)	0.31 g
cobalt chloride (CoCl ₂ * 6H ₂ O)	0.00125 g
cupric sulfate (CuSO ₄ *5H ₂ O)	0.00125 g
zinc sulfate (ZnSO ₄ * 7H ₂ O)	0.43 g
molybdic acid, sodium salt, <i>dihydrate</i> (NaMoO ₄ * 2H ₂ O)	0.0125 g
manganese sulfate (MnSO ₄ *H ₂ O)	0.845 g
potassium iodide (KI)	0.0415 g
- ✓ Heat and stir until boiling and dry ingredients have completely dissolved
- ✓ Bring to final volume (500 ml) with ddH₂O
- ✓ Dispense into desired vessel and store at 2-4 °C or aliquot and store in freezer

MS³ vitamin stock solution (100x) (liquid) – 500 ml

- ✓ To a small volume of double distilled water (ddH₂O) add:

glycine (free base)	0.1 g
myo-inositol	5.0 g
nicotinic acid (free base)	0.025 g
pyridoxine HCl	0.025 g
thiamine HCl	0.005 g

- ✓ Stir until ingredients are well blended
- ✓ Bring to final volume (500 ml) with ddH₂O
- ✓ Dispense into desired vessel and store at 2-4 °C or aliquot and store in freezer

Humulus growth medium II (HGII) (solid) - 1000 ml

- ✓ To a small volume of double distilled water (ddH₂O) add:

potassium nitrate (KNO ₃)	1.9 g
ammonium nitrate (NH ₄ KNO ₃)	1.65 g
magnesium sulfate (MgSO ₄)	0.1807 g
calcium chloride, <i>dihydrate</i> (CaCl ₂)	0.332 g
potassium phosphate, <i>monobasic</i> (KH ₂ PO ₄)	0.17 g
iron stock ²	20.0 ml
MS micronutrients ^{2,3}	10.0 ml
MS vitamins ^{2,3}	10.0 ml
glucose	20.0 g
BA (6-benzylaminopurine)	0.1 mg

- ✓ Stir until well blended
- ✓ Bring to final volume (1000 ml) with ddH₂O
- ✓ Adjust pH to 5.7
- ✓ Add:

gellan gum (Phytigel™ ⁴)	1.75 g
agar (Sigma® ⁴ A7002)	3.5 g

- ✓ Heat and stir until boiling
- ✓ Dispense into stacked Magenta®⁵ GA7 culture vessels (40 ml/vessel)
- ✓ Autoclave

Ca-free MS+3% (w/v) Na-alginate medium (liquid) – 100 ml

- ✓ To a small volume of double distilled water (ddH₂O) add:

potassium nitrate (KNO ₃)	0.19 g
ammonium nitrate (NH ₄ KNO ₃)	0.165 g
magnesium sulfate (MgSO ₄)	0.018 g
potassium phosphate, <i>monobasic</i> (KH ₂ PO ₄) ¹	0.017 g
iron stock ¹	1.0 ml
MS micronutrients ^{1,3}	1.0 ml
MS vitamins ^{1,3}	1.0 ml
sucrose	17.1 g
BA (6-benzylaminopurine)	0.1 mg
GA ₃ (gibberellic acid)	0.01 mg

- ✓ Stir until dry ingredients are completely dissolved
- ✓ Bring to final volume (100 ml)
- ✓ Adjust pH to 5.0
- ✓ Add:

alginate sodium salt, 2% viscosity (Sigma⁴ A2158) 3.0 g

To prevent clumping, add the alginate slowly to rapidly stirring medium. A homogenizer with a propeller-type stirring blade works well for this step.

- ✓ Continue stirring until well blended and alginate is completely dissolved (~20 minutes)
- ✓ Dispense into desired vessels
- ✓ Autoclave

100 mM calcium chloride+MS encapsulation medium (liquid) – 1000 ml

- ✓ To a small volume of double distilled water (ddH₂O) add:

MS basal medium w/vitamins ³	4.43 g	(pre-packaged as M519 ⁶)
calcium chloride (CaCl ₂), <i>dihydrate</i>	14.7 g	
sucrose	30.0 g	

- ✓ Stir until dry ingredients are completely dissolved
- ✓ Bring to final volume (1000 ml) with ddH₂O
- ✓ Adjust pH to 5.0
- ✓ Heat and stir until well blended
- ✓ Dispense into desired vessels
- ✓ Autoclave

⁶ *Phytotechnology Laboratories, Shawnee Mission, KS*

0.75 M sucrose+MS medium (liquid) – 1000ml

- ✓ To a small volume of double distilled water (ddH₂O) add:

MS basal medium w/vitamins ³	4.43 g	(prepackaged as M519 ⁵)
sucrose	256.72 g	

- ✓ Stir until dry ingredients are completely dissolved
- ✓ Bring to final volume (1000 ml) with ddH₂O
- ✓ Adjust pH to 5.0
- ✓ Mix and heat until boiling
- ✓ Dispense into desired vessels
- ✓ Autoclave

Ca-free MS medium (liquid) - 1000ml

- ✓ To a small volume of double distilled water (ddH₂O) add:

potassium nitrate (KNO ₃)	1.9 g
ammonium nitrate (NH ₄ KNO ₃)	1.65 g
magnesium sulfate (MgSO ₄)	0.1807 g
calcium chloride, <i>dihydrate</i> (CaCl ₂)	0.332 g
potassium phosphate, <i>monobasic</i> (KH ₂ PO ₄)	0.17 g
iron stock ²	20.0 ml
MS micronutrients ^{2,3}	10.0 ml
MS vitamins ^{2,3}	10.0 ml

- ✓ Adjust pH to 5.0
- ✓ Stir until well blended
- ✓ Bring to final volume (1000 ml) with ddH₂O
- ✓ Dispense into desired vessels
- ✓ Autoclave

Humulus recovery medium (solid) – 1000 ml

- ✓ To a small volume of double distilled water (ddH₂O) add:

potassium nitrate (KNO ₃)	1.9 g
ammonium nitrate (NH ₄ KNO ₃)	1.65 g
magnesium sulfate (MgSO ₄)	0.1807 g
calcium chloride, <i>dihydrate</i> (CaCl ₂)	0.332 g
potassium phosphate, <i>monobasic</i> (KH ₂ PO ₄)	0.17 g
iron stock ¹	20.0 ml
MS micronutrients ^{1,3}	10.0 ml
MS vitamins ^{1,3}	10.0 ml
glucose	20.0 g
GA ₃ (gibberellic acid)	0.05 mg
BA (6-benzylaminopurine)	0.5 mg

- ✓ Stir until dry ingredients are completely dissolved
- ✓ Bring to final volume (1000 ml) with ddH₂O
- ✓ Adjust pH to 5.0
- ✓ Add:

gellan gum (Phytigel™ ⁴)	1.75 g
agar (Sigma® ⁴ A7002)	3.5 g

- ✓ Heat and stir until well blended
- ✓ Autoclave
- ✓ In laminar flow hood, dispense slightly cooled liquid into Petri dishes