

## Attachment A. I. Column Chromatography, Question 2

Technique	TLC	Column		GC	HPLC
<b>Separation Technique</b>	Adsorption, partition or a combination of both effects	Adsorption and Partition		Partition between gas phase and stationary phase	Partition, adsorption or ion-exchange
<b>Phases</b>	<i>Stationary</i> Adsorbent – thin uniform layer or dry, finely powdered material such as silica gel or cellulose applied to a glass, plastic, or metal sheet or plate. <i>Mobile</i> Suitable solvent system	<i>Column Adsorption (CAC)</i> <i>Stationary</i> Adsorbent – activated alumina or silica gel as a dry solid or as a slurry <i>Mobile</i> Suitable solvent	<i>Column Partition (CPC)</i> <i>Stationary</i> Solvent adsorbed on a solid support <i>Mobile</i> Suitable solvent	<i>Stationary</i> Solid or immobilized liquid stationary phase. Liquid phases are found in packed or capillary columns. <i>Mobile</i> Gaseous mobile phase	<i>Stationary</i> Solid or immobilized liquid stationary phase <i>Mobile</i> Liquid mobile phase
<b>Equipment Needed</b>	Glass plates, storage rack, adsorbent, spreader, developing chamber, template, micropipette, sprayer and ultraviolet light source.	Chromatographic tube, delivery tube to control the flow rates of solvent and a tamping rod		Carrier gas source, injection port or auto-injectors, column, heated oven compartment, detector, and data handling system	Mobile Solvent, Solvent Delivery System (Pump), Injector (Autoinjector), column, detector, and data handling system
<b>Ease of Use</b>	User Friendly	User Friendly		Moderate to complex depending on Instrumentation	Moderate to complex depending on detectors and data handling system
<b>Accuracy</b>	Used for semi-quantitative or quantitative estimation	Used for quantitative analysis with the aid of titrimetric or spectrophotometric determinative step		Reliable quantitative results are obtainable especially with internal standards and using auto-injectors or auto-samplers	Reliable quantitative results are obtainable especially with internal standards and auto-injectors or auto-samplers.
<b>Sampling</b>	Apply the Test and	CAC-Compounds are		Compound of interest is to	Compounds are dissolved in

<b>Techniques</b>	Standard Solution as directed in the individual monograph and allow drying	dissolved in a small amount of solvent and added to the top of the column. CPC-A solution of the sample in a small volume of the mobile phase is added to the top of the column or a solution of the sample in a small volume of the immobile phase is mixed with the solid support and transferred to the column	be volatile and thermally stable when heated. The test mixture either in a solution or as a gas may be injected directly into the column	a suitable solvent. This technique allows for thermally unstable and non-volatile compounds to be chromatographed.
<b>Automation</b>	Not normally Multi spotting equipment is available.	None	Auto-injectors, auto-samplers	Auto-injectors, auto-samplers
<b>Cost</b>	Moderately Inexpensive (<\$500)	Inexpensive (<\$100)	Expensive (>15K +)	Expensive (>20K +)