

# Identification of CC10 using LSAB Kit in Formalin-Fixed, Paraffin Embedded Rodent Tissue

## Reagents:

[1X Automation Buffer](#)  
[3% Hydrogen Peroxide](#)  
[Antibody Diluent](#)  
[Citrate Buffer](#)  
[DAB Chromagen](#)  
[Hematoxylin](#)

## Antibody Information

### Avidin Biotin Blocking Kit

Vector Laboratories, Inc.  
Burlingame, CA 94010  
[www.vectorlabs.com](http://www.vectorlabs.com)  
1-800-227-6666  
Catalog #SP-2001

### Primary antibody: Goat anti-CC10 (T-18) Polyclonal Antibody

Santa Cruz Biotechnology, Inc.  
Santa Cruz, CA 95060  
[www.scbt.com](http://www.scbt.com)  
1-800-457-3801  
Catalog #sc-9772  
Concentration: 200ug/ml

### Dako Protein Block Serum-Free

DAKO Corporation  
Carpinteria CA 93013  
[www.dakousa.com](http://www.dakousa.com)  
Code no. X0909

### DAKO LSAB+ System HRP

DAKO Corporation  
Carpinteria CA 93013  
[www.dakousa.com](http://www.dakousa.com)  
Code No. K06901

## Staining Procedure

-Positive Control Tissue: Mouse lung; bronchiolar airways. Identifies clara cells

-Stain Localization: cytoplasm

Deparaffinize and hydrate slides through the following solutions.

Xylene	2 times	5 minutes
100% EtOH	2 times	3 minutes
95% EtOH	2 times	3 minutes
1X Automation Buffer	2 times	5 minutes

1. Quench endogenous peroxidase by placing slides in 3% hydrogen peroxide for 15 minutes.

2. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

3. Unmasking Techniques using the decloaker.

Add 500 ml D/W to the pan of the decloaker.

Place full rack of slides in 200 ml of 1X citrate buffer and place in the decloaker.

Decloak for 5 minutes. Pressure \_\_\_\_\_

Depressurize for 10 minutes.

Remove pan top and cool for 10 min. Temp \_\_\_\_\_

Rinse in D/W, 2x for 3 min each

4. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes.

5. Use Dako Serum-free Protein Block for 10 mins

Lot# \_\_\_\_\_ Exp. Date \_\_\_\_\_

**DO NOT RINSE SLIDES. CONTINUE TO AVIDIN-BIOTIN BLOCK.**

6. Apply Avidin/Biotin block

Lot# \_\_\_\_\_ Exp. Date \_\_\_\_\_ New Kit: yes / no

Apply avidin block - 15 min at RT.

Quick rinse in 1X AB.

Apply biotin block - 15 min at RT.

**DO NOT RINSE SLIDES WITH BUFFER BEFORE ADDING PRIMARY ANTIBODY.**

7. Apply primary antibody (Goat anti-CC-10) at 1:100 dilution and incubate for 30 mins.  
Lot#\_\_\_\_\_ Aliquoted yes / no Date Aliquoted\_\_\_\_\_

For negative control slides, normalize the protein concentration of normal goat serum to the protein concentration of the primary antibody (anti-CC10) and use this to make the 1:100 dilution. Apply normal goat serum to the slides and incubate for one hour.  
Lot#\_\_\_\_\_ Reconstituted Date\_\_\_\_\_

8. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

9. Apply Link (yellow bottle) from LSAB+ Kit for 15 mins  
Lot#\_\_\_\_\_ Exp\_\_\_\_\_

10. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

11. Apply label (red bottle) from LSAB+ kit and incubate for 15 minutes.

12. Rinse slides in 2 changes of 1X Automation Buffer for 5 minutes each.

13. Apply liquid Dako DAB Chromagen for 6 minutes in the dark.  
(Add 1 drop of DAB per ml of substrate)  
Exp.Date\_\_\_\_\_ New Kit: yes / no Lot#\_\_\_\_\_

14. Rinse in tap water 3 minutes.

15. Counterstain with Modified Harris Hematoxylin for 30 seconds.

16. Rinse in tap water until water is clear.

17. Place slides in 1X Automation Buffer for 1 minute with gentle agitation to blue slides.

18. Dehydrate through the following solutions.

95% Ethanol	1 change	3 minutes
100% EtOH	3 changes	3 minutes
Xylene	2 changes	5 minutes

19. Coverslip.  
updated 10/25/04