

510(k) SUMMARY

This Summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92

The assigned 510(k) number is: K073355

A. Introduction:

According to the requirements of 21 CFR 807.92 the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

B. Submitter's information

Name: Thermo Fisher Scientific Oy
Address: Ratastie 2
P.O. Box 100
FIN-01621 Vantaa
Finland
Phone: +358 (9) 329 100 tel
Fax: +358 (9) 3291 0500 fax
Contact person: Päivi Sormunen, Vice President of QRC
Date of Preparation: November 23rd , 2007

C. Device name

Proprietary name: Calcium, codes 981367 and 981772
Common name: Calcium
Classification: Clinical Chemistry
Class: II
Product Code: CJY

Proprietary name: sCal, code 981831
Common Name: Calibrator, Multi-Analyte Mixture
Classification: Clinical Chemistry
Class: II
Product Code: JIX

Proprietary name: Nortrol, code 981043
Common Name: Multi-analyte Controls (Assayed and unassayed)
Classification: Clinical Chemistry
Class: I
Product Code: JJY

Proprietary name: Abtrol, code 981044
Common Name: Multi-analyte Controls (Assayed and unassayed)
Classification: Clinical Chemistry
Class: I
Product Code: JJY

D. Intended Use

Calcium

For *in vitro* diagnostic use in the quantitative determination of the calcium concentration in human serum or plasma on T60 instruments.

sCal, code 981831

For *in vitro* diagnostic use on T60 instrument. sCal is used as a multicalibrator for quantitative measurements using methods defined by Thermo Fisher Scientific Oy.

Nortrol

For *in vitro* diagnostic use for quantitative testing on T60 instrument. Nortrol is a control serum to monitor trueness and precision of the analytes listed in the separate Nortrol value sheet. The given values are valid for T60 Clinical Chemistry Instruments using methods defined by Thermo Fisher Scientific Oy.

Abtrol

For *in vitro* diagnostic use for quantitative testing on T60 instrument. Abtrol is a control serum to monitor trueness and precision of the analytes listed in the separate Abtrol value sheet. The given values are valid for T60 Clinical Chemistry Instruments using methods defined by Thermo Fisher Scientific Oy.

E. Indications for use

The Calcium test system is intended for *in vitro* diagnostic use in the quantitative determination of the calcium concentration in human serum or plasma on T60 instruments. Calcium measurements are used in the diagnosis and treatment of parathyroid disease, a variety of bone diseases, chronic renal disease and tetany.

For sCal Calibrator, Nortrol and Abtrol see intended use.

F. Substantial Equivalence

Bayer Corporation, model Bayer ADVIA 2400 Chemistry System.

Bayer Corporation item: Bayer ADVIA Calcium assay.

G. Substantial equivalence -similarities

Calcium is substantially equivalent to other devices legally marketed in United States. We claim equivalence to the Bayer ADVIA Calcium assay (K991576).

The following table compares the Calcium with the predicate device
Table 1

Attribute	<u>New device #1</u>	<u>Predicate device #1</u>
Intended Use	For <i>in vitro</i> diagnostic use in the quantitative determination of the calcium concentration in human serum or plasma on T60 instruments.	For <i>in vitro</i> diagnostic use in the quantitative determination of calcium in human serum, plasma (lithium heparin), and urine on the ADVIA Chemistry systems. Such measurements are used in the diagnosis and treatment of parathyroid disease, a variety of bone diseases, chronic renal failure, and tetany.
Indication for Use	The Calcium test system is intended for quantitative <i>in vitro</i> diagnostic measurement of calcium concentration in human serum or plasma. Calcium measurements are used in the diagnosis and treatment of parathyroid disease, a variety of bone diseases, chronic renal disease and tetany	See intended use.
Assay Protocol	Calcium ions form a highly coloured complex with Arsenazo III at neutral pH. The amount of the complex is measured at 660 nm.	Calcium ions form a violet complex with <i>o</i> -cresolphthalein complexone in an alkaline medium. The reaction is measured at 545/658 nm.
Traceability/Standardization	The value of Calcium has been assigned by using NIST SRM 909b as a primary reference	The ADVIA CA method is traceable to a NIST atomic absorption reference method, which uses reference materials from the National Institute of Standards and Technology (NIST), via patient sample correlation.
Sample Type	Serum, plasma (Li-heparin)	Serum, plasma (Li-heparin) and urine
Reagent Storage	Reagents in unopened vials are stable at 2...25°C until the expiration date printed on the label. Keep away from sunlight.	For all systems, unopened reagents are stable until the expiration date printed on the product label when stored at 15° – 25°C. Do not freeze reagents.

Expected Values	Serum/ plasma 8.6 - 10.3 mg/dl (2.15 - 2.57 mmol/l)	Serum/Plasma: 8.3 – 10.6 mg/dL (2.08 – 2.65 mmol/L) Urine: 100 – 300 mg/day (2.50 – 7.50 mmol/day)
Instrument	T60 and DPC T60i, DPC T60i Kusti	ADVIA [®] 2400 Chemistry system.
Measuring Range	Serum/ plasma 2.8 - 16.0 mg/dl (0.70 - 4.00 mmol/l)	Serum/ Plasma: 1.0 – 15.0 mg/dL (0.25 – 3.75 mmol/L) Urine: 1.0 – 30.0 mg/dL (0.25 – 7.50 mmol/L)
Precision	<p>Serum:</p> <p>Within run</p> <p>Level 4.0 mg/ dL SD= 0.04 CV(%)= 1.0</p> <p>Level 8.4 mg/ dL SD= 0.07 CV(%)= 0.8</p> <p>Level 11.9 mg/ dL SD= 0.08 CV(%)= 0.7</p> <p>Between run</p> <p>Level 4.0 mg/ dL SD= 0.02 CV(%)= 0.6</p> <p>Level 8.4 mg/ dL SD= 0.03 CV(%)= 0.4</p> <p>Level 11.9 mg/ dL SD= 0.08 CV(%)= 0.7</p> <p>Total</p> <p>Level 4.0 mg/ dL SD= 0.06 CV(%)= 1.6</p> <p>Level 8.4 mg/ dL SD= 0.12 CV(%)= 1.5</p> <p>Level 11.9 mg/ dL SD= 0.18 CV(%)= 1.5</p>	<p>Serum:</p> <p>Within run</p> <p>Level 6.2 mg/dL SD= 0.06 CV(%)= 1.0</p> <p>Level 8.5 mg/dL SD= 0.17 CV(%)= 2.0</p> <p>Level 10.9 mg/dL SD= 0.18 CV(%)= 1.6</p> <p>Total</p> <p>Level 6.2 mg/dL SD= 0.12 CV(%)= 2.0</p> <p>Level 8.5 mg/dL SD= 0.21 CV(%)= 2.4</p> <p>Level 10.9 mg/dL SD= 0.23 CV(%)= 2.1</p> <p>Urine:</p> <p>Within run</p> <p>Level 5.8 mg/dL SD= 0.10 CV(%)= 1.7</p> <p>Level 12.7 mg/dL SD= 0.08 CV(%)= 0.6</p>

		<p>Total Level 5.8 mg/dL SD= 0.12 CV(%)= 2.1 Level 12.7 mg/dL SD= 0.64 CV(%)= 5.1</p>
Method Comparison	<p>Serum (Comparison to Bayer ADVIA 2400): $y = 1.04x - 0.002$ R = 0.994 range from 0.44 to 18.25 mg/dL N = 112</p>	<p>Serum (comparison to ADVIA 1650): $y = 1.01x + 0.27$ R = 0.988 Range: 4.6 – 12.1 mg/dL N = 242</p> <p>Serum (comparison to reference method): $y = 1.00x - 0.56$ R = 0.996 Range: 2.5 – 13.9 mg/dL N = 48</p> <p>Urine (comparison to ADVIA 1650): $y = 0.98x + 0.23$ R = 0.999 Range: 1.0 – 15.0 mg/dL N = 64</p>
Limitations	<p>Lipemia: No interference found up to 1000 mg/dL (10 g/l) of Intralipid.</p> <p>Hemolysate: No interference found up to 1000 mg/dl (10 g/l) of hemoglobin</p> <p>Bilirubin, conjugated: No interference found up to 58 mg/dL (1000 µmol/l) of conjugated bilirubin.</p> <p>Bilirubin, unconjugated: No interference found up to 58 mg/dL (1000 µmol/l) of unconjugated bilirubin.</p>	<p>Lipemia (from Intralipid): No significant interference found up to 625 mg/dl of Intralipid.</p> <p>Hemolysate: No significant interference found up to 525 mg/dl of hemoglobin.</p> <p>Bilirubin: No significant interference found up to 30 mg/dl.</p>



FEB 28 2008

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Thermo Fisher Scientific Oy
c/o Päivi Sormunen
Vice President of Industrial Solutions & QRC
Clinical Diagnostics Finland
Ratastie 2, P.O. Box 100
Fin-01621 Vantaa
Finland

Re: k073355
Trade Name: Calcium, sCal, Nortrol, Abtrol
Regulation Number: 21 CFR 862.1145
Regulation Name: Calcium Test System
Regulatory Class: Class II
Product Codes: CJY, JIX, JJY
Dated: January 24, 2008
Received: January 28, 2008

Dear Päivi Sormunen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0490. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address at <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Jean M. Cooper, M.S., D.V.M.

Director

Division of Chemistry and Toxicology

Office of *In Vitro* Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K073355

Device Name: Calcium, sCal, Nortrol, Abtrol

The Calcium test system is intended for quantitative *in vitro* diagnostic measurement of calcium concentration in human serum or plasma. Calcium measurements are used in the diagnosis and treatment of parathyroid disease, a variety of bone diseases, chronic renal disease and tetany

For *in vitro* diagnostic use on T60 instrument. sCal is used as a multicalibrator for quantitative measurements using methods defined by Thermo Fisher Scientific Oy

For *in vitro* diagnostic use for quantitative testing on T60 instrument. Nortrol is a control serum to monitor trueness and precision of the analytes listed in the separate Nortrol value sheet. The given values are valid for T60 Clinical Chemistry Instruments using methods defined by Thermo Fisher Scientific Oy.

For *in vitro* diagnostic use for quantitative testing on T60 instrument. Abtrol is a control serum to monitor trueness and precision of the analytes listed in the separate Abtrol value sheet. The given values are valid for T60 Clinical Chemistry Instruments using methods defined by Thermo Fisher Scientific Oy.


Prescription Use (21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use (21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)



Division Sign-Off
Office of In Vitro Diagnostic Device
Evaluation and Safety

510(k) K0 73355