

# Making RNA probes for *in situ* hybridization

## Reagents

**Nuclease free water**  
Sigma W4502

**NTPs**  
Roche order #1277057

**DIG NTP mix**  
Roche order #1277073

**Biotin NTP mix**  
Roche order #1685597

**Fluorescein NTP mix**  
Roche order #1685619

**DNP-11-UTP**  
Perkin Elmer NEL555001EA

**Template**  
via PCR, purified and resuspended in nuclease free water

**RNA polymerases**  
T3 - Ambion order #2062  
T7 - Ambion order #2073  
SP6 - Ambion order #2082

**RNAse inhibitor**  
SUPERase-In Ambion #2694

**RNAse free DNase I**  
Boehringer Mannheim 0776785

**RNAse free EtOH**

**RNAse free 3M NaOAc pH5.2**  
Dissolve NaOAc into ddH<sub>2</sub>O and adjust pH with HCl

## Make DNA templates via PCR

### Day 1

It's preferable to start with constructs that contain RNA polymerase start sites (T3, T7, or SP6) which allow use of primers outside the start-point. This gives the enzyme some docking space, improving transcription efficiency. If this is not possible, use extended primers with RNA-polymerase start sites in their 5' tails.

Purify PCR products with QiaQuick PCR purification kit and resuspend in nuclease free water.

## Make DNP-NTP mixture (if needed)

### Day 1

	Amount	Final concentration
<b>Sigma H2O</b>	5.7 µl	
<b>ATP 100mM</b>	2 µl	10 mM
<b>CTP 100mM</b>	2 µl	10 mM
<b>GTP 100mM</b>	2 µl	10 mM
<b>UTP 100mM</b>	1.3 µl	6.5 mM
<b>DNP-11-UTP 10mM</b>	7.0 µl	3.5 mM

## In vitro transcription reaction

### Day 1

	10 µl reaction	20 µl reaction
<b>Sigma H2O to final vol</b>		
<b>10X Buffer</b>	1 µl	2 µl
<b>10X NTP mix</b>	1 µl	2 µl
<b>RNAse inhibitor</b>	0.7 µl	1.3 µl
<b>template DNA</b>	>150ng	>300ng
<b>RNA polymerase</b>	1 µl	2 µl

*Heat 10X Buffer to dissolve salts before using.*

Incubate at 37°C for 2-4 hours.

## Eliminate template

### Day 1

Add 1  $\mu$ l DNaseI (RNase free) to each reaction.

Incubate for 15 min at 37°C. *This is a time sensitive step.*

---

## Precipitate probes

### Day 1 (immediately after eliminating the template)

Bring reaction to 50 $\mu$ l with RNase free water.

Add 1/10th volume 3M NaOAc pH5.2 (5 $\mu$ l)

Add 2.5 volumes of 100% EtOH (125 $\mu$ l)

Store at -20°C for at least 2 hours.

Spin 20 min at 4°C at 13200 rpm.

Remove supernatant.

*Pellet for DIG/ biotin reactions will be white, DNP reactions will be yellow and FITC reactions will be orange.*

Add 250-400 $\mu$ l 70% EtOH and spin 20 min at 4°C at 13200 rpm.

Remove supernatant and repeat wash with 100% EtOH.

For a 20 $\mu$ l reaction with a visible pellet, it is usually sufficient to dissolve the pellet directly in 100 $\mu$ l 50% formamide.

If necessary: resuspend probes in 50  $\mu$ l of RNase free water. Take 5  $\mu$ l for quantitation. Then add 45  $\mu$ l 100% formamide and store at -20°C.

For 10 $\mu$ l reaction : with visible pellet, it is sufficient to dissolve them directly in 50 $\mu$ l 50% formamide in RNase free water.

*Properly dissolved probes that have sufficiently high concentrations typically foam a lot. It will take a couple of minutes to dissolve a probe.*

---