



Section 3:

Prokaryotic Sample and Array Processing



Section 3






Contents

Section 3 *Prokaryotic Sample and Array Processing*

Chapter 1	<i>Direct Labeling of Enriched E. coli mRNA</i>	3.1.3
Chapter 2	<i>Preparation of Control Spike Transcripts for GeneChip E. coli Genome Array</i>	3.2.3
Chapter 3	<i>E. coli Target Hybridization</i>	3.3.3
Chapter 4	<i>GeneChip E. coli Genome Array: Washing, Staining, and Scanning</i>	3.4.3




AFFYMETRIX®



E. coli

Section 3, Chapter 1

Section 3, Chapter 1





Direct Labeling of Enriched *E. coli* mRNA

Introduction	3.1.4
Reagents and Materials Required	3.1.5
Reagent Preparation	3.1.7
Total RNA Isolation.	3.1.9
mRNA Enrichment Procedure	3.1.10
Step 1: cDNA Synthesis	3.1.10
Step 2: rRNA Digestion.	3.1.11
Step 3: cDNA Digestion	3.1.11
RNA Fragmentation and Labeling Reaction	3.1.12
Step 1: RNA Fragmentation	3.1.12
Step 2: RNA 5'-Thiolation	3.1.12
Step 3: Biotin Addition	3.1.13
Gel-Shift Assay	3.1.14

This Chapter Contains:

- Instructions for enriching, fragmenting, and labeling RNA target from *E. coli* cells.

After completing the procedures described in this chapter, the labeled and fragmented target is hybridized to the GeneChip *E. coli* Genome Array, as described in Section 3, Chapter 3.

Introduction

This chapter describes the assay procedures recommended for use with GeneChip *E. coli* Genome Array. These procedures differ from those described under the previous section primarily because, unlike their eukaryotic counterpart, the majority of *E. coli* messenger RNAs do not contain poly-A tails. To enrich for mRNA species, a procedure designed to remove 16S and 23S ribosomal RNAs, which constitute approximately 90% of the total RNA population, is utilized. The resulting RNA population is then directly labeled with biotin.

mRNA Enrichment Procedure

The enrichment procedure is a series of enzymatic steps that specifically eliminate the 16S and 23S rRNA species in the total *E. coli* RNA. Reverse transcriptase and primers specific to 16S and 23S rRNA are used to synthesize complementary DNAs. Then rRNA is removed enzymatically by treatment with RNase H, which specifically digests RNA within an RNA:DNA hybrid. The cDNA molecules are then removed by DNase I digestion and the enriched mRNA is purified on QIAGEN RNeasy columns.

RNA Fragmentation and Labeling Procedure

The direct labeling of RNA consists of the following steps:

- 1.** The RNA is fragmented by heat and ion-mediated hydrolysis.
- 2.** The 5'-end RNA termini are enzymatically modified by T4 polynucleotide kinase with γ -S-ATP.
- 3.** Biotin is conjugated to 5'-ends of the RNA. After purification of the product, the efficiency of the labeling procedure can be assessed using a gel-shift assay described in *Gel-Shift Assay* on page 3.1.14.

Reagents and Materials Required

The following reagents and materials are recommendations and have been tested and evaluated by Affymetrix scientists. For supplier phone numbers in the U.S. and Europe, please refer to the Supplier and Reagent Reference List, Appendix A, of this manual. Information and part numbers listed are based on U.S. catalog information. Additional reagents needed for the complete analysis are listed in the appropriate chapters. Appendix A contains a master list of all reagents used in this manual.

Total RNA Isolation

- MasterPure™ RNA Purification Kit, Epicentre Technologies, P/N MCR85102
- Isopropanol
- 75% Ethanol

mRNA Enrichment Procedure

- MMLV Reverse Transcriptase, New England BioLabs, P/N M0253L
- 10X MMLV Reverse Transcriptase Buffer, New England BioLabs, P/N M0253L (contains DTT and is also supplied with MMLV Reverse Transcriptase)
- 16S rRNA Primers

Name	Sequence
16S1514	5'-CCTAC GGTGA CCTTG TT-3'
16S889	5'-TTAAC CTTGC GGCCG TACTC-3'
16S541	5'-TCCGA TTAAC GCTTG CACCC-3'

- 23s rRNA Primers

Name	Sequence
23S2878	5'-CCTCA CGGTT CATTG GT-3'
23SE _{co} 2064	5'-CTATA GTAAA GGTTC ACGGG-3'
23SE _{co} 1595	5'-CCTGT GTCGG TTTGG GGT-3'
23S1022	5'-TCCCA CATCG TTTCC CAC-3'
23S539	5'-CCATT ATACA AAAGG TAC-3'

✓ Note

All primers should be purchased in the 1 μmole amount and purified by HPLC. Oligonucleotides purchased from Operon have been used successfully in the procedure described in this chapter.

- dATP, dCTP, dGTP, dTTP, Amersham Pharmacia Biotech, P/N 27-2035-01
- SUPERase•In™, Ambion, P/N 2696
- Ribonuclease H (RNase H), *E. coli*, Epicentre Technologies, P/N R0601K
- Deoxyribonuclease I (DNase I), Amersham Pharmacia Biotech, P/N 27-0514-01
- 0.5 M EDTA, pH 8.0, Invitrogen Life Technologies, P/N 15575-038
- RNeasy Mini Kit, QIAGEN, P/N 74104
- β-Mercaptoethanol
- 96% to 100% Ethanol
- Nuclease-free Water, Ambion, P/N 9930

RNA Fragmentation and Labeling Reaction

- T4 Polynucleotide Kinase and 10X NEBuffer, New England BioLabs, P/N 201L
- γ -S-ATP (20 μ moles), Roche Molecular Biochemical, P/N 1162306
- MOPS, Sigma-Aldrich, P/N M3183
- PEO-Iodoacetyl-Biotin (50 mg), Pierce Chemical, P/N 21334ZZ
- 3 M Sodium Acetate (NaOAc), pH 5.2, Sigma-Aldrich, P/N S 7899
- 96% to 100% Ethanol
- RNA/DNA Mini Column Kit, QIAGEN, P/N 14123
- β -Mercaptoethanol
- Isopropanol
- 70% Ethanol
- Glycogen (20 mg/mL), Roche Molecular Biochemical, P/N 901393

Gel-Shift Assay

- Novex XCell SureLock™ Mini-Cell, Invitrogen, P/N EI9001
- 4-20% TBE Gel, 1.0 mm, 12 well, Invitrogen, P/N EC62252
- Sucrose Gel Loading Dye, 5X, Amresco, P/N E-274
- 10X TBE Running Buffer
- SYBR Gold, Molecular Probes, P/N S-11494
- 10 bp and 100 bp DNA ladder, Invitrogen Life Technologies, P/N 10821-015 and 15628-019, respectively
- ImmunoPure NeutrAvidin, Pierce Chemical, P/N 31000ZZ
- 1M Tris pH 7.0, Ambion, P/N 9850G

Miscellaneous Supplies

- Hybridization Oven 640, Affymetrix, P/N 800139
- Sterile, RNase-free, microcentrifuge tubes, 1.5 mL, USA Scientific, P/N 1415-2600 (or equivalent)
- Micropipettors, (P-2, P-20, P-200, P-1000), Rainin Pipetman or equivalent
- Sterile-barrier pipette tips and non-barrier pipette tips
- Thermocycler
- Refrigerated microcentrifuge
- Spectrophotometer and quartz micro cuvettes
- Gel imaging system with appropriate filter for SYBR Green II or Gold.
- PCR tubes
- 0.2 μ m filters for liquid sterilization
- Tough Spots, Label Dots, USA Scientific, P/N 9902 (optional)

Reagent Preparation

mRNA Enrichment Procedure

rRNA Removal Primer Stock

Individual primers are maintained at a stock concentration of 100 μM . Each primer is diluted 10-fold in the final primer mix to a concentration of 10 μM . Individual stocks and the primer mix are stored at -20°C .

25 mM dNTP Mix

Add 50 μL of each dNTP stock solution (100 mM each) to a sterile microcentrifuge tube. Mix thoroughly. Store at -20°C .

RNA Fragmentation and Labeling Reaction

500 mM MOPS, pH 7.5

For 100 mL:

10.5 g MOPS

80 mL Distilled water

Mix and adjust pH to 7.5 with 1 M NaOH. Add Distilled water to 100 mL.

Filter through 0.2 μM filter. Store at room temperature.

Gel-Shift Assay

2 mg/mL NeutrAvidin

Resuspend 10 mg NeutrAvidin in 5 mL solution containing 50 mM Tris-Cl, pH 7.6, 200 mM NaCl. Store at 4°C .

Total RNA Isolation

Prior to the enrichment process, total RNA is isolated. For *E. coli*, Affymetrix has successfully used the Epicentre MasterPure RNA Purification Kit. Typical yields from a 10 mL log-phase culture range from 80 to 400 μg of total RNA. After purification, the RNA concentration is determined by absorbance at 260 nm on a spectrophotometer (1 absorbance unit = 40 $\mu\text{g}/\text{mL}$ RNA). The A_{260}/A_{280} ratio should be approximately 2.0, with ranges between 1.8 to 2.1 considered acceptable. We recommend checking the quality of the RNA by running it on an agarose gel prior to starting the assay. The figure below shows 1 μg samples from three acceptable RNA preparations. The 23S and 16S rRNA bands should be clear without any obvious smearing patterns.




Figure 3.1.1
Typical RNA preparations from *E. coli*

mRNA Enrichment Procedure

The following protocol starts with 100 µg of total RNA. Incubations are performed in a thermocycler. Each sample is equally split between four 200-µL PCR tubes. It takes approximately 5 hours to obtain enriched mRNA.

✓ Note

The enrichment procedure involves three enzymatic steps, and the amount of enzymes used and incubation time are critical to the success of the assay. Please exercise precautions and follow standard laboratory procedures when handling RNA samples.

Step 1: cDNA Synthesis

1. Prepare the following mixture in each of the four PCR tubes.

✓ Note

*Optionally, 1 µL of control transcript mix (lys, phe, dap, thr, and trp from *B. subtilis*) can be added to total RNA prior to the enrichment procedure for a final concentration of 10 pM of each control transcript to be applied to the arrays. Please refer to Section 3, Chapter 2 for the detailed protocol to prepare control sense RNA.*

*The detection limit of the assay is estimated to be around 5 pM. Assuming complete recovery of spike transcripts, the final hybridization mix contains each transcript at 10 pM concentration that is slightly above detection limit. Detection of these controls on the GeneChip *E. coli* Genome Array serves as indicators of the labeling efficiency. Alternatively, various control transcripts can also be spiked in at variable concentrations to demonstrate the dynamic range of the assay.*

Table 3.1.1

Primer Mixture for mRNA Enrichment

Components	Volume or Amount	Final Concentration
Total RNA	25.0 µg	0.83 µg/µL
10 µM rRNA Removal Primer Stock	7 µL	1.75 µM
500 pM Control Transcript (optional)	1 µL	
Nuclease-free DI H ₂ O	Up to 40.0 µL	—
Total Volume Added	40 µL	

2. Heat RNA and primer mixture to 70°C for 5 minutes and then cool to 4°C.
3. Add the following components to each of the four tubes.

Table 3.1.2

Reverse Transcription Components

Components	Volume	Final Concentration
10X MMLV RT Buffer	10.0 µL	1X
25.0 mM dNTP Mix	2.0 µL	0.5 mM
20 U/µL SUPERase•In	3.0 µL	0.6 U/µL
200 U/µL MMLV RT	2.5 µL	5 U/µL
Nuclease-free H ₂ O	42.5 µL	—
Total Volume Added	60 µL	

4. Incubate the reaction at 42°C for 25 minutes, then at 45°C for 20 minutes. Cool to 4°C.
5. Inactivate the enzyme at 65°C for 5 minutes, then hold at 4°C. Proceed immediately to *Step 2: rRNA Digestion*, below.

Step 2: rRNA Digestion

1. Add the following to each of the four tubes.

Table 3.1.3
RNase Digestion of rRNAs

Components	Volume	Final Concentration
10.0 U/μL RNase H	4.0 μL	0.4 U/μL
20 U/μL SUPERase•In	3.0 μL	0.6 U/μL
Total Volume Added	7 μL	

2. Incubate the reaction at 37°C for 25 minutes, then hold at 4°C. Proceed immediately to *Step 3: cDNA Digestion* below.

Step 3: cDNA Digestion

1. Add the following to each of the four tubes.

Table 3.1.4
Digestion of cDNA

Components	Volume	Final Concentration
5.0 U/μL DNase I*	4.0 μL	0.2 U/μL

*Dilute 10 U/μL DNase I to 5 U/μL with nuclease-free H₂O.

✓ Note

The concentration of DNase I is critical. Excessive digestion may cause degradation of mRNA. Lot-to-lot variation of enzyme activity may occur. Therefore, titrate DNase concentrations when necessary.

2. Incubate the reaction at 37°C for 20 minutes.
3. Inactivate the enzyme by adding 3 μL of 500 mM EDTA to each tube for a final concentration of 10.0 mM.

➔ IMPORTANT

The concentration and incubation time for DNase I are critical. Under-digestion is preferable to over-digestion.

4. Pool reaction product from all four tubes and clean up the enriched mRNA with QIAGEN RNeasy Mini column.
5. Quantify the enriched mRNA preparation by 260 nm absorbance. Typical yields for the procedure are 20 to 40 μg of RNA (1.0 A₂₆₀ unit = 40 μg/mL single strand RNA).

✓ Note

The enriched mRNA is stored at -20°C until ready for use in the subsequent fragmentation and labeling reaction.

RNA Fragmentation and Labeling Reaction

As in the previous procedure, the fragmentation and labeling reactions are done in PCR tubes in a thermocycler. A maximum of 20 µg of RNA per tube is used in the fragmentation step. Incomplete fragmentation may occur if excess RNA is used. Split samples into multiple tubes if the yield of RNA from the enrichment step is greater than 20 µg. Following the 5' thiolation and biotin addition reaction, the target is purified with ethanol precipitation or RNA/DNA Mini Columns. The reactions take about 6 hours to complete.



Note

The reaction of PEO-Iodoacetyl-Biotin is highly pH-dependent. MOPS is used as the buffer because of its inability to react with the iodoacetyl group under the reaction conditions. Other commonly used laboratory buffers containing amino groups such as Tris should not be used because of their reactivity with the iodoacetyl moiety.

Step 1: RNA Fragmentation

1. Prepare the following mixture.

Table 3.1.5
RNA Fragmentation Mix

Components	Volume or Amount	Final Concentration
10X NEBuffer for T4 Polynucleotide Kinase	10.0 µL	1.1X
enriched mRNA	up to 20.0 µg	—
DI H ₂ O	up to 88.0 µL total volume	—
Final Volume	88 µL	

2. Incubate the reaction at 95°C for 30 minutes.
3. Cool to 4°C.

Step 2: RNA 5'-Thiolation

1. Prepare the following mixture.

Table 3.1.6
RNA Thiolation Mix

Components	Volume	Final Concentration
Fragmented RNA (from Step 1)	88.0 µL	—
5 mM γ-S-ATP	2.0 µL	0.1 mM
10 U/µL T4 Polynucleotide Kinase	10.0 µL	1 U/µL
Final Volume	100.0 µL	

2. Incubate the reaction at 37°C for 50 minutes.
3. Inactivate the reaction by heating at 65°C for 10 minutes and then cool to 4°C.

4. Remove excess γ -S-ATP by ethanol precipitation. If you have started with multiple tubes, combine all samples in one sterile microcentrifuge tube. Add 1/10 volume of 3 M sodium acetate, pH 5.2, and 2.5 volumes of ethanol. Leave on ice for 15 minutes.
5. Spin at 14,000 rpm at 4°C for 30 minutes to pellet the RNA.
6. Resuspend the RNA pellet in 90 μ L of DI H₂O.

Step 3: Biotin Addition

1. Prepare the following mixture.

Table 3.1.7
RNA Labeling Mix

Components	Volume	Final Concentration
500 mM MOPS, pH 7.5	6.0 μ L	30 mM
Fragmented thiolated RNA (from Step 2)	90.0 μ L	—
50 mM PEO-Iodoacetyl-Biotin	4.0 μ L	2 mM
Final Volume	100.0 μL	

2. Incubate the reaction at 37°C for one hour.
3. Cool to 4°C.
4. Remove unincorporated biotin label using the QIAGEN RNA/DNA Mini Columns.

✓ Note

For increased RNA recovery, use one RNA/DNA column and 5.4 mL Buffer QRV2 for every 10.0 μ g RNA. It is recommended to add 50 μ g of glycogen as carrier during the precipitation step.

Alternatively, clean up the labeling product by ethanol precipitation instead of Qiagen columns (add 50 μ g of glycogen as carrier, 1/10 volume of 3 M sodium acetate and 2.5 volume of ethanol to samples), followed by twice washing the pellets with 750 μ L of 70% ethanol.

5. Dissolve the pelleted RNA in 20 to 30 μ L of nuclease-free water.
6. Quantify product by 260 nm absorbance. Typical yields for the procedure are 2 to 4 μ g of RNA. The labeled RNA should be stored at -20°C until ready for gel analysis or hybridization, as described in Section 3, Chapter 3.

✓ Note

The process may be stopped during ethanol precipitation following the RNA/DNA column purification.

Gel-Shift Assay

After purification of the target, the efficiency of the labeling procedure can be assessed using the following procedure. This quality control protocol prevents hybridizing poorly labeled target onto the probe array. The addition of biotin residues is monitored in a gel-shift assay where the fragments are incubated with avidin prior to electrophoresis. Biotin-containing residues are retarded or shifted during the electrophoresis due to avidin binding. The nucleic acids are then detected by staining. Affymetrix routinely obtains approximately 70% labeling efficiency as shown in the gel photograph (Figure 3.1.2). The procedure takes approximately 90 minutes to complete.

✓ Note

The absence of a shift pattern indicates poor biotin labeling. The problem should be addressed before proceeding to the hybridization step.




Figure 3.1.2
Gel-shift assay for monitoring *E. coli* target labeling efficiency


1. Prepare a NeutrAvidin solution of 2 mg/mL or higher. 50 mM Tris, pH 7.0 can be used to dilute the NeutrAvidin solution.
2. Place a 4-20% TBE gel into the gel holder and load system with 1X TBE Buffer.
3. For each sample to be tested, remove two 150 to 200 ng aliquots of fragmented and biotinylated sample to fresh tubes.
4. Add 5 μ L of 2 mg/mL NeutrAvidin to each tube.
5. Mix and incubate at room temperature for 5 minutes.
6. Add loading dye to all samples to a final concentration of 1X loading dye.
7. Prepare 10 bp and 100 bp DNA ladders (1 μ L ladder + 7 μ L water + 2 μ L loading dye for each lane).
8. Carefully load samples and two ladders on gel. (Each gel well can hold a maximum of 20 μ L.)

9. Run the gel at 150 volts until the front dye (red) almost reaches the bottom. The electrophoresis takes approximately 1 hour.
10. While gel is running, prepare at least 100 mL of a 1X solution of SYBR Green II or Gold for staining.

✓ Note

SYBR Green II and Gold are light sensitive. Therefore use caution and shield the staining solution from light. Prepare a new batch of stain at least once a week.

11. After the gel is complete, break open cartridge and stain the gel in 1X SYBR Green II or Gold for 10 minutes.
12. Place the gel on the UV light box and produce an image following standard procedure. Be sure to use the appropriate filter for SYBR Green II or Gold.




AFFYMETRIX®



E. coli

Section 3, Chapter 2

Section 3, Chapter 2





Preparation of Control Spike Transcripts for GeneChip® E. coli Genome Array

Overview	3.2.4
Reagents and Materials Required	3.2.5
Bacterial Plasmid DNA Preparation	3.2.6
Linearization of Plasmid DNA Preparation	3.2.7
Purification of Linearized Plasmid DNA	3.2.7
<i>In Vitro</i> Transcription (IVT) to Produce Control Sense Transcripts	3.2.8
Preparing the Control Transcript Mix	3.2.8

This Chapter Contains:

- Detailed steps for producing full-length control spike sense RNA.


After completing the procedures described in this chapter, the control sense transcripts can be added to purified *E. coli* RNA samples prior to enrichment and labeling procedure as described in Section 3, Chapter 1.

Overview

This chapter describes protocols used to generate sense RNA controls from *B. subtilis* genes. These control transcripts can be spiked into *E. coli* total RNA used for target preparation at a predetermined concentration to monitor labeling, hybridization, and staining efficiency.

To be used as control for assay performance, the GeneChip® *E. coli* Genome Array contains probe sets with sequences of *dap*, *thr*, *phe*, *lys*, and *trp* genes from *B. subtilis*. These genes have been cloned into Stratagene pBluescript as an *Xho I* to *Not I* insert, 5' to 3', respectively (see Section 2, Chapter 2, *Controls for Eukaryotic Arrays*).

pGIBS-lys	ATCC 87482
pGIBS-phe	ATCC 87483
pGIBS-thr	ATCC 87484
pGIBS-trp	ATCC 87485
pGIBS-dap	ATCC 87486



These clones can be digested with the *Not I* restriction enzyme to produce linear template DNA for the subsequent *in vitro* transcription (IVT) to produce sense strand RNA by T3 RNA polymerase as control molecules.

Bacteria containing these recombinant plasmids can be obtained from the American Type Culture Collection (ATCC).

Reagents and Materials Required

The following reagents and materials are recommendations and have been tested and evaluated by Affymetrix scientists. For supplier phone numbers in the U.S. and Europe, please refer to the Supplier and Reagent Reference List, Appendix A, of this manual. Information and part numbers listed are based on U.S. catalog information. Additional reagents needed for the complete analysis are listed in the appropriate chapters. Appendix A contains a master list of all reagents used in this manual.

- Expression Control Clones, American Type Culture Collection (ATCC)
 - pGIBS-lys ATCC 87482
 - pGIBS-phe ATCC 87483
 - pGIBS-thr ATCC 87484
 - pGIBS-trp ATCC 87485
 - pGIBS-dap ATCC 87486
- *Not I* restriction Endonuclease, New England BioLabs, P/N R0189S
- Phase Lock Gel, Brinkmann Instruments, P/N 955 15 415
- Phenol/chloroform/isoamyl alcohol, Ambion, P/N 9732
- MEGAscript™ T3 Kit, Ambion, P/N 1338

Miscellaneous Reagents

- 3 M NaAcetate (NaOAc)
- Absolute Ethanol
- 80% Ethanol
- RNeasy Mini Kit, QIAGEN, P/N 74104

Bacterial Plasmid DNA Preparation

- 1.** Grow *E. coli* bacterial cultures containing recombinant plasmids according to established protocols (a minimum 50 mL of culture volume is recommended).
- 2.** Prepare plasmid DNA from overnight cultures using standard procedures or commercial kits.
Affymetrix has obtained reliable results using QIAGEN Plasmid Kits for plasmid DNA isolation.

Linearization of Plasmid DNA Preparation

1. In a 50 μ L reaction volume, digest 10 μ g of plasmid with the restriction enzyme, *NotI*, according to the enzyme manufacturer's recommendations.
2. Analyze 50 ng of the uncut and linearized plasmid by gel electrophoresis on a 1% agarose gel. Complete digestion of the plasmid is required for IVT. Repeat restriction enzyme digestion, if necessary.

Purification of Linearized Plasmid DNA

Purify the linearized plasmid from restriction enzymes and potential RNase contaminants before proceeding to IVT using a Phase Lock Gel-phenol/chloroform extraction procedure.

Phase Lock Gels (PLG) form an inert, sealed barrier between the aqueous and organic phases of phenol-chloroform extractions. The solid barrier allows more complete recovery of the sample (aqueous phase) and minimizes interface contamination of the sample. PLG's are sold as premeasured aliquots in 1.5 mL tubes to which sample and phenol chloroform are directly added.

1. Pellet the Phase Lock Gel (1.5 mL tube with PLG I-heavy) in a microcentrifuge at $\geq 12,000 \times g$ for 20 seconds.
2. Dilute the linearized plasmid to final volume of 150 μ L with TE and add equal volume of (25:24:1) Phenol:chloroform:isoamyl alcohol (saturated with 10 mM Tris-HCl pH8.0/1 mM EDTA). Vortex.
3. Transfer the mix to the PLG tube and microcentrifuge at $\geq 12,000 \times g$ for 2 minutes.
4. Transfer the top aqueous phase to a new 1.5 mL tube.
5. Add 0.1 volumes (15 μ L) of 3 M NaOAc and 2.5 volumes (375 μ L) of absolute ethanol to the samples. Vortex.
6. Immediately centrifuge at $\geq 12,000 \times g$ in a microcentrifuge at room temperature for 20 minutes.
7. Carefully remove supernatant.
8. Wash pellet with 0.5 mL of 80% ethanol, then centrifuge at $\geq 12,000 \times g$ at room temperature for 5 minutes.
9. Remove the supernatant very carefully and air dry the pellet.
10. Resuspend DNA pellet in 15 μ L of RNase-free water.
11. Quantify the DNA by absorbance at 260 nm (50 μ g/mL of DNA for 1 absorbance unit at 260 nm).

✓ Note

The quality of DNA template can be monitored by the A_{260}/A_{280} ratio, which should be between 1.8 and 2.0 for pure DNA.

In Vitro Transcription (IVT) to Produce Control Sense Transcripts

Use MEGAscript™ T3 High Yield Transcription Kit for the IVT reaction.

1. To make up the reaction mix, follow the procedures in the instruction manual provided by Ambion.



Note

No tracer is involved in this assay.

2. Incubate the reaction for 4 hours at 37°C.
3. Cleanup the reaction product with RNeasy Mini column.
4. Quantify the transcript by absorbance at 260 nm (40 µg/mL RNA = 1 absorbance unit at 260 nm).



Note

It is recommended to examine the quality and integrity of the IVT product on an agarose gel.



IMPORTANT

Aliquot and freeze the IVT transcripts at -80°C. Avoid repeated freeze / thaw cycles.

Preparing the Control Transcript Mix

1. Prepare stock solutions for each of the five transcripts separately at 2.5 nM for each transcript.
Use the following table to calculate the amount of transcript needed to prepare 2.5 nM stock.

Table 3.2.1
Conversions for Preparing 2.5 nM Control Transcript Mix

Control RNA	Size (kb)	Molecular Weight	pMoles / µg
Lys	1	330,000	3.03
Phe	1.32	435,600	2.30
Dap	1.82	607,200	1.65
Thr	1.98	653,400	1.53
Trp	2.5 kb	825,000	1.21

2. Mix equal volume of all five transcripts so that final concentration for each is 500 pM.
3. Apply 4 µL of the transcript mix with each 100 µg of total RNA prior to the enrichment procedure as described in Chapter 1 of this section. Final concentration applied on the array for the control transcripts would be 10 pM, assuming 100% recovery.




Note

Since the hybridization mix has a volume of 200 µL, 10 pM would be the final concentration for the transcripts assuming 100% recovery throughout the enrichment and labeling procedure. Alternatively, different concentrations of transcript stock can be prepared to generate "variable" concentrations for different transcripts to monitor the dynamic range of the assay.



IMPORTANT

Aliquot and freeze the IVT transcripts at -80°C. Avoid repeated freeze / thaw cycles.




AFFYMETRIX®



E. coli

Section 3, Chapter 3

Section 3, Chapter 3





E. coli Target Hybridization

Reagents and Materials Required	3.3.5
Reagent Preparation	3.3.6
<i>E. coli</i> Target Hybridization	3.3.7

This Chapter Contains:

- Detailed steps for preparing the hybridization mix containing labeled target.
- Instructions for hybridizing the target mix to a GeneChip *E.coli* Genome Array.

After completing the procedures described in this chapter, the hybridized probe array is then ready for washing, staining, and scanning as detailed in Section 3, Chapter 4.

Reagents and Materials Required

The following reagents and materials are recommendations and have been tested and evaluated by Affymetrix scientists. For supplier phone numbers in the U.S. and Europe, please refer to the Supplier and Reagent Reference List, Appendix A, of this manual. Information and part numbers listed are based on U.S. catalog information. Additional reagents needed for the complete analysis are listed in the appropriate chapters. Appendix A contains a master list of all reagents used in this manual.

- Water, Molecular Biology Grade, BioWhittaker Molecular Applications / Cambrex, P/N 51200
- Acetylated Bovine Serum Albumin (BSA) solution (50 mg/mL), Invitrogen Life Technologies, P/N 15561-020
- Herring Sperm DNA, Promega Corporation, P/N D1811
- Micropure Separator, Millipore, P/N 42512 (optional)
- Control Oligo B2, 3 nM, Affymetrix, P/N 900301 (can be ordered separately)
- 5 M NaCl, RNase-free, DNase-free, Ambion, P/N 9760G
- MES Free Acid Monohydrate SigmaUltra, Sigma-Aldrich, P/N M5287
- MES Sodium Salt, Sigma-Aldrich, P/N M5057
- EDTA Disodium Salt, 0.5 M solution (100 mL), Sigma-Aldrich, P/N E7889

Miscellaneous Reagents

- Tough Spots, Label Dots, USA Scientific, P/N 9185 (optional)
- Tween-20, 10%, Pierce Chemical, P/N 28320

Miscellaneous Supplies

- Hybridization Oven 640, Affymetrix, P/N 800139
- Sterile, RNase-free, microcentrifuge tubes, 1.5 mL, USA Scientific, P/N 1415-2600 (or equivalent)
- Micropipettors, (P-2, P-20, P-200, P-1000), Rainin Pipetman or equivalent
- Sterile-barrier pipette tips and non-barrier pipette tips

Reagent Preparation

12X MES Stock

(1.22 M MES, 0.89 M [Na⁺])

For 1000 mL:

70.4 g MES-free acid monohydrate

193.3 g MES Sodium Salt

800 mL of Molecular Biology Grade water

Mix and adjust volume to 1000 mL.

The pH should be between 6.5 and 6.7. Filter through a 0.2 µm filter.



Do not autoclave, store at 2°C to 8°C, and shield from light. Discard solution if yellow.

2X Hybridization Buffer

(Final 1X concentration is 100 mM MES, 1 M [Na⁺], 20 mM EDTA, 0.01% Tween 20)

For 50 mL:

8.3 mL of 12X MES Stock

17.7 mL of 5 M NaCl

4.0 mL of 0.5 M EDTA

0.1 mL of 10% Tween 20

19.9 mL of water

Store at 2°C to 8°C, and shield from light

E. coli Target Hybridization

After determining that the fragmented RNA is labeled with biotin, prepare the hybridization solution mix. The recommended amount of RNA per assay is 1.5 to 4.0 μg . The solution is stable for approximately 6 to 8 hours at 4°C. The following protocol can be used for freshly prepared or frozen hybridization cocktail. Re-use of prokaryotic sample is not recommended, since the samples are end-labeled rather than internally labeled.

1. Prepare the following hybridization solution mix.

Table 3.3.1
Hybridization Solution Mix

Components	Volume	Final Concentration
2X MES Hybridization Buffer	100.0 μL	1X
3 nM Control Oligo B2	3.3 μL	50 pM
10 mg/mL Herring Sperm DNA	2.0 μL	0.1 mg/mL
50 mg/mL BSA	2.0 μL	0.5 mg/mL
Fragmented Labeled RNA	1.5 to 4.0 μg	—
Molecular Biology Grade Water	to a final volume of 200 μL	—
Final Volume	200 μL	

2. Equilibrate probe array to room temperature immediately before use.

✓ Note


It is important to allow the arrays to normalize to room temperature completely. Specifically, if the rubber septa are not equilibrated to room temperature, they may be prone to cracking which can lead to leaks.

3. Add the hybridization solution mix to the probe array.

✓ Note

It is necessary to use two pipette tips when filling the probe array cartridge: one for filling, and the second to allow venting of air from the hybridization chamber. After the addition of hybridization cocktail, the septa may be covered with Tough Spots to prevent evaporation.

4. Place probe array in the hybridization oven set at 45°C.
5. Avoid stress to the motor; load probe arrays in a balanced configuration around rotisserie axis. Rotate at 60 rpm.
6. Hybridize for 16 hours.
During the latter part of the 16-hour hybridization, proceed to Section 3, Chapter 4 to prepare reagents required immediately after completion of hybridization.




AFFYMETRIX®



E. coli

Section 3, Chapter 4

Section 3, Chapter 4





GeneChip *E. coli* Genome Array: Washing, Staining, and Scanning

Reagents and Materials Required	3.4.5
Reagent Preparation	3.4.6
Experiment and Fluidics Station Setup	3.4.7
Step 1: Defining File Locations	3.4.7
Step 2: Entering Experiment Information.	3.4.7
Step 3: Preparing the Fluidics Station.	3.4.8
Probe Array Wash and Stain	3.4.9
Probe Array Scan	3.4.13
Shutting Down the Fluidics Station	3.4.14
Customizing the Protocol	3.4.15

This Chapter Contains:

- Instructions for using the Fluidics Station 400 to automate the washing and staining of GeneChip *E.coli* Genome Arrays.
- Instructions for scanning probe arrays using the GeneArray Scanner.

After completing the procedures described in this chapter, the scanned probe array image (.dat file) is ready for analysis, as explained in the enclosed *GeneChip Expression Analysis: Data Analysis Fundamentals* booklet (P/N 701190).

Reagents and Materials Required

The following reagents and materials are recommendations and have been tested and evaluated by Affymetrix scientists. For supplier phone numbers in the U.S. and Europe, please refer to the Supplier and Reagent Reference List, Appendix A, of this manual. Information and part numbers listed are based on U.S. catalog information. Additional reagents needed for the complete analysis are listed in the appropriate chapters. Appendix A contains a master list of all reagents used in this manual.

- Water, Molecular Biology Grade, BioWhittaker Molecular Applications / Cambrex, P/N 51200
- Distilled water, Invitrogen Life Technologies, P/N 15230-147
- Acetylated Bovine Serum Albumin (BSA) solution (50 mg/mL), Invitrogen Life Technologies, P/N 15561-020
- R-Phycoerythrin Streptavidin, Molecular Probes, P/N S-866
- 5 M NaCl RNase-free, DNase-free, Ambion, P/N 9760G
- PBS, pH 7.2, Invitrogen Life Technologies, P/N 20012-027
- 20X SSPE (3 M NaCl, 0.2 M NaH₂PO₄, 0.02 M EDTA), BioWhittaker Molecular Applications / Cambrex, P/N 51214
- Goat IgG, Reagent Grade, Sigma-Aldrich, P/N I 5256
- Anti-streptavidin antibody (goat), biotinylated, Vector Laboratories, P/N BA-0500
- 10% surfact-Amps20 (Tween-20), Pierce Chemical, P/N 28320
- Bleach (5.25% Sodium Hypochlorite), VWR Scientific, P/N 21899-504 (or equivalent)
- ImmunoPure Streptavidin, Pierce Chemical, P/N 21125

Miscellaneous Supplies

- Sterile, RNase-free, microcentrifuge tubes, 1.5 mL, USA Scientific, P/N 1415-2600 (or equivalent)
- Micropipettors, (P-2, P-20, P-200, P-1000), Rainin Pipetman (or equivalent)
- Sterile-barrier pipette tips and non-barrier pipette tips
- Tygon Tubing, 0.04" inner diameter, Cole-Palmer, P/N H-06418-04

Reagent Preparation

Wash A: Non-Stringent Wash Buffer

(6X SSPE, 0.01% Tween 20)

For 1000 mL:

300 mL of 20X SSPE

1.0 mL of 10% Tween-20

699 mL of water

Filter through a 0.2 μ m filter

Wash B: Stringent Wash Buffer

(100 mM MES, 0.1 M [Na⁺], 0.01% Tween 20)

For 1000 mL:

83.3 mL of 12X MES Stock Buffer (see Section 3, Chapter 3 for reagent preparation)

5.2 mL of 5 M NaCl

1.0 mL of 10% Tween 20

910.5 mL of water

Filter through a 0.2 μ m filter

Store at 2°C to 8°C and shield from light

2X Stain Buffer

(Final 1X concentration: 100 mM MES, 1 M [Na⁺], 0.05% Tween 20)

For 250 mL:

41.7 mL 12X MES Stock Buffer (see Section 3, Chapter 3 for reagent preparation)

92.5 mL 5 M NaCl

2.5 mL 10% Tween 20

113.3 mL water

Filter through a 0.2 μ m filter

Store at 2°C to 8°C and shield from light

10 mg/mL Goat IgG Stock

Resuspend 50 mg in 5 mL PBS

Store at 4°C

1 mg/mL Streptavidin Stock

Resuspend 5 mg in 5 mL of PBS

Store at 4°C

Experiment and Fluidics Station Setup

Step 1: Defining File Locations

Before working with Affymetrix Microarray Suite it is important to define where the program stores and looks for files.

1. Launch Microarray Suite from the workstation and select **Tools** → **Defaults** → **File Locations** from the menu bar.

The File Locations window displays the locations of the following files:

- Probe Information (library files, mask files)
- Fluidics Protocols (fluidics station scripts)
- Experiment Data (.exp, .dat, .cel, and .chp files are all saved to location selected here)

2. Verify that all three file locations are set correctly and click **OK**.

Contact Affymetrix Technical Support if you have any questions regarding this procedure.

Step 2: Entering Experiment Information

To wash, stain, and scan a probe array, an experiment must first be defined in Microarray Suite.

1. Select **Run** → **Experiment Info** from the menu bar. Alternatively, click the New Experiment icon on the tool bar.
⇒ The Experiment Information dialog box appears allowing the experiment name to be defined along with several other parameters, such as probe array type, sample description, and comments.
2. Type in the **Experiment Name**.
3. In the **Probe Array Type** box, click the arrow and select **Ecoli** from the drop-down list.

Experiment name and probe array type are required. Complete as much of the other information as desired. The protocol information at the bottom of the dialog box is exported to the experiment information dialog box after the hybridization and scan are complete.

4. Save the experiment by selecting **Save**.

The name of the experiment is used by Microarray Suite to access the probe array type and data for the sample while it is being processed. Data files generated for the sample are automatically labeled to correspond to the experiment name. Microarray Suite automatically fills in the **Protocol** section of this dialog box with information on array processing from the fluidics station.

5. Close the Experiment Information dialog box.

Step 3: Preparing the Fluidics Station

The Fluidics Station 400 is used to wash and stain the probe arrays. It is operated using Microarray Suite.

Setting Up the Fluidics Station

1. Turn on the Fluidics Station using the toggle switch on the lower left side of the machine.
2. Select **Run** → **Fluidics** from the menu bar.
 - ⇒ The Fluidics Station dialog box appears with a drop-down list for selecting the experiment name for each of the fluidics station modules. A second drop-down list is accessed for choosing the Protocol for each of the four fluidics station modules.

**Note**

Refer to the Fluidics Station 400 User's Guide for instructions on connecting and addressing multiple fluidics stations.

Priming the Fluidics Station

Priming ensures that the lines of the fluidics station are filled with the appropriate buffers and the fluidics station is ready for running fluidics station protocols.

Priming should be done:

- when the fluidics station is first started.
 - when wash solutions are changed.
 - before washing, if a shutdown has been performed.
 - if the LCD window instructs the user to prime.
1. To prime the fluidics station, select **Protocol** in the Fluidics Station dialog box.
 2. Choose **Prime** for the respective modules in the Protocol drop-down list.
 3. Change the intake buffer reservoir A to **Non-Stringent Wash Buffer** and intake buffer reservoir B to **Stringent Wash Buffer**.
 4. Click **Run** for each module to begin priming.

Probe Array Wash and Stain

Following hybridization, the wash and stain procedures are carried out by the fluidics station using the **ProkGE-WS2** fluidics script. The procedure takes approximately 90 minutes to complete. The use of streptavidin in the first part of the stain procedure enhances the overall signal.

Preparing the Staining Reagents

1. Prepare the following stain and wash solutions the day of the procedure. The solutions are stable for approximately 6 to 8 hours at 4°C. Volumes given are sufficient for one probe array.

Streptavidin Phycoerythrin (SAPE) should be stored in the dark at 4°C, either foil-wrapped or kept in an amber tube. Remove SAPE from refrigerator and tap the tube to mix well before preparing stain solution. Do not freeze SAPE. Always prepare the SAPE stain solution immediately before use.

Table 3.4.1
Streptavidin Solution Mix

Components	Volume	Final Concentration
2X MES Stain Buffer	300.0 μ L	1X
50 mg/mL BSA	24.0 μ L	2 mg/mL
1 mg/mL Streptavidin	6.0 μ L	10 μ g/mL
DI H ₂ O	270.0 μ L	—
Total Volume	600 μL	

Table 3.4.2
Antibody Solution Mix

Components	Volume	Final Concentration
2X MES Stain Buffer	300.0 μ L	1X
50 mg/mL BSA	24.0 μ L	2 mg/mL
10 mg/mL Normal Goat IgG	6.0 μ L	0.1 mg/mL
0.5 mg/mL Biotin Anti-streptavidin	6.0 μ L	5 μ g/mL
DI H ₂ O	264.0 μ L	—
Total Volume	600 μL	

Table 3.4.3
SAPE Solution Mix

Components	Volume	Final Concentration
2X MES Stain Buffer	300.0 μ L	1X
50 mg/mL BSA	24.0 μ L	2 mg/mL
1 mg/mL Streptavidin Phycoerythrin	6.0 μ L	10 μ g/mL
DI H ₂ O	270.0 μ L	—
Total Volume	600 μL	

Table 3.4.4
Fluidics Protocols - Antibody Amplification for *E. coli* Targets

Standard Format ProkGE-WS2	
Post Hyb Wash #1	10 cycles of 2 mixes/cycle with Wash Buffer A at 25°C
Post Hyb Wash #2	4 cycles of 15 mixes/cycle with Wash Buffer B at 45°C
Stain	Stain the probe array for 10 minutes in Streptavidin Solution. Mix at 25°C
Post Stain Wash	10 cycles of 4 mixes/cycle with Wash Buffer A at 30°C
2nd Stain	Stain the probe array for 10 minutes in antibody solution. Mix at 25°C
3rd Stain	Stain the probe array for 10 minutes in SAPE Solution at 25°C
Final Wash	15 cycles of 4 mixes/cycle with Wash Buffer A at 30°C. The holding temperature is 25°C

- Wash Buffer A = non-stringent wash buffer
- Wash Buffer B = stringent wash buffer

2. In the Fluidics Station dialog box on the workstation, select the correct experiment name from the drop-down **Experiment** list. The probe array type will appear automatically.
3. In the **Protocol** drop-down list, select the **ProkGE-WS2** protocol to control the washing and staining of the probe array format being used.

**Note**

Three-stain protocols require the user to replace stain solutions as directed by the LCD window during staining steps.

4. Choose **Run** in the Fluidics Station dialog box to begin the washing and staining. Follow the instructions in the LCD window on the fluidics station.
If you are unfamiliar with inserting and removing probe arrays from the fluidics station modules, please refer to the *Fluidics Station 400 User's Guide*, *Fluidics Station 400 Video In-Service CD* (P/N 900374), or *Quick Reference Card* (P/N 08-0072).
5. Insert the appropriate probe array into the designated module of the fluidics station while the probe array lever is in the **EJECT** position. When finished, verify that the probe array lever is returned to the **ENGAGE** position.
6. Remove any microcentrifuge tube remaining in the sample holder of the fluidics station module(s) being used.
7. Place a microcentrifuge tube containing 600 μ L streptavidin solution into the sample holder, making sure that the metal sampling needle is in the tube with its tip near the bottom.
⇒ The Fluidics Station dialog box and the LCD window display the status of the washing and staining as they progress.
8. When the LCD window indicates, replace the microcentrifuge tube containing the streptavidin stain with a microcentrifuge tube containing 600 μ L antibody stain solution into the sample holder, making sure that the metal sampling needle is in the tube with its tip near the bottom.
9. When the LCD window indicates, replace the microcentrifuge tube containing antibody solution with the microcentrifuge tube containing the 600 μ L of SAPE solution.
10. When the protocol is complete, the LCD window displays the message **EJECT CARTRIDGE**.
11. Remove microcentrifuge tube containing stain and replace with an empty microcentrifuge tube.
12. Remove the probe arrays from the fluidics station modules by first moving the probe array holder lever to the **EJECT** position.
13. Check the probe array window for large bubbles or air pockets.
 - If bubbles are present, refer to **Table 3.4.5** on page 3.4.12.
 - If the probe array has no large bubbles, it is ready to scan on the GeneArray Scanner. **ENGAGE** wash block and proceed to *Probe Array Scan* on page 3.4.13.

If you do not scan the arrays right away, keep the probe arrays at 4°C and in the dark until ready for scanning.

If there are no more samples to hybridize, shut down the fluidics station following the procedure outlined in the section, *Shutting Down the Fluidics Station* on page 3.4.14.



Note

For proper cleaning and maintenance of the fluidics station, including the bleach protocol, refer to Section 4, Fluidics Station Maintenance Procedures.

Table 3.4.5

If bubbles are present

Return the probe array to the probe array holder. Latch the probe array holder by gently pushing it up until a light click is heard. Engage the washblock by firmly pushing up on the probe array lever to the **ENGAGE** position. The fluidics station will drain the probe array and then fill it with a fresh volume of the last wash buffer used. When it is finished, if the LCD window displays **EJECT CARTRIDGE** again, remove the probe array and inspect it again for bubbles. If no bubbles are present, it is ready to scan. Proceed to *Probe Array Scan* on page 3.4.13. If several attempts to fill the probe array without bubbles are unsuccessful, the array should be filled with **Wash A (non-stringent buffer)** manually, using a micropipette. Excessive washing will result in a loss of signal intensity.

Probe Array Scan

The scanner is also controlled by Affymetrix Microarray Suite. The probe array is scanned after the wash protocols are complete. Make sure laser is warmed up prior to scanning by turning the laser on at least 15 minutes before use. If probe array was stored at 4°C, warm to room temperature before scanning. Refer to the Microarray Suite online help and the appropriate scanner user's manual for more information on scanning.

If necessary, clean the glass surface of probe array with a non-abrasive towel or tissue before scanning. **Do not use alcohol to clean glass.**



Note

The scanner uses an argon-ion laser and is equipped with a safety interlock system. Defeating the interlock system may result in exposure to hazardous laser light.

1. Select **Run** → **Scanner** from the menu bar. Alternatively, click the Start Scan icon in the tool bar.
 - ⇒ The Scanner dialog box appears with a drop-down list of experiments that have not been run.
2. Select the experiment name that corresponds to the probe array to be scanned. A previously run experiment can also be selected by using the **Include Scanned Experiments** option box. After selecting this option, previously scanned experiments appear in the drop-down list.
3. By default, after selecting the experiment the number [2] is displayed in the **Number of Scans** box to perform the recommended 2X image scan.
4. Once the experiment has been selected, click the **Start** button.
 - ⇒ A dialog box prompts you to load a sample into the scanner.
5. Click the **Options** button to check for the correct pixel value and wavelength of the laser beam.
 - Pixel value = 3 μm
 - Wavelength = 570 nm
6. Open the sample door on the scanner and insert the probe array into the holder. Do not force the probe array into the holder. Close the sample door of the scanner.
7. Click **OK** in the Start Scanner dialog box.
 - ⇒ The scanner begins scanning the probe array and acquiring data. When **Scan in Progress** is Selected from the **View** menu, the probe array image appears on the screen as the scan progresses.

Shutting Down the Fluidics Station

1. After removing a probe array from the probe array holder, the LCD window displays the message **ENGAGE WASHBLOCK**.
2. Engage the washblock by firmly pushing up on the probe array lever to the **ENGAGE** position.
The fluidics station will automatically perform a Cleanout procedure. The LCD window will indicate the progress of the Cleanout procedure.
3. When the fluidics station LCD window indicates **REMOVE VIAL**, the Cleanout procedure is complete.
4. Remove the sample microcentrifuge tube from the sample holder.
5. If no other hybridizations are to be performed, place wash lines into a bottle filled with deionized water.
6. Select **Shutdown** for all modules from the drop-down **Protocol** list in the Fluidics Station dialog box. Click the **Run** button for all modules.
The Shutdown protocol is critical to instrument reliability. Refer to the *Fluidics Station 400 User's Guide* for more information.
7. After Shutdown protocol is complete, flip the ON/OFF switch of the fluidics station to the OFF position.

▶ IMPORTANT

To maintain the cleanliness of the fluidics station and obtain the highest quality image and data possible, a weekly bleach protocol and a monthly decontamination protocol are highly recommended. Please refer to Section 4, Fluidics Station Maintenance Procedures for further detail.

Customizing the Protocol

There may be times when the fluidics protocols need to be modified. Modification of protocols must be done before downloading the protocol to the fluidics station. Protocol changes will not affect runs in progress. For more specific instructions, refer to the Microarray Suite online help.

1. Select **Tools** → **Edit Protocol** from the menu bar.
2. In the Edit Protocol dialog box under Protocol Name, click the arrow to open a list of protocols. Click the protocol to be changed.
 - ⇒ The name of the protocol is displayed in the **Protocol Name** text box. The conditions for that protocol are displayed on the right side of the Edit Protocol dialog box.
3. Select the item to be changed and input the new parameters as needed, keeping parameters within the ranges shown below in **Table 3.4.6**.

Table 3.4.6
Valid Ranges for Wash/Stain Parameters

Parameter	Valid Range
Wash Temperature for A1, B, A2, or A3 (°C)	15 to 50
Number of Wash Cycles for A1, B, A2, or A3	0 to 99
Mixes / Wash cycle for A1, B, A2, or A3	15 to 50
Stain Time (seconds)	0 to 86,399
Stain Temperature (°C)	15 to 50
Holding Temperature (°C)	15 to 50

- Wash A1 corresponds to Post Hyb wash #1 in Table 3.4.4.
- Wash B corresponds to Post Hyb wash #2 in Table 3.4.4.
- Wash A2 corresponds to Post Stain Wash in Table 3.4.4.
- Wash A3 corresponds to Final Wash in Table 3.4.4.

4. To return to the default values for the protocol selected, click the **Defaults** button.
5. Once all the protocol conditions are modified as desired, change the name of the edited protocol in the **Protocol Name** box.

! CAUTION

*If the protocol is saved without entering a new **Protocol Name**, the original protocol parameters will be overwritten.*

6. Click **Save**, then close the dialog box.

Enter **0** (zero) for hybridization time if hybridization step is not required. Likewise, enter **0** (zero) for the stain time if staining is not required. Enter **0** (zero) for the number of wash cycles if a wash solution is not required.

