Appendix B2

Protocol for CHO Cells + hAR + Luciferase Assay

(Provided by Dr. Anne Marie Vinggaard, Institute of Food Safety and Toxicology, Danish Veterinary and Food Administration, Soborg, Denmark) [This page intentionally left blank]

Protocol for CHO Cells + hAR + Luciferase Reporter Gene Assay

Passage of cells

- **Day 2** 9a.m.: Cells are approx. 90 % confluent. Cells are transfected with cDNA and FuGene 6 (Roche) according to the table. Cells are incubated for 5 hours (14 a.m.) Turn the plates on paper towel and add media +/- hormones and chemicals (p.3). Cells are incubated for 20 hours (11a.m)
- **Day 3** Remove media by turning the plates on paper towel.

MgCl₂ is added to lysis buffer. All wells are added 20 μ l lysis buffer and incubated for 15 min on a shaker. Prepare luciferin solution containing (2 ml luciferin/ATP is added 2 ml lysis buffer with Mg Cl₂). Protect from light. Measure luciferase activity on the BioOrbit Galaxy Luminometer directly in the plates by injection of 40 μ l luciferin solution per well. The chemiluminiscense generated from each well is measured over a 1 sec interval after an incubation time of 2 sec.

Transfection scheme

The expression vector pSVAR0, AR13 and the MMTV-LUC reporter plasmid were both provided by Dr. Albert Brinkmann, Erasmus University, Rotterdam.

Optimum AR reporter gene assay conditions:

For 200 wells:	
DMEM/F12	: 940 µl
Fugene	$: 60 \mu l$ $: 15 \mu g$ $5 \mu l / well$
DNA (totally)	: 15 µg ^J

 $\begin{array}{ll} \underline{75 \text{ ng DNA per well}}\\ \overline{FuGene (ul) / DMEM} + \overline{FuGene (ul)} = 0.06\\ DNA \mu g / \overline{Fugene \mu l} = 0.25\\ psvAR0, MMTV-Luc = 1:100\\ AR13, MMTV-Luc = 2:100 \text{ (for cytotoxicity determination)} \end{array}$

	µg DNA	Plate 1 - 6	Plate x
			(cytotoxicity)
DMEM/F12 without serum		3290 µl	600 µl
FuGene		210 µl	38.3 µl
psvAR0 0.1 µg/µl (batch no.)	0.53 µg	5.3 µl	AR13(0.082μg/μl) 3.11 μl
MMTV-Luc 2.252 μg/μl (batch no.)	51.97 μg	23.1 µl	6.087 μl

DMEM/F12 without serum is added to a 15 ml plastic vial. FuGene is added without touching the walls of the vial. Gently mix and incubate for 5 min at room temperature. cDNA is added to another 15 ml vial. The diluted FuGene solution is added drop-wise to the cDNA. Gently mixing. The solution incubates for 15 min at room temperature. $8*250 \ \mu$ l is added to a column in a microtiter plate. 5 μ l is added to the each well containing the cells using a 100 μ l 8-channel pipette. Be sure that the cDNA is distributed well (DNA should lie as grain of sand in the media instantly or the day after).

Overview of plates

Don't use row A and H

50 μ l compound + 50 μ l media is added according to the scheme

Plate 1 Row B: Row C: Row D: Row E: Row F: Row G:	Final conc.: Compound x 0.01 nM R1881 0.025, 0.05, 0.10 0.20, 0.39, 0.78 1.56, 3.13; 6.25 12.5, 25, 50 0.01 nM R1881	0.02 nM R1881 x 12 wells + media uM x 4 wells + 0.02 nM R1881
Plate 2 Row B: Row C: Row D: Row E: Row F: Row G:	0; 0.001; 0.0023 nM 0.01; 0.023; 0,1 nM 0,23; 1,0; 2,3 nM 0; 1; 5 nM 10; 50; 100 nM 500; 1000; 5000nM	R1881 x 4 wells + media R1881 x 4 wells + media R1881 x 4 wells + media OHF x 4 wells + 0.02 nM R1881 OHF x 4 wells + 0.02 nM R1881 OHF x 4 wells + 0.02 nM R1881

Solvents

Hydroxyflutamide and R1881: Stock solutions in freezer no.

Positive antagonism control: Hydroxyflutamide (2 x conc.)

20 nM: $5 \ \mu l \ 10 \ \mu M + 2.5 \ ml media$ 10 nM: $5 \ \mu l \ 5 \ \mu M + 2.5 \ ml media$ 2 nM: $5 \ \mu l \ 1 \ \mu M + 2.5 \ ml media$	
---	--

Positive agonism control: R1881 (2 x conc.)

20 nM:	$5 \mu l 10 \mu M + 2.5 ml media - (not to be added)$
4.6 nM:	500 ul 20 nM + 1.665 ml media
2 nM:	5 μ l 1 μ M + 2.5 ml media
0.46 nM	500 ul 2 nM + 1.665 ml media
0.2 nM:	5 μ l 0.1 μ M + 2.5 ml media
0.02 nM:	40 µl 0.01 µM + 20 ml media
0.046 nM	500 ul 0.2 nM + 1.665 ml media
0.002 nM:	5 μ l 0.001 μ M + 2.5 ml media
0.0046 nM:	500 ul 0.02 nM + 1.665 ml media
0 nM:	$5 \mu l EtOH + 2.5 ml media.$
	-

Compound x	(2 x conc.)
100 µM:	$25 \ \mu l \ 20 \ mM + 5 \ ml \ media$
50.0 µM:	$1 \text{ ml } 100 \mu\text{M} + 1 \text{ ml media}$
25.0 µM:	1 ml 50 μ M + 1 ml media
12.5 µM:	1 ml 25 μ M + 1 ml media
6.25 µM:	$1 \text{ ml } 12.5 \mu\text{M} + 1 \text{ ml media}$
3.13 µM:	$1 \text{ ml } 6.25 \mu\text{M} + 1 \text{ ml media}$
1.56 µM:	$1 \text{ ml} 3.13 \mu\text{M} + 1 \text{ ml} \text{ media}$
0.78 µM:	$1 \text{ ml } 1.56 \mu\text{M} + 1 \text{ ml media}$
0.39 µM:	$1 \text{ ml } 0.78 \mu\text{M} + 1 \text{ ml media}$
0.20 µM:	$1 \text{ ml } 0.39 \mu\text{M} + 1 \text{ ml media}$
0.10 µM:	$1 \text{ ml } 0.20 \mu\text{M} + 1 \text{ ml media}$
0.05 µM:	$1 \text{ ml } 0.10 \text{\mu}\text{M} + 1 \text{ ml media}$

Media, buffers, compounds etc.

Lysis buffer:	25 mM Trisphosphate pH 7.8 (adjusted with phosphoric acid) 15 % glycerol 1 % Triton X-100 1 mM DTT Stored in sterile 50 ml vials at -20°C Before use: add 8 mM MgCl ₂ (8 μl 1M MgCl ₂ /ml buffer)	
Requirement:	x * 10 ml lysis buffer + x * 4.5 ml luciferin/ATP (Freezer no.)	
Requirement of media:	DMEM/F12 + 10 % DCC x ml + x ml for counting DMEM/F12 + 10 % FBS for a 25cm ² flask DMEM/F12 + 1% PSF	xx ml xx ml xx ml
Compounds:	Comp.1 Mw. xx g/mol Supplier:Lot no.: Purity: Stocksolution of 20 mM (x mg to x ml EtOH) Date: Person: Remarks:	

Ethanol: Merck pro analysis UN 1170, K 27773283-020

[This page intentionally left blank]