Attachment F

Region I Short Sheets and CLP Information Sheets

For hardcopy of Attachment F contact:

Steve Stodola, U.S. EPA Region I TEL: 781-860-4634 EMAIL: stodola.steve@epamail.epa.gov

For electronic copies of related documents see:

http://www.epa.gov/superfund/oerr/aoc/

(Click on CLP Products and Services)

USEPA CONTRACT LABORATORY PROGRAM SHORT SHEETS

AVAILABLE SHORT SHEETS

TITLE	DOCUMENT OR NUMBER		SHORT SHEET REVISION NUMBER
USEPA Contract Laboratory Program Statement of Work for Organic Analysis Multi-Media, Multi-Concentration	2/88		1.0
USEPA Contract Laboratory Program Statement of Work for Organic Analysis Multi-Media, Multi-Concentration	OLM01.0	1.0	
USEPA Contract Laboratory Program Statement of Work for Organic Analysis Multi-Media, Multi-Concentration	OLM01.9	3.0	
USEPA Contract Laboratory Program Statement of Work for Organic Analysis Multi-Media, Multi-Concentration	OLM02.1	2.0	
USEPA Contract Laboratory Program Statement of Work for Organic Analysis Multi-Media, High-Concentration	9/88		1.0
USEPA Contract Laboratory Program Statement of Work for Inorganic Analysis Multi-Media, Multi-Concentration	7/88		1.0
USEPA Contract Laboratory Program Statement of Work for Inorganic Analysis Multi-Media, Multi-Concentration	ILM01.0	1.0	
USEPA Contract Laboratory Program Statement of Work for Inorganic Analysis Multi-Media, Multi-Concentration	ILM02.1	2.0	
USEPA Contract Laboratory Program Statement of Work for Inorganic Analysis Multi-Media, Multi-Concentration	ILM03.0	2.0	
USEPA Contract Laboratory Program Statement of Work for Inorganic Analysis Multi-Media, High-Concentration	IHC01.2	1.0	

AVAILABLE SHORT SHEETS

TITLE	DOCUMENT I OR NUMBER		SHORT SHEET REVISION NUMBER
USEPA Contract Laboratory Program Statement of Work for Analysis of Polychlorinated Dibenzo-p-Dioxins (PCDD) and PolychlorinatedDibenzofurans (PCDF) Multi-Media, Multi-Concentration	DFLM01.2		3.0
Superfund Analytical Methods for Low Concentration Water for Organics Analysis	10/92		2.0
USEPA Contract Laboratory Program Water Quality Parameters in Multi-Concentration Water (WQP)	6/93		2.0
USEPA Contract Laboratory Program Statement of Work for Inorganic Analysis Multi-Media, Multi-Concentration	ILM04.0	0.0	
USEPA Contract Laboratory Program Statement of Work for Organic Analysis Multi-Media, Multi-Concentration	OLM03.2	0.0	

SOW SHORT SHEET GUIDELINES

The following is a description of the categories included in the Statement of Work (SOW) Short Sheets.

- ! Title: Name of the document as it appears on the cover page of the document.
- ! Document Number: Number of the document as it appears on the cover of the document. Not all SOWs have associated document numbers. A SOW may be referenced by either the document date or the document number, depending on which one is applicable.
- ! Document Date: Date the document was initially issued according to the cover page of the document. Not all SOWs have associated document dates. A SOW may be referenced by either the document date or the document number, depending on which one is applicable.
- ! Effective Dates: Range of dates for which the Regions can submit samples for analysis under a particular SOW. If contracts have not been awarded yet, then the expected award date is listed.
- ! Concentration: Range of sample concentrations for which the SOW is applicable, such as low to medium, high, or > 20 mg/kg.
- ! Data Turnaround: Number of days the laboratory has to submit the complete data package after sample receipt.
- ! Matrices: Sample matrices for which the SOW is applicable, such as aqueous, soil, sediment, multi-phase, etc.
- ! Significant Features: Information about the SOW which distinguishes it from other SOWs. This section highlights critical items such as holding times, concentrations and matrices which may be different.
- ! Revisions/Modifications: Revisions from the previous SOW which may significantly affect data useability.
- ! Recommended Uses: Explanation of appropriate Superfund activities for which the SOW may be utilized.
- ! Analytes/CRQLs: References Attachment I which lists the parameters included in the analysis and their respective CRQLs or CRDLs. The aqueous and soil CRQLs and CRDLs are listed. "Notes" are provided at the bottom of each Attachment to document deviations from the CRQLs and CRDLs listed.

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS

MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	Not Applicable
DOCUMENT DATE:	February 1988
EFFECTIVE DATES:	January 20, 1989 through September 10, 1991
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	35 Days
MATRICES:	Aqueous/Soil/Sediment*

SIGNIFICANT FEATURES

- ! The parameters include volatile, semivolatile, and pesticide/PCB compounds.
- ! Volatiles and semivolatiles are analyzed by GC/MS; pesticides/PCBs are analyzed by GC/ECD.
- ! Major Tentatively Identified Compounds (TICs) are reported for GC/MS analyses.
- ! Second column confirmation by GC/ECD is required for all pesticides/PCBs. Pesticides/PCBs which are identific concentrations above 10 ng/uL are confirmed by GC/MS analysis.

REVISIONS/MODIFICATIONS

The 9/88 and 4/89 revisions to the 2/88 SOW do not significantly affect data useability.

RECOMMENDED USES

This Routine Analytical Services (RAS) method is recommended for broad spectrum analysis to define the na of potential site contamination during SSI, LSI, and RI/FS activities. This method is suitable when a thirty five day to results is adequate. It is recommended for samples from known or suspected hazardous waste sites where potential co be present at significant risk levels.

Sediment samples with high moisture content should be solicited as RAS + SAS (Special Analytical Service) achieve the CRQLs.

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in .

ATTACHMENT I (page 1 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION FEBRUARY, 1988

TARGET COMPOUND LIST - VOLATILES

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Chloromethane	10	10
Bromomethane	10	10
Vinyl Chloride	10	10
Chloroethane	10	10
Methylene Chloride	5	5
Acetone	10	10
Carbon Disulfide	5	5
1,1-Dichloroethene	5	5
1,1-Dichloroethane	5	5
1,2-Dichloroethene(total)	5	5
Chloroform	5	5
1,2-Dichloroethane	5	5
2-Butanone	10	10
1,1,1-Trichloroethane	5	5
Carbon Tetrachloride	5	5
Vinyl Acetate	10	10
Bromodichloromethane	5	5

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
1,2-Dichloropropane	5	5
cis-1,3-Dichloropropene	5	5
Trichloroethene	5	5
Dibromochloromethane	5	5
1,1,2-Trichloroethane	5	5
Benzene	5	5
trans-1,3-Dichloropropene	5	5
Bromoform	5	5
4-Methyl-2-pentanone	10	10
2-Hexanone	10	10
Tetrachloroethene	5	5
Toluene	5	5
1,1,2,2-Tetrachloroethene	5	5
Chlorobenzene	5	5
Ethyl Benzene	5	5
Styrene	5	5
Total Xylenes	5	5

- ! THE SAMPLE-SPECIFIC CRQLs FOR SOIL SAMPLES WILL BE ADJUSTED FOR PERCENT MOISTURE AND WILL BE HIGHEI ABOVE.
- ! MEDIUM LEVEL SOIL CRQL = 125 x AQUEOUS CRQL REPORTED IN UG/KG.

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION FEBRUARY, 1988

TARGET COMPOUND LIST - SEMIVOLATILES

	TARGET	COMPOUND L.
Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Phenol	10	330
bis(2-Chloroethyl)ether	10	330
2-Chlorophenol	10	330
1,3-Dichlorobenzene	10	330
1,4-Dichlorobenzene	10	330
Benzyl alcohol	10	330
1,2-Dichlorobenzene	10	330
2-Methylphenol	10	330
bis(2-Chloroisopropyl) ether	10	330
4-Methylphenol	10	330
N-Nitroso-di-n-dipropylamine	10	330
Hexachloroethane	10	330
Nitrobenzene	10	330
Isophorone	10	330
2-Nitrophenol	10	330
2,4-Dimethylphenol	10	330
Benzoic acid	50	1600
bis(2-Chloroethoxy) methane	10	330
2,4-Dichlorophenol	10	330
1,2,4-Trichlorobenzene	10	330
Naphthalene	10	330
4-Chloroaniline	10	330
Hexachlorobutadiene	10	330
4-Chloro-3-methylphenol (para-chloro-meta-cresol)	10	330
2-Methylnaphthalene	10	330
Hexachlorocyclopentadiene	10	330
2,4,6-Trichlorophenol	10	330
2,4,5-Trichlorophenol	50	1600
2-Chloronaphthalene	10	330
2-Nitroaniline	50	1600
Dimethylphthalate	10	330
Acenaphthylene	10	330
2,6-Dinitrotoluene	10	330

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
3-Nitroaniline	50	1600
Acenaphthene	10	330
2,4-Dinitrophenol	50	1600
4-Nitrophenol	50	1600
Dibenzofuran	10	330
2,4-Dinitrotoluene	10	330
Diethylphthalate	10	330
4-Chlorophenyl-phenyl ether	10	330
Fluorene	10	330
4-Nitroaniline	50	1600
4,6-Dinitro-2-methylphenol	50	1600
N-nitrosodiphenylamine	10	330
4-Bromophenyl-phenylether	10	330
Hexachlorobenzene	10	330
Pentachlorophenol	50	1600
Phenanthrene	10	330
Anthracene	10	330
Di-n-butylphthalate	10	330
Fluoranthene	10	330
Pyrene	10	330
Butylbenzylphthalate	10	330
3,3'-Dichlorobenzidine	20	660
Benzo(a)anthracene	10	330
Chrysene	10	330
bis(2-Ethylhexyl)phthalate	10	330
Di-n-octylphthalate	10	330
Benzo(b)fluoranthene	10	330
Benzo(k)fluoranthene	10	330
Benzo(a)pyrene	10	330
Indeno(1,2,3-cd)pyrene	10	330
Dibenz(a,h)anthracene	10	330
Benzo(g,h,i)perylene	10	330

NOTE:

HIGHER TH

[!] THE SAMPLE-SPECIFIC CRQLs FOR SOIL SAMPLES WILL BE ADJUSTED FOR PERCENT MOISTURE AND WILL BE

! MEDIUM LEVEL SOIL CRQL = 1980 x AQUEOUS CRQL REPORTED IN UG/KG. ATTACHMENT I (page 3 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION FEBRUARY, 1988

TARGET COMPOUND LIST - PESTICIDES/PCBs

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
alpha-BHC	0.05	8.0
beta-BHC	0.05	8.0
delta-BHC	0.05	8.0
gamma-BHC(Lindane)	0.05	8.0
Heptachlor	0.05	8.0
Aldrin	0.05	8.0
Heptachlor epoxide	0.05	8.0
Endosulfan I	0.05	8.0
Dieldrin	0.10	16.0
4,4'-DDE	0.10	16.0
Endrin	0.10	16.0
Endosulfan II	0.10	16.0
4,4'-DDD	0.10	16.0
Endosulfan sulfate	0.10	16.0

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
4,4'-DDT	0.10	16.0
Methoxychlor	0.5	80.0
Endrin ketone	0.10	16.0
alpha-Chlordane	0.5	80.0
gamma-Chlordane	0.5	80.0
Toxaphene	1.0	160.0
Aroclor-1016	0.5	80.0
Aroclor-1221	0.5	80.0
Aroclor-1232	0.5	80.0
Aroclor-1242	0.5	80.0
Aroclor-1248	0.5	80.0
Aroclor-1254	1.0	160.0
Aroclor-1260	1.0	160.0

- ! THE SAMPLE-SPECIFIC CRQLs FOR SOIL SAMPLES WILL BE ADJUSTED FOR PERCENT MOISTURE AND WILL BE HIGHER 1 ABOVE.
- ! MEDIUM LEVEL SOIL CRQL = 15 x LOW LEVEL SOIL CRQL REPORTED IN UG/KG.

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS

MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	OLM01.0
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	September 28, 1990 through February 1994
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	14 Days or 35 Days
MATRICES:	Aqueous/Soil/Sediment *

SIGNIFICANT FEATURES

- ! The parameters include volatile, semivolatile, and pesticide/PCB compounds.
- ! Volatiles and semivolatiles are analyzed by GC/MS, pesticides/PCBs are analyzed by GC/ECD.
- ! Major Tentatively Identified Compounds (TICs) are reported for GC/MS analyses.
- ! Second column confirmation by GC/ECD is required for all pesticides/PCBs. Pesticides/PCBs which are identified by GC/EC above 10 ng/uL are confirmed by GC/MS analysis.

REVISIONS/MODIFICATIONS

The following is a list of the significant changes from the 2/88 SOW that are incorporated in the OLM01.0 SOW:

- ! Selected volatile CRQLs have been raised; pesticide/PCB low soil CRQLs have been lowered; and selected pesticide/PCB aqu changed.
- ! Target Compound List (TCL) changes include the elimination of vinyl acetate from the volatile TCL; the elimination of benzy acid from the semivolatile TCL; the addition of carbazole to the semivolatile TCL; and the addition of endrin aldehyde to the p semivolatile TCL compound bis(2-chloroisopropyl)ether was renamed 2,2'-oxybis(1-chloropropane).
- ! A new method for analysis of pesticides/PCBs is used. Changes include the use of wide bore capillary columns; new surrogate techniques.
- ! Pesticide/PCB quantitation is performed using both the primary and secondary columns. The lower value is reported by the lal

The only significant change in the OLM01.1 (December, 1990) and OLM01.1.1 (February, 1991) revisions to the OLM01.0 lowering of selected semivolatile CRQLs. The significant changes in the OLM01.1 through OLM01.7 revisions to the OLM01.0 SC selected semivolatile CRQLs and options for either a 14 day or 35 day data turnaround.

RECOMMENDED USES

This Routine Analytical Services (RAS) method is recommended for broad spectrum analysis to define the nature and extent contamination during SSI, LSI, and RI/FS activities. This method is suitable when a fourteen day or thirty five day turnaround for r recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at signific

* Sediment samples with high moisture content should be solicited as RAS + SAS (Special Analytical Service) in order to acl

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in Attachment 1

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM01.0

TARGET COMPOUND LIST - VOLATILES

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Chloromethane	10	10
Bromomethane	10	10
Vinyl Chloride	10	10
Chloroethane	10	10
Methylene Chloride	10*	10*
Acetone	10	10
Carbon Disulfide	10*	10*
1,1-Dichloroethene	10*	10*
1,1-Dichloroethane	10*	10*
1,2-Dichloroethene (total)	10*	10*
Chloroform	10*	10*
1,2-Dichloroethane	10*	10*
2-Butanone	10	10
1,1,1-Trichloroethane	10*	10*
Carbon Tetrachloride	10*	10*
Bromodichloromethane	10*	10*
1,2-Dichloropropane	10*	10*

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
cis-1,3-Dichloropropene	10*	10*
Trichloroethene	10*	10*
Dibromochloromethane	10*	10*
1,1,2-Trichloroethane	10*	10*
Benzene	10*	10*
trans-1,3-Dichloropropene	10*	10*
Bromoform	10*	10*
4-Methyl-2-pentanone	10	10
2-Hexanone	10	10
Tetrachloroethene	10*	10*
Toluene	10*	10*
1,1,2,2-Tetrachloroethene	10*	10*
Chlorobenzene	10*	10*
Ethyl Benzene	10*	10*
Styrene	10*	10*
Xylenes (total)	10*	10*

- ! THE SAMPLE-SPECIFIC CRQLs FOR SOIL SAMPLES WILL BE ADJUSTED FOR PERCENT MOISTURE AND WILL BE HIGHER TH
- ! MEDIUM LEVEL SOIL CRQL = 120 x AQUEOUS CRQL REPORTED IN UG/KG.

^{*} CRQLs previously 5 ug/L and 5 ug/Kg in 2/88 SOW.

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM01.0

TARGET COMPOUND LIST - SEMIVOLATILES

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Phenol	10	330
bis(2-Chloroethyl)ether	10	330
2-Chlorophenol	10	330
1,3-Dichlorobenzene	10	330
1,4-Dichlorobenzene	10	330
1,2-Dichlorobenzene	10	330
2-Methylphenol	10	330
2,2'-oxybis(1-Chloropropane)	10	330
4-Methylphenol	10	330
N-Nitroso-di-n-dipropylamine	10	330
Hexachloroethane	10	330
Nitrobenzene	10	330
Isophorone	10	330
2-Nitrophenol	10	330
2,4-Dimethylphenol	10	330
bis(2-Chloroethoxy) methane	10	330
2,4-Dichlorophenol	10	330
1,2,4-Trichlorobenzene	10	330
Naphthalene	10	330
4-Chloroaniline	10	330
Hexachlorobutadiene	10	330
4-Chloro-3-methylphenol	10	330
2-Methylnaphthalene	10	330
Hexachlorocyclopentadiene	10	330
2,4,6-Trichlorophenol	10	330
2,4,5-Trichlorophenol	25*	800*
2-Chloronaphthalene	10	330
2-Nitroaniline	25*	800*
Dimethylphthalate	10	330
Acenaphthylene	10	330
2,6-Dinitrotoluene	10	330
3-Nitroaniline * CRQLs previously 50 ug/L and 1600	25*	800* W ** CRQLs p

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Acenaphthene	10	330
2,4-Dinitrophenol	25*	800*
4-Nitrophenol	25*	800*
Dibenzofuran	10	330
2,4-Dinitrotoluene	10	330
Diethylphthalate	10	330
4-Chlorophenyl-phenyl ether	10	330
Fluorene	10	330
4-Nitroaniline	25*	800*
4,6-Dinitro-2-methylphenol	25*	800*
N-nitrosodiphenylamine	10	330
4-Bromophenyl-phenylether	10	330
Hexachlorobenzene	10	330
Pentachlorophenol	25*	800*
Phenanthrene	10	330
Anthracene	10	330
Carbazole	10	330
Di-n-butylphthalate	10	330
Fluoranthene	10	330
Pyrene	10	330
Butylbenzylphthalate	10	330
3,3'-Dichlorobenzidine	10**	330**
Benzo(a)anthracene	10	330
Chrysene	10	330
bis(2-Ethylhexyl)phthalate	10	330
Di-n-octylphthalate	10	330
Benzo(b)fluoranthene	10	330
Benzo(k)fluoranthene	10	330
Benzo(a)pyrene	10	330
Indeno(1,2,3-cd)pyrene	10	330
Dibenz(a,h)anthracene	10	330
Benzo(g,h,i)perylene	10	330

^{*} CRQLs previously 50 ug/L and 1600 ug/Kg in 2/88 SOW ** CRQLs previously 20 ug/L and 660 ug/Kg in 2/88 SOW.

[!] THE SAMPLE-SPECIFIC CRQLs FOR SOIL SAMPLES WILL BE ADJUSTED FOR PERCENT MOISTURE AND WILL BE HIGHER THAN THOSE LISTED ABOVE.

! MEDIUM LEVEL SOIL CRQL = 1000 x AQUEOUS CRQL REPORTED IN UG/KG.
ATTACHMENT I (page 3 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM01.0

TARGET COMPOUND LIST - PESTICIDES/PCBs

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil** CRQL (ug/Kg,ppb)
alpha-BHC	0.05	1.7
beta-BHC	0.05	1.7
delta-BHC	0.05	1.7
gamma-BHC (Lindane)	0.05	1.7
Heptachlor	0.05	1.7
Aldrin	0.05	1.7
Heptachlor epoxide	0.05	1.7
Endosulfan I	0.05	1.7
Dieldrin	0.10	3.3
4,4'-DDE	0.10	3.3
Endrin	0.10	3.3
Endosulfan II	0.10	3.3
4,4'-DDD	0.10	3.3
Endosulfan sulfate	0.10	3.3

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil** CRQL (ug/Kg,ppb)
4,4'-DDT	0.10	3.3
Methoxychlor	0.5	17.0
Endrin ketone	0.10	3.3
Endrin aldehyde	0.10	3.3
alpha-Chlordane	0.05*	1.7
gamma-Chlordane	0.05*	1.7
Toxaphene	5.0*	170.0
Aroclor-1016	1.0*	33.0
Aroclor-1221	2.0*	67.0
Aroclor-1232	1.0*	33.0
Aroclor-1242	1.0*	33.0
Aroclor-1248	1.0*	33.0
Aroclor-1254	1.0	33.0
Aroclor-1260	1.0	33.0

Aqueous CRQLs changed from 2/88 SOW to the following:

* Aqueous CRQLs (ug/L) - alpha- and gamma- Chlordane from 0.5 to 0.05;

Toxaphene from 1.0 to 5.0;

Aroclors-1016, 1232, 1242, and 1248 from 0.5 to 1.0;

Aroclor-1221 from 0.5 to 2.0.

All low soil CRQLs changed from 2/88 SOW to the following:

** Low Soil CRQLs (ug/Kg) -alpha-BHC through Endosulfan I from 8.0 to 1.7;

Dieldrin through 4,4'-DDT and Endrin ketone from 16.0 to 3.3;

Methoxychlor from 80.0 to 17.0;

alpha- and gamma-Chlordane from 80.0 to 1.7;

Toxaphene from 160.0 to 170.0;

Aroclor-1016, 1232, 1242, and 1248 from 80.0 to 33.0;

Aroclor-1221 from 80.0 to 67.0;

Aroclor-1254 and 1260 from 160.0 to 33.0.

NOTE:

! THE SAMPLE-SPECIFIC CRQLs FOR SOIL SAMPLES WILL BE ADJUSTED FOR PERCENT MOISTURE AND WILL BE HIGHER ABOVE.

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	OLM01.9
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	July 1993 through February 1995
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	14 Days or 35 Days
MATRICES:	Aqueous/Soil/Sediment*

SIGNIFICANT FEATURES

- ! The parameters include volatile, semivolatile, and pesticide/PCB compounds.
- ! Volatiles and semivolatiles are analyzed by GC/MS, pesticides/PCBs are analyzed by GC/ECD.
- ! Major Tentatively Identified Compounds (TICs) are reported for GC/MS analyses.
- ! Second column confirmation by GC/ECD is required for all pesticides/PCBs. The lower of the two concentrations detected on both columns is reported. Pesticides/PCBs which are detected at concentrations above 10 ng/uL are confirmed by GC/M

REVISIONS/MODIFICATIONS

The only significant revisions to the OLM01.0 SOW in the OLM01.8 revision were changes in the format and content of the Agency Standard diskette deliverable.

The following is a list of the significant changes from the OLM01.8 revision that are incorporated in the OLM01.9 rev

- ! MS/MSD analysis is not required for SDGs containing only equipment/trip blanks or PE samples.
- ! Specific instructions are given regarding resolution of problems with reduced sample volume and MS/MSD sample

RECOMMENDED USES

This Contract Laboratory Program (CLP) method is recommended for broad spectrum analysis to define the nature and extent of potential site contamination during Site Investigation (SI) and Remedial Investigation/Feasibility Study (RI/FS) activities and to verify that Remedial Design/Remedial Action (RD/RA) activities comply with pre-determined clean-up standards. This method is suitable when a 14 day or 35 day tumaround for results is adequate. It is recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at significant risk levels.

* This method is not applicable to sediment samples with high moisture content.

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in Attacl

ATTACHMENT 1 (page 1 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM01.9

TARGET COMPOUND LIST - VOLATILES

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Chloromethane	10	10
Bromomethane	10	10
Vinyl chloride	10	10
Chloroethane	10	10
Methylene chloride	10	10
Acetone	10	10
Carbon disulfide	10	10
1,1-Dichloroethene	10	10
1,1-Dichloroethane	10	10
1,2-Dichloroethene (total)	10	10
Chloroform	10	10
1,2-Dichloroethane	10	10
2-Butanone	10	10
1,1,1-Trichloroethane	10	10
Carbon tetrachloride	10	10
Bromodichloromethane	10	10
1,2-Dichloropropane	10	10

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg/ppb)
cis-1,3-Dichloropropene	10	10
Trichloroethene	10	10
Dibromochloromethane	10	10
1,1,2-Trichloroethane	10	10
Benzene	10	10
trans-1,3-Dichloropropene	10	10
Bromoform	10	10
4-Methyl-2-pentanone	10	10
2-Hexanone	10	10
Tetrachloroethene	10	10
1,1,2,2-Tetrachloroethane	10	10
Toluene	10	10
Chlorobenzene	10	10
Ethylbenzene	10	10
Styrene	10	10
Xylenes (Total)	10	10

- ! The sample-specific CRQLs for soil samples will be adjusted for percent moisture and will be higher than those listed above.
- ! Medium level soil CRQL = 120 x aqueous CRQL reported in ug/kg.

ATTACHMENT I (page 2 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM01.9

TARGET COMPOUND LIST - SEMIVOLATILES

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Phenol	10	330
bis(2-chloroethyl)ether	10	330
2-Chlorophenol	10	330
1,3-Dichlorobenzene	10	330
1,4-Dichlorobenzene	10	330
1,2-Dichlorobenzene	10	330
2-Methylphenol	10	330
2,2'-oxybis(1-chloropropane)	10	330
4-Methylphenol	10	330
N-Nitroso-di-n-propylamine	10	330
Hexachloroethane	10	330
Nitrobenzene	10	330
Isophorone	10	330
2-Nitrophenol	10	330
2,4-Dimethylphenol	10	330
bis-(2-Chloroethoxy)methane	10	330
2,4-Dichlorophenol	10	330
1,2,4-Trichlorobenzene	10	330
Naphthalene	10	330
4-Chloroaniline	10	330
Hexachlorobutadiene	10	330
4-Chloro-3-methylphenol	10	330
2-Methylnaphthalene	10	330
Hexachlorocyclopentadiene	10	330
2,4,6-Trichlorophenol	10	330
2,4,5-Trichlorophenol	25	800
2-Chloronaphthalene	10	330
2-Nitroaniline	25	800
Dimethylphthalate	10	330
Acenaphthylene	10	330
2,6-Dinitrotoluene	10	330
3-Nitroaniline	25	800

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Acenaphthene	10	330
2,4-Dinitrophenol	25	800
4-Nitrophenol	25	800
Dibenzofuran	10	330
2,4-Dinitrotoluene	10	330
Diethylphthalate	10	330
4-Chlorophenyl-phenylether	10	330
Fluorene	10	330
4-Nitroaniline	25	800
4,6-Dinitro-2-methylphenol	25	800
N-Nitrosodiphenylamine	10	330
4-Bromophenyl-phenylether	10	330
Hexachlorobenzene	10	330
Pentachlorophenol	25	800
Phenanthrene	10	330
Anthracene	10	330
Carbazole	10	330
Di-n-butylphthalate	10	330
Fluoranthene	10	330
Pyrene	10	330
Butylbenzylphthalate	10	330
3,3'-Dichlorobenzidine	10	330
Benzo(a)anthracene	10	330
Chrysene	10	330
bis(2-ethylhexyl)phthalate	10	330
Di-n-octylphthalate	10	330
Benzo(b)fluoranthene	10	330
Benzo(k)fluoranthene	10	330
Benzo(a)pyrene	10	330
Indeno(1,2,3-cd)pyrene	10	330
Dibenzo(a,h)anthracene	10	330
Benzo(g,h,i)perylene	10	330

- ! The sample-specific CRQLs for soil samples will be adjusted for percent moisture and will be higher than those listed above.
- ! Medium level soil CRQL = 1000 x aqueous CRQL in ug/Kg.

ATTACHMENT 1 (page 3 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM01.9

TARGET COMPOUND LIST - PESTICIDE/PCB

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
alpha-BHC	0.05	1.7
beta-BHC	0.05	1.7
delta-BHC	0.05	1.7
gamma-BHC(Lindane)	0.05	1.7
Heptachlor	0.05	1.7
Aldrin	0.05	1.7
Heptachlor Epoxide	0.05	1.7
Endosulfan I	0.05	1.7
Dieldrin	0.10	3.3
4,4'-DDE	0.10	3.3
Endrin	0.10	3.3
Endosulfan II	0.10	3.3
4,4'-DDD	0.10	3.3
Endosulfan sulfate	0.10	3.3

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/L,ppb)
4,4'-DDT	0.10	3.3
Methoxychlor	0.5	17.0
Endrin ketone	0.10	3.3
Endrin aldehyde	0.10	3.3
alpha-Chlordane	0.05	1.7
gamma-Chlordane	0.05	1.7
Toxaphene	5.0	170.0
Aroclor-1016	1.0	33.0
Aroclor-1221	2.0	67.0
Aroclor-1232	1.0	33.0
Aroclor-1242	1.0	33.0
Aroclor-1248	1.0	33.0
Aroclor-1254	1.0	33.0
Aroclor-1260	1.0	33.0

NOTE:

! The sample-specific CRQLs for soil samples will be adjusted for percent moisture and will be higher than those listed above.

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	OLM02.1
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	No contracts have been awarded
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	14 Days or 35 Days
MATRICES:	Aqueous/Soil/Sediment*

SIGNIFICANT FEATURES

- ! The parameters include volatile, semivolatile, and pesticide/PCB compounds.
- ! Volatiles and semivolatiles are analyzed by GC/MS, pesticides/PCBs are analyzed by GC/ECD.
- ! Major Tentatively Identified Compounds (TICs) are reported for GC/MS analyses.
- ! Second column confirmation by GC/ECD is required for all pesticides/PCBs. The lower of the two concentrations detected on both columns is reported. Pesticides/PCBs which are detected at concentrations above 10 ng/uL are confirmed by GC/MS analysis.

REVISIONS/MODIFICATIONS

The following is a list of the significant changes from the OLM01.9 SOW that are incorporated in the OLM02.0 SOW (including revision Companies).

- ! For volatiles analysis, if gaseous compounds 1) fail to exhibit narrow symmetrical peak shape, 2) are not separated from the solvent front, or 3) are not resolved greater than 90% from each other; then a subambient oven controller must be used and the initial temperature must be $\leq 10^{\circ}$ C.
- ! Background subtraction must be performed utilizing a spectrum obtained no greater than 20 scans prior to the elution of BFB or DFTPP.
- ! The final column temperature must be held for three minutes following the elution of the last BNA target compound.
- ! Samples which contain alkane series in the TICs will be evaluated using the mass chromatograms of m/z 43, 57, and 71 and the alkane series will be reported as one TIC along with a total estimated concentration.
- ! Within 72 hours of detecting a multi-component pesticide/PCB in a field sample, a standard must be analyzed.
- ! The number of TICs for volatiles have been raised from 10 to 30 and for semivolatiles have been raised from 20 to 30.

RECOMMENDED USES

This Contract Laboratory Program (CLP) method is recommended for broad spectrum analysis to define the nature and extent of potential site contamination during Site Investigation (SI) and Remedial Investigation/Feasibility Study (RI/FS) activities and to verify that Remedial Design/Remedial Action (RD/RA) activities comply with the predetermined clean-up standards. This method is suitable when a 14 day or 35 day turnaround for results is adequate. It is recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at significant risk levels.

* This method is not applicable to sediment samples with high moisture content.

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in Attachment 1.

ATTACHMENT 1 (page 1 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM02.1

TARGET COMPOUND LIST - VOLATILES

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Chloromethane	10	10
Bromomethane	10	10
Vinyl chloride	10	10
Chloroethane	10	10
Methylene chloride	10	10
Acetone	10	10
Carbon disulfide	10	10
1,1-Dichloroethene	10	10
1,1-Dichloroethane	10	10
1,2-Dichloroethene (total)	10	10
Chloroform	10	10
1,2-Dichloroethane	10	10
2-Butanone	10	10
1,1,1-Trichloroethane	10	10
Carbon tetrachloride	10	10
Bromodichloromethane	10	10
1,2-Dichloropropane	10	10

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg/ppb)
cis-1,3-Dichloropropene	10	10
Trichloroethene	10	10
Dibromochloromethane	10	10
1,1,2-Trichloroethane	10	10
Benzene	10	10
trans-1,3-Dichloropropene	10	10
Bromoform	10	10
4-Methyl-2-pentanone	10	10
2-Hexanone	10	10
Tetrachloroethene	10	10
1,1,2,2-Tetrachloroethane	10	10
Toluene	10	10
Chlorobenzene	10	10
Ethylbenzene	10	10
Styrene	10	10
Xylenes (Total)	10	10

- ! The sample-specific CRQLs for soil samples will be adjusted for percent moisture and will be higher than those listed above.
- ! Medium level soil CRQL = 120 x aqueous CRQL reported in ug/kg.

ATTACHMENT I (page 2 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM02.1

TARGET COMPOUND LIST - SEMIVOLATILES

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
Phenol	10	330
bis(2-chloroethyl)ether	10	330
2-Chlorophenol	10	330
1,3-Dichlorobenzene	10	330
1,4-Dichlorobenzene	10	330
1,2-Dichlorobenzene	10	330
2-Methylphenol	10	330
2,2'-oxybis(1-chloropropane)	10	330
4-Methylphenol	10	330
N-Nitroso-di-n-propylamine	10	330
Hexachloroethane	10	330
Nitrobenzene	10	330
Isophorone	10	330
2-Nitrophenol	10	330
2,4-Dimethylphenol	10	330
bis-(2-Chloroethoxy)methane	10	330
2,4-Dichlorophenol	10	330
1,2,4-Trichlorobenzene	10	330
Naphthalene	10	330
4-Chloroaniline	10	330
Hexachlorobutadiene	10	330
4-Chloro-3-methylphenol	10	330
2-Methylnaphthalene	10	330
Hexachlorocyclopentadiene	10	330
2,4,6-Trichlorophenol	10	330
2,4,5-Trichlorophenol	25	800
2-Chloronaphthalene	10	330
2-Nitroaniline	25	800
Dimethylphthalate	10	330
Acenaphthylene	10	330
2,6-Dinitrotoluene	10	330
3-Nitroaniline	25	800

Compound	Aqueous CRQL	Low Soil CRQL
	(ug/L,ppb)	(ug/Kg,ppb)
Acenaphthene	10	330
2,4-Dinitrophenol	25	800
4-Nitrophenol	25	800
Dibenzofuran	10	330
2,4-Dinitrotoluene	10	330
Diethylphthalate	10	330
4-Chlorophenyl-phenylether	10	330
Fluorene	10	330
4-Nitroaniline	25	800
4,6-Dinitro-2-methylphenol	25	800
N-Nitrosodiphenylamine	10	330
4-Bromophenyl-phenylether	10	330
Hexachlorobenzene	10	330
Pentachlorophenol	25	800
Phenanthrene	10	330
Anthracene	10	330
Carbazole	10	330
Di-n-butylphthalate	10	330
Fluoranthene	10	330
Pyrene	10	330
Butylbenzylphthalate	10	330
3,3'-Dichlorobenzidine	10	330
Benzo(a)anthracene	10	330
Chrysene	10	330
bis(2-ethylhexyl)phthalate	10	330
Di-n-octylphthalate	10	330
Benzo(b)fluoranthene	10	330
Benzo(k)fluoranthene	10	330
Benzo(a)pyrene	10	330
Indeno(1,2,3-cd)pyrene	10	330
Dibenzo(a,h)anthracene	10	330
Benzo(g,h,i)perylene	10	330

- ! The sample-specific CRQLs for soil samples will be adjusted for percent moisture and will be higher than those listed above.
- ! Medium level soil CRQL = 1000 x aqueous CRQL in ug/Kg.

ATTACHMENT 1 (page 3 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM02.1

TARGET COMPOUND LIST - PESTICIDE/PCB

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/Kg,ppb)
alpha-BHC	0.05	1.7
beta-BHC	0.05	1.7
delta-BHC	0.05	1.7
gamma-BHC(Lindane)	0.05	1.7
Heptachlor	0.05	1.7
Aldrin	0.05	1.7
Heptachlor Epoxide	0.05	1.7
Endosulfan I	0.05	1.7
Dieldrin	0.10	3.3
4,4'-DDE	0.10	3.3
Endrin	0.10	3.3
Endosulfan II	0.10	3.3
4,4'-DDD	0.10	3.3
Endosulfan sulfate	0.10	3.3

Compound	Aqueous CRQL (ug/L,ppb)	Low Soil CRQL (ug/L,ppb)
4,4'-DDT	0.10	3.3
Methoxychlor	0.5	17.0
Endrin ketone	0.10	3.3
Endrin aldehyde	0.10	3.3
alpha-Chlordane	0.05	1.7
gamma-Chlordane	0.05	1.7
Toxaphene	5.0	170.0
Aroclor-1016	1.0	33.0
Aroclor-1221	2.0	67.0
Aroclor-1232	1.0	33.0
Aroclor-1242	1.0	33.0
Aroclor-1248	1.0	33.0
Aroclor-1254	1.0	33.0
Aroclor-1260	1.0	33.0

NOTE:

! The sample-specific CRQLs for soil samples will be adjusted for percent moisture and will be higher than those listed above.

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS

MULTI-MEDIA, HIGH CONCENTRATION

DOCUMENT NUMBER:	Not Applicable
DOCUMENT DATE:	September 1988
EFFECTIVE DATES:	June 7, 1989 through December 26, 1991
CONCENTRATION:	High: Greater than 20 ppm
DATA TURNAROUND:	35 Days
MATRICES:	Liquid/Solid/Multi-phase

SIGNIFICANT FEATURES

- ! No holding times are designated for high concentration samples.
- ! The analyses are suitable for highly contaminated samples (>20 mg/Kg).
- ! The analyses are acceptable for liquid, solid, or multiphase samples. Multi-phase samples are separated into water miscible liquid, water immiscible liquid, or solid phases. Each phase is analyzed separately.
- ! Volatile, extractable (semivolatiles and pesticides), and multicomponent extractable (Aroclors and Toxaphene) cor
- ! Volatiles and extractables are analyzed by GC/MS; Aroclors and Toxaphene are analyzed by GC/ECD.
- ! Second column confirmation by GC/ECD is required for Aroclors and Toxaphene.
- ! Major Tentatively Identified Compounds (TICs) are reported for GC/MS analyses.

REVISIONS/MODIFICATIONS

The 1/89 and 4/89 revisions to the 9/88 SOW do not significantly affect data useability.

RECOMMENDED USES

This Routine Analytical Services (RAS) method is recommended for pre-remedial, remedial, or removal projects where high concentrations of organic contaminants (greater than 20 mg/Kg) are suspected and a thirty five day tumaround for results is adequate. It is recommended for samples obtained from drummed material, waste pits or lagoons, piles of waste, tanker trucks, on site tanks, or apparent contaminated soil areas. The waste material may be industrial process waste, byproducts, raw materials, intermediates and contaminated products. Samples may be spent oil, spent solvents, paint wastes, metal treatment wastes, and polymer formulations.

The method is suitable for solids, liquids, or multiphase samples, a phase being either water miscible liquid, water immiscible liquid, or solid. Various methods of phase separation may be utilized depending on the number and types of phases in a sample

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, HIGH CONCENTRATION SEPTEMBER, 1988

ATTACHMENT I (page 2 of 3)

TARGET COMPOUND LIST - VOLATILES

Compound	CRQL (mg/Kg,ppm)
Chloromethane	5.0
Bromomethane	5.0
Vinyl Chloride	5.0
Chloroethane	5.0
Methylene Chloride	2.5
Acetone	5.0
Carbon Disulfide	2.5
1,1-Dichloroethene	2.5
1,1-Dichloroethane	2.5
1,2-Dichloroethene (total)	2.5
Chloroform	2.5
1,2-Dichloroethane	2.5
2-Butanone	5.0
1,1,1-Trichloroethane	2.5
Carbon Tetrachloride	2.5
Vinyl Acetate	5.0
Bromodichloromethane	2.5

Compound	CRQL (mg/Kg,ppm)
1,2-Dichloropropane	2.5
cis-1,3-Dichloropropene	2.5
Trichloroethene	2.5
Dibromochloromethane	2.5
1,1,2-Trichloroethane	2.5
Benzene	2.5
trans-1,3-Dichloropropene	2.5
Bromoform	2.5
4-Methyl-2-pentanone	5.0
2-Hexanone	5.0
Tetrachloroethene	2.5
1,1,2,2-Tetrachloroethane	2.5
Toluene	2.5
Chlorobenzene	2.5
Ethylbenzene	2.5
Styrene	2.5
Xylene (Total)	2.5

- ! ALL CRQLs ARE BASED ON WET WEIGHT AND APPLY TO SOLID AND LIQUID SAMPLES.
- ! RESULTS FOR BOTH SOLID AND LIQUID SAMPLES ARE REPORTED AS MG/KG, WET WEIGHT.

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, HIGH CONCENTRATION SEPTEMBER, 1988

TARGET COMPOUND LIST - EXTRACTABLES

Compound	CRQL (mg/Kg,ppm)
Phenol	20
bis(2-Chloroethyl)ether	20
2-Chlorophenol	20
1,3-Dichlorobenzene	20
1,4-Dichlorobenzene	20
Benzyl alcohol	20
1,2-Dichlorobenzene	20
2-Methylphenol	20
bis(2-Chloroisopropyl) ether	20
4-Methylphenol	20
N-Nitroso-di-n-dipropylamine	20
Hexachloroethane	20
Nitrobenzene	20
Isophorone	20
2-Nitrophenol	20
2,4-Dimethylphenol	20
Benzoic acid	100
bis(2-Chloroethoxy) methane	20
2,4-Dichlorophenol	20
1,2,4-Trichlorobenzene	20
Naphthalene	20
4-Chloroaniline	20
Hexachlorobutadiene	20
4-Chloro-3-methylphenol (para-chloro-meta-cresol)	20
2-Methylnaphthalene	20
Hexachlorocyclopentadiene	20
2,4,6-Trichlorophenol	20
2,4,5-Trichlorophenol	100
2-Chloronaphthalene	20
2-Nitroaniline	100
Dimethylphthalate	20

Compound	CRQL (mg/Kg,ppm)
Acenaphthylene	20
2,6-Dinitrotoluene	20
3-Nitroaniline	100
Acenaphthene	20
2,4-Dinitrophenol	100
4-Nitrophenol	100
Dibenzofuran	20
2,4-Dinitrotoluene	20
Diethylphthalate	20
4-Chlorophenyl-phenylether	20
Fluorene	20
4-Nitroaniline	100
4,6-Dinitro-2-methylphenol	100
N-nitrosodiphenylamine	20
4-Bromophenyl-phenylether	20
alpha-BHC	20
Hexachlorobenzene	20
beta-BHC	20
Pentachlorophenol	100
gamma-BHC (Lindane)	20
Phenanthrene	20
Anthracene	20
delta-BHC	20
Heptachlor	20
Aldrin	20
Di-n-butylphthalate	20
Fluoranthene	20
Heptachlor epoxide	20
Monochlorobiphenyl	100
Dichlorobiphenyl	100
Trichlorobiphenyl	100

- ! ALL CRQLs ARE BASED ON WET WEIGHT AND APPLY TO SOLID AND LIQUID SAMPLES. ! RESULTS FOR BOTH SOLID AND LIQUID SAMPLES ARE REPORTED AS MG/KG, WET WEIGHT.

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, HIGH CONCENTRATION SEPTEMBER, 1988

TARGET COMPOUND LIST - EXTRACTABLES (CONTINUED)

Compound	CRQL (mg/Kg,ppm)
Tetrachlorobiphenyl	100
Pyrene	20
gamma-Chlordane	20
Endosulfan I	20
alpha-Chlordane	20
4,4,'-DDE	20
Dieldrin	20
Hexachlorobiphenyl	100
Pentachlorobiphenyl	100
Endrin	20
Endosulfan II	20
4,4'-DDD	20
Heptachlorobiphenyl	100
Butylbenzylphthalate	20
Endosulfan sulfate	20
4,4'-DDT	20

Compound	CRQL (mg/Kg,ppm)
Endrin ketone	20
Benzo(a)anthracene	20
Methoxychlor	20
Chrysene	20
Octachlorobiphenyl	200
3,3'-Dichlorobenzidine	40
bis(2-Ethylhexyl)phthalate	20
Nonachlorobiphenyl	200
Decachlorobiphenyl	200
Di-n-octylphthalate	20
Benzo(b)fluoranthene	20
Benzo(k)fluoranthene	20
Benzo(a)pyrene	20
Indeno(1,2,3-cd)pyrene	20
Dibenz(a,h)anthracene	20
Benzo(g,h,i)perylene	20

${\sf TARGE} \underline{{\sf T}} \, {\sf COMPOUND} \, {\sf LIST} \, {\sf -MULTICOMPONENT} \, {\sf EXTRACTABLES}$

Compound	CRQL (mg/Kg,ppm)
Toxaphene	50
Aroclor 1016	10
Aroclor 1221	10
Aroclor 1232	10
Aroclor 1242	10
Aroclor 1248	10
Aroclor 1254	10
Aroclor 1260	10

- ! ALL CRQLs ARE BASED ON WET WEIGHT AND APPLY TO SOLID AND LIQUID SAMPLES. ! RESULTS FOR BOTH SOLID AND LIQUID SAMPLES ARE REPORTED AS MG/KG, WET WEIGHT.

TITLE: USEPA CONTRACT LABORATORY PROGRAM

STATEMENT OF WORK FOR INORGANIC ANALYSIS

MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	Not Applicable
DOCUMENT DATE:	July 1988
EFFECTIVE DATES:	May 23, 1989 through March 29, 1992
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	35 Days
MATRICES:	Aqueous/Soil/Sediment *

SIGNIFICANT FEATURES

! The analyses are suitable for aqueous, soil or sediment samples at low to medium concentration levels.

REVISIONS/MODIFICATIONS

2/89 - Method for Total Cyanide (CN) Analysis by Midi Distillation Method 335.2 CLP-M was added.

6/89 Revisions to the 7/88 SOW do not significantly affect data useability.

RECOMMENDED USES

This Routine Analytical Services (RAS) method is recommended for broad spectrum analysis to define the nature and extent of potential site contamination during SSI, LSI, and RI/FS activities. This method is suitable when a thirty five day turnaround for results is adequate. It is recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at significant risk levels.

* Sediment samples with high moisture content should be solicited as RAS + SAS (Special Analytical Service) in order to achieve the CRDLs.

ANALYTES/CRDLs

The parameters included in the analysis and the Contract Required Detection Limits (CRDLs) are listed in Attachm

ATTACHMENT I (page 1 of 1)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR INORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION JULY, 1988

TARGET ANALYTE LIST

Analyte	AQUEOUS CRDL (ug/L,ppb)	SOIL CRDL (mg/Kg,ppm)
Aluminum	200	40
Antimony	60	12
Arsenic	10	2
Barium	200	40
Beryllium	5	1
Cadmium	5	1
Calcium	5000	1000
Chromium	10	2
Cobalt	50	10
Copper	25	5
Iron	100	20
Lead	3	0.6

Analyte	AQUEOUS CRDL (ug/L,ppb)	SOIL CRDL (mg/Kg,ppm)
Magnesium	5000	1000
Manganese	15	3
Mercury	0.2	0.1
Nickel	40	8
Potassium	5000	1000
Selenium	5	1
Silver	10	2
Sodium	5000	1000
Thallium	10	2
Vanadium	50	10
Zinc	20	4
Cyanide	10	2

NOTE:

! THE SAMPLES PECIFIC CRDLS FOR SOIL SAMPLES WILL BE ADJUSTED FOR PERCENT MOISTURE AND WILL BE HIGHER THAN THOSE LISTED ABOVE.

TITLE: USEPA CONTRACT LABORATORY PROGRAM

STATEMENT OF WORK FOR INORGANIC ANALYSIS

MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	ILM01.0
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	September 7, 1990 through September 26, 1993
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	35 Days
MATRICES:	Aqueous/Soil/Sediment *

SIGNIFICANT FEATURES

- ! The analyses are suitable for aqueous, soil or sediment samples at low to medium concentration levels.
- ! This Statement of Work includes the midi distillation for cyanide analysis and the microwave digestion for GFAA and ICP analyses. These two sample preparation procedures require less sample volume than the traditional Statement of Work sample preparation procedures.

REVISIONS/MODIFICATIONS

None to Date

RECOMMENDED USES

This Routine Analytical Services (RAS) method is recommended for broad spectrum analysis to define the nature and extent of potential site contamination during SSI, LSI, and RI/FS activities. This method is suitable when a thirty five day turnaround for results is adequate. It is recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at significant risk levels.

* Sediment samples with high moisture content should be solicited as RAS + SAS (Special Analytical Service) in order to achieve the CRDLs.

ANALYTES/CRDLs

The parameters included in the analysis and the Contract Required Detection Limits (CRDLs) are listed in Att

ATTACHMENT I (page 1 of 1)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR INORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION ILM01.0

TARGET ANALYTE LIST

Analyte	Aqueous CRDL (ug/L,ppb)	Soil CRDL (mg/Kg,ppm)
Aluminum	200	40
Antimony	60	12
Arsenic	10	2
Barium	200	40
Beryllium	5	1
Cadmium	5	1
Calcium	5000	1000
Chromium	10	2
Cobalt	50	10
Copper	25	5
Iron	100	20
Lead	3	0.6

Analyte	Aqueous CRDL (ug/L,ppb)	Soil CRDL (mg/Kg,ppm)
Magnesium	5000	1000
Manganese	15	3
Mercury	0.2	0.1
Nickel	40	8
Potassium	5000	1000
Selenium	5	1
Silver	10	2
Sodium	5000	1000
Thallium	10	2
Vanadium	50	10
Zinc	20	4
Cyanide	10	2

NOTE:

 $\hbox{!} \quad \hbox{THESAMPLE-SPECIFIC CRDLs FOR SOIL SAMPLES WILL BE ADJUSTED FOR PERCENT MOISTURE AND WILL BE HIGHER THAN } \\ \quad \hbox{THOSE LISTED ABOVE}.$

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR INORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	ILM02.1
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	March 1993 through October 1994
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	14 Days or 35 Days
MATRICES:	Aqueous/Soil/Sediment*

SIGNIFICANT FEATURES

! The analyses are suitable for aqueous, soil, or sediment samples at low to medium concentration levels.

REVISIONS/MODIFICATIONS

The following is a list of the significant changes from the ILM01.0 SOW that are incorporated in the ILM02.1 SOW:

- ! Specific analysis instructions are presented when samples are received for dissolved metals analysis.
- ! Requirements for contract reports/deliverables distribution are included for 14 day turnaround contracts.
- ! SOPs are now required to be distributed by the laboratories to EPA-NEIC.
- ! Microwave digestion for soil/sediment samples is not appropriate for quantitative recovery of antimony.

RECOMMENDED USES

This Contract Laboratory Program (CLP) method is recommended for broad spectrum analysis to define the nature and extent of potential site contamination during Site Investigation (SI) and Remedial Investigation/Feasibility Study (RI/FS) activities and to verify that Remedial Design/Remedial Action (RD/RA) activities comply with pre-determined clean-up standards. This method is suitable when a 14 day or 35 day tumaround for results is adequate. It is recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at significant risk levels.

* This method is not applicable to sediment samples with high moisture content.

ANALYTES/CRDLs

The parameters included in the analysis and the Contract Required Detection Limits (CRDLs) are listed in Attachn

ATTACHMENT 1

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR INORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION ILM02.1

TARGET ANALYTE LIST

Analyte	Aqueous CRDL (ug/L,ppb)	Soil CRDL (mg/Kg,ppm)
Aluminum	200	40
Antimony	60	12
Arsenic	10	2
Barium	200	40
Beryllium	5	1
Cadmium	5	1
Calcium	5000	1000
Chromium	10	2
Cobalt	50	10
Copper	25	5
Iron	100	20
Lead	3	0.6

Analyte	Aqueous CRDL (ug/L,ppb)	Soil CRDL (mg/Kg,ppm)
Magnesium	5000	1000
Manganese	15	3
Mercury	0.2	0.1
Nickel	40	8
Potassium	5000	1000
Selenium	5	1
Silver	10	2
Sodium	5000	1000
Thallium	10	2
Vanadium	50	10
Zinc	20	4
Cyanide	10	2

NOTE:

! The sample-specific CRDLs for soil samples will be adjusted for percent moisture and will be higher than those listed a

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR INORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	ILM03.0
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	May 1993 through March 1996
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	14 Days or 35 Days
MATRICES:	Aqueous/Soil/Sediment*

SIGNIFICANT FEATURES

! The analyses are suitable for aqueous, soil, or sediment samples at low to medium concentration levels.

REVISIONS/MODIFICATIONS

The following is a list of the significant changes from the ILM02.1 SOW that are incorporated in the ILM03.0 SOW:

- ! An analytical spike and an aqueous Laboratory Control Sample (LCS) are not required when analyzing field samples fo
- ! Terminology was added to require that sample coolers be returned to the appropriate sampling office within 14 days following shipment receipt.
- ! Additional instrumentation requirements were added for greater than 500 samples per month capacity.
- ! For cyanide water analysis, the LCS requirement was changed from "not required" to "using the distilled ICV as the LC

RECOMMENDED USES

This Contract Laboratory Program (CLP) method is recommended for broad spectrum analysis to define the nature and extent of potential site contamination during Site Investigation (SI) and Remedial Investigation/Feasibility Study (RI/FS) activities and to verify that Remedial Design/Remedial Action (RD/RA) activities comply with pre-determined clean-up standards. This method is suitable when a 14 day or 35 day turnaround for results is adequate. It is recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at significant risk levels.

* This method is not applicable to sediment samples with high moisture content.

ANALYTES/CRDLs

The parameters included in the analysis and the Contract Required Detection Limits (CRDLs) are listed in Attachment

ATTACHMENT 1

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR INORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION ILM03.0

TARGET ANALYTE LIST

Analyte	Aqueous CRDL (ug/L,ppb)	Soil CRDL (mg/Kg,ppm)
Aluminum	200	40
Antimony	60	12
Arsenic	10	2
Barium	200	40
Beryllium	5	1
Cadmium	5	1
Calcium	5000	1000
Chromium	10	2
Cobalt	50	10
Copper	25	5
Iron	100	20
Lead	3	0.6

Analyte	Aqueous CRDL (ug/L,ppb)	Soil CRDL (mg/Kg,ppm)
Magnesium	5000	1000
Manganese	15	3
Mercury	0.2	0.1
Nickel	40	8
Potassium	5000	1000
Selenium	5	1
Silver	10	2
Sodium	5000	1000
Thallium	10	2
Vanadium	50	10
Zinc	20	4
Cyanide	10	2

NOTE:

! The sample-specific CRDLs for soil samples will be adjusted for percent moisture and will be higher than those listed a

TITLE: USEPA CONTRACT LABORATORY PROGRAM

STATEMENT OF WORK FOR INORGANIC ANALYSIS

MULTI-MEDIA, HIGH CONCENTRATION

DOCUMENT NUMBER:	IHC01.2
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	May 15, 1991 through November 30, 1993
CONCENTRATION:	High
DATA TURNAROUND:	35 Days
MATRICES:	Liquid/Solid/Multi-phase

SIGNIFICANT FEATURES

- ! The analyses are suitable for highly contaminated samples.
- ! The analyses are acceptable for liquid, solid, or multiphase samples. Multi-phase samples are separated into water miscible liquid, water immiscible liquid, or solid phases. Each phase is analyzed separately.
- ! The analyses include conductivity and pH; potassium is not included.

REVISIONS/MODIFICATIONS

The IHC01.1 and IHC01.2 revisions to the IHC01.0 SOW do not significantly affect data useability.

RECOMMENDED USES

This Routine Analytical Services (RAS) method is recommended for pre-remedial, or removal projects where high concentrations of inorganic contaminants are suspected and a thirty five day turnaround for results is adequate. It is recommended for samples obtained from drummed material, waste pits or lagoons, piles of waste, tanker trucks, on site tanks, or apparent contaminated soil areas. The waste material may be industrial process waste, byproducts, raw materials, intermediates and contaminated products. Samples may be spent oil, spent solvents, paint wastes, metal treatment wastes, and polymer formulations.

The method is suitable for solids, liquids, or multiphase samples, a phase being either water miscible liquid, water immiscible liquid, or solid. A phase separation step is applied prior to digestion. Each phase is analyzed and reported as a separat

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in .

ATTACHMENT I (page 1 of 1)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR INORGANIC ANALYSIS MULTI-MEDIA, HIGH CONCENTRATION IHC01.2

TARGET ANALYTE LIST

Analyte	CRQL (mg/Kg, ppm)
Aluminum	80
Antimony	20
Arsenic	5
Barium	80
Beryllium	5
Cadmium	10
Calcium	80
Chromium	10
Cobalt	20
Copper	40
Iron	20
Lead	10
Magnesium	80

Analyte	CRQL (mg/Kg, ppm)
Manganese	10
Mercury	0.3
Nickel	20
Selenium	5
Silver	10
Sodium	80
Thallium	20
Vanadium	20
Zinc	10
Cyanide	1.5
pН	NA
Conductivity	3.0 (µmhos/cm)

- ! ALL CRQLs ARE BASED ON WET WEIGHT AND APPLY TO SOLID AND LIQUID SAMPLES.
- ! RESULTS FOR BOTH SOLID AND LIQUID SAMPLES ARE REPORTED AS MG/KG, WET WEIGHT.

TITLE: USEPA CONTRACT LABORATORY PROGRAM

STATEMENT OF WORK FOR ANALYSIS OF POLYCHLORINATED DIBENZO-p-DIOXINS (PCDD) AND POLYCHLORINATED DIBENZOFURANS (PCDF) MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	DFLMO1.2
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	May 1992 through May 1994
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	45 days
MATRICES:	Aqueous/Soil/Fly Ash/Chemical Waste

SIGNIFICANT FEATURES

- ! The parameters include 2,3,7,8-substituted Tetra-, Penta-, Hexa-, Hepta-, and Octachlorinated dibenzo-p-dioxins (PCDDs) and dibenz
- ! Total homologue concentrations are reported for a given level of chlorination (i.e. Total TCDD, Total PCDF, etc.).
- ! 2,3,7,8-TCDD toxic equivalence is determined using all 2,3,7,8-substituted isomers.
- ! The dioxins and furans are analyzed by High Resolution Gas Chromatography and Low Resolution Mass Spectrometry (HRGC/LRM
- ! Second column confirmation is required if the toxic equivalence is greater than or equal to 0.7ppb (soil or fly ash), 7ppt (aqueous) or
- ! Chemical waste includes oils, stillbottoms, oily sludge, wet fuel oil, and surface water heavily contaminated with oils. Soil samples i

REVISIONS/MODIFICATIONS

Revisions DFLMO1.1 (September 1991) and DFLMO1.2 (April 1992) do not significantly affect data useability.

RECOMMENDED USES

This Contract Laboratory Program (CLP) method is recommended for analysis of polychlorinated dioxins and furans to define the nat site contamination during Site Investigation (SI) and Remedial Investigation/Feasibility Study (RI/FS) activities and to verify that Remed (RD/RA) activities comply with pre-determined clean-up standards. This method is suitable when a forty five day turnaround for results recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at significant ri applicable to sediment samples with high moisture content.

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in Table 1.

TABLE 1

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ANALYSIS OF POLYCHLORINATED DIBENZO-p-DIOXINS (PCDD) AND POLYCHLORINATED DIBENZOFURANS (PCDF) MULTI-MEDIA, MULTI-CONCENTRATION DFLMO1.2

TARGET COMPOUND LIST

Compound	Aqueous CRQL (ng/L,ppt)	Low Soil CRQL ^A (ug/Kg,ppb)	Fly Ash CRQL ^A (ug/Kg,ppb)	Chemical Waste CRQL (ug/Kg,ppb)	TEFs ^B
2,3,7,8-TCDD	10	1.0	1.0	10	1.0
2,3,7,8-TCDF	10	1.0	1.0	10	0.10
1,2,3,7,8-PeCDF	25	2.5	2.5	25	0.05
1,2,3,7,8-PeCDD	25	2.5	2.5	25	0.50
2,3,4,7,8-PeCDF	25	2.5	2.5	25	0.50
1,2,3,4,7,8-HxCDF	25	2.5	2.5	25	0.10
1,2,3,6,7,8-HxCDF	25	2.5	2.5	25	0.10
1,2,3,4,7,8-HxCDD	25	2.5	2.5	25	0.10
1,2,3,6,7,8-HxCDD	25	2.5	2.5	25	0.10
1,2,3,7,8,9-HxCDD	25	2.5	2.5	25	0.10
2,3,4,6,7,8-HxCDF	25	2.5	2.5	25	0.10
1,2,3,7,8,9-HxCDF	25	2.5	2.5	25	0.10
1,2,3,4,6,7,8-HpCDF	25	2.5	2.5	25	0.01
1,2,3,4,6,7,8-HpCDD	25	2.5	2.5	25	0.01
1,2,3,4,7,8,9-HpCDF	25	2.5	2.5	25	0.01
OCDD	50	5.0	5.0	50	0.001
OCDF	50	5.0	5.0	50	0.001

Homologue Compounds	Number of Possible Isomers	Number of 2,3,7,8-Substituted Isomers
TOTAL TCDD	22	1
TOTAL TCDF	38	1
TOTAL PeCDD	14	1
TOTAL PeCDF	28	2
TOTAL HxCDD	10	3
TOTAL HXCDF	16	4
TOTAL HpCDD	2	1
TOTAL HpCDF	4	2
OCDD	1	1
OCDF	1	1

NOTE:

A. The sample-specific CRQLs for soils and fly ash samples will be adjusted for percent moisture and will be higher than those listed above

B. TEF = 2,3,7,8-TCDD Toxic Equivalence Factors (TEFs) for the PCDDs/PCDFs.

TITLE: SUPERFUND ANALYTICAL METHODS FOR LOW CONCENTRATION WATER FOR ORGANICS

DOCUMENT NUMBER:	Not Applicable
DOCUMENT DATE:	October 1992
EFFECTIVE DATES:	December 1992 through June 1994
CONCENTRATION:	Low
DATA TURNAROUND:	14 Days
MATRICES:	Aqueous

SIGNIFICANT FEATURES

- ! The parameters include volatile, semivolatile, and pesticide/PCB compounds.
- ! Volatiles and semivolatiles are analyzed by GC/MS; pesticides/PCBs are analyzed by GC/ECD.
- ! All parameters require significantly reduced CRQLs as compared to the OLM01.9 SOW.
- ! A 25mL aliquot of sample is purged for volatiles analysis.
- ! Analysis of a Laboratory Control Sample (LCS) for each parameter is required.
- ! Semivolatile and pesticide/PCB samples must be extracted by continuous liquid-liquid extraction procedures.

REVISIONS/MODIFICATIONS

The following is a list of the significant changes from the 6/91 Low Concentration SOW:

- ! The requirement for a diskette deliverable was removed.
- ! Technical acceptance criteria for the volatile LCS were established.
- ! Potential action against a laboratory for Performance Evaluation (PE) scores below 75% was added.

RECOMMENDED USES

This Contract Laboratory Program (CLP) method is recommended for broad spectrum analysis to define the nature and level organic* contamination in water supplies during Site Investigation (SI) and Remedial Investigation/Feasibility Study verify that Remedial Design/Remedial Action (RD/RA) activities comply with pre-determined clean-up standards. This method is made that the sources are expected to be from drinking water and well/ground water sources around Superfund sites. This method is suit turnaround for results is adequate.

* This method should not be used for low concentration volatile organics analyses in Region I when comparison to the dr Contaminant Levels (MCLs) is required. The Region I EPA 524.2 standard specifications should be utilized for this pu

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in Attachmet

ATTACHMENT 1 (page 1 of 3)

SUPERFUND ANALYTICAL METHODS FOR LOW CONCENTRATION WATER FOR ORGANICS ANALYSIS 10/92

TARGET COMPOUND LIST - VOLATILES

COMPOUND	CRQL (ug/L, ppb)
Chloromethane	1
Bromomethane	1
Vinyl chloride	1
Chloroethane	1
Methylene chloride	2
Acetone	5
Carbon disulfide	1
1,1-Dichloroethene	1
1,1-Dichloroethane	1
cis-1,2-Dichloroethene	1
trans-1,2-Dichloroethene	1
Chloroform	1
1,2-Dichloroethane	1
2-Butanone	5
Bromochloromethane	1
1,1,1-Trichloroethane	1
Carbon tetrachloride	1
Bromodichloromethane	1
1,2-Dichloropropane	1
cis-1,3-Dichloropropene	1

COMPOUND	CRQL (ug/L, ppb)
Trichloroethene	1
Dibromochloromethane	1
1,1,2-Trichloroethane	1
Benzene	1
trans-1,3-Dichloropropene	1
Bromoform	1
4-Methyl-2-pentanone	5
2-Hexanone	5
Tetrachloroethene	1
1,1,2,2-Tetrachloroethane	1
1,2-Dibromoethane	1
Toluene	1
Chlorobenzene	1
Ethylbenzene	1
Styrene	1
Xylenes (Total)	1
1,3-Dichlorobenzene	1
1,4-Dichlorobenzene	1
1,2-Dichlorobenzene	1
1,2-Dibromo-3-chloropropane	1

ATTACHMENT I (page 2 of 3)

SUPERFUND ANALYTICAL METHODS FOR LOW CONCENTRATION WATER FOR ORGANICS ANALYSIS 10/92

TARGET COMPOUND LIST - SEMIVOLATILES

COMPOUND	CRQL (ug/L, ppb)
Phenol	5
bis-(2-chloroethyl)ether	5
2-Chlorophenol	5
2-Methylphenol	5
2,2'-oxybis(1-chloropropane)	5
4-Methylphenol	5
N-Nitroso-di-n-propylamine	5
Hexachloroethane	5
Nitrobenzene	5
Isophorone	5
2-Nitrophenol	5
2,4-Dimethylphenol	5
bis-(2-Chloroethoxy)methane	5
2,4-Dichlorophenol	5
1,2,4-Trichlorobenzene	5
Naphthalene	5
4-Chloroaniline	5
Hexachlorobutadiene	5
4-Chloro-3-methylphenol	5
2-Methylnaphthalene	5
Hexachlorocyclopentadiene	5
2,4,6-Trichlorophenol	5
2,4,5-Trichlorophenol	20
2-Chloronaphthalene	5
2-Nitroaniline	20
Dimethylphthalate	5
Acenaphthylene	5
2,6-Dinitrotoluene	5
3-Nitroaniline	20
Acenaphthene	5

COMPOUND	CRQL (ug/L, ppb)
2,4-Dinitrophenol	20
4-Nitrophenol	20
Dibenzofuran	5
2,4-Dinitrotoluene	5
Diethylphthalate	5
4-Chlorophenyl-phenylether	5
Fluorene	5
4-Nitroaniline	20
4,6-Dinitro-2-methylphenol	20
N-Nitrosodipenylamine	5
4-Bromophenyl-phenylether	5
Hexachlorobenzene	5
Pentachlorophenol	20
Phenanthrene	5
Anthracene	5
Di-n-butylphthalate	5
Fluoranthene	5
Pyrene	5
Butylbenzylphthalate	5
3,3'-Dichlorobenzidine	5
Benzo(a)anthracene	5
Chrysene	5
bis(2-ethylhexyl)phthalate	5
Di-n-octylphthalate	5
Benzo(b)fluoranthene	5
Benzo(k)fluoranthene	5
Benzo(a)pyrene	5
Indeno(1,2,3-cd)pyrene	5
Dibenzo(a,h)anthracene	5
Benzo(g,h,i)perylene	5

ATTACHMENT 1 (page 3 of 3)

SUPERFUND ANALYTICAL METHODS FOR LOW CONCENTRATION WATER FOR ORGANICS ANALYSIS 10/92

TARGET COMPOUND LIST - PESTICIDE/PCB

PARAMETER	CRQL (ug/L, ppb)
alpha-BHC	0.01
beta-BHC	0.01
delta-BHC	0.01
gamma-BHC (Lindane)	0.01
Heptachlor	0.01
Aldrin	0.01
Heptachlor epoxide	0.01
Endosulfan I	0.01
Dieldrin	0.02
4,4'-DDE	0.02
Endrin	0.02
Endosulfan II	0.02
4,4'-DDD	0.02
Endosulfan sulfate	0.02

PARAMETER	CRQL (ug/L, ppb)
4,4'-DDT	0.02
Methoxychlor	0.10
Endrin ketone	0.02
Endrin aldehyde	0.02
alpha-Chlordane	0.01
gama-Chlordane	0.01
Toxaphene	1.0
Aroclor-1016	0.2
Aroclor-1221	0.4
Aroclor-1232	0.2
Aroclor-1242	0.2
Aroclor-1248	0.2
Aroclor-1254	0.2
Aroclor-1260	0.2

TITLE: USEPA CONTRACT LABORATORY PROGRAM WATER QUALITY PARAMETERS (WQP) IN MULTI-CONCENTRATION WATER

DOCUMENT NUMBER:	Not Applicable
DOCUMENT DATE:	June 1993
EFFECTIVE DATES:	July 1993 through June 1994
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	14 Days or 35 Days
MATRICES:	Aqueous

SIGNIFICANT FEATURES

- ! The parameters include Alkalinity, Ammonia-Nitrogen, Total Organic Carbon (TOC), Chemical Oxygen Demand (COD), Chloride, Nitrate/Nitrite, Total Phosphorous, Total Dissolved Solids (TDS), Total Suspended Solids (TSS), and Sulfate.
- ! Ion chromatography may be used in place of conventional methods for the determination of Chloride, Nitrate/Nitrite, Phosphorous, and Sulfate.
- ! The laboratory IDL for a parameter may exceed the CRDL if the sample concentrations are greater than 5xIDL.

REVISIONS/MODIFICATIONS

The following is a list of the significant changes from the 2/93 version that are incorporated in the 6/93 SOW:

- ! The CRDL for Total Dissolved Solids was elevated from 4,000 ug/L to 10,000 ug/L.
- ! The analysis procedure for Total Organic Carbon (TOC) was changed from following a step-by-step procedure for calibrating and standardizing the TOC analyzer to performing these functions according to the manufacturers specifications.

RECOMMENDED USES

These Contract Laboratory Program (CLP) methods are intended for use only with aqueous samples. They are recommended for analysis of selected water quality parameters to define the nature and extent of potential site contamination during Site Investigation (SI) and Remedial Investigation/Feasibility Study (RI/FS) activities and to verify that Remedial Design/Remedial Action (RD/RA) activities comply with pre-determined clean-up standards. These methods may also be used for monitoring the wastewater treatment processes of pretreatment plants and Publicly Owned Treatment Works (POTWs). This method is suitable when a 14 day or 35 day tumaround for results is adequate. It is recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at significant risk levels.

ANALYTES/CRDLs

The parameters included in the analysis and the Contract Required Detection Limits (CRDLs) are listed in Attachment

ATTACHMENT 1

USEPA CONTRACT LABORATORY PROGRAM WATER QUALITY PARAMETERS (WQP) IN MULTI-CONCENTRATION WATER 6/93

PARAMETER LIST - WATER QUALITY PARAMETERS

PARAMETER	CRDL (ug/L, ppb)
Alkalinity	2000
Ammonia-Nitrogen	1000
Total Organic Carbon (TOC)	100
Chemical Oxygen Demand (COD)	3000
Chloride	2000
Nitrate/Nitrite	100
Total Phosphorous	100
Total Dissolved Solids (TDS)	10000
Total Suspended Solids (TSS)	4000
Sulfate	2000

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR INORGANICS ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	ILM04.0
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	August 1995 through May 1999
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	14 Days or 35 Days
MATRICES:	Aqueous/Soil/Sediment*

SIGNIFICANT FEATURES

- ! Analyses are suitable for aqueous, soil or sediment samples at low to medium concentration levels.
- ! Metals except mercury are analyzed by furnace AA, flame AA, and/or ICP-AES; mercury by cold vapor AA; and cyanide by spe distillation.

RECOMMENDED USES

This Contract Laboratory Program (CLP) method is recommended for broad spectrum analysis to define the nature and extent of puring Site Investigation (SI) and Remedial Investigation/Feasibility Study (RI/FS) activities and to verify that Remedial Design/Remactivities comply with pre-determined clean-up standards. This method is suitable when a 14 day or 35 day turnaround for results is a for samples from known or suspected hazardous waste sites where potential contamination may be present at significant risk levels.

* This method is not applicable to soil/sediment samples with high moisture content.

REVISIONS/MODIFICATIONS

The following is a list of significant changes from the SOW ILM03.0 that are incorporated in the SOW ILM04.0:

- ! New procedure: The lab must measure the sample shipping cooler temperature at time of sample receipt using the cooler temperat present, located in the cooler. The lab must contact CLASS if the cooler's temperature exceeds 10°C.
- ! If dissolved metals are required by the EPA Regional offices and there are no instructions on the Traffic Report, the lab must dige dissolved metals.
- ! If elements (e.g., As, Pb, Se, Tl) traditionally analyzed by GFAA are analyzed by ICP, the spiking levels for Furnace AA must be are ≤CRDL).
- ! Additional analysis frequency requirement for ICP CRDL and ICS standards: also analyzed at a frequency of not greater than 20 a analysis run.
- For analytes with CRDLs $\leq 10 \,\mu\text{g/L}$, ICP ICSA results must be within $\pm 2x$ CRDL of the true value; otherwise, these analytes must alternate method (e.g., GFAA) for samples analyzed since the last compliant ICSA.
- ! For analytes with CRDLs ≤ 10 μg/L, the ICP ICSA results must be reported from an undiluted sample analysis.
- ! Independent ICP Interference Check Samples: Elements As, Sb, Se, Tl were added (0.1, 0.6, 0.05, and 0.1 mg/L, respectively); Aş lowered (0.2 and 0.05 mg/L, respectively).
- ! Alternate methods for catastrophic ICP failure are no longer included.
- ! Cyanide spiking concentration for aqueous and soil matrix spikes: $100 \,\mu\text{g/L}$ in the distillate (i.e., the final sample solution prepare the amount of sample used.
- ! Mercury clarification of the CRDL standard requirement in the calibration curve for manual cold vapor AA; a linear regression ε automated cold vapor AA analysis.

ANALYTES/CRDLS

The parameters included in the analysis and the Contract Required Detection Limits (CRDLs) are listed in Attachment I.

ATTACHMENT I

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR INORGANICS ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION ILM04.0

TARGET ANALYTE LIST

Analyte	Aqueous CRDL (ug/L, ppb)	Soil CRDL (mg/kg, ppm)
Aluminum	200	40
Antimony	60	12
Arsenic	10	2
Barium	200	40
Beryllium	5	1
Cadmium	5	1
Calcium	5000	1000
Chromium	10	2
Cobalt	50	10
Copper	25	5
Iron	100	20
Lead	3	0.6

Analyte	Aqueous CRDL (ug/L, ppb)	Soil CRDL (mg/kg, ppm)
Magnesium	5000	1000
Manganese	15	3
Mercury	0.2	0.1
Nickel	40	8
Potassium	5000	1000
Selenium	5	1
Silver	10	2
Sodium	5000	1000
Thallium	10	2
Vanadium	50	10
Zinc	20	4
Cyanide	10	2.5

NOTE:

- ! The CRDLs for soils are based on the following: for all metals (except mercury) 1 gram sample, 200 mL final digestate volume mercury 0.2 gram sample, 100 mL final sample solution volume, 100% solids; and for cyanide 1 gram sample, 250 mL final d distillate taken for manual spectrophotometric determination, and 100% solids.
- ! The sample-specific CRDLs for soil/sediment samples will be adjusted for percent solids and will be higher than those listed abov

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION

DOCUMENT NUMBER:	OLM03.2
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	July 1994 through February 1999
CONCENTRATION:	Low to Medium
DATA TURNAROUND:	14 Days or 35 Days
MATRICES:	Aqueous/Soil/Sediment*

SIGNIFICANT FEATURES

The parameters include volatile, semivolatile, and pesticide/PCB compounds.

Volatiles and semivolatiles are analyzed by GC/MS; pesticides/PCBs are analyzed by GC/ECD.

Major Tentatively Identified Compounds (TICs) are reported for GC/MS analyses.

All pesticide/PCB analyses require second column confirmation by GC/ECD. Pesticide/PCB compounds are confirmed by GC concentration to be detected by GC/MS. (Concentrations in the sample extract at or above 10 ng/µL for pesticides, 50 ng/µL fc Toxaphene should enable the lab to confirm by GC/MS analysis.)

RECOMMENDED USES

This Contract Laboratory Program (CLP) method is recommended for broad spectrum analysis to define the nature and extent of contamination during Site Investigation (SI) and Remedial Investigation/Feasibility Study (RI/FS) activities and to verify that Reme (RD/RA) activities comply with the pre-determined clean-up standards. This method is suitable when a 14 day or 35 day turnaroun recommended for samples from known or suspected hazardous waste sites where potential contamination may be present at signific

This method is not applicable to soil/sediment samples with high moisture content.

REVISIONS/MODIFICATIONS

The following is a list of significant changes from the SOW OLM01.9 SOW that are incorporated in the SOW OLM03.1 (inclu

Volatiles

- Medium level soil/sediment extract and purge solution preparation procedures were revised for the method blank, MS/MSD, an
- Minimum sample amount for low level soil/sediment analysis was lowered to 0.5 g.

 Minimum initial and continuing calibration Relative Response Factor (RRF) for 1,1,2,2-Tetrachloroethane was lowered to 0.30 In addition to % recovery criteria, System Monitoring Compounds are evaluated based on Relative Retention Time (RRT) crite
- The concentration of Methylene Chloride in method, storage, and instrument blanks must be less than 2.5 times the CRQL.

Semivolatiles

- Low level soil CRQLs were changed from 800 ug/kg to 830 ug/kg for the following compounds: 2,4,5-Trichlorophenol, 2-Nitro Dinitrophenol, 4-Nitrophenol, 4-Nitroaniline, 4,6-Dinitro-2-methylphenol, and Pentachlorophenol.

 The minimum RRF for initial and continuing calibrations for Acenaphthylene and Acenaphthene was changed to 0.900.
- The RRT of each surrogate must be within + 0.06 RRT units of its RRT in the continuing calibration standard. The sample(s) 1 RRT criteria are not met.
- The GPC blank must contain target compounds at less than the CRQL, except phthalate esters (less than 5 times the CRQL).

Volatiles and Semivolatiles

- GC/MS Performance Check: Background subtraction must be performed using a single scan acquired no more than 20 scans pr (VOA) or DFTPP (SVOA).
- The number of volatile organic TICs was raised from 10 to 30; the number of semivolatile organic TICs was raised from 20 to and summed separately (not included in the 30). For both fractions, the data collection window is specified as 30 sec before the min after the last target analyte.
- The non-TCL library search is performed using NIST/EPA/NIH (1992 or later) and/or Wiley (1991 or later) or equivalent.

Pesticides/PCBs

- Packed columns must not be used for analysis.
- Initial and continuing calibration PEM Resolution Criteria are $\geq 90\%$; %D criteria are $\geq -25\%$ and $\leq 25\%$. Initial calibration %RSD criteria are $\leq 20\%$ except for alpha-BHC and delta-BHC which are $\leq 25\%$. PEMs, INDAs, INDBs and instrument blanks must meet initial calibration criteria.

- Surrogate advisory limits are 30-150%. Method blank surrogates must meet criteria.
- Detailed Sulfur Cleanup Blank preparation and QC acceptance/corrective action procedures were added.

Semivolatiles and Pesticides/PCBs

- Continuous L/L extraction procedures with and without hydrophobic membranes are provided for sample preparation.
- GPC calibration solution concentrations changed for bis(2-Ethylhexyl)phthalate (0.5 mg/mL) and Methoxychlor (0.1 mg/mL).

All Fractions

- The lab must measure the sample shipping cooler temperature at time of sample receipt using the cooler temperature indicator the cooler. The lab must contact CLASS if the cooler's temperature exceeds 10°C.
- Specifications for analysis of multiphase samples are included.
- Detailed requirements were added for sample collection and preservation, standard analysis and documentation, and corrective standards, mixture standards, and blanks.
- The MS/MSD are analyzed and reported at the same dilution as the least dilute sample for which the original sample results are

The following are significant changes from the OLM03.0 or OLM03.1 revision that are incorporated in the OLM03.2 revision:

- Semivolatiles: Limits number of searched alkanes to 20 suspected alkane peaks of greatest apparent concentration. Alkanes ar
- Semivolatiles: Revised for clarification if Internal Standard (IS) recoveries were outside criteria for a sample used for the MS sample extract is required only if IS recovery criteria were met in both the MS and MSD.
- Pesticide/PCB: The pesticide GPC calibration check solution contains Endrin and Dieldrin at 0.5 ug/mL for a 2 mL GPC loop; contains Aroclor 1016 and Aroclor 1260 in methylene chloride at 0.5 ug/mL for a 2 mL GPC loop.
- All fractions: Ampulated standard solution extracts in glass vials may be used until the manufacturer's expiration date, or 2 yrs manufacturer's expiration date is provided.

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in Attachment I.

2 of 2 Revision 0, 7/22/96

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM03.2

TARGET COMPOUND LIST - VOLATILES

Compound	Aqueous CRQL ug/L, ppb	Low Soil CRQL ug/kg, ppb	Medium Soil CRQL ug/kg, ppb
Chloromethane	10	10	1200
Bromomethane	10	10	1200
Vinyl Chloride	10	10	1200
Chloroethane	10	10	1200
Methylene Chloride	10	10	1200
Acetone	10	10	1200
Carbon Disulfide	10	10	1200
1,1-Dichloroethene	10	10	1200
1,1-Dichloroethane	10	10	1200
1,2-Dichloroethene (total)	10	10	1200
Chloroform	10	10	1200
1,2-Dichloroethane	10	10	1200
2-Butanone	10	10	1200
1,1,1-Trichloroethane	10	10	1200
Carbon Tetrachloride	10	10	1200
Bromodichloromethane	10	10	1200
1,2-Dichloropropane	10	10	1200
cis-1,3-Dichloropropene	10	10	1200
Trichloroethene	10	10	1200
Dibromochloromethane	10	10	1200
1,1,2-Trichloroethane	10	10	1200
Benzene	10	10	1200
trans-1,3-Dichloropropene	10	10	1200
Bromoform	10	10	1200
4-Methyl-2-Pentanone	10	10	1200
2-Hexanone	10	10	1200
Tetrachloroethene	10	10	1200
1,1,2,2-Tetrachloroethane	10	10	1200
Toluene	10	10	1200
Chlorobenzene	10	10	1200
Ethylbenzene	10	10	1200
Styrene	10	10	1200

Xylenes (Total)	10	10	1200

NOTE:

! The sample-specific CRQLs for soil/sediment samples will be adjusted for percent moisture and will be higher than those listed above.

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM03.2

TARGET COMPOUND LIST - SEMIVOLATILES

Compound	Aqueous CRQL ug/L, ppb	Low Soil CRQL ug/kg, ppb	Medium Soil CRQL ug/kg, ppb
Phenol	10	330	10000
bis(2-Chloroethyl)ether	10	330	10000
2-Chlorophenol	10	330	10000
1,3-Dichlorobenzene	10	330	10000
1,4-Dichlorobenzene	10	330	10000
1,2-Dichlorobenzene	10	330	10000
2-Methylphenol	10	330	10000
2,2'-oxybis(1-Chloropropane)	10	330	10000
4-Methylphenol	10	330	10000
N-Nitroso-di-n-propylamine	10	330	10000
Hexachloroethane	10	330	10000
Nitrobenzene	10	330	10000
Isophorone	10	330	10000
2-Nitrophenol	10	330	10000
2,4-Dimethylphenol	10	330	10000
bis(2-Chloroethoxy)methane	10	330	10000
2,4-Dichlorophenol	10	330	10000
1,2,4-Trichlorobenzene	10	330	10000
Naphthalene	10	330	10000
4-Chloroaniline	10	330	10000
Hexachlorobutadiene	10	330	10000
4-Chloro-3-methylphenol	10	330	10000
2-Methylnaphthalene	10	330	10000
Hexachlorocyclopentadiene	10	330	10000
2,4,6-Trichlorophenol	10	330	10000
2,4,5-Trichlorophenol	25	830*	25000
2-Chloronaphthalene	10	330	10000
2-Nitroaniline	25	830*	25000
Dimethylphthalate	10	330	10000
Acenaphthylene	10	330	10000
2,6-Dinitrotoluene	10	330	10000
3-Nitroaniline	25	830*	25000

NOTE:

- The sample-specific CRQLs for soil/sediment samples will be adjusted for percent moisture and will be higher than those listed above. Low level soil CRQLs were changed from 800 μ g/kg in SOW OLM01.9 to 830 μ g/kg in SOW OLM03.2.

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM03.2

TARGET COMPOUND LIST - SEMIVOLATILES

Compound	Aqueous CRQL ug/L, ppb	Low Soil CRQL ug/kg, ppb	Medium Soil CRQL ug/kg, ppb
Acenaphthene	10	330	10000
2,4-Dinitrophenol	25	830*	25000
4-Nitrophenol	25	830*	25000
Dibenzofuran	10	330	10000
2,4-Dinitrotoluene	10	330	10000
Diethylphthalate	10	330	10000
4-Chlorophenyl-phenylether	10	330	10000
Fluorene	10	330	10000
4-Nitroaniline	25	830*	25000
4,6-Dinitro-2-methylphenol	25	830*	25000
N-Nitrosodiphenylamine	10	330	10000
4-Bromophenyl-phenylether	10	330	10000
Hexachlorobenzene	10	330	10000
Pentachlorophenol	25	830*	25000
Phenanthrene	10	330	10000
Anthracene	10	330	10000
Carbazole	10	330	10000
Di-n-butylphthalate	10	330	10000
Fluoranthene	10	330	10000
Pyrene	10	330	10000
Butylbenzylphthalate	10	330	10000
3,3'-Dichlorobenzidine	10	330	10000
Benzo(a)anthracene	10	330	10000
Chrysene	10	330	10000
bis(2-Ethylhexyl)phthalate	10	330	10000
Di-n-octylphthalate	10	330	10000
Benzo(b)fluoranthene	10	330	10000
Benzo(k)fluoranthene	10	330	10000
Benzo(a)pyrene	10	330	10000
Indeno(1,2,3-cd)pyrene	10	330	10000
Dibenzo(a,h)anthracene	10	330	10000

Benzo(g,h,i)perylene	10	330	10000
----------------------	----	-----	-------

NOTE:

The sample-specific CRQLs for soil/sediment samples will be adjusted for percent moisture and will be higher than those listed above. Low level soil CRQLs were changed from 800 μ g/kg in SOW OLM01.9 to 830 μ g/kg in SOW OLM03.2.

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANIC ANALYSIS MULTI-MEDIA, MULTI-CONCENTRATION OLM03.2

TARGET COMPOUND LIST - PESTICIDES/PCBS

Compound	Aqueous CRQL ug/L, ppb	Soil CRQL ug/kg, ppb
alpha-BHC	0.050	1.7
beta-BHC	0.050	1.7
delta-BHC	0.050	1.7
gamma-BHC (Lindane)	0.050	1.7
Heptachlor	0.050	1.7
Aldrin	0.050	1.7
Heptachlor epoxide	0.050	1.7
Endosulfan I	0.050	1.7
Dieldrin	0.10	3.3
4,4'-DDE	0.10	3.3
Endrin	0.10	3.3
Endosulfan II	0.10	3.3
4,4'-DDD	0.10	3.3
Endosulfan sulfate	0.10	3.3
4,4'-DDT	0.10	3.3
Methoxychlor	0.50	17
Endrin ketone	0.10	3.3
Endrin aldehyde	0.10	3.3
alpha-Chlordane	0.050	1.7
gamma-Chlordane	0.050	1.7
Toxaphene	5.0	170
Aroclor-1016	1.0	33
Aroclor-1221	2.0	67
Aroclor-1232	1.0	33
Aroclor-1242	1.0	33
Aroclor-1248	1.0	33
Aroclor-1254	1.0	33
Aroclor-1260	1.0	33

NOTE:

The sample-specific CRQLs for soil/sediment samples will be adjusted for percent moisture and will be higher than those listed above.

TITLE: USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANICS ANALYSIS LOW CONCENTRATION WATER

DOCUMENT NUMBER:	OLC02.1
DOCUMENT DATE:	Not Applicable
EFFECTIVE DATES:	September 1996 through September 2000
CONCENTRATION:	Low
DATA TURNAROUND:	14 Days
MATRICES:	Aqueous*

SIGNIFICANT FEATURES

- ! The parameters include volatile, semivolatile, and pesticide/PCB compounds.
- ! Volatiles and semivolatiles are analyzed by GC/MS; pesticides/PCBs are analyzed by GC/ECD.
- ! All parameters have significantly reduced CRQLs as compared to SOW OLM03.2.
- ! A 25 ml aliquot of sample is purged for volatiles analysis.
- ! Semivolatile samples are extracted by continuous liquid-liquid extraction procedures; pesticide/PCB samples are extracted by separatory funnel or continuous liquid-liquid extraction procedures.
- ! Analysis of a Laboratory Control Sample (LCS) for each parameter is required.
- ! Analysis of an MS/MSD duplicate pair is <u>not</u> required.

REVISIONS/MODIFICATIONS

The following is a list of the significant changes from the 10/92 Low Concentration SOW that are incorporated in SOW OLC02.1:

- ! The lab must measure the sample shipping cooler temperature at time of sample receipt using the cooler temperature indicator blank, if present, located in the cooler. The lab must contact CLASS if the cooler's temperature indicator exceeds 10°C.
- ! The term Laboratory Evaluation Sample (LES) replaced Performance Evaluation Sample.
- ! Compound 1,2,4-Trichlorobenzene was removed from the semivolatile Target Compound List (TCL) and added to the volatile TC
- ! The number of volatile TICs was raised from 10 to 30; the number of semivolatile TICs was raised from 20 to 30. Up to 10 alkanes, which are not part of the 30 semivolatile TICs, are searched.
- ! The separatory funnel extraction method for pesticide/PCB samples was added.
- ! Control limits for recovery of volatile, semivolatile, and pesticide/PCB Laboratory Control Sample (LCS) compounds were modifi

RECOMMENDED USES

This Contract Laboratory Program (CLP) method is recommended for broad spectrum analysis to define the nature and extent of potential low level organic contamination in water supplies during Site Investigation (SI) and Remedial Investigation/Feasibility Study (RI/FS) activities and to verify that Remedial Design/Remedial Action (RD/RA) activities comply with pre-determined clean-up standards. This method attains lower detection limits than SOW OLM03.2 and can aid in the determination of low level contamination in public drinking water supplies. The majority of samples are expected to be from drinking water and well/ground water sources around Superfund sites. This method is suitable when a 14 day turnaround for results is adequate.

* This method may not be applicable for analysis of aqueous low concentration volatile organic samples when project DQOs require comparison of sample results to the drinking water Maximum Contaminant Levels (MCLs).

ANALYTES/CRQLs

The parameters included in the analysis and the Contract Required Quantitation Limits (CRQLs) are listed in Attachment 1.

ATTACHMENT 1 (Page 1 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANICS ANALYSIS LOW CONCENTRATION WATER OLC02.1

TARGET COMPOUND LIST - VOLATILES

COMPOUND	CRQL (ug/L, ppb)
Chloromethane	1
Bromomethane	1
Vinyl chloride	1
Chloroethane	1
Methylene chloride	2
Acetone	5
Carbon disulfide	1
1,1-Dichloroethene	1
1,1-Dichloroethane	1
cis-1,2-Dichloroethene	1
trans-1,2-Dichloroethene	1
Chloroform	1
1,2-Dichloroethane	1
2-Butanone	5
Bromochloromethane	1
1,1,1-Trichloroethane	1
Carbon Tetrachloride	1
Bromodichloromethane	1
1,2-Dichloropropane	1
cis-1,3-Dichloropropene	1
Trichloroethene	1

COMPOUND	CRQL (ug/L, ppb)
Dibromochloromethane	1
1,1,2-Trichloroethane	1
Benzene	1
trans-1,3-Dichloropropene	1
Bromoform	1
4-Methyl-2-pentanone	5
2-Hexanone	5
Tetrachloroethene	1
1,1,2,2-Tetrachloroethane	1
1,2-Dibromoethane	1
Toluene	1
Chlorobenzene	1
Ethylbenzene	1
Styrene	1
Xylenes (Total)	1
1,3-Dichlorobenzene	1
1,4-Dichlorobenzene	1
1,2-Dichlorobenzene	1
1,2-Dibromo-3-chloropropane	1
1,2,4-Trichlorobenzene	1

ATTACHMENT 1 (Page 2 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANICS ANALYSIS LOW CONCENTRATION WATER OLC02.1

TARGET COMPOUND LIST - SEMIVOLATILES

COMPOUND	CRQL (ug/L, ppb)
Phenol	5
bis(2-Chloroethyl)ether	5
2-Chlorophenol	5
2-Methylphenol	5
2,2'-oxybis(1-Chloropropane)	5
4-Methylphenol	5
N-Nitroso-di-n-propylamine	5
Hexachloroethane	5
Nitrobenzene	5
Isophorone	5
2-Nitrophenol	5
2,4-Dimethylphenol	5
bis(2-Chloroethoxy)methane	5
2,4-Dichlorophenol	5
Naphthalene	5
4-Chloroaniline	5
Hexachlorobutadiene	5
4-Chloro-3-methylphenol	5
2-Methylnaphthalene	5
Hexachlorocyclopentadiene	5
2,4,6-Trichlorophenol	5
2,4,5-Trichlorophenol	20
2-Chloronaphthalene	5
2-Nitroaniline	20
Dimethylphthalate	5
Acenaphthylene	5
2,6-Dinitrotoluene	5
3-Nitroaniline	20

COMPOUND	CRQL (ug/L, ppb)
4-Nitrophenol	20
Dibenzofuran	5
2,4-Dinitrotoluene	5
Diethylphthalate	5
4-Chlorophenyl-phenylether	5
Fluorene	5
4-Nitroaniline	20
4,6-Dinitro-2-methylphenol	20
N-Nitrosodiphenylamine	5
4-Bromophenyl-phenylether	5
Hexachlorobenzene	5
Pentachlorophenol	20
Phenanthrene	5
Anthracene	5
Di-n-butylphthalate	5
Fluoranthene	5
Pyrene	5
Butylbenzylphthalate	5
3,3'-Dichlorobenzidine	5
Benzo(a)anthracene	5
Chrysene	5
bis(2-Ethylhexyl)phthalate	5
Di-n-octylphthalate	5
Benzo(b)fluoranthene	5
Benzo(k)fluoranthene	5
Benzo(a)pyrene	5
Indeno(1,2,3-cd)pyrene	5
Dibenzo(a,h)anthracene	5

Acenaphthene	5
2,4-Dinitrophenol	20

Benzo(g,h,i)perylene	5
----------------------	---

ATTACHMENT 1 (Page 3 of 3)

USEPA CONTRACT LABORATORY PROGRAM STATEMENT OF WORK FOR ORGANICS ANALYSIS LOW CONCENTRATION WATER OLC02.1

TARGET COMPOUND LIST - PESTICIDES/PCBS

COMPOUND	CRQL (ug/L, ppb)
alpha-BHC	0.01
beta-BHC	0.01
delta-BHC	0.01
gamma-BHC (Lindane)	0.01
Heptachlor	0.01
Aldrin	0.01
Heptachlor epoxide	0.01
Endosulfan I	0.01
Dieldrin	0.02
4,4'-DDE	0.02
Endrin	0.02
Endosulfan II	0.02
4,4'-DDD	0.02
Endosulfan sulfate	0.02

COMPOUND	CRQL (ug/L,ppb)
4,4'-DDT	0.02
Methoxychlor	0.10
Endrin ketone	0.02
Endrin aldehyde	0.02
alpha-Chlordane	0.01
gamma-Chlordane	0.01
Toxaphene	1.0
Aroclor-1016	0.20
Aroclor-1221	0.40
Aroclor-1232	0.20
Aroclor-1242	0.20
Aroclor-1248	0.20
Aroclor-1254	0.20
Aroclor-1260	0.20



Analytical Services Available Through Superfund's Analytical Operations Branch/Contract Laboratory Program (CLP)

e Analytical Operations Branch of the Hazardous Site Evaluation Division has 18 analytical services available to rect of variety of needs for Superfund decision-making. The following table shows types of analyses available, otement of work and fact sheet numbers, dates, turnaround times, lower limit of required quantitation, and uses, nitations, average costs, and morthly laboratory capacity. For additional information about a specific service, cfuding target analytes, consult the method-specific fact sheet or call your Regional Sample Control Coordinator or LP Technical Project Officer, or call the Analytical Operations Branch on (703) 603-8870.

Type of Analysis (State) near of Work No., Feet Sheet Normber or Date)	Data Turnaround Timelai	Lower Limit of Registed Quentitetion*	August 199
Full Organics. Voletiles, Semivolatiles, Pesticides in Water and Soil IOUMO.S, 9240.0-08FS)	35 days 14 days	WATER: Volation: 10 µg/L Samheolecius: 10-25 µg/L Anadecides: 0.05-0.10 µg/L, accept methoxychlor 0.5 µg/L: Aradiare 1-2 µg/L, taxaphene 5 µg/L : SORS: Volation: Low cone, 10 µg/Kg; madium cone, 1,200 µg/Kg	USES: To determine extent of organice contamination; assess potential for the to homen health and environment; determine appropriate clean-up action; determine when remedial actions are complete; determine absence of organic contaminants. LIMITATIONS: Main survey method to determine general organics concentration. Use for other than low concentration ground and drinking water complete. For suspected high hexard samples, use high occupants on
		Semiralstins: Low cone. 320-800 μg/Kg; medium cone. 10.000-25,000 μg/Kg; medium cone. 10.000-25,000 μg/Kg. Participae: 1.7-3.3 μg/Kg. except methoxychics: 17 μg/Kg: Arcolors 33-67 μg/Kg, toxaphens 170 μg/Kg.	organic method. AVE COST: 1 790 per sample (35-day) all fractional 1.356 per sample (14-day) all fractional 2.366 per sample (14-day) all fractional 2.366 per sample per month (35-day) 700 samples per month (14-day)
.aw concentration /oistiles. Samivolatiles, Pastucides in Water (SOW Rev. 19/92, 9240.0- (DFS) (available in SAS)	14 days	WATER: Volentes: 1 μg/L (kerones 5 μg/L) Sembolectes: 5-20 μg/L Protectes: 0.01-0.02 μg/L except methoxychior 0, 1 μg/L: Aroclors 0.2-0.4 μg/L, toxaphene 1.0 μg/L	USES: To determine extent of well/ground water contem- instent; essees potential for risk to human health; determine appropriate destrup actions; determine when remedial actions are completed; determine absence of organic conteminants.
		LIMITATIONS: Use only for drinking/ground water type samples when concentrations for undiluted complex are not expected to exceed the upper limit of the cultivation curves as follows: 25 µg/L for volatiles, 30 µg/L for samivofatiles, and 0.32 µg/L for pesticidas except for Aroclars 6.4 µg/L and recepted 16 µg/L.	
	;		AVE COST: \$875 per sample (all fractions) \$135 per sample (VOA only)
			CAPACTY: 50 samples per month (all tractions) 100 samples per month (VGA only)
Okins Furans in Water, ly Asn. Sca. and Waxte 4: DFLM01 1, 09240.0- 7F5)	45 daye	<u>WATER</u> : 0.01-0.05 μg:L <u>FLY ASM</u> : 1-5 μg·Kg <u>SC</u> IL: 1-5 μg·Kg	USES: To determine extent of diskinifusen contamination: **seese potential for risk to burnen health and environment; determine appropriate clean-up actions, determine when remedial actions are complete; determine absence of diskinifuran contaminants.
	İ	<u>WASTE</u> . 10-50 بوم (q	<u>LIMITATIONS</u> : the only for specified matrices, analytes, and concentrations.
i			AVE COST: 5635 per sample
			CAPACITY: 50 samples per month

The required quantitation level is a range of quantitation levels for analytics in a service. See the last sheet for analytic and technique-specific information.

Type of Analysis (Statement of Work No., Fact Sheet Number of Dete)	©ota *Turneround *Time(a)	Lower Limit of Required Quartitetion	UseelLimitations/Gosts/Capacity
Paped Turnsround Dicking in Solids, Asphalt, Water, Air and Wipe (SOW Rev. 11/92)	T6 haurs (stactroxic) 40 hears (stactronic, sir) 7 days (hardcopy)	SOLIDS: 0.3 µgrKg ASPHALT: 0.7 µgrKg YGPE, WATER, AIR: 1.0 ng/sample	USES: To address dioxin situations requiring a quack answer or on-site feedback; identify critical semples for confirmatory analyses: optimite analysis conditions for confirmatory analyses. LIMITATIONS: Use only for semples likely to contain low fevel dioxin and when dota are needed quickly AVE COST: 4166 per sample CAPACITY: 3,500 samples per month
Quick Turnerbund Organics in Weser, Sail/Solid, Wipe (QTM SOW Ray, 2/93, Ovart Fact Sheet 7/93) Availsbie in SASI	48 hours* lelectronic) 7 days (hardcopy) * validated data	WATER: PAH: 20 µg/L Phenoles: 50 µg/L Anodors: 1-2 µg/L (toxophen= 5µg/L) Perceides: 0.1 µg/L Volatios: 20 µg/L \$0 L/50L D: PAHs: 330 µg/Kg Phenole: 830 µg/Kg Arodors: 17-33 µg/Kg (toxophene 83 µg/Kg) Perceides: 1.7 µg/Kg Volatios: 40 µg/Kg	USES. To address organics situations requiring a quick answer or on-site feedback; to direct sampling efforts, monitoring well piacament; selection of acreen intervals; to monitor descript and treatments for offectiveness. LIMITATION: Service is most affective when contaminants of concern are known or suspected and can be tootised on analyte(s) or inactional of concern. Use with the full organics service. AVE COST: \$250 per fractional analysis CAPACITY: 90 fractional analyses per day 2,160 fractional analyses per month.
Inorganics. Total Metals, Discolved Metals. Cyanide in Water and Soil (ILM03.0.9240.0-09F5)	35 days 14 days	WATER: Total metals; Disserved metals: 3-60 μg/L Hg 0.2 μg/L IA200 μg/L IGa, K, Mg. Na 5,000 μg/L Cyanida: 10 μg/L	USES: To determine extern of morganice contamination: assesse potential for risk in human health and anvironment: determine appropriate clean-up action. determine when remedial actions are complete; determine absence of inorganic contaminants. LIMITATIONS Main survey method to determine general inorganics concentration. Use for other than low concentration ground and drinking water samples for suspected high hazard samples, use high concentration inorganic method AVE COST. SIDS per sample (35-day) \$1.65 per sample (114-day) CAPACITY 4,600 samples per month (35-day) 900 samples per month (14-day)
Taxicity Characteristic Leasting Procedure (TCLP SCW Rev. 6/93) Javailable in SAS)	35 days	ORGANICS Valuation, Seminoplesian: 50 μg/L Ipentechlorophenol 325 μg/Li Perticides: 1 0 μg/L (g-BHC 0.5 μg/L). Inaxaphene 50 μg/L: METALS 100:200 μg/L (Pd. Hg 40:60 μg/L, Ba 4:000 μg/L)	USES: Colormonation of toxicity characteristic for coupesal of waste. LIMITATION Committation in its are nightly mair and dependent. AVE COST: Organica: 3675 per sample fall focusors inorganics: \$135 per sample. CAPACITY: Organics: 50 samples per month inorganics. 50 samples per month.

ι.

Type of Analysis (State ment of Work No., Fact State Number or Date)	Osta Tjameroland Tjanela)	Knew Limit of Required Openitudion	Gave Limitations/Costs)Capacity
Multi-Concentration Weter Quality Perameters ISOW Rev. 7/93, Draft Fact Sheet 6/931 [available in SAS)	14 daya	PARAMETER ALK, CL SO4 2,000 μg/L TOC, AR, TP 100 μg/L MHZ 1,000 μg/L; CO0 3,000 μg/L 7SS 4,000 μg/L; TOS 10,000 μg/L	USES: To determine water quality: determine appropriate—clean-up actions: determine when remedial actions are complete; determine presence or absence of weter quality-parameters. LIMITATIONS: Use when water samples contain low or medium lavels of up to 10 water quality parameters. Quantitation limits are highly matrix-dependent. AVE COST: \$215 per sample fall parameters! CAPACITY: 250 samples per month.
High Concentration Total Matala and Cyanida in Liquid, Solid, and Multiphase Samples (HCO1, 2, 9240,0-18FS)	3≤ daya	CIQUID, SCLID, AND MULTIPHASE: Fold charals: 5,000-80,000 µg/Kg (for Hg 300 µg/Kg) Cyardor: 1,500 µg/Kg	USES: To determine extent of high fevel contemination; determine drum contemination; seeses potential for risk to human health and environment; determine appropriate clean-up actions; determine presence or obsence of high fevels of inorganic contemination. LIMITATION: Use for determination of other than low or medium concentration inorganic conteminants. AVE COST: \$95 per sample CAPACITY: 100 samples per month:
Hight Concentration Organiza In Water - Immisciplie Liquids and Solide ISOW Rev. 9:88 and 4:89] (avadable in SAS)	35 daya	LIG <u>UIDS and SDLEDS:</u> Volation: 2.5-5.0 mg/Kg Extractables: 20:200 mg/Kg American: 10 mg/Kg Toxapheres: 50 mg/Kg	USES: To detections extent at high level contermination determine drum contamination; seekes potential for total human health and environment; determine appropriate clean-up actions, determine presence or obsence of high levels of organic conteminents. UNITATIONS: Use for determination of other than low or inedium concentration organic contaminents. ANY COST \$650 (all fractions)
Low Concentration Total Matella, Gyanida, Total Nitrogen, Fluorida in Water (7/92 Scatt), 9240.0-11FSI (avadagle in SAS)	14 day∎	WATER: Fotal metals: 1-10 μgA (Hg B 2 μg/L) (As, Ni, Zn 20μgA; Fa 100 μg/L) (Ca. Mg. Na 500 μg/L, K /50 μg/L) Cysnich: 10 μg/L Fot Aic. 100 μg/L Fluoride: 200 μg/L	USES: To determine extent of well/ground water contents attent parameters for rick in human health, determine appropriate clean-up actions: determine when remedial actions are complete, determine absence or storganic conteminants. LIMITATIONS: Use only for drinking/ground water type samples when concentrations for uncitured earnies are not expected to exceed the upper from of the calibration curves: 200 pg/L for Mg. 200-400 pg/L for openide, 1 mg/L for total nitrogen, 1000 mg/L for fluoride, and for the following analytex using flame AA - 7 mg/L for Ca. 0.5 mg/L for Mg. 2 mg/L for K, and 1 mg/L for Na. For analytes not described above, no upper calibration limit is specified. AVE COST: S215 per sample
			CAPACITY: Project-apacific requests

Type of Analysis (State) ment of Work No., Fect Sheet Number or Date)	Date Vineround Timels)	Lower Limit at Required Quartitution	Unser(Limitational Costs)Capacity
Analysis of Ambient Aid ISOW Rev. VCAA01.0, VTAA01.0.5VAA01.0 end MAA01.0, 9240.0- 15FSI levelable in SAS)	35 daye	ORGANICS. Volation: centerer 2-5 ppbv: Tenta 2-48 ng on column Sommirolistier. Perticitor: 31-193 ng/m² NORGANICS: Yotal metals: 1-109 ng/m²	USES: To determine extent of contamination; assess posential for risk in human health and environment; determine appropriate clean-up action; determine when remedial ections are complete; determine presence or absence of air contamination. LIMITATEONS: (Voluntes! For contater, insufficient flow rates and sampling times may prevent collection of sufficient quantity of sample; for Tenas, contain analytes are not continued. AVE COST: \$700 per sample (organics) 4375 per sample (metals)
e.			CAPACITY Project-specific requests
General Radiochemical Analytical Services Protocal in Water and Seil (GRASP SOW Ray, 6/1/92) (available in SASI	. 45 d∎ya	WATER Gross Alpha: 3 pCi/L Gross Beta: 2 pCi/L Metals: 0.05 pCi/L (for Ra, Sr-90 1 pCi/L: Ce-137 10 pCi/L; T 1.000 pCi/L SCIL/SEDIMENT Metals: 0.05 pCi/g (for Ra, Sr-90 0.1 pCi/g)	USES: To determine presence or absence of sadiosectivity; assess potential for risk in human health and environment; determine appropriate clean-up action; determine when remedial educate are complete LIMITATION: Main survey method to determine general radioactivity AVE COST: \$1,300 per xemple (all analytes) CAPACITY: Project-specific requests

.



Analytical Services Available Through Superfund's Analytical Operations Branch/Contract Laboratory Program (CLP)

T Unalytical Operations Branch of the Hazardous Site Evaluation Division has 16 analytical services available to test a Variety of needs for Superfund decision-making. The following table shows types of analyses available, fatement of work and fact sheet numbers, dates, turnaround times, lower limit of required quantitation, and uses, mitations, average costs, and monthly laboratory capacity. For additional information about a specific service, tluding target analytes, consult the method-specific fact sheet or call your Regional Sample Control Coordinator or Technical Project Officer, or call the Analytical Operations Branch on (703) 603-8870.

pe of Analysis (State- ent of Work No., Fact est Number of Date)	Date Turneround Time(a)	Lower Limit of Required Countrieston	Less/Lindertores/Codes/Capacity
il Organica. Valatiles, mivolaules, Pasticides Water and Soit LMQ.S., 9240,0-08FSI vailable in RAS	35 daya 14 daya	WATER: Volenker: 10 pg/L ServiceSer: 10.25 pg/L Particides: 0.05-0.10 pg/L, except methoxychlor 0.5 pg/L: Asociore 1-2 pg/L, toxephene 5 pg/L SOILS: Volenker: Low conc. 10 pg/Kg; medium conc. 1.200 pg/Kg SensholetRer: Low conc. 330-800 pg/Kg; medium conc. 10.000-25,000 pg/Kg: medium conc. 10.000-25,000 pg/Kg Particides: 1.7-3.3 pg/Kg, except methoxychlor 17 pg/Kg; Aradore 33-67 pg/Kg, toxephene 170 pg/Kg	USES: To determine extent of organics contamination; asses potential for risk to human health and envisonment; determine appropriate clean-up action; determine when remodal ections are complete; determine absence of organic contaminants. LIMITATIONS: Main survey method to determine general organics concentration. Use for other then low concentration ground and drinking water eamples. For suspected high hazard semples, use high concentration organic method. AVE COST: \$ 790 per eample [35-day; ell fractional \$1,355 per sample (14-day; all fractional \$2,400 complete per month [35-day]
w concentration :lasses, Semivolatiles, esticides in Water (SOW v, 10/92, 9240.0- :FS) sigble as MQ Multi- ant SAS	14 days	WATER: Valuation: 1 μg/L katones 5 μg/L Semination: 5-20 μg/L Pertialder: 0.01-0.02 μg/L except mathoxyclalor 0.1 μg/L; Aradions 0.2-0.4 μg/L, texephene 1.0 μg/L	DSES: To determine extent of well/ground water contentination; seases potential for risk to human health; determine appropriate clean-up actions; determine when remedial actions are completed; determine absence of organic conteminants. LIMITATIONS: Use only for drinking/ground werer type samples when concentration for undiluted samples are not expected to exceed the upper limit of the calibration curves as follows: 25 µg/L for volations, 80 µg/L for semivoration, and 0.32 µg/L for perticides except for Aroctors 5.4 µg/L and toxaphene 16 µg/L. AVE COST: 1875 per semple [all fractions] 135 per semple [VOA only] CAPACITY: 50 semples per month [all fractions] 100 samples per month (VOA only)
opine/Furene in Water, * Ash. Soil. and Waste #LM01,3, 09240.0- *F\$)	45 days	<u>WATEM</u> : 0 01-0.05 µg/L <u>FLY 48</u> H: 1-5 µg/Kg <u>SCIL</u> : 1-5 µg/Kg <u>WASTE</u> : 10-50 µg/Kg	USES: To determine extent of Gozinifuran contemination; passes potential for risk to human health and environment; determine appropriate mean-up accorns: determine when remedial actions are complete; determine absence of didainifurant contempanants. LIMITATIONS: Use only for specified metrices, analytes, and concentrations. AVE COST: 9635 per semple

⁾ The required quantitation level is a range of quantitation levels for analysis in a service. See the fact cheet for analysis and technique-specific information.

Type of Analysis (Statement of Work No Foot Sheet Humber or Date)	Data Turnground Time(a)	Lower Limit of Required Guentication	Uses/Limitations/Control/Egal-City
Repid Turnéround Dissins en Solice, Asphalt, Weter, Air and Wipe (SOW Rev. 11/92) Available in RAS	16 hours (electronic) 48 hours (electronic) 8ir] 7 days [hardcopy)	<u>şατίθ\$</u> : 0.3 μg/Kg <u>ASPHALT</u> : 0.7 μg/Kg <u>wnPE, WATER, AIR</u> : 1,0 ag/sample	USES: To address dioxin situations requiring a quick snawer or on site teachack; identify critical samples tooffrmatory enalyses; optimize enalysis conditions to confirmatory enalyses. LIMITATIONS: Use only for samples likely to contain level dioxin and when oats are needed quickly. AVE COST: 4186 per sample CAPACTLY: 3.600 samples per month
Cauck Turnsround Organics in Weter, Schl'Solid, Wipe (QTM SQW Rev. 2/93, Draft Feet Sheet 7/93) Available se HQ Multi- Client SAS	48 hours* (electronic) 7 days (hardcopy) * validated date	γνατερ: PANs: 20 μg/L Phenole: 50 μg/L Arouloid: 1-2 μg/L (toxephene Sμg/L) - Persoider: 0.1 μg/L Volution: 20 μg/L SDIL/SOLIQ: PANs: 330 μg/Kg Phenole: 37-33 μg/Kg (toxephene 83 μg/Kg) Persoider: 17-33 μg/Kg (toxephene 83 μg/Kg) Persoider: 17-37 μg/Kg Volution: 40 μg/Kg	USES: To address argenizes situations requiring a quic enswer or on-site leedback; direct exampling efforts; locate exampling areas; monitor the placement of well effect ecreen intervels; nunctor cleanups and treatmeter effectiveness. Can be used to *citor PRP activisand provide feedback at public meetings. LIMITATION: Service is most effective when contements of concern are known or suspected and can be focused on analytetal or tractional of concern. Uses the full organics service AVE COST: 1250 per tractional analysis CAPACITY: 90 tractional analyses per day 2,160 frectional analyses per month
Inorganice, Total Macels, Dissolved Metals. Cygnide in Weter and Soli (ILM03.0.9240.0-(-#\$) Available in RAS	35 days 14 days	WATER: Texal metals, Dissolved metals: 3-69 µg/L [Hg 0.2 µg/L) [Al 200 µg/L); [Ca, K. Mg. Na 5,000 µg/L] Cysnide: 10 µg/L SOIL: Total metals: 600-10 ⁴ µg/Kg (Hg 100 µg/Kg) Cysnide: 1,000 µg/Kg	USES: To determine extent of inorganics contaminate seases parential for rick in human health and environment; determine appropriate clean-up ection; determine when remedial actions are complete; determine when remedial actions are complete; determine observe of inorganic contaminants. LIMITATIONS: Mass survey method to determine gazinesgances concentration. Use for other man low concentration ground and drinking water sampless. For examples the determine gazinesgances concentration ground and drinking water sampless. For examples, use high concentration in the concentration and the concentration are made to the concentration and the concentration are made to the concentration and the concentration are provided to the concentration and the concentration are concentrated to the concent
Toxicity Characteristic Leaching Procedure (TCLP SOW Rov. 5/93) Available on HC Multi- Client SAS	35 daye	ORGANICS Volacios: Semivolatios: 50 μg/L Ipentecniorophenol and 2,4,5- suchlarophenol 125 μg/L) Perpicion: 1.0 μg/L Ig/BHC 0.5 μg/L! Incaphane 50 μg/L! METALS 100-200 μg/L (Pd. Hg/40-60 μg/L, Ba/4,000 μg/L)	USES: Determination of toxicity characteristic for displaying the second of weater and toxicity characteristic for displaying the second of weater and the second of the s

pe of Analysis (State - series) Wesk No., Feet News Number of Date!	Dete Turnercentd Timele:	Lower Emit of Required Quantitetion.	Uses Limitations)Costs/Capacity
lulti-Concentration Jarer Quality Parameters TOW Rev. 7/93, Dreft act Sheet 8/93) vailable as HQ Multi- cient SAS	14 days	<u>PARAMETER</u> ALC CL 504 2.000 μg/L TOC, ACL TP 100 μg/L NN3 1,000 μg/L; COD 3,000 μg/L TSS 4,000 μg/L; TDS 10,000 μg/L	USES: To determine water quality; determine appropriets clean-up actions: determine when remedial actions are complete; determine presence of tossade of water quality parameters. LIMITATIONS: Use when water samples contain low or medium levels of up to 10 water quality parameters. Quantization limits are highly matrix-dependent. AVE COST: \$215 per sample (all parameters) CAPACITY: 250 samples per month
igh Concentration Total state and Cyanide in quid, Solld, and luttiphaza Samples HC01.2, 3240.0-16FS) vailable in RAS	35 dey4	LIQUID SOLID, AND MULTIPHASE: Total matein: 5,000-80,000 μg/Kg Har Hg 200 μg/Kg! Connide: 1,500 μg/Kg	USES: To determine extent of high level contamination; determine drum contamination; assess potential for risk to human health and environment; determine appropriate clean-up actions; determine presence or absence of high levels of inorganic contaminants; LIMITATION: Use for determination of other than low or medium concentration inorganic contaminants. AVE COST: \$95 per temple CAPACITY: 100 semples per month
igh Concentration rgaries in Weter - smissible Uquids and osds (SDW Rev. 9/88 nd 4/89) valiable in SAS	⊋S daya	LIQUIDS and SCRIDS: Volation: 2.5-5.0 mg/Kg Entractables: 20-200 mg/Kg Aro-dors: 10 mg/Kg Zazapharas: 50 mg/Kg	USES: To determine extent of high level contemination: determine drum contemination: serves potential for risk to human health and environment; determine appropriate chan-up actions; determine presents or absence of high levels of organic conteminants. LIMITATIONS: Use first determination of other than low or medium contemination organic conteminants. AVE COST \$650 (all fractions)
ew Concentration Total letais. Cyanida. Total Progen, Physrida in Pater (7/92 Draft. 240 G-11F51 Pailable as HQ Multi- Cans SAS	14 days	<u>WATER</u> : Total anotate: 1-10 μg/L (Mg 0.2 μg/L) (As, N), In 20μg/L; Fe 100 μg/L) (Ca, Mg, Na 500 μg/L; K 750 μg/L) Cyanide: 10 μg/L Tot Aic. 100 μg/L Stanide: 200 μg/L	USES: To determine essent of well-ground water contemination: essees potential for risk in human health: determine appropriate clean-up accorat; determine when remedial actions are complete: determine absence of inorganic conteminants. LINITATIONS: Use only for denking/ground water type samples when concentrations for unpillured samples are not expected to exceed the upper limit of the coultration curves: 200 µg/L for Hg, 200-409 µg/L for cyanide. 1 mg/L for total nitrogen, 1000 mg/L for the following snalytes using flame AA · 7 mg/L for Cs. D.5 mg/L for Mg, 2 mg/L for X, and 1 mg/L for Ns. For ensystes not described above, no upper calibration limit is specified. AVE COST: 1215 per sample. GAPACITY: Project epecific requests

Type of Analysis (State- nent of Work No., Foot Sheet Number or Date)	Data Turnqround Turnq(a)	Lower Limit of Required Quantitation	Uses/Liteitetions/Ceets/Capacity
Analysis of Ambient Air SQW Rev. VCAA01.0. /TAA01.0.SVAA01.0 and MAA01.0. 9240.0- ISFS:	35 days	ORGANICS: Valenties: canacter 2-5 ppbs; Tenas 2-46 ng on column Servivolente, Participies: 37-183 ng/m² INORGANICS: Yould rootals: 1-109 ng/m²	USES: To determine extent of contamination: assess potential for risk in human health and environment; determine spropriate clean-up action; determine when remedial actions are complete; determine presence or absence of sir contamination. LIMITATIONS: (Volation) For careater, insufficient flow rates and sampling times may prevent collection of sufficient quantity of sample; for Tenax, certain analytes are not captured. AVE COST: \$700 per sample (organical \$375 per sample (metals))
General Radiochemical Analytical Services - Protocol in Water and Soil CRASP SOW Rev. 5/1/921 Qualleble in SAS	45 daye	WATER Gross Alpha: 3 pCIAL Gross Beta: 2 pCIAL Metals: 0.05 pCiA, [for Re. 5s-90 1 pCiA, Cs-137 10 pCIA; T 1.000 pCiA) SOIL:SEDIMENT Metals: 0.05 pCi/g (for Re. 9s-90 0.1 pCi/g)	USES: To determine presence or obserce of radioactivity: assess potential for risk in human health and environment; determine appropriate clean-up action; determine when remedial actions are complete LIMITATION: Main survey method to determine general redioactivity AVE COST: \$1,000 per sample (48 analyses) CAPACITY: Project-specific requests

.