NOV 8 2002

K011471

510(k) Notification 43S11 Polygram 98 pH Testing System (incl. bile)

Page 1 of 5

510(k) SUMMARY

as required per 807.92(c)

1. Submitters Name, Address:

Medtronic Functional Diagnostics A/S Tonsbakken 16-18 DK-2740 SKOVLUNDE Tel: + 45 44 57 95 02 Fax: + 45 44 57 90 10 Contact person for this submission: Tove Kjaer Date submission was prepared: May 10, 2001

2. Trade Name, Common Name and Classification Name:

A. Trade Name: Polygram 98 pH Testing System (incl. bile)

B. Common Name, Classification Name, Class and Regulation Number:

Common Name	Classificat ion Number	Class	Regulation Number
Bilitec 2000	78 FFX	П	21 CFR 876.1725
Optical Fiber Probe	78 FFX	П	21 CFR 876.1725
Polygram 98 pH Testing Application	78 FFX	Π	21 CFR 876.1725

3. Predicate Device Identification:

The functionality and intended use of the Polygram 98 pH Testing System (incl. bile) is equivalent to Medtronic Functional Diagnostics A/S's Polygram 98 pH Testing System (K 981733).

4. Device Description:

The system is an ambulatory system for use in evaluating reflux disorders in the gastrointestinal tract. The system measures bile online using sensors in the patient. The data is captured and recorded. The data is uploaded from the Bilitec 2000 by use of the Polygram 98 pH Testing Application software for later display, analysis and reporting. In its daily use, a trained technician and/or physician are the main user of the system.

The Main tasks when performing a reflux testing procedure:

- Prepare equipment including calibration
- Enter patient demographic information
- perform procedure and obtain relevant data
- Review, analysis and post procedure activities
- Create and print a report

5. Intended Use:

The Polygram 98 pH Testing System is intended to record, store, view and analyze esophageal and gastric refluxate data to diagnose reflux disorders. The Polygram 98 pH Testing System can also be used to locate the position of the proximal Lower Esophageal Sphincter (LES) manometrically, to assist in the accurate placement of the pH catheter.

6. Table of Device Similarities and differences to predicate device

Manufacturer	Medtronic Synectics AB	Medtronic Functional Diagnostics A/S	
510(k) number	Predicate Device Polygram 98' pH Testing System, i.e Digitrapper pH Polygram '98, pH Testing Application - K 981733	Modified DevicePolygram 98 pH Testing Systemv2.2 (incl. Bile) , i.eBilitec 2000Polygram 98, pH Testing Application v2.2 (including Bile)	Polygram 98 software has been updated in order to facilitate upload and analysis of data recorded by the Bilitec 2000 device

K011471

-

510(k) Notification 43S11 Polygram 98 pH Testing System (incl. bile)

Resolution

Dimension

monitoring

Display

Weight

On-line

Better than 0.01pH

display. 115 x 140 x 50 mm

Via own Screen

300 gr.

119 by 73 pixel graphic

Page 3 of 5

General:	<u>Predicate devics:</u> - all	<u>Modified Device</u> - pH System	Explanation of the differences compared to the Predicate devices
Intended Use / Indication of Use	The pH System is intended to record, store, view and analyze esophageal and gastric pH data to diagnose reflux disorders. The pH System can also be used to locate the position of the proximal Lower Esophageal Sphincter (LES) manometrically, to assist in the accurate positioning of the pH catheter	The pH System is intended to record, store, view and analyze esophageal and gastric refluxate to diagnose reflux disorders. The pH System can also be used to locate the position of the proximal Lower Esophageal Sphincter (LES) manometrically, to assist in the accurate positioning of the pH catheter	Same
Intended Populations	Pediatric to Adults	Same	
Sterilization	Accessories are not supplied sterile, manufacturer label the accessories with cleaning instructions.	Same	
Biocompatibility	Catheters are the only part that come into contact with the patients.	Same	
Technical Features:	Predicate devices - Digitrapper pH	Modified Device - Bilitec 2000	Explanation of the differences compared To the Predicate devices
Number of	4 channels pH, 1	2 channels Bile	
Channels	pressure		
Sampling rate	1/4 Hz	1 Hz with 4, 8 or 16sec averaging time	
Memory Size	2Mb Flash Data RAM	128 Kbytes flash memory	
Event markers Built in memory back up	3 Now has Flash ram (no backup required). Clock has 14 days backup.	1 event marker Has Flash ram.	Flash RAM is non-volatile therefor requires no 'backup'. The Bilitec has no clock and no back-up of date/time information is then required
Power supply	2 x 1.5V AA alkaline	4 x 1,5V AA size alkaline batteries	
Measuring range	0-9 pH	Absorbance from 0 to 1 unit	The spectrophotometric measurement of the absorbance of Bile (bilirubin) is measured in the range of 0 to 1 (0-100%)
Recording time	24 Hours	24 Hours	
	20mA	30 mA in normal mode and 16	The set-up mode and the spectrophotometric measurement
Current consumption		mA on average in Sampling/idle mode	requires different levels of current consumption
	IrDA Communication (Infra Red) Better than 0.01nH		requires different levels of current

<= 0,004 absorbance unit

200mm x 100mm x 35mm

characters

450 gr.

Same

LCD consists of 4 rows with 16

510(k) Notification 43S11 Polygram 98 pH Testing System (incl. bile)

Page 4 of 5

Features:

•

Predicate Device - Polygram '98, pH Testing Application

<u>Modified Device</u> <u>- Polygram '98, pH</u> <u>Testing Application</u> <u>v2.2 (including Bile)</u> Explanation of the differences compared to the Predicate devices

Signals to	рН	PH and Bile	
analyze			
User	Menu selections,	Same	
commands	keyboard combinations,		
	screen "buttons"		
Calculated parameters	 Maximum, Minimum Duration of period Number of acid refluxes Number of long acid refluxes Longest acid reflux Total time pH below 4 Symptom index Symptom Association Probability Number of alkaline shifts Longest alkaline shifts Longest alkaline shift Total tome pH 	 Same for pH, plus for Bile analysis: Maximum, Minimum Duration of period Number of bile refluxes Number of long bile refluxes Longest bile reflux Total time bile above 0,14 Fraction time bile above 0,14 Symptom index Symptom Association Probability 	Calculated parameters for bile are similar to the pH analysis
Scoring, Normals	 above 8 14. Fraction time pH above 8 15. DeMeester & Johnson (adult) 16. Boix-Ochoa (pediatric) 17. Infant normals percentile graph 	Same for pH. NA for the Bile analysis	No scoring system exists for Bile analysis
	(ESPGAN normals)		
Reports	Signal tracings and reports. Optional selections only.	Same	
Patient database	Relational database with logical patient- recording relations	Same	-
Additional data	User definable additional patient/recording parameters	Same	
User help system	Online help system with descriptions of procedures	Same	
Signal review method	Time - tracing based	Same	
Recording control	Real time monitoring of signals	Same	
Recording configuration	A template is used for each type of recording. User definable.	Same	

7. Assessment of non-clinical performance data for equivalence:

Verifications results show that the enhanced system performs as its predicate system.

8. Assessment of clinical performance data for equivalence:

Clinical trials have not been performed. This new system does not raise any new safety or performance issues.

9. Biocompatability:

The Optical Fiber probe has been tested for biocompatibility.

10. Sterilization:

Not applicable

11. Standards and Guidances:

The Polygram 98 pH testing System (incl bile) conforms to the following voluntary and mandatory standards:

• EN 60601-1, Medical equipment

The following guidances were followed:

- DRAERD Premarket Notification 510(k) Screening checklist, RRG Rev. 3/14/95
- ODE Guidance for the Content of Premarket Submission for Medical Device Containing Software Draft Document

DEPARTMENT OF HEALTH & HUMAN SERVICES



Public Health Service

NOV 8 2002

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Mr. Toni Kennet Jørgensen Regulatory Affairs Specialist Medtronic Functional Diagnostics A/S Tonsbakken 16-18 DK-2740 Skovlunde DENMARK Re: K011471 Trade/Device Name: Polygram '98 pH Testing System (incl Bile), Bilitec™ 2000 and Optical Fiber Probe Regulation Number: 21 CFR §876.1725 Regulation Name: Gastrointestinal motility monitoring system Regulatory Class: II Product Code: 78 FFX Dated: August 9, 2002 Received: August 12, 2002

Dear Mr. Jørgensen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html.

Sincerely yours,

Mancy C Grogdon

Nancy C. Błogdon Director, Division of Reproductive, Abdominal, and Radiological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

Indication for Use Statement

Page 1 of 1

510(k) Number (if known): <u>K011471</u>

Device Name: Polygram 98 pH Testing System (incl. bile)

Indications for Use:

The Polygram 98 pH Testing System is intended to record, store, view and analyze esophageal and gastric refluxate data to diagnose reflux disorders. The Polygram 98 pH Testing System can also be used to locate the position of the proximal Lower Esophageal Sphincter (LES) manometrically, to assist in the accurate placement of the pH catheter.

MRI Compatibility Statement:

The Ambulatory pH System is not compatible for use in a MRI magnetic field.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence o	f CDRH. Office of [Device Evaluation (ODE)
Prescription Use (Per 21 CFR 801.109)	OR	Over-The-Counter Use (Optional Format 1-2-96)
(Division Sign Division of Re and Radiologi 510(k) Numbe	eproductive, Abdomir	nal. 11471