PCB Congener Analysis of XAD-2 Resins and GFF Filters Using GC/ECD

Eric Crecelius and Lisa Lefkovitz Pacific Northwest National Laboratory Battelle Marine Sciences Laboratory 1529 West Sequim Bay Road Sequim, WA 98382

Standard Operating Procedure MSL-M-093-00

June 1994

PCB Congener Analysis of XAD-2 Resins and GFF Filters Using GC/ECD

1.0 Application and Scope

- 1.1 This SOP is applicable to the analysis of GFF filter and XAD-2 resin extracts prepared according to SOPs MSL-M-091 and MSL-M-092 for polychlorinated biphenyls (PCB) and trans-nonachlor by capillary gas chromatography with ⁶³Ni electron-capture detection.
- 1.2 This procedure provides typical gas chromatography (GC) conditions for the detection of trace levels of PCBs and trans-nonachlor, methods for identifying the analytes, and methods for analyte quantitation using the internal standard method. Attachment 1 lists the most frequently analyzed compounds. However, this list may be amended to meet requirements of specific projects.

2.0 Definitions

The following terms and acronyms may be associated with this procedure:

| ECD | Electron capture detector or detection |
|-----------|---|
| GC | Gas chromatography |
| PCB | Polychlorinated biphenyl |
| RF | Response factor |
| RRF | Relative response factor; response factor of analyte normalized to the response |
| | factor of the internal standard. |
| RSD | Relative standard deviation (%) |
| RT | Retention time |
| IS | Internal standard - compound(s) added just prior to analysis on instrument. |
| Surrogate | Compound(s) added prior to extraction to assess efficiency of method. |

In addition, it should be noted that the numbering scheme used for PCB congeners, e.g. PCB 3, is that used by Ballschmiter and Zell.

3.0 Responsible Staff

Project Manager: A Scientist responsible for 1) administration of the project; 2) providing project specific quality control requirements to the laboratory; 3) defending the data in a Quality Assurance Audit; and 4) reporting results to client.

Laboratory Supervisor: A Technical Specialist or Scientist having expertise in the principles involved with this procedure and in the use of the GC. Responsible for 1) ensuring that analysts are trained in operation of the GC; 2) appropriate quality control samples are included with the sample analysis to monitor precision and accuracy of the analysis; 3) checking the analysts' work to ensure that data are collected and interpreted correctly; 4) making decisions regarding problems with the analysis or deviations from the SOP; 5) defending the data in a Quality Assurance Audit; and 6) reporting results to project manager or client.

Analyst: A Technician, Technical Specialist, or Scientist assigned to conduct analyses using this procedure. Responsible for 1) understanding the proper use and maintenance of the GC;
2) recording information regarding instrument use and maintenance in the appropriate log books;
3) analyzing the appropriate number of quality assurance samples for each batch of samples analyzed;
4) tabulating all sample and QC data and reviewing the quality of the data based on QC guidelines presented in this SOP and any other project-specific QC guidelines;
5) reporting results to the Project Manager; and 6) defending the data during an audit.

Quality Assurance Representative: A qualified staff member assigned to the Quality Assurance Unit. Responsible for monitoring the project activities and conducting Quality Assurance Audits to ensure that 1) analysts have conducted the analysis according to the SOP and that deviations from the SOP have been noted in the appropriate log book or in the project files; 2) instrument use and maintenance records are kept correctly; and 3) data have been reported and presented accurately.

4.0 Procedures

The GC must be maintained and operated as described in Battelle SOP No. MSL-M-075.

4.1 GC Preparation

The primary quantitation GC column is a 60 meter x 0.25 mm (i.d.) fused silica capillary column coated with a 5% phenyl-, 95% methyl-polysiloxane film of 0.25 μ m thickness. (J&W Scientific, Inc., 60 meter DB-5 or equivalent). The confirmation column is a 60 m x 0.25 mm (i.d.) fused silica capillary column coated with a 14% cyanopropylphenyl-86% methyl polysiloxane film of 0.25 m thickness (J&W Scientific, Inc. 60 meter DB-1701 or equivalent). Both columns are installed in a single splitless injection port using a 2-hole ferrule.

4.2 Sample Collection, Preservation, and Handling

To conduct this analysis, the analyst should receive the samples as solvent extracts reduced to an appropriate volume (as specified in SOP MSL-M-091 or MSL-M-092). Holding times to be followed are those specified in the project specific QAPjP; normally the holding time for extracts is 40 days from date of extraction. If holding times have been exceeded, the Project Manager should be notified immediately. Refer to project-specific plans or protocols for sample collection, preservation, and handling methods.

4.3 Sample Specifications

Sample preparation methods may vary depending on the sample matrix and project needs; refer to project-specific protocols. Normally, this method is used to analyze extracts prepared according to Battelle SOPs MSL-M-091 and MSL-M-092. Samples and standards for analysis using this SOP should be prepared in hexane.

Table 1: Suggested Instrument Conditions for PCB and Chlorinated Pesticide Analysis

| Injection port type | S | Splitless | |
|----------------------------|--------|------------|--|
| Injection port temperature | 2 | 250°C | |
| Detector temperature | 3 | 300°C | |
| Initial temperature | 5 | 50°C | |
| Initial hold | 1 | .5 min | |
| Ramp 1 rate | 1 | 0°C/min | |
| Final temp 1 | 1 | l05°C | |
| Final time 1 | 0 |).0 min | |
| Ramp 2 rate | 1 | °C/min | |
| Final temp 2 | 2 | 225°C | |
| Ramp 3 rate | 1 | 0°C/min | |
| Final temp 3 | 2 | 280°C | |
| Final time 3 | 2 | 20 min | |
| Carrier gas | H | Hydrogen | |
| Carrier gas velocity | 6 | 50 cm/sec | |
| ECD make-up gas | N | Nitrogen | |
| Make-up gas flow | 3 | 30 cc/min | |
| Split vent flow | 2 | 215 cc/min | |
| Split vent on-time | 1 | 1.5 min | |
| Injection volume | 5 | 5 μL | |
| Injection speed | 1 | 0 μL/sec | |
| Hot needle time 0.0 | 07 min | | |
| Needle residence time | 0 |).2 min | |
| | | | |

Primary Column (DB-5) Conditions:

4.4 Analyte Identification

Prior to sample analysis, the elution order of the analytes of interest must be determined by analyzing the analytes individually or in combination with other analytes having known or predetermined retention times. The retention times of the analytes have all been verified on both the quantitation (DB-5) and confirmation columns (DB-1701) specified above under the GC conditions listed in Table 1, and are listed in Attachments 1 and 2. The elution order of the congeners was determined from data provided by another laboratory using similar GC conditions (Mullins et al. 1985, 1994 personal communications). Additional information on analyte identification is discussed in Section 4.6.1.

4.5 Instrument Calibration

Before the sample is injected into the GC, the detectors must be calibrated to determine the response of the detector to the analytes of interest. Demonstration of linearity of detector response is required before sample analysis. Calibration checks must be analyzed at a minimum frequency

of once every 10 samples during sample analysis.

- 4.5.1 *Initial Calibration:* Prior to initiating any sample analyses initial calibration is performed by analyzing calibration standards at a total of seven levels spanning the concentration range of 9.1 ng to 1830 ng as ?Total PCBs". Not all the congeners respond or are clearly resolved from other congeners at all these levels. On the DB-5 quantitation column, eleven congeners (BZ #s 12, 13, 100, 147 & 124, 134, 130, 129, 167, 173, and 189) are only determinable at three of these concentration levels. The calibration curves for these consist of the three concentrations plus zero. All other congeners are calibrated using a minimum of four concentrations. Calibration standards include the appropriate surrogates (Surr.) and internal standards (IS) at concentrations identical to their concentrations in the samples (see MSL-M-091 and MSL-M-092). An initial calibration must also be run if any GC conditions have changed significantly or if continuing calibration acceptance criteria (CCAC) are not met. Initial calibration acceptance criteria are a correlation coefficient (r²) of 0.95 for a minimum 4 point curve using a quadratic or linear fit.
- 4.5.2 *Continuing Calibration (CCAC):* A mid-level calibration solution is analyzed as a calibration check minimally every 10 samples while samples are being analyzed. All sample analyses must be bracketed by two calibration check standards that meet calibration criteria. Continuing calibration acceptance criteria for the primary quantitation column are a 75-125% recovery relative to the total PCB value expected.
 - 4.5.2.1 *Performance Criteria*: Additional performance criteria that will be evaluated within each CCAC include: a recovery of 50-150% for PCBs 6 & 205 (which represent small peaks), and a recovery of 75-125% for PCBs 101, 185, 44, & 180 (which represent average and large peaks). To ensure proper identification, the retention time of the internal standard reference peaks, PCB 30 and PCB 204, cannot shift by more than ± 0.4 min (see Section 4.6.1).
 - 4.5.2.2 *Internal Standard Criteria:* To ensure that the internal standards are not interfered with the area or height ratios between the internal standards PCB 30 and PCB 204 in the samples are monitored. If the area or height ratios observed in the samples differ from those observed during initial calibration by more than $\pm 15\%$ relative percent difference (RPD), all congeners must be quantitated relative to the uncontaminated internal standard and the data flagged accordingly.

4.6 Sample Analysis Procedure

Samples are analyzed under the same analytical conditions as the calibration standards. Samples must be bracketed by acceptable calibrations (see acceptance criteria in Section 4.5.2). Criteria for accepting peaks as analytes of interest are explained in Sections 4.6.1 through 4.6.4.

4.6.1 *Retention Time Windows:* Analytes are identified by the data system by setting allowable time windows for reference peak and analyte peak retention times. Initial retention times are set during the initial calibration of the instrument. From that point forward, the system utilizes retention time windows in which to "look" for peaks of interest and reference peaks. The reference peak windows are set at ± 0.4 min. and the analyte peak windows are set at ± 0.10 min. Since the internal standard reference peaks are large and clearly resolved

using a larger window for recognizing them is justified. The data system then uses the

actual retention time of the internal standard reference peaks to adjust the position of the analyte peak windows, which are much narrower, to compensate for any minor chromatographic drift. One internal standard reference peak shall also be designated an RRT reference peak to further assist in proper peak identification through the use of relative retention times.

- 4.6.2 *Second Column Confirmation*. Since each injection is split between the quantitation column (DB-5) and confirmation column (DB-1701) all sample results will be confirmed by whether or not a peak is observed at the appropriate retention time on both columns. The retention time criteria described in Section 4.6.1 apply to the confirmation column also. Trans-nonachlor can only be quantitated on the DB-1701 since it is interfered with by PCB 99 on the DB-5. As a result of this, second column confirmation of trans-nonachlor is not possible except as PCB 99 and trans-nonachlor on the DB-5. When trans-nonachlor is confirmed in this manner the PCB 99 value must be flagged accordingly.
- 4.6.3 *Resolution*. Resolution on the primary (DB-5) column should be sufficient to separate congeners #17 and #18 into two peaks with a valley less than half the height of PCB #17. If this cannot be accomplished a new column must be installed or the instrument conditions optimized further.
- 4.6.4 *Minimum Area or Height*: Peaks with a signal-to-noise ratio of three or less should be regarded as not detected unless otherwise noted in a specific project plan and/or documented by project management.

5.0 Data Analysis and Reporting

5.1 Data Recording

Data quantitation and calculations will be performed on personal computers using commercial software such as Varian Star Chromatography Software (Version 4.0 or higher), Microsoft Excel (Version 4.0 or higher) or database software (Access Version 2.0). All transfers of data to forms and data reductions (e.g., concentration calculations, means, standard deviations) will be checked by the analyst and approved by the Laboratory Supervisor. Hard copies of GC printouts of calibrations and sample data and spreadsheet reports will be kept in the GC/ECD files. A copy of the summary sheets and extraction logs will be placed in the appropriate project file in the Chemistry Group Central Files. Hard copies of chromatograms from each sample and all calibrations will be kept in the GC/ECD files unless otherwise noted in a specific project plan. All GC data files will be archived on magnetic tape.

5.2 Sample Quantitation

The internal standard method is used for quantitation. PCB 30 is used to quantify all PCB congeners with retention times up to and including PCB 110, all other congeners are quantitated vs. PCB 204. PCB 30 is used to quantitate the surrogates PCB 14 and PCB 65 and PCB 204 is used to quantitate the surrogate PCB 166. In addition, the results reported are corrected for the recovery of the surrogates PCB 65 and PCB 166; the surrogate PCB 14 is not used because it is oftentimes interfered with. The recovery of PCB 65 is used to correct the congeners quantitated vs.

PCB 30 and the recovery of PCB 166 is used to correct the congeners quantitated vs. PCB 204. If the ratio of internal standard areas exceeds the criteria set in Section 4.5.2.2, all congeners must be quantitated relative to the uncontaminated internal standard and the data must be flagged accordingly.

The concentration of a specific analyte in a sample is calculated by the Varian Star chromatography workstation. The system uses the concentration amounts for the individual analytes entered into the peak table and the results from the calibration standards to generate the coefficients of a polynomial curve fit that is, in turn, used to calculate the analytical concentrations. This result is then adjusted for internal standard and surrogate recoveries in an Excel spreadsheet:

Star result = $Amt_{(star)}$

Where:

Amt_(star) is obtained by solving:

 $Area(or \ height)_{(unk)} = A * (Amt_{(star)}^{2} + B * Amt_{(star)} + C)$

Where:

A, B & C = coefficients of the polynomial equation $Amt_{(star)} = amount of compound present in extract unadjusted for internal or surrogate$ stds $<math>Area(or height)_{(unk)} = area or height of the peak for the selected compound found in the$ analysis run.

At this point the result is not corrected by any internal standards. The results are then imported to an Excel spreadsheet using various star macros and adjusted for internal and surrogate standard recoveries.

 $Result = Amt_{(ca)} / (Sample Volume)$

Where:

 $\begin{aligned} Rec_{(istd)} &= Amt_{(istd)} / Amt \ expected_{(istd)} \\ Rec_{(surr)} &= Amt_{(surr)} / Amt \ expected_{(surr)}) / Rec_{(istd)} \\ Amt_{(adj)} &= Amt_{(star)} / Rec_{(istd)} \\ Amt_{(ca)} &= Amt_{(adj)} / Rec_{(surr)} \end{aligned}$

Where:

 $Amt_{(ca)} = Amount of compound present in analysis run$ $Amt_{(istd)} = Amount of internal standard present in extract, calculated by Varian Star$ $Amt expected_{(istd)} = Amount of internal standard added to extract$ $Amt_{(surr)} = Amount of surrogate standard present in extract, calculated by Varian Star$ $Amt expected_{(surr)} = Amount of surrogate standard added to extract$

5.3 Confirmation Data

The identification of particular PCB congeners is confirmed in all the samples by analysis on a column of dissimilar polarity (DB-1701). Confirmation is strictly by retention time only. The same criteria for retention time windows as described for the quantitation column (Section 4.6.1) are applied to the confirmation column.

5.4 Surrogate and Spike Recovery Calculation

Surrogate and Spike recoveries are calculated from the quantitation column. Calculations are as follows:

Surrogate % Recovery = $\underline{Q}_d \times 100$ Q_a

 $Q_d = Quantity determined by analysis$ $Q_a = Quantity added$

 $Matrix Spike Recovery = (\underline{SSR-SR}) \times 100$ \underline{ESR}

 $SSR = Spike \ sample \ result$ $SR = Sample \ result$ $ESR = Expected \ sample \ result = Q_a/Divisor, \ where \ Divisor = amount \ sample \ analyzed$ The relative percent difference (RPD) between replicates is calculated as follows:

 $RPD = \frac{|SR - SDR|}{1/2} x 100$ $\frac{1/2}{(SR + SDR)}$ SDR = Sample Duplicate Result

6.0 Quality Control

Some quality control considerations associated with this SOP are described in the individual sections to which they apply. The following additional quality control criteria apply unless otherwise specified in a project specific QAPjP:

| Field and Method Blanks | ≤10 ng Total PCB |
|-------------------------|---|
| Surrogate Recoveries | 40-140% (excluding PCB 14) |
| Matrix Spike Recoveries | 50 to 150% Total PCB & 70% of analytes within 50-150% range |
| Blank Spike Recoveries | 70-130% Total PCB, & 70% of analytes within 60-140% range |
| Sample Duplicate RPD | \leq 50% Total PCB and \leq 100% for analytes \geq 5X MDL |

7.0 Safety

All analysts following this procedure should be aware of routine laboratory safety concerns, including the following:

- 1. Protective clothing and eyeglasses should be worn when appropriate
- 2. Proper care must be exercised when using syringes
- 3. Certain areas of the GC system are heated. Avoid bodily contact with these areas and use care in handling flammable solvents in and around the GC system.

8.0 Training

All analysts following this procedure will be directly supervised by the Principal Investigator, qualified analyst, or laboratory supervisor until they have demonstrated to the satisfaction of the supervisor that they are capable of operating the GC independently. At a minimum, the analyst trainee should be competent in operation and maintenance of the GC (SOP No. MSL-M-075). The analyst trainee should also be able to analyze and quantify a multi-point calibration and quantitate a sample of known concentration (e.g., a reference material or matrix spike) within established control limits. Documentation of training will be recorded on training assignment and on-the-job training forms from SOP MSL-A-006. Records of this training will be kept by the laboratory Quality Assurance Representative.

9.0 References

- 9.1 MSL-A-006 Marine Sciences Laboratory Training.
- 9.2 MSL-M-075 Routine GC Maintenance.
- 9.3 MSL-M-091 Extraction and Cleanup of Resin Cartridges for Polychlorinated Biphenyls and trans Nonachalor.
- 9.4 MSL-M-092 Extraction and Cleanup of Glass Fiber Filters for Polychlorinated Biphenyls and Nonachlor.
- 9.5 Quality Assurance Plan Green Bay Mass Balance Study "Cleaning Methods for XAD-2 Resin and Filters," U.S. Environmental Protection Agency (EPA). 1986.
- 9.6 Analytical Quality Assurance Project Plan (QAPjP) for the EPA Lake Michigan PCB Mass Balance Study, DRAFT, dated October 25, 1994.
- 9.7 K Ballschmiter and Zell, "Analysis of Polychlorinated Biphenyls by Capillary Gas chromatography," Fresenius Z. Analytical Chemistry, #302 pp. 20-31 (1980)

ł

Attachment 1: Mixed Congener Standard on DB-5

(Page 1 of 5)



1

Attachment 1: Mixed Congener Standard on DB-5

(Page 2 of 5)



•

Ł

Attachment 1: Mixed Congener Standard on DB-5 (Page 3 of 5)

| Title : 7 pt pp-376 GLNPO june using height for cal Run File : C:\STAR\MODULE16\728JN4\7284013.RUN Method File : C:\STAR\MODULE16\728JN4\7284014.MTH Sample ID : PP-376C 183 | | | | | | | | |
|---|---|--|--|--|--|--|--|--|
| Injection Date: 26-MAY-95 4:51 PM | Calculation Date: 22-JUN-95 2:56 PM | | | | | | | |
| Operator : TJF Workstation: Instrument : Varian Star #1 Channel : A = 5 | Detector Type: ADCB (10 Volts) Bus Address : 16 Sample Rate : 5.00 Hz Run Time : 152.503 min | | | | | | | |

********** Star Chromatography Software ********** Version 4.01 ************

Run Mode : Analysis Peak Measurement: Peak Height Calculation Type: External Standard

| | | | Ret. | Time | | | Width | |
|------|-----------|-----------|--------|--------|----------|----------|-------|--------|
| Peak | Peak | Result | Time | Offset | Height | Sep. | 1/2 | Status |
| No. | Name | (NG/L) | (min) | (min) | (counts) | Code | (sec) | Codes |
| | | | | | | | | |
| 1 | 1 | 12.452705 | 18.578 | 0.006 | 264 | VB | 6.2 | |
| 2 | 3 | 7.748489 | 22.887 | 0.005 | 154 | BV | 6.1 | |
| 3 | 4+10 | 3.673396 | 25.507 | 0.009 | 257 | BB | 6.6 | |
| 4 | 7+9 | 1.111898 | 28.924 | 0.030 | 437 | BB | 11.9 | |
| S | 6 | 1.810442 | 30.360 | 0.010 | 269 | BB | 7.5 | |
| 6 | 8+5 | 14.648691 | 31.133 | 0.006 | 1047 | BB | 7.3 | |
| 7 | 14 (SURR) | 19.450569 | 33.155 | 0.001 | 1303 | BB | 7.4 | |
| 8 | 19 | 0.278382 | 33.893 | 0.003 | 92 | BB | 6.3 | |
| 9 | 30(ISTD) | 10.014711 | 35.261 | 0.041 | 3154 | BB | 7.3 | R |
| 10 | 12 | 0.075297 | 36.707 | 0.033 | 12 | BV | 1.8 | |
| 11 | 13 | 0.178523 | 36.975 | -0.013 | 92 | VB | 10.1 | |
| 12 | 18 | 3.918240 | 37.545 | 0.007 | 831 | BV | 7.0 | |
| 13 | 15+17 | 3.817946 | 37.760 | 0.008 | 501 | VB | 7.8 | |
| 14 | 24 | 0.052955 | 39.012 | 0.005 | .100 | BV | 7.1 | |
| 15 | 27 | 0.200533 | 39.119 | -0.011 | 94 | VB | 8.5 | |
| 16 | 16 | 2.078020 | 40.253 | 0.005 | 478 | BV | 7.3 | |
| 17 | 32 | 2.000064 | 40.409 | 0.008 | 430 | VB | 8.7 | |
| 18 | 29 | 0.051388 | 42.527 | 0.000 | 16 | VB | 5.6 | |
| 19 | 26 | 0.683344 | 43.543 | 0.019 | 174 | BV | 7.4 | |
| 20 | 25 | 0.299107 | 43.835 | 0.014 | 99 | VB | 8.1 | |
| 21 | 31+28 | 10.189798 | 44.834 | 0.014 | 1700 | VB | 12.0 | |
| 22 | 21 | 0.031160 | 46.056 | 0.023 | 10 | BV | 9.7 | |
| 23 | 33 | 3.577584 | 46.341 | 0.021 | 829 | vv | 7.8 | |
| 24 | 53 | 0.672404 | 46.558 | 0.002 | 272 | VB | 9.7 | |
| 25 | 51 | 0.179196 | 47.259 | -0.000 | 82 | BV | 6.9 | |
| 26 | 22 | 2.948125 | 47.523 | 0.020 | 637 | vv | 7.8 | |
| 27 | 45 | 0.937357 | 48.246 | 0.009 | 272 | VB | 6.8 | |
| 28 | 46 | 0.377848 | 49.529 | 0.004 | 95 | BB | 6.5 | |
| 29 | 52 | 4.857748 | 50.845 | 0.008 | 941 | BV | 7.5 | |
| 30 | 43 | 0.266150 | 51.156 | 0.001 | 73 | vv | 8.0 | |
| 31 | 49 | 2.580539 | 51.451 | 0.012 | 828 | vv | 7.4 | |
| 32 | 48+47 | 1.112766 | 51.893 | 0.006 | 505 | vv | 0.0 | |
| 33 | 65(SURR) | 4.909873 | 52.106 | 0.003 | 2146 | VB | 7.5 | |
| 34 | 44 | 4.793020 | 53.978 | 0.013 | 1277 | BV | 7.7 | |
| 35 | 42+37 | 1.452151 | 54.423 | 0.009 | 580 | VR | 9.1 | |
| 36 | 41+71 | 2.459767 | 55.776 | 0.009 | 751 | BU | 8.2 | |
| 37 | 64 | 1.899905 | 55.929 | 0.012 | 974 | VR VR | 10.2 | |
| 38 | 40 | 0.956243 | 57.084 | 0.004 | 345 | BV | 7.5 | |
| | | | | | - 15 | | | |

ŧ

.

Attachment 1: Mixed Congener Standard on DB-5 (Page 4 of 5)

| 39 | 103(ISTD) | 9.980792 | 57.584 | 0.039 | 2275 | VB | 7.6 R |
|----|-------------|-----------|---------|--------|-------------|-----------|-------|
| 40 | 100 | 0.102683 | 58.510 | 0.003 | 33 | 88 | 9.2 |
| 41 | 63 | 0.205672 | 59.395 | 0.015 | 68 | BV | 8.2 |
| 42 | 74 | 2.026874 | 59.972 | 0.021 | 767 | vv | 7.7 |
| 43 | 70+76 | 3.643651 | 60.712 | 0.018 | 1401 | vv | 8.0 |
| 44 | 66 | 5.517148 | 61.166 | 0.017 | 1361 | ٧V | 9.2 |
| 15 | 05 | 2 298095 | 61 326 | -0.000 | 719 | VB | 11.5 |
| 16 | 0) | 0 527131 | 67 313 | 0 005 | 164 | BV | 7.2 |
| 40 | 51 51 | 2 021600 | 67 706 | 0.000 | 1497 | RV | R.2 |
| 47 | 00-04 | 3.021030 | 64 600 | 0.023 | 101 | vv | 7 7 |
| 40 | 32704 | 1.050027 | 64 077 | 0.004 | 43 | VR | 6 5 |
| 49 | 09 | 1.0000110 | 04.3/3 | 0.012 | 635 | 70 | 8.5 |
| 50 | 101 | 1.932/18 | 03.314 | 0.007 | 923 | ממ ממ | 7.6 |
| 51 | 99 | 0.758792 | 00.2/4 | 0.008 | 3201 | פס | 7.0 |
| 52 | 83 | 0.143149 | 67.975 | 0.004 | 57 | DD | 7.0 |
| 53 | 97 | 0.568289 | 68.821 | 0.006 | 327 | BV | 7.9 |
| 54 | 81 | 0.184326 | 69.423 | 0.010 | 81 | VV | 1.9 |
| 55 | 87 | 1.073545 | 69.704 | 0.013 | 551 | VV | 8.3 |
| 56 | 85 | 0.709927 | 70.280 | 0.005 | 328 | vv | 8.0 |
| 57 | 136 | 0.778711 | 70.625 | 0.005 | 210 | VB | 7.4 |
| 58 | 110+77 | 2.019769 | 71.348 | 0.013 | 911 | BB | 8.4 |
| 59 | 82 | 0.469813 | 72.937 | 0.006 | 273 | BB | 7.5 |
| 60 | 151 | 1.803968 | 73.559 | 0.008 | 801 | BB | 8.0 |
| 61 | 135+144 | 0.948432 | 74.298 | -0.000 | 333 | BV | 8.1 |
| 62 | 147+124 | 0.395914 | 74.639 | 0.028 | 23 | VB | 7.6 |
| 63 | 107 | 0.124099 | 75.149 | 0.013 | 51 | BV | 7.6 |
| 64 | 123+149 | 3.043263 | 75.576 | 0.008 | 1224 | VV | 8.1 |
| 65 | 118 | 1.255572 | 75.882 | 0.020 | 597 | VB | 8.3 |
| 66 | 134 | 0.068224 | 77.027 | 0.006 | 40 | BB | 6.4 |
| 67 | 114+131 | 0 141561 | 77 481 | 0.020 | 55 | BV | 8.3 |
| 20 | 146 | 0 400207 | 70 000 | 0 016 | 204 | RR | 7.8 |
| 60 | 13341634106 | 1 502065 | 70.550 | 0 013 | 1858 | RR | 9.5 |
| 70 | 13271337103 | 1 041740 | 01 006 | 0.010 | 1285 | 22 | 8 1 |
| 70 | 141 | 1.041300 | 07.002 | 0.005 | 1205 | 100 | 7 9 |
| 11 | 13/+1/0 | 0.2/1340 | 03.030 | 0.005 | 230 | 5V 17D | 7.7 |
| 12 | 130 | 0.0/1001 | 01 220 | 0.015 | 1102 | 0.0 | 10.6 |
| 15 | 103+138 | 2.926401 | 84.228 | 0.009 | 1103 | BV | 10.0 |
| 74 | 158 | 0.252452 | 84.586 | 0.010 | 160 | VB | 8.8 |
| 75 | 129 | 0.015677 | 85.365 | -0.003 | 4.3 | VP | 17.4 |
| 76 | 178 | 1.193907 | 85.885 | 0.006 | 518 | PV | 8.1 |
| 77 | 166(SURR) | 4.808515 | 86.215 | 0.000 | 3252 | VV | 8.4 |
| 78 | 175 | 0.194360 | 86.738 | 0.002 | 100 | TS | 0.0 |
| 79 | 187+182 | 4.078599 | 87.260 | 0.001 | 2366 | VB | 8.4 |
| 80 | 183 | 1.905587 | 88.047 | 0.002 | 1101 | BV | 8.1 |
| 81 | 128 | 0.094070 | 88.478 | 0.023 | 126 | VB | 9.5 |
| 82 | 167 | 0.041623 | 89.243 | 0.014 | 15 | BV | 6.6 |
| 83 | 185 | 0.494288 | 89.556 | 0.007 | 605 | VB | 8.0 |
| 84 | 174 | 3.529921 | 90.910 | -0.001 | 1802 | BB | 8.6 |
| 85 | 177 | 1.909970 | 91.710 | 0.005 | 934 | BB | 8.0 |
| 86 | 202+171 | 0.871405 | 92.411 | 0.005 | 731 | BV | 8.0 |
| 87 | 156 | 0.065165 | 92.687 | 0.011 | 82 | VB | 8.2 |
| 88 | 173 | 0.034467 | 93.124 | -0.000 | 35 | BB | 7.0 |
| 89 | 157+200 | 0.439246 | 93.709 | -0.002 | 326 | BV | 8.4 |
| 90 | 204(1570) | 9 931970 | 93.938 | 0.050 | 5983 | VA | 8.3 R |
| 91 | 172+197 | 0 595701 | 94 772 | 0.004 | 185 | 89 | 8.7 |
| 07 | 180 | 6 010375 | 95 740 | 0.004 | 105 1761 | 00 01 | 9 C |
| 07 | 101 | 0.310323 | 06 120 | 0.012 | 1404 111 | 101 | 0.0 |
| 30 | 101 | 0.444703 | 101.00 | 0.003 | 323 | V 13 | 0.4 |
| 74 | 100 | 0.110430 | 30.//0 | 0.002 | 70 | 88 | 0.7 |
| 32 | 133 | 0.430443 | 97.250 | -0.007 | 886 | 88 | 8.1 |
| 96 | 1/0+190 | 1.942022 | 100.228 | 0.006 | 1556 | BV | 9.1 |
| 97 | DBC (SURR) | 10.149161 | 100.727 | ~0.013 | 5274 | VB | 8.4 |
| 98 | 198 | 0.125405 | 101.285 | -0.002 | 179 | TS | 0.0 |

Attachment 1: Mixed Congener Standard on DB-5 (Page 5 of 5)

| 99 | 201 | 4.806856 | 101.905 | -0.002 | 2272 | BB | 8.3 | |
|-----|---------|------------|---------|--------|--------|----|------|--|
| 100 | 203 | 2.407035 | 102.656 | -0.012 | 1611 # | BV | 8.3 | |
| 101 | 196 | 1,977705 | 102.750 | -0.013 | 1323 | VB | 10.8 | |
| 102 | 189 | 0.039184 | 105.148 | -0.020 | 26 | BB | 8.1 | |
| 103 | 208+195 | 0.896954 | 107.209 | 0.005 | 1064 | BB | 7.9 | |
| 104 | 207 | 0.093248 | 108.476 | -0.015 | 133 | BB | 7.4 | |
| 105 | 194 | 2.084478 | 110.864 | -0.002 | 2088 | PB | 8.2 | |
| 106 | 205 | 0.112808 | 111.588 | -0.003 | 112 | BB | 7.7 | |
| 107 | 206 | 0.777460 | 116.694 | -0.007 | 849 | BB | 7.9 | |
| 108 | 209 | 0.013712 | 121.396 | 0.000 | 48 | BB | 7.2 | |
| | | | | | | | | |
| | Totals: | 249.217190 | | 0.826 | 82548 | | | |

Status Codes:

ł

R - Reference peak

Total Unidentified Counts : 5777 counts

Detected Peaks: 196 Rejected Peaks: 0 Identified Peaks: 108

Amount Standard: N/A Multiplier: 1.000000 Divisor: 1.000000

Baseline Offset: -11 microVolts

Noise (used): 140 microVolts - monitored before this run

Rack: 1 Vial: 12 Injection Number: 1 Injection Volume: 5.0 ul

Original Notes:

MULLIN CAL STD W SINGLE LEVEL SURROGATES FOR SURROGATE INTERNAL STD CALIBRATION

.

Appended Notes:

Attachment 2. Mixed Congener Standard on DB-1701

(Page 1 of 5)



Attachment 2: Mixed Congener Standard on DB-1701

(Page 2 of 5)



Attachment 2: Mixed Congener Standard on DB-1701 (Page 3 of 5)

.

Title : june cal DB-1701 5 ul inj using height for cal Run File : C:\STAR\MODULE16\728JN4\7284013.RUN Method File :- C:\STAR\MODULE16\728JN4\7284014.MTH ----Sample ID : PP-376C 183 Injection Date: 26-MAY-95 4:51 PM Calculation Date: 22-JUN-95 3:05 PM Operator : TJF Detector Type: ADCB (10 Volts) Bus Address : 16 Sample Rate : 5.00 Hz Workstation: Instrument : Varian Star #1 Channel : B = 01Run Time : 152.503 min

.

Run Mode : Analysis Peak Measurement: Peak Height Calculation Type: External Standard

| | | | Ret. | Time | | | Width | |
|------|--------------|-----------|--------|--------|----------|-----------|-------|--------|
| Peak | Peak | Result | Time | Offset | Reight | Sen | 1/2 | Statue |
| No. | Name | (NG/L) | (min) | /minl | (counte) | Code | 1000) | Codee |
| | ***** | | | | | | (sec) | Codes |
| 1 | 10 | 0.176037 | 45.261 | 0.000 | 65 | RA | 6.2 | |
| 2 | 4 | 3.018902 | 45.844 | 0.014 | 50 | RR | 6.1 | |
| 3 | 7 | 0.265841 | 49.117 | 0.018 | 95 | BB | 7.8 | • |
| 4 | 9 | 0.335552 | 49.281 | 0.010 | 57 | BB | 9.7 | |
| 5 | 6 | 1.816117 | 52.143 | 0.011 | 174 | BB | 6.8 | |
| 6 | 6 | 14.159795 | 53.171 | 0.005 | 593 | BV | 7.3 | |
| _ 1 | 5 | 0.278926 | 53.431 | -0.002 | 77 | VB | 7.2 | |
| 8 | 14SURR | 5.036743 | 54.919 | 0.005 | 841 | BB | 7.5 | |
| 9 | 30ISTD | 10.156128 | 56.083 | 0.039 | 2043 | BP | 7.2 | R |
| 10 | 19 | 0.281359 | 57.367 | 0.004 | 59 | BB | 7.2 | |
| 11 | 12+13+17+18 | 6.391715 | 60.958 | 0.010 | 643 | BB | 8.7 | |
| 12 | 24 | 0.052619 | 62.559 | 0.017 | 49 | BB | 7.0 | |
| 13 | 15 | 1.490831 | 62.831 | 0.014 | 54 | BB | 8.2 | |
| - 14 | 27 | 0.203341 | 63.221 | 0.006 | 51 | BB | 7.3 | |
| 15 | 32 | 1.854917 | 64.585 | 0.004 | 223 | BB | 7.5 | |
| 16 | 16+29 | 2.099022 | 65.381 | 0.011 | 273 | BB | 7.6 | |
| 17 | 26+25 | 1.021916 | 67.942 | 0.008 | 128 | BB | 13.1 | |
| 18 | 28+31 | 10.125882 | 69.381 | 0.022 | 1314 | 8B | 9.0 | |
| 19 | 21+33+53 | 4.102441 | 71.875 | 0.018 | 542 | BB | 8.1 | |
| 20 | 51 | 0.023595 | 72.245 | -0.019 | 10 | BB | 11.4 | |
| 21 | 22 | 2.925104 | 73.676 | 0.015 | 402 | BB | 7.8 | |
| - 22 | 45 | 0.928908 | 74.149 | 0.016 | 156 | 88 | 7.6 | |
| 23 | 43+46+52 | 4.576473 | 76.102 | 0.011 | 476 | BB | 8.1 | |
| 24 | 65+47+49 | 3.705014 | 76.477 | 0.004 | 1378 | 88 | 9.5 | |
| - 25 | 48 | 0.847172 | 76.761 | 0.003 | 172 | BB | 8.0 | |
| 26 | 44 | 4.734031 | 80.575 | 0.016 | 803 | BV | 8.2 | |
| - 27 | 42 | 1.447829 | 80.786 | 0.009 | 296 | VB | 9.7 | |
| 28 | 103ISTD | 9.717182 | 82.020 | -0.002 | 1354 | BB | 7.7 | |
| 29 | 41+71+100+64 | 3.625286 | 82.314 | 0.016 | 534 | BB | 11.7 | |
| 30 | 40+63 | 1.256547 | 85.098 | 0.004 | 260 | BB | 7.5 | |
| - 11 | 74 | 2.065444 | 85.858 | 0.014 | 494 | BB | 7.9 | |
| 32 | 70 | 3.283874 | 87.639 | 0.017 | 864 | 89 | 8.0 | |
| 33 | 66+95 | 5.817988 | 87.991 | 0.017 | 911 | BB | 9.9 | |
| - 34 | 91 | 0.533801 | 88.717 | 0.009 | 90 | BB | 7.7 | |
| 35 | 60+101 | 5.634208 | 91.523 | 0.005 | 760 | BB | 8.7 | |
| 36 | 56+99 | 3.559681 | 92.018 | 0.014 | 572 | 8B | 9.8 | |
| - 37 | 89+84 | 1.446028 | 92.557 | 0.001 | 174 | BR | 7.9 | |
| 38 | 83 | 0.158861 | 95.090 | 0.004 | 46 | 8B | 8.1 | |
| | | | | | •• | | | |

- • .

Attachment 2: Mixed Congener Standard on DB-1701 (Page 4 of 5)

| 39 | 92+97 | 1.070954 | 95.977 | -0.002 | 245 | BB | 8.7 | | |
|--------|----------------------------------|------------|---------|----------|-------|----|-------|---|--|
| 40 | trans nonach | 0.006340 | 96.512 | -0.013 | 2 | BB | 2.6 | | |
| 41 | 87 | 1.052906 | 97.292 | 0.005 | 296 | BB | 7.8 | | |
| 42 | 85 | 0.695014 | 97.675 | 0.004 | 187 | BB | 7.8 | | |
| 43 | 136 | 0.726305 | 97.945 | 0.001 | 104 | BB | 7.9 | | |
| 44 | 110 | 2.009048 | 99.362 | 0.009 | 530 | BB | 8.0 | | |
| 45 | 151 | 1.796553 | 99.955 | 0.003 | 458 | BB | 7.6 | | |
| 46 | 144 | 0.630569 | 100.531 | 0.002 | 92 | BB | 7.2 | | |
| 47 | 77 | 0.228873 | 101.073 | -0.003 | 133 | BB | 8.4 | | |
| 48 | 82 | 0.464388 | 101.797 | 0.015 | 161 | BB | 7.6 | | |
| 49 | 135 | 0.156599 | 102.165 | 0.014 | 10 | BB | 8.0 | | |
| 50 | 149 | 3.115006 | 102.432 | 0.011 | 686 | BB | 8.1 | | |
| 51 | 107+123 | 0.071140 | 102.667 | 0.011 | 22 | BB | 11.3 | | |
| 52 | 118 | 1.281861 | 103.558 | 0.019 | 385 | BB | 7.8 | | |
| 53 | 114 | 0.356901 | 103.899 | -0.008 | 24 | ₿B | 8.1 | | |
| 54 | 134 | 0.076659 | 104.496 | 0.006 | 28 | BB | 8.3 | | |
| 55 | 131+137 | 0.035547 | 105.067 | 0.012 | 45 | BB | 7.7 | | |
| 56 | 146 | 0.403681 | 105.455 | 0.019 | 106 | BB | 7.9 | | |
| 57 | 153 | 2.840417 | 106.450 | 0.015 | 895 | BB | 8.3 | | |
| 58 | 132 | 0.684852 | 108.216 | 0.016 | 555 | BB | 8.3 | | |
| 59 | 141+176+105 | 2.702293 | 109.348 | -0.009 | 515 | BB | 14.4 | | |
| 60 | 129 | 0.014328 | 110.015 | 0.013 | 16 | BB | 7.1 | | |
| . 61 | 130 | 0.081114 | 111.250 | 0.007 | 20 | BB | 6.9 | | |
| 62 | 163+138+158+ | 4.417425 | 112.199 | 0.010 | 757 | BB | 10.6 | | |
| 63 | 182 | 0.697081 | 112.930 | 0.018 | 48 | BB | 8.5 | | |
| 64 | 187 | 3.900236 | 113.722 | 0.023 | 1328 | BV | 9.6 | | |
| 05 | 166515 | 4.921273 | 113.888 | -0.006 | 2041 | vv | 8.8 | | |
| 66 | 183 | 1.816817 | 114.360 | 0.019 | 569 | VB | 8.3 | | |
| 67 | 185 | 0.465152 | 116.771 | 0.018 | 302 | BB | 8.2 | | |
| 68 | 167 | 0.037056 | 117.391 | -0.003 | 9 | BB | 6.7 | | |
| 69 | 202+167+128 | 0.311932 | 117.624 | -0.006 | 191 | BB | 8.0 | | |
| 70 | 174+200+2041 | 11.663553 | 118.659 | 0.045 | 3932 | BB | 9.5 | R | |
| /1 | 177 | 1.840270 | 119.537 | 0.014 | 492 | BB | 8.3 | | |
| 12 | 171+197 | 0.363120 | 120.146 | 0.008 | 214 | BB | 8.2 | | |
| /3 | 1/3 | 0.037543 | 121.310 | 0.023 | 19 | BB | 8.2 | | |
| /4 | 156 | 0.070713 | 121.833 | 0.019 | 54 | BB | 7.5 | | |
| 15 | 1/2 | 0.576366 | 122.556 | 0.028 | 204 | BB | 7.6 | | |
| /0 | 157+180+199 | 7.206979 | 123.651 | 0.021 | 2487 | BB | 8.1 | | |
| | 193 | 0.435179 | 124.422 | 0.021 | 159 | BB | 8.2 | | |
| 78 | 191 | 0.124112 | 124.824 | 0.021 | 41 | BB | 7.6 | | |
| /9 | 198 | 0.068043 | 128.187 | 0.003 | 109 | BB | 0.0 | | |
| 80 | 1/0+190+201 | 0.939771 | 128.836 | 0.007 | 209 | BB | 6.4 | | |
| 81 | 196+203 | 4.480081 | 129.174 | -0.008 | 1660 | BB | 6.6 | | |
| 82 | 208 | 0.101339 | 130.523 | -0.013 | 206 | BB | 4.3 | | |
| 83 | 207 | 0.095886 | 131.078 | -0.015 | 148 | BB | 4.7 | | |
| 64 | 189 | 0.036086 | 131.360 | -0.016 | 47 | BB | 0.0 | | |
| 60 | 195 | 0.703326 | 131.705 | -0.016 | 1035 | BB | 3.5 | | |
| 00 | 194 | 1.901737 | 133.084 | -0.020 | 2564 | BB | 3.2 | | |
| 87 | 205 | 0.093067 | 133.394 | -0.020 | 146 | BB | 0.0 | | |
| 88 | 206 | 0.724105 | 134.676 | -0.020 | 995 | BB | 3.6 | | |
| 63 | 20,9 | 0.010731 | 135.533 | -0.019 | 55 | BB | 3.3 | | |
| | | | | | | | ~~~~~ | | |
| | TOCATA: | 181-095437 | | 0.651 | 42619 | | | | |
| State | us Codect | | | | | | | | |
| R - 1 | us coues: Reference | 4 | | | | | | | |
| 41 - 1 | wererence bea | ~ | | | | | | | |
| Tota | Total Unidentified Counter: 6147 | | | | | | | | |
| | | a counca . | 01 | a counts | | - | | | |
| | | | | | - | | | | |

Detected Peaks: 293

Rejected Peaks: 6

Identified Peaks: 89

Attachment 2: Mixed Congener Standrd on DB-1701 (Page 5 of 5)

 Amount Standard: N/A
 Multiplier: 1.000000
 Divisor: 1.000000

 Baseline Offset: -14 microVolts

 Noise (used): 50 microVolts - monitored before this run

 Rack: 1
 Vial: 12

 Injection Number: 1
 Injection Volume: 5.0 ul