

ANTIMICROBIAL RESISTANCE TESTING

Salmonella, generic *E. coli* and *Enterococcus* antimicrobial susceptibility testing are conducted using the Sensititre™ System (Trek Diagnostics, Inc., Westlake, Ohio) as per manufacturers directions. Antimicrobials include those used in both human and veterinary medicine and are configured in a 96 well custom made panel. National Committee for Clinical Standards (NCCLS) guidelines are followed throughout the testing procedure. Quality control strains (*Escherichia coli* ATCC 25922, *Pseudomonas aeruginosa* ATCC 27853, and *Enterococcus faecalis* 29212) are included in susceptibility testing on a weekly basis.

Campylobacter isolates are evaluated for susceptibility to 8 antimicrobials using the E-Test (AB Biodisk, Piscataway, NJ). Briefly, E-test strips are placed onto 150 mm Mueller Hinton plus 5% lysed horse blood plates (B-D Biosciences) previously swabbed with a suspension of the *Campylobacter* test isolate (density equivalent of 0.5 McFarland standard) prepared from a fresh culture. The plates are incubated for 48 h at 42° C under microaerobic conditions (5% O₂, 10% CO₂ and 85% N₂). Minimum inhibitory concentrations of the antimicrobials are determined by reading the numeric concentration (µg/ml) printed on the strip at the point of growth inhibition. Quality control strains (*Escherichia coli* ATCC 25922, *Staphylococcus aureus* ATCC 25923 and *Campylobacter jejuni* ATCC 33560) are included in susceptibility testing on a weekly basis.

References

- Barrett, T.J., H. Lior, J.H. Green, R. Khakhria, J.G. Wells, B.P. Bell, K.D. Greene, J. Lewis, and P.M. Griffin. 1994. Laboratory investigation of a multistate food-borne outbreak of *Escherichia coli* O157:H7 by using pulsed-field gel electrophoresis and phage typing. *J Clin Microbiol.* 32(12):3013-3017.
- Jackson, C. R., P. J. Fedorka-Cray, and J. B. Barrett. 2004. Use of a genus-and species-specific multiplex PCR for identification of enterococci. *J. Clin. Microbiol.* 42(8):3558-65
- Stern, N.J., B. Wojton, and K. Kwiatek. 1992. A differential-selective medium and dry ice-generated atmosphere for recovery of *Campylobacter jejuni*. *J. Food Prot.* 55:514-517.