

(B) For the assay to be valid:

(1) The positive control sera must give a result within one dilution of the previously determined titer.

(2) The negative control sera must be negative.

(3) The backtitration of the antigen must be 1:4 or 1:8.

(4) The RBC control must give tight, non-hemolyzed buttons.

(5) Sera controls (well A of each test sera) must not have non-specific agglutination or hemolysis. If negative, report as "negative with non-specific agglutination or non-specific hemolysis" or "unable to evaluate due to non-specific agglutination or hemolysis" or treat the serum to remove the non-specific agglutination and repeat the test. (See paragraph (e)(2)(v) of this section.)

(v) *Treatment to remove non-specific agglutination.* (A) *Purpose.* Treatment of serum to remove non-specific agglutination that is interfering with HI assays.

(B) *Specimen.* Serum.

(C) *Materials.* Homologous RBC's (chicken or turkey), 50 percent solution PBS, centrifuge, incubator, 4C (refrigerator).

(D) *Procedure.* (1) Prepare a 1:5 dilution of test serum by adding 50 µL of serum to 200 µL of PBS.

(2) Prepare a 50 percent solution of RBC's by adding equal volumes of packed RBC's to PBS. Mix well.

(3) Add 25 µL of 50 percent RBC solution to the serum dilutions.

(4) Vortex gently to mix.

(5) Incubate at 4 °C for 1 hour.

(6) Centrifuge to pellet the RBC's.

(7) Use the supernatant to perform the HI assay. Modify the dilution scheme in the assay to consider the initial 1:5 dilution prepared in the treatment. For the 1:5 dilution scheme, do not add PBS to row A. Add 50 µL of the 1:5 treated supernatant to row A. Serially dilute 25 µL from rows A through H. This prepares a serum dilution of 1:10 through 1:640 in rows B through H.

[49 FR 19803, May 10, 1984, as amended at 57 FR 57342, Dec. 4, 1992; 59 FR 12799, Mar. 18, 1994; 63 FR 3, Jan. 2, 1998; 67 FR 8469, Feb. 25, 2002]

§ 147.8 Procedures for preparing egg yolk samples for diagnostic tests.

The following testing provisions may be used for retaining the classification U.S. M. Gallisepticum Clean under § 145.23(c)(1)(ii)(C) and § 145.33(c)(1)(ii)(C), and for retaining the classification U.S. M. Synoviae Clean under § 145.23(e)(1)(ii)(b) and § 145.33(e)(1)(ii)(b) of this chapter.

(a) Under the supervision of an Authorized Agent or State Inspector, the eggs which are used in egg yolk testing must be selected from the premises where the breeding flock is located, must include a representative sample of 30 eggs collected from a single day's production from the flock, must be identified as to flock of origin and pen, and must be delivered to an authorized laboratory for preparation for diagnostic testing.

(b) The authorized laboratory must identify each egg as to the breeding flock and pen from which it originated, and maintain this identity through each of the following:

(1) Crack the egg on the round end with a blunt instrument.

(2) Place the contents of the egg in an open dish (or a receptacle to expose the yolk) and prick the yolk with a needle.

(3) Using a 1 ml syringe without a needle, aspirate 0.5 ml of egg yolk from the opening in the yolk.

(4) Dispense the yolk material in a tube. Aspirate and dispense 0.5 ml of PBS (phosphate-buffered saline) into the same tube, and place in a rack.

(5) After all the eggs are sampled, place the rack of tubes on a vortex shaker for 30 seconds.

(6) Centrifuge the samples at 2500 RPM (1000×g) for 30 minutes.

(7) Test the resultant supernatant for *M. gallisepticum* and *M. synoviae* by using test procedures specified for detecting IgG antibodies set forth for testing serum in § 147.7 (for these tests the resultant supernatant would be substituted for serum); except that a single 1:20 dilution hemagglutination inhibition (HI) test may be used as a screening test in accordance with the procedures set forth in § 147.7.

NOTE: For evaluating the test results of any egg yolk test, it should be remembered

§ 147.9

9 CFR Ch. I (1–1–05 Edition)

that a 1:2 dilution of the yolk in saline was made of the original specimen.

[50 FR 19900, May 13, 1985; 63 FR 3, Jan. 2, 1998]

§ 147.9 Standard test procedures for avian influenza.

(a) The agar gel immunodiffusion (AGID) test should be considered the basic screening test for antibodies to Type A influenza viruses. The AGID test is used to detect circulating antibodies to Type A influenza group-specific antigens, namely the ribonucleoprotein (RNP) and matrix (M) proteins. Therefore, this test will detect antibodies to all influenza A viruses, regardless of subtype. The AGID test can also be used as a group-specific test to identify isolates as Type A influenza viruses. The method used is similar to that described by Beard.⁶ The basis for the AGID test is the concurrent migration of antigen and antibodies toward each other through an agar gel matrix. When the antigen and specific antibodies come in contact, they combine to form a precipitate that is trapped in the gel matrix and produces a visible line. The precipitin line forms where the concentration of antigen and antibodies is optimum. Differences in the relative concentration of the antigen or antibodies will shift the location of the line towards the well with the lowest concentration or result in the absence of a precipitin line. Electrolyte concentration, pH, temperature, and other variables also affect precipitate formation.

(1) *Materials needed.*

- (i) Refrigerator (4 °C).
- (ii) Freezer (–20 °C).
- (iii) Incubator or airtight container for room temperature (approximately 25 °C) incubations.
- (iv) Autoclave.
- (v) Hot plate/stirrer and magnetic stir bar (optional).
- (vi) Vacuum pump.
- (vii) Microscope illuminator or other appropriate light source for viewing results.
- (viii) Immunodiffusion template cutter, seven-well pattern (a center well

surrounded by six evenly spaced wells). Wells are 5.3 mm in diameter and 2.4 mm apart.

(ix) Top loading balance (capable of measuring 0.1 gm differences).

(x) Pipetting device capable of delivering 50µl portions.

(xi) Common laboratory supplies and glassware—Erlenmeyer flasks, graduated cylinders, pipettes, 100 × 15 mm or 60 × 15 mm petri dishes, flexible vacuum tubing, side-arm flask (500 mL or larger), and a 12- or 14-gauge blunt-ended cannula.

(2) *Reagents needed.*

(i) Phosphate buffered saline (PBS), 0.01M, pH 7.2 (NVSL media #30054 or equivalent).

(ii) Agarose (Type II Medium grade, Sigma Chemical Co. Cat.# A-6877 or equivalent).

(iii) Avian influenza AGID antigen and positive control antiserum approved by the Department and the Official State Agency.

(iv) Strong positive, weak positive, and negative control antisera approved by the Department and the Official State Agency (negative control antisera optional).

(3) *Preparing the avian influenza AGID agar.* (i) Weigh 9 gm of agarose and 80 gm of NaCl and add to 1 liter of PBS (0.01 M, pH 7.2) in a 2 liter Erlenmeyer flask.

(ii) To mix the agar, either:

(A) Autoclave the mixture for 10 minutes and mix the contents by swirling after removing from the autoclave to ensure a homogeneous mixture of ingredients; or

(B) Dissolve the mixture by bringing to a boil on a hot plate using a magnetic stir bar to mix the contents in the flask while heating. After boiling, allow the agar to cool at room temperature (approximately 25 °C) for 10 to 15 minutes before dispensing into petri plates.

(iii) Agar can be dispensed into small quantities (daily working volumes) and stored in airtight containers at 4 °C for several weeks, and melted and dispensed into plates as needed.

NOTE: Do not use agar if microbial contamination or precipitate is observed.

(4) *Performing the AGID—(i) Detection of serum antibodies.* (A) Dispense 15 to 17 mL of melted agar into a 100 × 15 mm

⁶Beard, C.W. Demonstration of type-specific influenza antibody in mammalian and avian sera by immunodiffusion. Bull. Wld. Hlth. Org. 42:779–785. 1970.