MYNIAID

Sprint 1 Wireframes - Change Request

April 23, 2015



Design Tenets

Three design tenets will guide the design of MyNIAID


Scannable

Each view of MyNIAID should be scannable so that the user can quickly ascertain what the most important information on the page is.




Responsive

We are designing with responsiveness in mind, so that MyNIAID can be used seamlessly across devices.





Quick Actions

User should be able to quickly and easily take action on items.



MyNIAID — Web View





Web - MyNIAID - Compliance

NIAID-WF-352

User has clicked on Compliance filter.

USER NAM

User's name appears. On click, menu opens up for user to go to Subscriptions and Important Contacts.


FILTER BEHAVIOR

By default, all types, sources, and priorities are displayed. A user can choose to filter the list of items by clicking on a row within each facet. When a user clicks on a row to filter by, the row is highlighted. The user can only choose one within each bucket.

When filtering by type, the label updates to indicate the number and type of items that are being displayed.

The number next to the un-selected filter is removed, since none of those items are being displayed.


The numbers in adjacent facets are updated once a filter has been chosen.



Web - MyNIAID - Notifications


NIAID-WF-353

User has clicked on Notifications filter.



vorites App Catalog User Name ▼

Last updated 4:55pm C



You have not snoozed

any notification or

compliance items.


Web - MyNIAID - Snoozed -

None NIAID-WF-354

User has not snoozed any items and clicks on the Snoozed tab.

NO SNOOZES

No number is displayed if the user has not snoozed any items. A message is displayed that indicates that the user has not snoozed any items.



Search within...

Compliance

Notifications

Type

All

Source

All

Due Date

All

Overdue

Due Soon

Upcoming

Not Provided

EDRMS

Source Application 2

Source Application 3

Q

12

33

35

6

19

9

15


2

My Notifications

Favorites

App Catalog

User Name ▼



Web - MyNIAID - All Active - Sort Menu NIAID-WF-355

User has clicks on the Sort By menu and reveals their options for sorting.

TYPE

Order things should be displayed if by type:

• Type - Compliance, Notifications

Within each type, display in the following order:

- Due Date (closest to farthest)
- Source (alphabetical)Title (alphabetical)
- TRY IT: Click to view items sorted by priority

SOURCE

If user sorts by source, items appear alphabetical order by source. If there are multiple items from a single source, the items should be sorted in following order:

- Type (Compliance, Notifications)
- Due Date (closest to farthest)
- Title (alphabetical)


DUE DATE

`If user sorts by priority, items are displayed in the following order:

- Overdue
- Due Soon
- Upcoming
- Not Provided

If there are multiple items for a due date, the items should then be sorted by:

- Type (Compliance, Notifications)
- Source (alphabetical)
- Title (alphabetical)




Food Allergies have been lowered in the last 2 years

3

Web - MyNIAID - All Active - Sorted by Due Date NIAID-WF-356

User has sorted items by due date.

TRY IT: Click to filter by due date



App Catalog

Favorites

User Name ▼

MyNIAID

Web - MyNIAID - All Active - Filtered by Due Date

NIAID-WF-357

User has filtered items by due date.

ITEM COUNT

If user filters by due date or source, the item count updates to show X of X items are displayed.


ADJACENT FACETS

When a user filters in one facet, the counts in the other facets update. A user can filter across multiple facets at the same time (e.g. filter by type AND due date).

TRY IT: Click to add filters back in


DEFAULT


By default, all items within a filter type are displayed.



FILTER SELECTED

When user selects a facet, the state changes to indicate the facet that is selected. The facets that were not selected change state to indicate that they are not being displayed. The numbers should be displayed whether the facet is active or not. User can only select one at a time. If there are zero items in a facet, the item should be displayed with a 0 next to it.





Web - MyNIAID - All Active - Add Snooze - Selecting NIAID-WF-358

User clicked on snooze icon.


SNOOZE OPTIONS

User chooses a date to snooze until from a date picker. Important states:

- Dates in past Not selectable
- Dates in future Selectable
- Selected Date should be highlighted

User clicks arrows to move across months.

TRY IT: Click to snooze item until next week



Web - MyNIAID -All Active -Add Snooze - Snoozed Undo

NIAID-WF-359

User snoozed an item.

NUMBER OF ITEMS

After snoozing an item, the active items number is updated to reflect that an item has been snoozed. A number also appears next to the Snoozed tab item.

TRY IT: Click to dismiss message prior to it disappearing by itself

TRY IT: Click to undo snooze

SNOOZE MESSAGE


After a user chooses how long they want to snooze, the item disappears and a message appears in its place. The snoozed message indicates the date (e.g. Mar 13) that the item will be snoozed until.

UNDO

On click, undo removes the message and the item that was snoozed reappears in the Active list.

CLOSE


On click, close dismisses the message. The item remains in the Snoozed tab. If the user does not dismiss the message after 3 seconds, it will automatically disappear.



Web - MyNIAID - All Active - Snooze Added NIAID-WF-360

User dismissed snooze message or snooze message disappeared after 3 seconds.

TRY IT: Click to view snoozed items



MyNIAID

Web - MyNIAID - Snoozed

NIAID-WF-361

User Name ▼

App Catalog

Favorites

User clicked on Snoozed tab and there are items that have been snoozed.

NUMBER OF ITEMS

Number of snoozed items.


FACETS

Item numbers are updated to reflect number that are snoozed.

TRY IT: Click to edit snoozed item

ACTION BUTTON

If an item is snoozed, the snoozed icon is given a different style than the act and dismiss buttons.



MyNIAID

Web - MyNIAID - Snoozed - Edit Snooze NIAID-WF-362

User clicks on Snooze icon on an already snoozed item.

EDIT SNOOZE & UNSNOOZE


User Name ▼

App Catalog

Favorites

If an item has already been snoozed, when the user clicks on the snooze icon, the calendar appears. The calendar has the previously selected snooze date selected. Below the calendar, "Unsnooze" appears. On click of "Unsnooze," the item is no longer snoozed and moves back to the Active list of items.

TRY IT: Click to dismiss the snooze



MyNIAID

Web - MyNIAID - Snoozed - Undo Edit Snooze NIAID-WF-363

User dismissed snooze.

UPDATED ITEM NUMBERS

User Name ▼

App Catalog

Favorites

-- The item numbers update if a user un-snoozes an item.

TRY IT: Click to dismiss message

TRY IT: Click to undo dismissing an item

UN-SNOOZE MESSAGE


If the user dismisses the snooze, the snooze disappears from the Snoozed tab and moves to the Active tab. In its place, a message appears.

UNDO UN-SNOOZE

On click of UNDO, the item that was un-snoozed is added back to the Snoozed list.


CLOSE

The user can close the message. If the user does not close the message, it disappears after 3 seconds.



Web - MyNIAID - Snoozed - Snoozed Undone NIAID-WF-364

User dismissed snooze - message closed or disappeared.



Web - MyNIAID - Dismiss Item

NIAID-WF-365

User dismissed an item.

UPDATED ITEM NUMBERS

==The item numbers update if a user dismisses an item.

DISMISS MESSAGE

If the user dismisses an item, a message is displayed that indicates that the action was taken.

UNDO UN-SNOOZE


On click of UNDO, the item that was dismissed is added back to the

CLOSE

The user can close the message. If the user does not close the message, it disappears after 3 seconds.


TRY IT: Click to dismiss message

TRY IT: Click to undo dismissing an item



Web - MyNIAID - After Dismissing Item NIAID-WF-366

User dismissed item - message closed or disappeared.



Web - MyNIAID - Refresh Message NIAID-WF-367

User has not refreshed the page in ten minutes.

--_NO

If user clicks no, list of items is not refreshed.

YES

If user clicks yes, list of items is refreshed.

Thank you.

