# Hard X-ray tomography on the nanometer length scale

#### C. Rau

University of Illinois at Urbana-Champaign, Urbana, IL 61801,USA UNICAT/Advanced Photon Source, ANL, Argonne, IL 60439, USA

- Introduction
- Microscope at 34 ID-C
- Limits with current equipment
- Prospects for the future
  - A microscope with 50nm resolution Summary

## **Experimental set-up**

#### **Detector system Tomography stage**



-Tomography and in-line phase contrast imaging

-High quality stages Rotation Stage with air bearing (run out < 20nm)



Sample
 Translation stages
 Rotation stage

### Scientific application (coll. T. Terry, M. Harris, Purdue University)

#### Projection image



#### Reconstructed slice



- SiO<sub>2</sub> and ZrO<sub>2</sub> spheres
- Hollow spheres with high relative surface
- biomedical applications; chemical reactions
- Size distribution in function of the preparation procedure
- In-line phase contrast
- Features over several length scales; wall thickness
- More detailed study  $\rightarrow$  microscope

#### **Scientific application**

(coll. C.-P. Richter, Northwestern University)



http://www.iurc.montp.inserm.fr/cric51/audition/en glish/cochlea/fcochlea.htm

- -Better understanding of hearing process
- in function of frequency, diff. parts stimulated
- Dynamic study
- Again: features on different length scales
- Unique approach to the problem
- Challenging project

## **Magnified Imaging**

#### **Resolution mainly limited by detector**

#### How to obtain higher spatial resolution?

- lensless imaging (-> talk C. Jacobsen)
- imaging using hard X-ray optics
  - Scanning techniques with focused beam
  - Fan beam geometry
  - Full-field imaging

## **Magnified Imaging**

#### Hard X-ray optics Compound Refractive Lens



#### **Planar Lens**



#### **Fresnel-Zone Plate**



#### **Kirkpatrick-Baez Mirror**



## **Magnified Imaging**



#### Preface

## Experiments in the past with impact for current development

## Hard X-ray full-field microscope

#### **Experiments at ESRF ID 22**

#### **Monochromatic**



-beamline nicely adapted (lenses with f~1-2m) - $I_{pink}$ ~50xI<sub>mono</sub>;  $\Delta E/E=10^{-2}(pink);10^{-4}(mono)$ 

#### Parabolic Compound Refractive Lens (C.G. Schroer, B. Lengeler, RWTH Aachen)

#### Set of parabolic lenses (typ. 50-200)



Focal length f can be varied with energy E and number of lenses N

f=2 m with 50 lenses at 20 keV
effective aperture ~ 0.2 mm
diffr. limit 0.3 micron
transmission 3 %



R: 0.2 mm R<sub>0</sub>: 0.45 mm d: 0.01 mm material: Al

## X-ray magnified tomography

#### <u>Sample</u> Boron fiber with tungsten core



Combined absorption-contrast (tungsten) and phase-contrast (boron) object

## Tomography: phase contrast and pink beam







3.4x mono 3.4x pink

10x pink

Al parabolic CRL; E=19.7 keV

- no significant difference pink mono (3.4x)
- resolution  $10x\sim 1\mu m$
- study of mechanism of image formation
- beryllium lenses overcome actual limit

Ref.: C. Rau et al., NIM A, (2001).

T. Weitkamp et al., SPIE proc., (2002).

C. Schroer et al., SPIE proc., (2002).

### **KB-FZP** microscope

C. Rau, U. Neuhäusler, G. Schneider, O. Hignette, O. Leupold



<u>Condenser</u> (O. Hignette, G. Rostaing, ESRF) **KB multilayer mirror**: high reflectivity, matches Objective aperture

Objective(M. Panitz, IRP Göttingen):micro-FZP:high resolution, reasonable efficiency

#### **Kirkpatrik-Baez System** (O. Hignette, A. Rostaing, ESRF Grenoble)

Design of the system<sup>1</sup>



Ni-Be<sub>4</sub>C multilayer d=24 Å

**f**<sub>hor</sub>=**0.9 m f**<sub>vert</sub> =**0.6 m** 

**Reflectivity 33 %** 

<sup>1</sup> courtesy O. Hignette

#### Micro Fresnel-Zone plate

(M. Panitz, G. Schneider, M. Peuker, G. Schmahl: Uni Götttingen) Scheme<sup>1</sup>



<sup>1</sup>courtesy U. Neuhäusler

Au **65** μ**m 70 nm** 85 nm **Thickness** 500 nm Efficiency 2 % (10% at 13 keV at 4 keV) **2x10**<sup>4</sup> 4.7 cm f at 13keV

focal length f=2·r<sub>n</sub>·dr<sub>n</sub>/ $\lambda$ resolution  $\delta$ = 0.61·  $\lambda$  /N.A.=1.22·dr<sub>n</sub>

## **High-resolution imaging**

Imaging with sub-100nm lateral resolution with a combined KB<sup>1</sup>-FZP<sup>2</sup> microscope

KB system as a condenser: high flux
micro-FZP: high resolution 85 nm
exposure time: < 1min, E= 13 keV</li>

Significant improvement

- Nano-tomography
- Zernike phase-contrast
- <sup>1</sup> O. Hignette, ESRF
   <sup>2</sup> M. Panitz, IRP Göttingen



Why are these experiments important?

- Proof of principle
- Tomography possible
- Use of pink beam
- Magnified Imaging of phase contrast objects
- High resolution with FZPs



## *The Present* Realization of the microscope

## Microscope at UNICAT/APS

<u>UNICAT</u>: 4 different partners (UIUC, ORNL, NIST, UOP) 3 beamlines in sector 33&34 Methods: USAXS, XRD, Spectroscopy, ...

Tomography setup at UNICAT/APS

- Tomography
- KB-FZP microscope

Nano-Tomography available for large scientific community

Sub-100nm tomography

development Instrument – valuable Experiments

## Layout Microscope 34 ID

- Hard X-rays
  - Penetration depth , sample handling, depth of focus
- Resolution <100nm
  - Use KB-FZP system
- -Applications in Materials Science and Biomedicine

- Flexible design
- Evolutionary improvement of setup
- Technical development accompanied by scientific experiments

## Microscope at 34 ID

- Energy range 6-12 keV (currently 9keV)
  - Depth of focus of FZP: 50-100µm
  - $\rightarrow 100 \mu m$ -diameter samples
  - $\rightarrow$  fits to Pixel array of CCD
  - Elements like Cu, Zr, Au, Ag can be studied using absorption contrast
  - Biomaterials have to be studied by making use of phase contrast
  - Detector efficient at this energy

#### Microscope at 34 ID-C



-Experiments in "parasitic mode" (E-Station main user)
-> amount of beamtime

### Kirkpatrik-Baez System (P. Eng, CARS, University of Chicago/J. Qian, R. Conley,L. Assoufid, APS)

#### Photo of the system<sup>1</sup>



**Pt-coated Si crystal** 

**Reflectivity >80 %** (theoretical)

**Roughness**  $< 1 \text{\AA}$ 

<sup>1</sup> courtesy P. Eng, CARS, U. of Chicago

## Microscope at 34 ID



## Microscope at 34 ID

#### Profile of an GaAs-edge 20x X-ray magnification



Measured Resolution

165 nm 20x X-ray magnification135 nm 30x X-ray magnification... nm 40x X-ray magnification

#### **Phase II-IV**

*The Near and far future* Next developments and potential of microscope

## Microscope at 34 ID

#### Limits of current setup

Resolution limit:

 $\begin{array}{ll} \underline{Higher \ order \ of \ ZP} \\ Resolution & \delta = 1.22 dr_n / \mathbf{m} & m: order \\ Efficiency & \eta \sim 1/m^2 & m \ is \ impair \end{array}$ 

<u>Finer outer zones</u> dr<sub>n</sub> of 40nm

(Phase II)

Reaching 10nm resolution? (Phase III)

Zone plates with finer structure?

## Microscope at 34 ID

Even finer structures possible?

Problem aspect ratio: Etching techniques limited (> 1:20??) →Deposition techniques: sputtering

Way of fabrication: Circular lenses Linear lenses

Thin lens approximation valid?

## **Prospects for the future**

#### A microscope with 50nm resolution

- Depth of focus
- Field of View
- Intensity
- Contrast
- Mechanical Stability

## **Depth of focus**

Depth of Focus:

 $\sim \delta^2/\lambda$ 

DOF= $2(dr_n)^2/\lambda$ 

 $\rightarrow$  DOF 15-30µm for E = 6-12 keV

higher energies - contrast Only partially sharp imaging – O.K. for Tomography Smaller samples – handling, preparation Field of view of Detector

## **Field of view**

#### Field of view:

- 2048 Pixel with 12.5 nm size  $\rightarrow$  25  $\mu$ m field of view
- Larger Detector Array (e.g. 4k x 4k)
- Tomography with 360° rotation?
- Partial Tomography?

- More general (valid also for higher resolution): combine low and high resolution tomographic data high resolution around rotation center
- → fits also with DOF problem: only central part of sample can be imaged with high resolution

## Intensity

Intensity:

Why important?

Not only smaller but also thinner features (3D)  $\rightarrow$  contrast decrease $\rightarrow N_{Ph}^{2}$ 

with constant detector resolution and field of view:

Increase magn.  $\rightarrow$  illum. Detector Area<sup>2</sup> $\rightarrow$  N<sub>Ph</sub><sup>2</sup>

 $\delta \sim N_{Ph}^{4}$ 

## Microscope at 34 ID

Intensity: Use of Pink beam 50-100 higher Intensity:

-Condition for FZP (chromatic aberration):  $N < E/\Delta E = 10^2$  N: Number of Zones -but also 100 Zones minimum required for functioning of FZP

| For $dr_n = 40$ nm @ 91     | keV              |                                     |
|-----------------------------|------------------|-------------------------------------|
| $D = 4N dr_n$               | $D = 16 \ \mu m$ | D: diameter of Zone Plate           |
| And                         |                  | dr <sub>n</sub> : finest Zone width |
| $f = 4N (dr_n)^2 / \lambda$ | f = 5 mm         | N: number of Zones                  |
|                             |                  | f: focal length                     |

- In principle possible- Real quality and feasibility of such a FZP?

- Perhaps intermediate solution? N=500 $\rightarrow \Delta E/E=2*10^{-3}$ 

### Contrast

Contrast:

Number of counts necessary to visualize a given contrast:

Rose Signal-to-Noise Ratio  $\Delta SNR_{Rose} = C * Sqrt(N_q)$  C: contrast  $N_q$ : mean number of quanta/Pixel  $SNR_{Rose}$  should be at least 5 (for the table 10)

#### With actual equipment

| С  | N <sub>q</sub> | Field of View                              | t <sub>exp</sub> |
|----|----------------|--|------------------|
| 1% | 106            | 30x20 μm <sup>2</sup><br>(2400x1600 Pixel) | minutes          |

#### Space for improvements (pink beam, efficiencies,...)

Absorption/Transmission\* of different elements with d=50nm as a function of E



- Go to lower energies? May be, but

- Problem DOF

Considering "free standing sample" - Penetration depth
 → hard X-rays = phase contrast \*http://www-cxro.lbl.gov

#### **Phase Contrast - Refractive Index**

## For very thin samples ( $\Delta \Phi \ll 1$ ) and $\Delta \Phi = (2\pi/\lambda) \, \delta \, 1$ supposing pure phase contrast object $\Delta \Phi \sim C$



Sufficient contrast in interesting energy range

\*http://www-cxro.lbl.gov

## Microscope at 34 ID

Phase Contrast: Making use of Zernike Phase contrast method:

Basic Idea:

Transform a phase information into an amplitude information
Diffracted beam(1<sup>st</sup> order) phase shifted in respect to undiffracted beam(0<sup>th</sup> order).

-1<sup>st</sup> and 0<sup>th</sup> order spatially separated in back focal plane -Phase ring in back focal plane:

phase shift ( $\pi/2$  or  $3\pi/2$ ) and attenuation of 0<sup>th</sup> order -In image plane increased contrast can be achieved

## **Zernike Phase Contrast**



1<sup>st</sup> order 0<sup>th</sup> order Resulting

## **Zernike Phase Contrast**



## Summary

Development of sub-100 nm Microscope at 34 ID-C -KB-FZP optic -Zernike Phase Contrast

Current project for user-friendly scientific applications

In future 50nm imaging/tomography

- min.exposures/optimization

Ultimate resolution 10nm? Relation to other methods? Lensless imaging? let's listen to the next talk!

#### Acknowledgements

K. Peterson, P. Jemian (UNICAT)
U. Neuhäusler (U. Bielefeld), G. Schneider (BESSY)
S. Vogt, P. Illinski, C. Roehrig, E. Wrobel, B. Tieman (APS)
R. Divan, B. Lai (APS)
G. Falkenberg (Hasylab), I. Robinson (UIUC)

The UNICAT facility at the Advanced Photon Source (APS) is supported by the U.S. DOE under Award No. DEFG02-91ER45439, through the Frederick Seitz Materials Research Laboratory at the University of Illinois at Urbana-Champaign, the ORNL (U.S. DOE contract DE-AC05-00OR22725 with UT-Battelle LLC), the NIST (U.S. Department of Commerce) and UOP LLC. The APS is supported by the U.S. DOE, Basic Energy Sciences, Office of Science under contract No. W-31-109-ENG-38.