PHASE TRANSITIONS IN POLYMERIC AND MICELLAR SYSTEMS

1. INTRODUCTION

The most active Small-Angle Neutron Scattering (SANS) research areas are in "Polymers" and "Complex Fluids" which include micellar systems. Understanding the thermodynamics of mixing of polymers and/or complex fluids is important. This involves insight into the major phase transitions and the various structures that are formed.

Polymer mixtures comprise polymer solutions, polymer blends, copolymers as well as their mixtures. Micelles are formed when a surfactant, oil and water are mixed. Adding oil forms a microemulsion. These systems form structures in the nanometer length scale that can be quantified by the SANS technique.





Figure 1: Various polymeric and micellar systems. The small dots represent solvent molecules.

2. PHASE TRANSITION BASICS

Polymer mixtures phase separate either through cooling and are characterized by an Upper Critical Spinodal Temperature (UCST), or upon heating and are characterized by a Lower Critical Spinodal Temperature (LCST). Combinations of LCST and UCST are also possible.



PHASE SEPARATION ISLAND

Figure 2: Various possible phase diagrams.

Composition fluctuations are enhanced when phase transition lines are approached from the mixed (1-phase) region. This increases the SANS intensity thereby making this technique an appropriate probe for thermodynamic (phase transition) studies.

The "binodal" condition is the phase transition line between the homogeneously mixed 1phase region and the phase separated (i.e., demixed) 2-phase region in the temperaturecomposition phase diagram. The spinodal line is the deeper part of the demixed 2-phase region. The region between the binodal and the spinodal lines is the nucleation-andgrowth region while the region within the spinodal line is the spinodal region. Most water-soluble systems phase separate upon heating due to the softening of hydrogenbonding at elevated temperatures. Note that the binodal line is mapped out using light

scattering (cloud point measurements) whereas the spinodal line is measured using SANS.




Figure 3: Schematic representation of the phase diagram for an LCST polymer blend.

For block copolymers, the spinodal line is referred to as the Order-Disorder-Temperature (ODT). The ODT defines a region which comprises three main forms of microphase separation; these correspond to spherical, cylindrical and lamellar morphologies. The LCST line is referred to as the Lower Critical Ordering Temperature (LCOT) for copolymers. The copolymer composition is varied by increasing the length of one of blocks in the copolymer.




Figure 4: Representation of the various copolymer microphases for a Lower Critical Ordering Temperature (LCOT) system.

Surfactants are amphiphilic molecules which contain a hydrophilic part and a hydrophobic part. Surfactants can be ionic when the hydrophilic part of the molecule is charged. They can also be nonionic when the molecule is completely uncharged. When surfactants are mixed with water and are above a certain concentration and temperature, they form micelles in solution. The concentration threshold is called the Critical Micelle Concentration (CMC) and the temperature threshold is the Critical Micelle Temperature (CMT). One example of a nonionic surfactant is a Pluronic copolymer described next.

3. WHAT ARE PLURONICS?


Poly(ethylene oxide) (referred to as PEO) is the simplest water soluble polymer with a monomer structure –CH₂CH₂O–. Its neighbor in the homologous series poly(propylene oxide) (or PPO) does not dissolve in water at ambient temperature. Pluronics are triblock copolymers composed of one central PPO block connected to two PEO blocks. The PEO block dissolves well in aqueous media because it is mostly hydrophilic, while the PPO block does not dissolve because it is mostly hydrophobic at ambient temperature. This amphiphilic nature of pluronics molecules makes them form micelles at ambient temperature. At low temperature, this balance does not hold and both PEO and PPO blocks dissolve in water thereby remaining individual molecules. Pluronics are commercially available materials used in the cosmetics and pharmaceutical industries. The CMT and CMC vary depending on the block composition of the various pluronics which allows selection of the optimal surfactant for the desired application. The P85 pluronic considered here consists of 26 EO monomers in each of the outside blocks and

40 PO monomers in the middle block. P85 is referred to as $EO_{26}PO_{40}EO_{26}$. The molecular weight of P85 is around 4,600 g/mol.

4. SANS FROM P85 PLURONIC

The interest here is in micelle-formation as well as in the investigation of phase transitions that form various microstructures. A series of P85 samples of different weight fractions in D_2O (d-water) were prepared and measured by SANS.

In order to investigate micelle formation, SANS data were taken from 10 % by weight P85 in deuterated water (d-water) at various temperatures (from 10 °C to 60 °C). Deuterated water was used in order to enhance the neutron contrast. Micelles are well formed at 30 °C. When micelles are formed, the scattering is characterized by two features: a peak characteristic of inter-micelles interactions, and decay at high Q characterizing the tail of the single-particle form factor. In-between these two features, a second shoulder (around $Q = 0.15 \text{ Å}^{-1}$) can be observed. This shoulder is also characteristic of the single particle form factor (oscillation of the spherical Bessel function) and is more or less visible depending on the size polydispersity of the micelles.



10 % P85 in d-water

Figure 5: SANS data from the 10 % by weight P85 in D_2O at various temperatures. Data statistics are different for the two (low-Q and high-Q) instrument configurations. The overlap region for the two SANS configurations used is around 0.04 Å⁻¹.

In order to determine the Critical Micelle Temperature, the SANS intensity at a low Q value ($Q = 0.02 \text{ Å}^{-1}$ is chosen) is plotted with increasing temperature. The CMT is seen to be around 22 °C for the 10 % by weight P85/d-water sample. This characteristic temperature is observed as a jump in the low-Q intensity.





Figure 6: Sigmoid shape variation of the low-Q SANS intensity with temperature for the 10 % by weight P85/d-water sample.

SANS data are shown also for a 0.5 % by weight P85/d-water sample for a wider temperature range (10 °C to 100 °C). The low P85 weight fraction allows the observation of single-particle effect (form factor) without interference from the inter-particle interaction peak. The intensity is seen to increase with increasing temperature but the shape of the form factor changes from unimers (below the CMT) to spherical micelles, then to cylindrical micelles and finally to lamellar micelles. This is a rich phase behavior and structural transitions which are deduced based on standard clues in the SANS spectra. For example, scattering for unimers (low temperature) is characterized by a $1/Q^2$ variation in the high-Q (so-called Porod) region. Spherical and cylindrical micelles are characterized by a $1/Q^4$ variation in the Porod region (rod-like scattering). Lamellae are characterized by a $1/Q^2$ dependence in the low-Q region (scattering from a 2D object) and the formation of an inter-lamellar Bragg peak (referred to as the "long period" in crystalline polymers).



0.5 % P85 in d-water

Figure 7: SANS data for the 0.5 % P85 in D₂O for increasing temperature.

When the PPO central block of P85 becomes hydrophobic, spherical micelles form. Upon heating, hydrogen-bonding between the PEO outer blocks and D_2O starts softening yielding cylindrical micelles. Upon further heating, the PEO blocks become hydrophobic and lead to a phase transition producing the lamellar structure. In the lamellar phase, the triblock molecules combine to form multilayer structures rich in PEO or PPO; D_2O is mostly excluded and forms solvent-rich pockets.

Here also, a plot of the low-Q SANS intensity shows clearly the transition temperatures between the observed phases. Heating and cooling cycles show that these phase transitions are reversible.




Figure 8: Variation of the low-Q SANS intensity for increasing temperature for the 0.5 % by weight P85/d-water sample.

5. THE P85/D-WATER PHASE DIAGRAM

The P85/d-water phase diagram has been mapped out (Mortensen, 1996) using the SANS technique (among other techniques). The main phases (unimers, spherical micelles, cylindrical micelles and lamellar phase) can be observed at low P85 weight fraction. Other phases (ordered spherical micelles, hexagonal phase and a disordered phase) are not discussed here.




Figure 9: The P85/d-water phase diagram.

6. MICELLAR STRUCTURES MODELS

At low P85 weight fraction, unimers (dissolved copolymer chains), then spherical micelles, then cylindrical micelles, and finally lamellae form with increasing temperature. These structures are described here. If we assume that the micelles scattering length density is uniform, then the SANS cross section (scattering intensity) can be expressed as follows:

$$\frac{d\Sigma(Q)}{d\Omega} = \left(\frac{N}{V}\right) \Delta \rho^2 V_P^2 P(Q) S_I(Q).$$
(1)

Here N/V is the micelle number density, $\Delta \rho^2$ is the contrast factor, V_P is the micellar volume, P(Q) is the single micelle form factor and $S_I(Q)$ is the inter-micelles structure factor.

The form factor P(Q) for a Gaussian polymer coil of radius of gyration R_g is given by the so-called Debye function:

$$P(Q) = \frac{2}{Q^4 R_g^4} \left[exp(-Q^2 R_g^2) - 1 + Q^2 R_g^2 \right].$$
 (2)

For a sphere of radius R:

$$\mathbf{P}(\mathbf{Q}) = \left[\frac{3\mathbf{j}_1(\mathbf{Q}\mathbf{R})}{\mathbf{Q}\mathbf{R}}\right]^2 = \left[\frac{3}{\mathbf{Q}\mathbf{R}}\left(\frac{\sin(\mathbf{Q}\mathbf{R})}{(\mathbf{Q}\mathbf{R})^2} - \frac{\cos(\mathbf{Q}\mathbf{R})}{\mathbf{Q}\mathbf{R}}\right)\right]^2.$$
 (3)

The spherical Bessel function of first order $j_1(x)$ has been defined.

For a randomly oriented cylinder of length L and radius R:

$$P(Q) = \frac{1}{2} \int_{-1}^{1} d\mu \left[\frac{\sin(Q\mu L/2)}{Q\mu L/2} \right]^{2} \left[\frac{2J_{1}(Q\sqrt{1-\mu^{2}}R)}{Q\sqrt{1-\mu^{2}}R} \right]^{2}.$$
 (4)

A lamellar stack of length L and (extended) radius R can also be represented by the same form factor:

$$P(Q) = \frac{1}{2} \int_{-1}^{1} d\mu \left[\frac{\sin(Q\mu L/2)}{Q\mu L/2} \right]^{2} \left[\frac{2J_{1}(Q\sqrt{1-\mu^{2}}R)}{Q\sqrt{1-\mu^{2}}R} \right]^{2}.$$
 (5)

Here $J_1(x)$ is the cylindrical Bessel function of first order.

These models are used to perform nonlinear least squares fits to the SANS data.

7. SANS DATA ANALYSIS

Using the relevant IGOR SANS Data Analysis models, fits are performed to the 0.5 % $P85/D_2O$ SANS data at various temperatures. Since the P85 weight fraction was low, no inter-micelle interactions were considered; i.e., structure factor effects were neglected $(S_I(Q) = 1)$ for the unimers, spherical micelles and cylindrical micelles structures. An inter-lamellar structure factor is built into the lamellar structure model considered since lamellae are fundamentally different from other structures. The lamellar phase results from the breaking of hydrogen bonds making PEO blocks hydrophobic.

Since micelle formation is driven by the hydrophobic (PPO) part of the P85 triblock copolymer, core-shell micellar structures are considered in each case (spherical, cylindrical and lamellar). Models that include polydispersity and instrumental smearing effects were used.




Figure 10: Representation of core-shell sphere, cylinder and lamellar stack.

Using the smeared Debye function model, fits are performed to the 10 °C data. There are three varying parameters: an overall scale factor, a radius of gyration and a Q-independent (mostly incoherent scattering) background. The low-Q SANS upturn characteristic of some form of clustering has been excluded from the fit. The results are included here.

Fit to SmearedDebye,




Figure 11: SANS data for the 0.5 % P85/D₂O at 10 $^{\circ}$ C and fit to the smeared Debye model.

The smeared polydisperse core-shell sphere model is used to fit the 50 °C data representing spherical micelles. Fitting parameters include an overall scale factor, the sphere core radius and its polydispersity and the shell thickness, the scattering length densities for the core, the shell and the solvent regions and a constant background. The solvent scattering length density is kept constant equal to the value for D_2O . Note that since the PEO blocks are hydrated, the scale factor is not necessarily equal to the P85 volume fraction.

Fit to SmearedPolyCoreForm,

Data file: X0_5__p85_50C_txt

```
scale =
                     0.00883924 \pm 0.000177863
avg core rad (A) =
                     46.2927 \pm 0.178344
core polydisp (0,1) = 0.141087 \pm 0.00137233
shell thickness (A) = 23.7195 \pm 0.479928
SLD core (A-2) =
                     4.38452e-07 \pm 1.93293e-08
SLD shell (A-2) =
                     5.58624e-06 \pm 3.15427e-08
SLD solvent (A-2) = 6.3e-06 \pm 0
bkg(cm-1) =
                     0.0706587 \pm 0.000121782
                     502.904
chisq =
Npnts =
                     235
                                    Sqrt(X^2/N) = 1.46288
Fitted range = [0,234] = 0.003837 < Q < 0.4054
```




Figure 12: SANS data for the 0.5 % P85/D₂O at 50 $^{\circ}$ C and fit to the smeared polydisperse core-shell spherical micelles model.

The smeared polydisperse core-shell cylinder model is used to fit the 80 °C data. The fitting parameters include the scale factor, the mean cylinder core radius and its polydispersity, the cylinder core length, the radial shell thickness and the face shell thickness, the three scattering length densities (core, shell and solvent) and finally the constant background. The solvent scattering length density for D_2O is fixed during the fits. The numerical integration used to perform the orientational averaging slows down the fitting process. Typical fit results are included.

Fit to SmearedPolyCoShCylinder, 7

Data file: X0_5__p85_80C_txt

scale = 0.0126776 ± 0.627644 mean CORE radius (A) = 19.2283 ± 0.80576 radial polydispersity (sigma) = 0.356149 ± 0.0125716 CORE length (A) =5770.6 ± 604.38 radial shell thickness (A) = 28.9607 ± 0.790341 face shell thickness (A) = 1613.72 ± 316.611 SLD core $(A^{-2}) =$ $1.67113e-06 \pm 0.000114591$ SLD shell $(A^-2) =$ $3.74124e-06 \pm 6.33398e-05$ SLD solvent $(A^-2) =$ $6.3e-06 \pm 0$

incoh. bkg (cm⁻¹) = $0.0749588 \pm 0.000157727$ chisq = 1769.76Npnts = 235 Sqrt(X²/N) = 2.74425 Fitted range = [0,234] = 0.003837 < Q < 0.4054FitError = No Error FitQuitReason = No Error




Figure 13: SANS data for the 0.5 % P85/D₂O at 80 °C and fit to the smeared polydisperse core-shell cylindrical micelles model.


The lamellar structure is harder to model. A crude attempt is included here. Note that this model requires an inter-lamellar structure factor. For the sake of simplicity, details have been left out; some details can be found elsewhere (Hammouda, 2008).

This lamellar structure model assumes a stack of lamellae (think of this as a cylindrical stack of disks) with very large disk radius and alternating lamellae (PEO-PPO-PEO, then PEO-PPO-PEO, etc). PEO forms the core part and PEO forms the layer part of the lamellae. This model assumes a small number of lamellae per stack. This model is used to fit the $95 \,^{\circ}C$ data. Fit results are included.

Fit to SmearedStackedDiscs,

Data file: X0_5__p85_95C_txt

scale = 0.0036907 ± 0.363266 Disc Radius (A) = 9417.1 ± 0 Core Thickness (A) = 34.0089 ± 68.4543






Figure 14: SANS data for the 0.5 % P85/D₂O at 95 $^{\circ}$ C and fit to the smeared stacked lamellae model.

This lamellar stack model does not adequately fit the fully developed lamellar structure observed at $100 \,^{\circ}$ C.

REFERENCES

K. Mortensen, "Structural Studies of Aqueous Solutions of PEO-PPO-PEO Triblock Copolymers: Their Micellar Aggregates and Mesophases; A SANS Study" J. Phys.; Condensed Matter <u>8</u>, A13 (1996). B. Hammouda, "Probing Nanoscale Structures: the SANS Toolbox", (2008). The pdf file is available online at http://www.ncnr.nist.gov/staff/hammouda/the_SANS_toolbox.pdf. Check for instance chapters 35 and 44.