| Method 1622/1623 Spiking Suspension Enumeration Form: <i>Crypt</i> One) | osporidium Giardia (Circle |
|--|----------------------------|
| Internal spiking suspension enumeration ID number (note on Line 7 of bench sheet) | |
| Source of stock Cryptosporidium oocyst/ Giardia cyst suspension | |
| Stock Cryptosporidium oocyst/ Giardia cyst suspension lot number | |
| Date that stock Cryptosporidium oocyst/ Giardia cyst suspension was received | |
| Date of this spiking suspension enumeration | |
| Analyst performing enumeration | |
| | |
| Well number (or hemacytometer chamber number if not using well-slide enumeration) | Oocysts/cysts counted |
| 1 | |
| 2 | |
| 3 | |
| 4 | |
| 5 | |
| 6 | |
| 7 | |
| 8 | |
| 9 | |
| 10 | |
| | |
| Time enumeration completed (samples must be spiked within 24 hours of enumeration) | |
| | |
| Mean of the 10 counts | |
| Standard deviation of the 10 counts | |
| Relative standard deviation | |
| Volume used per well (in μL) | |
| | |
| | |

Laboratory ICR ID (if applicable):

Laboratory name:

| Comments: | | |
|-----------|--|--|
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

July 16, 1998, Revision - Draft