

Prosthetics Remove IFCAP Item Description From 2319

Release Notes

Patch RMPR*3*99

Version 3.0

November 2004

Department of Veterans Affairs VISTA Health System Design and Development (HSD&D)

Table of Contents

View Prosthetics Item Transactions (Read Only)	Patch RMPR*3*99 Overview NIF IFCAP Changes in the 2319 (Roll and Scroll)	
GUI 2319 Button on DOR Menu – Appliance Transactions (Tab 4)		
	• • • • • • • • • • • • • • • • • • • •	
GUI 2319 Button on DOR Menu - Home Oxygen (Tab 8)	GUI 2319 Button on DOR Menu - Home Oxygen (Tab 8)	
	endix A	
Appendix A	Glossary	

Remove IFCAP Item Description from 2319 Release Notes Patch RMPR*3*99 Overview

Introduction

Patch RMPR*3*99 (Remove IFCAP Item Description from 2319) was developed for Prosthetics to provide an interim solution for the continuation of the National Item File (NIF) data cleansing effort.

On the national level, the IFCAP Item file is pulled from each facility's database and merged into a VISN database. Then the VISN's NIF team reviews the Items and edits, deletes and adds descriptions (which are considered to be the "data cleansing" process).

If the IFCAP Item (at the national level) was a Prosthetic Item that was changed, added or deleted, then the Prosthetic Patient Treatment Record (2319) was historically being changed to match the data cleansing changes.

Patch RMPR*3*99 prevents the Prosthetic 2319 from historically being changed when data cleansing happens. This patch protects the **PROSTHETICS ITEM DESCRIPTIONS** in our Prosthetic 2319 from being overridden when the NIF make changes to the **IFCAP DESCRIPTIONS**.

Patch RMPR*3*99 changes

Patch RMPR*3*99 affects multiple **Prosthetics Menu** options. All three types of the Prosthetics 2319 have changes including:

- 1) Roll and Scroll
- 2) Read Only Roll and Scroll
- 3) GUI 2319 from the Delayed Order Report (DOR) application.

This patch will change the Prosthetics 2319 *Appliance Transactions* and *Home Oxygen Items*, switching the **IFCAP Item Description** to the **HCPCS Description** display.

Note to IRM

There is no new .EXE file to install. All the changes are M routines and a new cross reference in File 660. The **IFCAP Item Description** will remain in the database until phase 2 of the NIF implementation.

This is the menu option name that IRM references: **RMPR VIEW 2319 READ ONLY.**

Note: Please read these Release Notes for further information.

Patch RMPR*3*99 Overview, Continued

No HCPCS entered prior to 1997

Since HCPCS codes were not used prior to 1997, you will continue to see the **IFCAP Item Description** that was selected at that time. This is considered to be protected data and will not change.

Even though the **Appliance Transaction (Screen 4)** column heading displays **HCPCS** (in the Roll and Scroll), the information will continue to display the original **IFCAP Item Description** for those transactions entered prior to 1997.

Even though the GUI version of the **Appliance Transaction (Tab 4)** displays **Item** as the column header, it is actually the **IFCAP Item** for those transactions entered prior to 1997. The column heading was **not** changed to **HCPCS**.

NIF IFCAP Changes in the 2319 (Roll and Scroll)

2319 Roll and Scroll changes

Access the **Prosthetics Official's Menu** option, **Display/Print Patient 2319 (23)** to view the **Appliance Transactions (Screen 4)** or **Home Oxygen (Screen 8)**. Below are the sample screens for the Roll and Scroll 2319 that were affected by the NIF IFCAP changes.

Display/Print Patient 2319 Menu (23)

```
Select Prosthetic Official's Menu Option: DD <Enter> Display/Print
23
      Display/Print Patient 2319
ΤI
      Transaction Inquiry
ΙP
      Print All Prosthetic Items
ΒI
      Print Prosthetic Billings for MAS
IH
      Item History
SE
      Search for Recalled Item
SI
      Site Parameter Inquiry
VI
      Vendor Inquiry
Select Display/Print Option: 23 <Enter> Display/Print Patient 2319
```

Select the patient

```
Select PROSTHETIC PATIENT:
                              PROSPATIENT, ONE <Enter>
                                                        LOMA
LINDA VAMC 6-10-31
                    70%
000040707
            YES
                           SC VETERAN
                                           LL/SD/
MODULE III - GIM PCP:
PROSPROVIDER, ONE CX **ADVANCE DIRECTIVE - NO **
DEVICE: HOME// <Enter> VIRTUAL CONNECTION
                                   Right Margin: 80// <Enter>
Current Disability Codes are:
AO/DIS
               ALL OTHER S/C
                              S/C
```

Select the Appliance Transactions (Screen 4)

```
Select one of the following:
                  PATIENT DEMOGRAPHICS
                  CLINIC ENROLLMENTS/CORRESPONDENCE
         3
                  ENTITLEMENT INFORMATION
         4
                  APPLIANCE TRANSACTIONS
         5
                  AUTO ADAPTIVE INFORMATION
                  CRITICAL COMMENTS
                   ADD/EDIT DISABILITY CODE
                   HOME OXYGEN ITEMS
Enter 10-2319 screen to VIEW (1-8),'^' to EXIT, or 'return' to
continue : 4 <Enter> APPLIANCE TRANSACTIONS
                   SSN: 000-04-0707 DOB: JUN 10,1931 CLAIM# SS
HHUU,UXKHUS G
```

NIF IFCAP Changes in the 2319 (Roll and Scroll), Continued

HCPCS Column

Notice the **HCPCS** column heading below is highlighted. This is the column of data that has been changed with Patch RMPR*3*99 to the **HCPCS Description** from the **IFCAP Item Description** (which was the **Item** column heading).

New HCPCS Column Title

Date Qty HCPCS Type Vendor Sta Serial Delivery Date Tot Cost 1. 05/29/02 1 AUTO BLOOD I VENDOR, ONE 605 05/29/02 42.01
2. 07/26/00 1 VASCULAR G 1 VENDOR, ONE 605 05/29/02 42.01
End of Appliance/Repair records for this veteran!
+=Turned-In *=Historical Data I=Initial X=Repair S=Spare R=Replacement
Enter 1-2 to show full entry, '^' to exit or `return` to
continue. 1 <enter></enter>
PROSPATIENT, ONE SSN: 000-00-0001 LOMA LINDA VAMC DOB: 99-99 -1931
APPLIANCE/REPAIR LINE ITEM DETAIL <4-1>
TYPE OF FORM: STOCK ISSUE INITIATOR: PROSPROVIDER, THREE
DATE: MAY 29, 2002
DELIVER TO:
TYPE TRANS: INITIAL ISSUE QTY: 1 SOURCE: COMMERCIAL
VENDOR: PROSVENDOR, ONE
VENDOR PHONE: 800-555-5555
28690 N. Any DRIVE
LAKE FORREST, ILLINOIS 60045
DELIVERY DATE: MAY 29, 2002
TOTAL COST: \$42.01 OBL:
REMARKS: [If you type a specific note here, it appears on
Screen 4 under the transaction with notessee below.]
DISABILITY SERVED: SC/OP
ITEM DESCRIPTION: BLOOD PRESSURE UNIT DIGITAL
APPLIANCE: BLOOD PRESSURE UNIT DIGITAL
PSAS HCPCS: A4670 AUTO BLOOD PRESSURE MONITOR
ICD-9 Code: 401.9 HYPERTENSION NOS
CPT MODIFIER:
DESCRIPTION: AUTO BLOOD PRESSURE MONITOR
EXTENDED DESCRIPTION:
AUTOMATIC BLOOD PRESSURE MONITOR
Enter RETURN to continue or '^' to exit:

The Item Description is new!!

Remarks

If you enter a specific note in the **Remarks** field, it will be easier to locate it to identify it for future Item needs. The note displays on Screen 4 below the transaction.

Item Description / Appliance prompt

The **IFCAP Item Description** field display above is new. This is the field you should review to determine which **Prosthetic Item** was provided to the patient. You should no longer review the **Appliance** field. When the **Item Description** and the **Appliance** field are identical (as shown above), then the "new note" that is described on the next page will not appear. See next page.

NIF IFCAP Changes in the 2319 (Roll and Scroll), Continued

Item **Description NEW Note**

With Patch RMPR*3*99, whenever an **IFCAP Item Description** has been changed, edited or deleted, you will see a new note under this field as follows:

*** See Above for Original Item Description ***

Due to the data cleansing, the **Appliance** field (which was the original Item selected) has a new description and the original description is now the **Item Description** field.

After Patch Install - Item **Description** with NOTE

TYPE OF FORM: STOCK ISSUE INITIATOR: PROSPROVIDER, ONE

DATE: MAY 29, 2002

DELIVER TO:

TYPE TRANS: INITIAL ISSUE QTY: 1 SOURCE: COMMERCIAL

VENDOR: PROSVENDOR, ONE VENDOR PHONE: 800-555-5555

28690 N. Any DRIVE

LAKE FORREST, ILLINOIS 60045

DELIVERY DATE: MAY 29, 2002

TOTAL COST: \$42.01 OBL:

REMARKS:

DISABILITY SERVED: SC/OP

ITEM DESCRIPTION: BLOOD PRESSURE UNIT DIGITAL *** See Above For Original Item Description ***

APPLIANCE: BLOOD PRESSURE CUFF DIGITAL

PSAS HCPCS: A4670 AUTO BLOOD PRESSURE MONITOR

ICD-9 Code: 401.9 HYPERTENSION NOS

CPT MODIFIER:

DESCRIPTION: AUTO BLOOD PRESSURE MONITOR

EXTENDED DESCRIPTION:

AUTOMATIC BLOOD PRESSURE MONITOR

Enter RETURN to continue or '^' to exit:

To interpret the example above, the original **Item Description** was "Blood Pressure Unit Digital." Due to data cleansing, the description was changed to "Blood Pressure Cuff Digital."

The Appliance Field was changed to the new description: "Blood Pressure Cuff Digital" and the original or saved description appears in the new field: **Item** Description.



5

View Prosthetics Item Transactions (Read Only)

"Read Only" patch change

The change with Patch RMPR*3*99 that was previously explained also affects the Roll and Scroll version of Prosthetics and provides the following new (Read only) menu option:

View Prosthetics Item Transactions

Below is a sample screen.

View 2319 (Read Only)

```
Select OPTION NAME: VIEW PROSTHETICS ITEM TRANSACT <Enter> RMPR VIEW 2319
READ ONLY
View Prosthetics Item Transactions
SITE: SUPPORT ISC// <Enter>
                                                        578A5
Select PROSTHETIC PATIENT: PROSPATIENT, ONE <Enter>
                                                       SUPPORT ISC
                                                                       4-5-61
000111111
          SC VETERAN
  YES
PATIENT IS DECEASED. DATE OF DEATH WAS OCT 21,1993@08:00
Would you Like to continue Processing this Patient? NO// YES <Enter>
DEVICE: HOME// <Enter> TELNET Right Margin: 80// <Enter>
     Select one of the following:
                   ITEM TRANSACTIONS
                   HOME OXYGEN ITEMS
Enter DATA screen to VIEW (Item Transactions or Home Oxygen),'^' to EXIT, or
'return' to continue :
```


GUI 2319 Button on DOR Menu – Appliance Transactions (Tab 4)

Delayed Order Reports (DOR)

The **GUI 2319 Appliance Transactions (Tab 4)** window from the Delayed Order Report (DOR) below shows the heading: **Item** column; however this information is now the **HCPCS Description** with Patch RMPR*3*99.

The column heading **Item** will not be changed to **HCPCS**.

Appliance Transactions (Tab 4)



Next Step

Click the **View Detail** button at the bottom of this window to display the **Appliance Transaction Detail** window and proceed to the next page.


GUI 2319 Button on DOR Menu – Appliance Transactions (Tab 4), Continued

Appliance Field

The **Appliance** field listed in the **Transaction Detail** window below is protected from changes to IFCAP (from the NIF group's data cleansing).

The example below, "Wheelchair 9000XT 18-inch" in the **Appliance** field will remain the same with Patch RMPR*3*99.

Note: The **Appliance** field below will always display the originally selected **IFCAP Item Description**. The new "**Item Description**: field in Roll and Scroll and the new note "*See Above*…" is not necessary in the GUI 2319. This is because the GUI 2319 **IFCAP Item Description** will always be displayed and protected from changes.



Close button


To close the **Transaction Detail** window and return to the **View 2319** window, click the **Close** button.

GUI 2319 Button on DOR Menu - Home Oxygen (Tab 8)

Home Oxygen (Tab 8)

The 2319 in the Delayed Order Report (DOR) shown below of the **Home Oxygen** (Tab 8) window shows the **Item** column which is the **HCPCS Description** with Patch RMPR*3*99.

Home Oxygen window



View Detail button


Click the **View Detail** button at the bottom of this window to display the **Home Oxygen Detail** and proceed to the next page.

GUI 2319 Button on DOR Menu - Home Oxygen (Tab 8),

Continued

Home Oxygen (Tab 8) Transaction Detail Notice the **Appliance** field on the **Transaction Detail** window from the **Home Oxygen** (**Tab 8**) remains the same, even when changes are made to the NIF IFCAP Item Description.

Note: The **Appliance** field below will always display the originally selected **IFCAP Item Description**. The new "**Item Description**: field in Roll and Scroll and the new note "*See Above*…" is not necessary in the GUI 2319. This is because the GUI 2319 **IFCAP Item Description** will always be displayed and protected from changes.



Appendix A

Glossary

DOR	Delayed Order Report – Prosthetics application.
GUI	Graphical User Interface. This is a windows-based application with buttons and tabs instead of prompts in a Roll and Scroll environment for an application. Microsoft Windows is a GUI environment.
HCPCS	Healthcare Common Procedure Coding System – formerly Health Care Financing Administration Common Procedures Coding System - is distributed by Centers for Medicare & Medicaid Services (CMS) and the Health and Human Services (HHS) for all other substances, equipment, supplies, or other items used in Health Care Services. These items include, but are not limited to: medical supplies; orthotic and prosthetic devices; and durable medical equipment.
NIF	The National Item File (NIF) database is used to standardize naming practices so that facilities will be able to identify, track and evaluate purchases.
	The NIF is also a group that works on this database file.
VistA	Veterans Health Information Systems and Technology Architecture, formerly known as Decentralized Hospital Computer Program, encompasses the complete information environment at VA medical facilities. It consists of hardware, software packages, and comprehensive support for system-wide and station specific, clinical and administrative automation needs.